## Trypsine-EDTA

**REF**      **CEZTDA00-0U**  
**CEZTDA01-0U**

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7, avenue de Scandinavie  
91940 Les Ulis  
FRANCE

## 15. Bibliographie

- Balk, S.D., Lestourgeon, D. and Mitchell, R.S. 5-methyltetrahydrofolic acid, 5-formyltetrahydrofolic acid (folinic acid), and folate acid requirements of normal and Rous sarcoma virus-infected chicken fibroblasts. *Cancer Res.*, 1978, 38, 3966-3968.
- Carrel, A. and Burrows, M.T. Cultivation of adult tissues and organs outside of the body. *J. Am. Med. Ass.*, 1910, 55, 1379-1381.
- Eagle, H. The specific amino acid requirements of a human carcinoma cell (strain HeLa) in tissue culture. *J. Exp. Med.*, 1955, 102, 37-48.
- Eagle, H. The minimum vitamins requirements of the L and HeLa cells in tissue culture, the production of specific vitamin deficiencies and their cure. *J. Exp. Med.*, 1955, 102, 595-600.
- Eagle, H. Nutrition needs of mammalian cells in tissue culture. *Science*, 1955, 122, 501-504.
- Eagle, W.R. Production of malignancy in vitro-IV the mouse fibroblast cultures and changes seen in the living cells. *J. Nat. Cancer. Inst.*, 1943, 4, 165-212.
- Fischer, A., Astrup, P., Ehrens-Vård, G. and Ohlenschlager, V. Growth of animal tissue cells in artificial media. *Proc. Soc. Exp. Biol. Med.*, 1948, 67, 40-46.
- Gey, G.O. An improved technique for massive tissue culture. *Am. J. Cancer*, 1933, 17, 752-756.
- Graham, A.F. and Siminovitch, L. Prolifération of monkey kidney cells in rotating cultures. *Proc. Soc. Exp. Biol. Med.*, 1955, 89, 326-327.
- Griffith, J.B. and Pirt, S.J. The uptake of amino acids by mouse cells (strain LS) during growth in batch culture and chemostat culture : the influence of cell growth rate. *Proc. Roy. Soc. Serie B*, 1967, 168, 421-438.
- Ham, R.G. Clonal growth of mammalian cells in chemically defined medium. *Proc Natl. Acad. Sci.*, 1965, 53, 288-293.

# Fiche Technique Trypsine-EDTA

---

- Me Keehan, W.L., Mc Keehan, K.A., Hanunon, S.L. and Ham, R.G. Improved media for clonal growth of human diploid fibroblasts at low concentration of serum protein. *In vitro*, 1977, 13, 399-416.
- Monard; D., Solomon, F., Rentsch, M. and Gysin, R. Glia-induced morphological differentiation in neuroblastoma cells. *Proc. Nat. Acad. Sci.*, 1973, 70, 1894-1897.
- Morgan, J.F., Morton, H.J. and Parker, R.C. Nutrition of animal cells in tissue culture. I. Initial studies on a synthetic medium. *Proc. Soc. Exp. Biol. Med.*, 1950, 73, 1-8.
- Moscona, A. Cell suspensions from organ rudiments of chick embryos. *Exp. Cell Res.*, 1952, 3, 535-539.
- Reitzer, J.L., Wice, B.M. and Kennell, D. Evidence that glutamine, not sugar, is a major energy source of cultured HeLa cells. *J. Biol. Chem.*, 1979, 254, 2669-2676.
- Rous, P. and Jones, F.A. A method of obtaining suspensions of living cells from the fixed tissues and from the plating out of individual cells. *J. Exp. Med.* 1916, 23, 549-555.
- Waymouth, C. Osmolality of mammalian blood and media for culture of mammalian cells. *In vitro*, 1970, 6, 109-127.
- F. Devreker, K. Hardy, M. Van den Bergh, A.S. Vannin, S. Emiliani and Y. Englert Amino acids promote human blastocyst development in vitro .*Human Reproduction* Vol.16, No.4 pp. 749-756, 2006.
- Mitalipov SM1, White KL, Farrar VR, Morrey J, Reed WA. Development of nuclear transfer and parthenogenic rabbit embryos activated with inositol 1,4,5-Trisphosphate. *Biol Reprod.* 1999 Apr;60(4):821-7.

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---

## 1. General Information

This is the most widely used protease in cell culture. Isolated from pig and bovine pancreas, it catalyzes the cleavage reaction at the lysine or arginine level. Its optimal pH of action is 8.0; it does not require the presence of specific ions or cofactors. It is insensitive to the presence of EDTA which is used in combination for cell detachment. The action of trypsin is inhibited by serum and divalent ion.

## 2. Intended use

Trypsin is used for the preparation of mammal cells and tissues and is in combination with EDTA for the detachment of cells.

The product is intended for in-vitro use, not for human therapy or veterinary applications.

Trypsin-EDTA should be used by qualified medical laboratory personnel.

# Instructions for use Trypsin EDTA

---

## 3. Symbols

|                                                                                     |                                        |
|-------------------------------------------------------------------------------------|----------------------------------------|
|    | Reference                              |
|    | Batch number                           |
|    | Temperature limit                      |
|    | Expiration date                        |
|    | Manufacturer                           |
|  | CE marked product                      |
|  | In vitro diagnostic medical device     |
|  | Instructions for use                   |
|  | Warning, read the instructions for use |
|  | Date of manufacture                    |
|  | Do not use if the packaging is damaged |

## 4. Packaging

The packaging is 100ml for all references. Different compositions and concentrations are available.

# Instructions for use Trypsin EDTA

---

| Description                         | Packaging | Ref.        |
|-------------------------------------|-----------|-------------|
| Trypsin-EDTA (Versene) 1X, Sterile  | 100 mL    | CEZTDA00-0U |
| Trypsin-EDTA (Versene) 10X, Sterile | 100 mL    | CEZTDA01-0U |

## 5. Characteristics

| Description                         | Composition                                        | Ref.        |
|-------------------------------------|----------------------------------------------------|-------------|
| Trypsin-EDTA (Versene) 1X, Sterile  | 0,5g/l trypsin and 0,2 g/l EDTA without phenol red | CEZTDA00-0U |
| Trypsin-EDTA (Versene) 10X, Sterile | 5g/l trypsin and 2g/l EDTA without phenol red      | CEZTDA01-0U |

## 6. Conservation and Storage

Store the Trypsin-EDTA solution at -15°C/-22°C until the expiration date indicated on the label.

After opening the bottle, store the digestive solution at -15°C/-22°C for 1 month.

## 7. Shipping

The delivery is made in dry ice.

## 8. Required material non provided

Depending on the application, not supplied materials may be required (pipets, flasks, micropipettes...).

## 9. Warnings and precautions

The solution must be taken up in balanced saline solution without Ca or Mg for 10X concentrations.

The user should check in advance whether the use with another reagent affects or not the performance of the device.

The device should be handled with care by the user in order to maintain the verified microbiological status of the product.

# Instructions for use Trypsin EDTA

---

Do not use if the packaging is damaged.

## 10. Protocole

|   |                                                                                   |
|---|-----------------------------------------------------------------------------------|
| 1 | 10X dilution of Trypsin                                                           |
| 2 | Rinse the cells with a little diluted Trypsin                                     |
| 3 | Cover the cells completely with Trypsin previously heated to 37°C in a water bath |
| 4 | Incubate the mixture at 37°C until the cell mat is detached.                      |

## 11. Quality control

The Liquid saline solutions undergo quality control tests and microbiological controls. Trypsin solutions are controlled for mycoplasma.

### **Physical and chemical controls:**

pH and osmolarity as well as density are measured by standard equipment and standards. And according to standardized procedures.

### **Microbiological controls:**

Bacterial and fungal sterility controls are performed according to the protocol. The samples to be tested are previously incubated at two temperatures (20°C and 35°C) for 14 days. At the end of the incubation period, a direct sub-culture is performed.

The culture media used are:

- for aerobic germs BTCS
- for anaerobic germs Thioglycolate

In parallel, controls are sown to check if the media are able to allow the growth of a small number of organisms and that there is no inhibiting effect of the medium. Signs of bacterial growth are checked at regular intervals.

### **Biological controls:**

The absence of cytotoxicity of the saline solutions is verified on MRC5 cells grown in monolayer and known to be sensitive to cytotoxic elements. After several washes, the cells maintained with the saline solution to be tested are observed for 48 hours to ensure the absence of signs of cytotoxicity.

## 12. Waste disposal

Eliminate all waste according to local legislation.

---

# Instructions for use Trypsin EDTA

---

## 13. Incident report

Any serious incident occurring in connection with the device shall be notified to EUROBIO SCIENTIFIC and to the competent authority of the Member State in which the user and/or patient is established.

## 14. Technical assistance

For assistance with our products, please contact our technical support.

Eurobio Scientific customer service can be reached by e-mail at [adv@eurobio-scientific.com](mailto:adv@eurobio-scientific.com) or by phone at +33 (0)1.69.79.64.80

7, avenue de Scandinavie  
91940 Les Ulis  
FRANCE

## 15. Bibliography

- Balk, S.D., Lestourgeon, D. and Mitchell, R.S. 5-methyltetrahydrofolic acid, 5-formyltetrahydrofolic acid (folinic acid), and folate acid requirements of normal and Rous sarcoma virus-infected chicken fibroblasts. *Cancer Res.*, 1978, 38, 3966-3968.
- Carrel, A. and Burrows, M.T. Cultivation of adult tissues and organs outside of the body. *J. Am. Med. Ass.*, 1910, 55, 1379-1381.
- Eagle, H. The specific amino acid requirements of a human carcinoma cell (strain HeLa) in tissue culture. *J. Exp. Med.*, 1955, 102, 37-48.
- Eagle, H. The minimum vitamins requirements of the L and HeLa cells in tissue culture, the production of specific vitamin deficiencies and their cure. *J. Exp. Med.*, 1955, 102, 595-600.
- Eagle, H. Nutrition needs of mammalian cells in tissue culture. *Science*, 1955, 122, 501-504.
- Eagle, W.R. Production of malignancy in vitro-IV the mouse fibroblast cultures and changes seen in the living cells. *J. Nat. Cancer. Inst.*, 1943, 4, 165-212.
- Fischer, A., Astrup, P., Ehrens-Vård, G. and Ohlenschlager, V. Growth of animal tissue cells in artificial media. *Proc. Soc. Exp. Biol. Med.*, 1948, 67, 40-46.
- Gey, G.O. An improved technique for massive tissue culture. *Am. J. Cancer*, 1933, 17, 752-756.
- Graham, A.F. and Siminovitch, L. Prolifération of monkey kidney cells in rotating cultures. *Proc. Soc. Exp. Biol. Med.*, 1955, 89, 326-327.

# Instructions for use Trypsin EDTA

---

- Griffith, J.B. and Pirt, S.J. The uptake of amino acids by mouse cells (strain LS) during growth in batch culture and chemostat culture : the influence of cell growth rate. Proc. Roy. Soc. Serie B, 1967, 168, 421-438.
- Ham, R.G. Clonal growth of mammalian cells in chemically defined medium. Proc Natl. Acad. Sci., 1965, 53, 288-293.
- Me Keehan, W.L., Mc Keehan, K.A., Hanunon, S.L. and Ham, R.G. Improved media for clonal growth of human diploid fibroblasts at low concentration of serum protein. In vitro, 1977, 13, 399-416.
- Monard; D., Solomon, F., Rentsch, M. and Gysin, R. Glia-induced morphological differentiation in neuroblastoma cells. Proc. Nat. Acad. Sci., 1973, 70, 1894-1897.
- Morgan, J.F., Morton, H.J. and Parker , R.C. Nutrition of animal cells in tissue culture. I. Initial studies on a synthetic medium. Proc. Soc. Exp. Biol. Med., 1950, 73, 1-8.
- Moscona, A. Cell suspensions from organ rudiments of chick embryos. Exp. Cell Res., 1952, 3, 535-539.
- Reitzer, J.L., Wice, B.M. and Kennell, D. Evidence that glutamine, not sugar, is a major energy source of cultured HeLa cells.. J. Biol. Chem., 1979, 254, 2669-2676.
- Rous, P. and Jones, F.A. A method of obtaining suspensions of living cells from the fixed tissues and from the plating out of individual cells. J. Exp. Med. 1916, 23, 549-555.
- Waymouth, C. Osmolality of mammalian blood and media for culture of mammalian cells. In vitro, 1970, 6, 109-127.
- F. Devreker, K. Hardy, M. Van den Bergh, A.S. Vannin, S. Emiliani and Y. Englert Amino acids promote human blastocyst development in vitro .Human Reproduction Vol.16, No.4 pp. 749-756, 2006.
- Mitalipov SM1, White KL, Farrar VR, Morrey J, Reed WA. Development of nuclear transfer and parthenogenic rabbit embryos activated with inositol 1,4,5-Trisphosphate. Biol Reprod. 1999 Apr;60(4):821-7.