Fiche technique MILIEU DE SEPARATION DES LYMPHOCYTES $d = 1.078$

MILIEU DE SEPARATION DES LYMPHOCYTES $d = 1.078$

REF

CMSMSL01-0U
CMSMSL01-01
CMSMSL01-1A

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7, avenue de Scandinavie
91940 Les Ulis
FRANCE

15. Bibliographie

- ❑ Balk, S.D., Lestourgeon, D. and Mitchell, R.S. 5-methyltetrahydrofolic acid, 5-formyltetrahydrofolic acid (folinic acid), and folie acid requirements of normal and Rous sarcoma virus-infected chicken fibroblasts. *Cancer Res.*, 1978, 38, 3966-3968.
- ❑ Carrel, A. and Burrows, M.T. Cultivation of adult tissues and organs outside of the body. *J. Am. Med. Ass.*, 1910, 55, 1379-1381.
- ❑ Eagle, H. The specific amino acid requirements of a human carcinoma cell (strain HeLa) in tissue culture. *J. Exp. Med.*, 1955, 102, 37-48.
- ❑ Eagle, H. The minimum vitamins requirements of the L and HeLa cells in tissue culture, the production of specific vitamin deficiencies and their cure. *J. Exp. Med.*, 1955, 102, 595-600.
- ❑ Eagle, H. Nutrition needs of mammalian cells in tissue culture. *Science*, 1955, 122, 501-504.
- ❑ Eagle, W.R. Production of malignancy in vitro-IV the mouse fibroblast cultures and changes seen in the living cells. *J. Nat. Cancer. Inst.*, 1943, 4, 165-212.
- ❑ Fischer, A., Astrup, P., Ehrens-Vård, G. and Ohlenschlager, V. Growth of animal tissue cells in artificial media. *Proc. Soc. Exp. Biol. Med.*, 1948, 67, 40-46.
- ❑ Gey, G.O. An improved technique for massive tissue culture. *Ai'n. J. Cancer*, 1933, 17, 752-756.
- ❑ Graham, A.F. and Siminovitch, L. Prolifération of monkey kidney cells in rotating cultures. *Proc. Soc. Exp. Biol. Med.*, 1955, 89, 326-327.
- ❑ Griffith, J.B. and Pirt, S.J. The uptake of amino acids by mouse cells (strain LS) during growth in batch culture and chemostat culture : the influence of cell growth rate. *Proc. Roy. Soc. Serie B*, 1967, 168, 421-438.
- ❑ Ham, R.G. Clonal growth of mammalian cells in chemically defined medium. *Proc Natl. Acad. Sci.*, 1965, 53, 288-293.
- ❑ Me Keehan, W.L., Mc Keehan, K.A., Hanunon, S.L. and Ham, R.G. Improved media for clonal growth of human diploid fibroblasts at low concentration of serum protein. *In vitro*, 1977, 13, 399-416.
- ❑ Monard; D., Solomon, F., Rentsch, M. and Gysin, R. Glia-induced morphological différentiation in neuroblastoma cells. *Proc. Nat. Acad. Sci.*, 1973, 70, 1894-1897.

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Instructions for use Lymphocytes separation medium

d = 1.078

1. General information

The Lymphocyte Separation Medium (LSM), a gradient combining Ficoll® and Diatrizoate (d = 1.078+/- 0.002), allows the recovery of lymphocytes contained in blood by simple centrifugation. After preparation of the Diatrizoate buffer and neutralization by caustic soda, Ficoll® is added, the medium is then sterilized by filtration.

2. Intended use

The Lymphocyte Separation Medium (LSM) is a product intended for general laboratory use. This medium commercialised by Eurobio Scientific is particularly used before in vitro diagnostic examination procedures on human samples.

The Lymphocyte Separation Medium (LSM) is therefore suitable for use In Vitro diagnostic procedures.

Lymphocyte Separation Medium (LSM) is a medical device for in vitro diagnostic procedures and cannot be recycled.

The Lymphocyte Separation Medium (LSM) must be used by qualified medical laboratory personnel.

Instructions for use Lymphocytes separation medium

d = 1.078

3. Symbols

Reference

Batch number

Temperature limit

Expiration date

Manufacturer

Date of manufacture

CE marked product

In vitro diagnostic medical device

Instructions for use

Warning, read the instructions for use

Do not use if the packaging is damaged

Keep out of direct sunlight

Instructions for use Lymphocytes separation medium d = 1.078

4. Packaging

Description	Packaging	Référence
Lymphocytes separation medium d = 1,078	100 ml	CMSMSL01-0U
Lymphocytes separation medium d = 1,078	500 ml	CMSMSL01-01
Lymphocytes separation medium d = 1,078 in pierceable flask	150 ml	CMSMSL01-1A

5. Characteristics

Lymphocytes separation medium d = 1,078+/- 0,002 is a medium whose sterile state is claimed.

Formulation :

Components g/l	CMSMSL01
EDTA Na2	0.0530
DIATRIZOIC ACID	94.100
NaOH	5.8892
FICOLL 400®	56.000

6. Conservation and storage

Store the Separation lymphocytes medium at +2°C/+30°C until the expiration date indicated on the label.

After opening the bottle, store the medium at à +2°C/+30°C for 1 month.

7. Shipping

The delivery can be made at room temperature. Indeed, a temporary transit at room temperature does not alter product characteristics.

8. Required material non provided

Depending on the application, materials not supplied may be required (pipets, flasks, micropipettes...).

9. Warnings and precautions

The device should be handled with care by the user in order to maintain the verified microbiological status of the product.

Instructions for use Lymphocytes separation medium

$d = 1.078$


Do not use if the packaging is damaged.

10. Protocole

The whole blood must first be defibrinated by shaking it for a few minutes (5 to 6 beads for 5 to 8 ml of whole blood) or collected in the presence of heparin.

- Dilute the blood sample by half with saline.
- In a 10 ml tube, sequentially add 3 ml of MSL then 4 ml of blood, without mixing.
- Centrifuge at room temperature for 30 to 40 minutes at 400g.
- Using a Pasteur pipette, collect the dense opalescent layer as shown in the diagram (erythrocytes and other leukocytes sediment at the bottom of the tube).
- Replace the isolated lymphocyte suspension in an equivalent volume of Hanks' saline.
- Centrifuge at room temperature for 10 minutes at 160g.
- Resuspend the pellet in an equivalent volume of Hanks' saline.
- Centrifuge at room temperature for 10 minutes at 160g.
- Retain the lymphocyte pellet.

Note: Some immature granulocytes may accompany the lymphocytes in uremic patients or those on immunosuppressive therapy.

Instructions for use Lymphocytes separation medium

d = 1.078

11. Quality control

Physical and chemical controls:

pH and osmolarity as well as density are measured by standard equipment and calibrators.

Microbiological controls:

Bacterial and fungal sterility controls are performed according to the protocol. The samples to be tested are previously incubated at two temperatures (20°C and 35°C) for 14 days. At the end of the incubation period, a direct sub-culture is performed.

The culture media used are:

- for aerobic germs BTCS
- for anaerobic germs Thioglycolate

In parallel, controls are sown to check that the media are able to allow the growth of a small number of organisms and that there is no inhibiting effect from the medium. Signs of bacterial growth are checked at regular intervals.

The endotoxin level is measured according to the European Pharmacopoeia.

The quality of the separation is controlled using horse blood.

12. Waste disposal

Eliminate all waste according to local legislation related to waste from healthcare activities with risks of infection.

13. Incident report

Any serious incident occurring in connection with the device shall be notified to EUROBIO SCIENTIFIC and to the competent authority of the Member State in which the user and/or patient is established.

14. Technical assistance

For product assistance, please contact our technical support.

Eurobio Scientific customer service can be reached via e-mail at adv@eurobio-scientific.com or via phone at +33 (0)1.69.79.64.80

Instructions for use Lymphocytes separation medium

d = 1.078

7, avenue de Scandinavie
91940 Les Ulis
FRANCE

15. Bibliography

- ❑ Balk, S.D., Lestourgeon, D. and Mitchell, R.S. 5-methyltetrahydrofolic acid, 5-formyltetrahydrofolic acid (folinic acid), and folate acid requirements of normal and Rous sarcoma virus-infected chicken fibroblasts. *Cancer Res.*, 1978, 38, 3966-3968.
- ❑ Carrel, A. and Burrows, M.T. Cultivation of adult tissues and organs outside of the body. *J. Am. Med. Ass.*, 1910, 55, 1379-1381.
- ❑ Eagle, H. The specific amino acid requirements of a human carcinoma cell (strain HeLa) in tissue culture. *J. Exp. Med.*, 1955, 102, 37-48.
- ❑ Eagle, H. The minimum vitamins requirements of the L and HeLa cells in tissue culture, the production of specific vitamin deficiencies and their cure. *J. Exp. Med.*, 1955, 102, 595-600.
- ❑ Eagle, H. Nutrition needs of mammalian cells in tissue culture. *Science*, 1955, 122, 501-504.
- ❑ Eagle, W.R. Production of malignancy in vitro-IV the mouse fibroblast cultures and changes seen in the living cells. *J. Nat. Cancer. Inst.*, 1943, 4, 165-212.
- ❑ Fischer, A., Astrup, P., Ehrens-Vård, G. and Ohlenschlager, V. Growth of animal tissue cells in artificial media. *Proc. Soc. Exp. Biol. Med.*, 1948, 67, 40-46.
- ❑ Gey, G.O. An improved technique for massive tissue culture. *Am. J. Cancer*, 1933, 17, 752-756.
- ❑ Graham, A.F. and Siminovitch, L. Prolifération of monkey kidney cells in rotating cultures. *Proc. Soc. Exp. Biol. Med.*, 1955, 89, 326-327.
- ❑ Griffith, J.B. and Pirt, S.J. The uptake of amino acids by mouse cells (strain LS) during growth in batch culture and chemostat culture : the influence of cell growth rate. *Proc. Roy. Soc. Serie B*, 1967, 168, 421-438.
- ❑ Ham, R.G. Clonal growth of mammalian cells in chemically defined medium. *Proc Natl. Acad. Sci.*, 1965, 53, 288-293.
- ❑ McKeehan, W.L., McKeehan, K.A., Hanunon, S.L. and Ham, R.G. Improved media for clonal growth of human diploid fibroblasts at low concentration of serum protein. *In vitro*, 1977, 13, 399-416.
- ❑ Monard, D., Solomon, F., Rentsch, M. and Gysin, R. Glia-induced morphological differentiation in neuroblastoma cells. *Proc. Nat. Acad. Sci.*, 1973, 70, 1894-1897.

Instructions for use Lymphocytes separation medium

d = 1.078

- ❑ Morgan, J.F., Morton, H.J. and Parker, R.C. Nutrition of animal cells in tissue culture. I. Initial studies on a synthetic medium. Proc. Soc. Exp. Biol. Med., 1950, 73, 1-8.
- ❑ Moscona, A. Cell suspensions from organ rudiments of chick embryos. Exp. Cell Res., 1952, 3, 535-539.
- ❑ Reitzer, J.L., Wice, B.M. and Kennell, D. Evidence that glutamine, not sugar, is a major energy source of cultured HeLa cells. J. Biol. Chem., 1979, 254, 2669-2676.
- ❑ Rous, P. and Jones, F.A. A method of obtaining suspensions of living cells from the fixed tissues and from the plating out of individual cells. J. Exp. Med. 1916, 23, 549-555.
- ❑ Waymouth, C. Osmolality of mammalian blood and media for culture of mammalian cells. In vitro, 1970, 6, 109-127.
- ❑ F. Devreker, K. Hardy, M. Van den Bergh, A.S. Vannin, S. Emiliani and Y. Englert Amino acids promote human blastocyst development in vitro. Human Reproduction Vol.16, No.4 pp. 749-756, 2006.
- ❑ Mitalipov SM1, White KL, Farrar VR, Morrey J, Reed WA. Development of nuclear transfer and parthenogenic rabbit embryos activated with inositol 1,4,5-Trisphosphate. Biol Reprod. 1999 Apr;60(4):821-7.