
EurobioPlex DermaXtract

For lysis of nail samples

Reference EBE-077 Version 1.02 – 10/16/2024

Instructions for use

Available on www.eurobio-scientific.com

Instructions for Use EurobioPlex DermaXtract

ENGLISH.....	1
---------------------	----------

FRANÇAIS.....	22
----------------------	-----------

Table of content

1. General information.....	3
2. Intended use.....	3
3. Symbols.....	4
4. Packaging.....	5
5. Principle of the test	5
6. Composition	5
7. Storage/stability.....	5
8. Materials required but not provided	5
9. Cautions and recommendation for users.....	6
10. Sample collection	6
11. Procedure	6
12. Waste Disposal.....	7
13. Incident report.....	7
14. Technical assistance.....	8

Instructions for Use EurobioPlex DermaXtract

1. General information

The buffers in this kit allow lysis of nail samples and the stabilization of the DNA obtained after lysis. The DNA obtained with this kit can be tested with the EurobioPlex DermaScreen (EBX-073) kit.

It is the responsibility of the end user to qualify these products for their specific application. These products are not qualified or intended for any animal or human therapeutic use.

2. Intended use

The DermaXtract kit is a two component reagent kit for extraction of PCR-ready genomic DNA from nail samples.

Both buffers are ready to use solutions. The Extraction buffer allows the lysis of the nails and the Stabilisation buffer guarantees long-term stability of the genomic DNA.

This reagent has been designed for use in combination with the EurobioPlex DermaScreen kit (EBX-073). Its use for other purposes must be validated by the user.

The DermaXtract is a medical device for in vitro diagnosis use.

The DermaXtract must be used by qualified medical laboratory personnel.

The DermaXtract is suitable for use in *In Vitro* diagnostic procedures.

Instructions for Use EurobioPlex DermaXtract

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Date of manufacture

CE marked product

In vitro diagnostic medical device

Instructions for use

Do not use if packaging is damaged

Caution

Instructions for Use EurobioPlex DermaXtract

4. Packaging

The EurobioPlex DermaXtract kit (EBE-077) contains two ready to use buffers (see table below).

Cap colour	Content of the kit	50 reactions	100 reactions	200 reactions	Reconstitution
Red	Extraction Buffer	4x 1500 µl	8x 1500 µl	16 x 1500 µL	Ready to use
Transparent	Stabilisation Buffer	4x 1500 µl	8x 1500 µl	16 x 1500 µL	Ready to use

5. Principle of the test

The Extraction buffer results in the lysis of nail samples and enables the extraction of the DNA from the nails. The stabilisation buffer is used to stabilize the DNA for long-term storage at -20°C. Both buffers are ready to use.

6. Composition

The kit does not contain hazardous components.

7. Storage/stability

All reagents from EBE-077 should be stored at room temperature (+15/+30°C) before opening and is stable until the date printed on the tube labels.

8. Materials required but not provided

The following materials and equipment are required but not provided

- Biological cabinet/UV cabinet
- Centrifuge for microtubes
- Heating block
- Vortex
- Adjustable micropipettes (2-20 µl, 20-200 µl, 200-1000 µl)

Instructions for Use EurobioPlex DermaXtract

- DNase-free and RNase-free filter tips for micropipettes
- Sterile microtubes
- Gloves (talc-free)
- Sterile DNase free 1.5 ml tubes with safety lock

9. Cautions and recommendation for users

Read these instructions carefully before starting the procedure.

- Do not use the product if the individual packaging is damaged.
- Upon receipt of the device, it is necessary to ensure the integrity of the packaging and containers that guarantee the quality of the reagents. Preliminary examination of the reagents must not reveal any precipitate, particles or leaks.
- This extraction should only be performed by qualified laboratory employees.
- The reagents are exclusively for in vitro diagnostic use.
- This experiment should be performed according to good laboratory practice.
- Use a labcoat and disposable gloves.
- Do not use this kit after the expiration date.
- Use DNase-free filter tips for micropipettes

10. Sample collection

- Nail samples can be stored at room temperature.
- Samples should be considered as potentially infectious material and are preferably prepared in a laminar flow cabinet or biosafety cabinet.
- Transfer nail part(s) to a 1.5 mL sterile tube with safety lock before starting with the DNA-extraction.

11. Procedure

DNA-extraction with the DermaXtract kit

- Add 50 µL* of Extraction buffer to the sterile 1.5 mL tube with safety lock and make sure the nail part is completely covered in the Extraction buffer.
- Add 5 µL of the Inhibition Control from the EBX-073 kit to the sample.
- Incubate 10 minutes at 100°C (do not exceed these 10 minutes!).
- Briefly spin down the sample(s).
- Add 50 µL of Stabilisation buffer and vortex briefly.

Instructions for Use EurobioPlex DermaXtract

- Incubate 1 min at room temperature.
- Centrifugate 1 min at 10000 rpm.
- Transfer the supernatant to new (labeled) tubes with screwcaps and avoid transfer of nail parts**.
- Store samples at -20°C when not immediately required for PCR testing.

**To cover larger nail samples, also 100 µL Extraction buffer can be used. This requires 10 µL Inhibition control and 100 µL stabilization buffer.*

***The transfer to a new tube is not essential as long as no nail parts are transferred to the PCR plate.*

Inhibition control

The Armored DNA Inhibition Control is part of the EurobioPlex DermaScreen PCR kit (EBX-073) and is a DNA sequence in a protective protein coat. It is used to differentiate true negative samples from negative samples as a result of a non-efficient DNA-extraction or PCR failure.

12.Waste Disposal

Dispose of all waste in accordance with DASRI legislation or to the local applicable regulations.

13.Incident report

Any serious incident involving the reagents is reported to Eurobio Scientific and to the competent authority of the Member state in which the user and/or patient is established.

14. Technical assistance

For assistance with our products, please contact our technical support.

Eurobio Scientific's customer service can be contacted by e-mail at adv@eurobio.scientific.com or by phone at +33 (0)1.69.79.64.80

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