EurobioPlex

Bordetella Pertussis and Parapertussis REAL-TIME PCR

For **qualitative** real-time PCR

EBX-001-UN

REF **EBX-001-48**
EBX-001-24

24/48/96 reactions

Version 6.03 of 13/11/2024

Validated on:

- CFX96™ Real Time PCR detection system (Biorad) with analysis on CFX Manager v 3.1 (Biorad)
- T-COR 8®-IVD (Tetcore Inc.) with T-COR 8 SmartCT™ (auto v1) software

Reagents compatible with the following real-time PCR systems:

Applied Biosystems Prism® 7500; qTOWER 2.2; Cepheid Smart Cyclor II®;
Qiagen Rotor-Gene 3000/6000; Stratagene MX3005P*; Biorad Chromo4/CFX 96; Roche LightCycler 480...

* ongoing

Storage conditions: keep all reagents between -15°C and -22°C until use
and after first use

Instructions for use

Available on request from info@eurobio-scientific.com

TABLE OF CONTENTS

ENGLISH.....	1
FRANÇAIS.....	23
DEUTSCH.....	46

INTENDED USE.....	3
INTRODUCTION	3
PRINCIPLE OF DETECTION	4
DESCRIPTION AND CONTENT OF THE KIT.....	4
STORAGE.....	6
CAUTIONS AND NOTES.....	6
COLLECTION OF SAMPLES, TRANSPORT AND STORAGE	6
PROCEDURE.....	7
VALIDATION OF THE EXPERIMENT**	10
DATA ANALYSIS AND INTERPRETATION	10
PERFORMANCES ANALYSIS.....	11
SPECIFICITIES OF THE REAL-TIME PCR T-COR 8®-IVD INSTRUMENT	14
BIBLIOGRAPHY.....	21
WASTE ELIMINATION	21
SYMBOLS	22

INTENDED USE

The *Bordetella pertussis/parapertussis* test uses real-time polymerase chain reaction (PCR) amplification and is designed for the qualitative detection of *Bordetella pertussis* and *Bordetella parapertussis* from a nucleic acid extract. This test is indicated to diagnose the occurrence of an infection by *Bordetella pertussis* or by *Bordetella parapertussis* in humans. Extracted DNA is the starting material for the Eurobioplex *Bordetella* kit.

The Eurobioplex *Bordetella pertussis/parapertussis* has been validated on the following type of samples:

- Nasopharyngeal aspirations

This amplification system has been validated on 107 nasopharyngeal aspiration samples and 6 quality control samples from a French National Reference Center.

INTRODUCTION

Whooping cough is a very contagious disease of the respiratory system caused by Gram negative bacteria, *Bordetella pertussis* and *Bordetella parapertussis*. The clinical symptoms are a chronic cough, which can be very strong and can be followed by whistling breezing and throwing up. For the most severe cases, most often in newborns, this can lead to hypoxia, irreversible lesions and death. In older children and adults, particularly the ones that have been vaccinated or contracted the disease, it is not as severe and is characterized by a chronic cough.

On the first part of the XXth century, the prevalence was very high with 1% of the population suffering each year, and in most cases affecting young children. In the 1940s, vaccination considerably decreased prevalence, but nevertheless the disease remained. Since the end of the 1990s, there was a dramatic increase in the number of cases in developed countries performing vaccination. In 2004, 25 827 cases were reported in the USA, the highest number since 1959. Similar observations were made in Canada and Europe. The World Health organization estimates 50 million cases each year in the world, with 350 000 deaths, with a prevalence in all group ages. In most cases, *Bordetella pertussis* is the causing agent. 3% to 35% are caused by *Bordetella parapertussis*, which is responsible of a mild form of the disease similar to whooping cough. Although the causes of the reappearance of *Bordetella pertussis* infections is not clear, experts agree overall that the immunity conferred by the vaccine declines after 7 to 10 years, explaining disease during teenage years. These individuals could transmit the disease to newborn within the family.

Different techniques, including cultures, serology, and PCR are available for diagnosis. Culture is the technical reference with a specificity of 100%, but with a lower sensitivity. Culture can have 56% sensitivity at the beginning of the disease but decreases with time and is lower for patients that had an antimicrobial treatment or prior vaccination. *B. pertussis* culture requires a specific medium and the growth takes 7 to 14 days, much too long for acute cases. Serology uses a comparison between serum of the patient in the acute phase and recovery and requires at least 4 weeks interval, preventing initial diagnosis.

Specific serology for IgG antitoxins have been developed but samples must be taken 2 weeks after start of the symptoms. Other serology tests exist on the market but have not been validated and cannot discriminate between a recent infection, old infection or vaccination.

PRINCIPLE OF DETECTION

The *Bordetella pertussis/parapertussis* Eurobioplex test uses real time polymerase chain reaction (PCR) amplification, that simultaneously in one reaction, detects *Bordetella pertussis* (target region IS481) and *Bordetella parapertussis* (target region IS1001) from a nucleic acid extract as well as a human cellular control (b actin). The specific probe for *B. pertussis* is labelled with the FAM fluorophore. The specific probe for *B. parapertussis* is labeled with the HEX fluorophore. The specific probe for the cellular control is labeled with the Cy5 fluorophore. All emit fluorescence when hydrolyzed during the elongation of the amplification product of the internal control. Fluorescence measurements over time is correlated to accumulation of amplified products.

DESCRIPTION AND CONTENT OF THE KIT

The Real time PCR Eurobioplex *Bordetella* kit is ready for use for the specific detection of *B. pertussis* and *B. parapertussis* and human cellular control. Fluorescence is emitted and individually recorded through optical measurements during the PCR. The detection of the amplified fragment is performed by a fluorimeter using the channels shown in Table 1 below.

The kit contains reagents and enzyme for the amplification of *B. pertussis* and *B. parapertussis* DNA and human cellular control (Table 1).

Table 1:

Target	Fluorophore	Excitation	Emission
<i>B. pertussis</i>	FAM	495 nm	515 nm
<i>B. parapertussis</i>	HEX	535 nm	555 nm
Human Cellular Control	Cy5	647 nm	667 nm

Equivalent channels on different PCR instruments:

- Channel **FAM** (ABI Systems, SmartCycler II, Chromo 4/CFX96, Systems Mx, T-COR8®-IVD), Channel 530 (LC 480), Channel Green (RotorGene)
- Channel **HEX** (Chromo 4/CFX96, Mx Systems, T-COR8®-IVD), Channel VIC (ABI Systems), Channel Alexa532 (SmartCycler II), Channel 580 (LC 480), Channel Yellow (RotorGene),
- Channel **Cy5** (ABI Systems, Chromo 4/CFX96, Mx Systems, T-COR8®-IVD), Channel 660 (LC 480), Channel Alexa647 (SmartCyclerII), Channel Red (RotorGene)

Note: On LC480 instrument II, apply color compensation for the following channels: FAM-HEX/VIC- (465-510, 533-580).

Table 2:

Cap color	Components	24 reactions	48 reactions	96 reactions
-----------	------------	--------------	--------------	--------------

Red	ENZYME	4 x 100µL	2 x 375µL	4 x 375µL
Transparent	OLIGOMIX	4 x 60µL	2 x 210µL	4 x 210µL
Yellow	Positive control (<i>Bordetella</i> +human) (CP)	4 x 40µL	2 x 40µL	4 x 40µL
Blue	Water=negative control(CN-H ₂ O)	1 x 1000 µl	1 x 1000 µl	1 x 1000 µl

Required material not provided:

- ◇ Biological Hood
- ◇ Real time PCR instrument
- ◇ Micro centrifuge
- ◇ Vortex
- ◇ Plates / tubes for Real time PCR
- ◇ Micropipettes
- ◇ DNase-free RNase free filter tips for micropipettes
- ◇ Sterile microtubes
- ◇ Powder-free gloves

STORAGE

All reagents must be stored between -15°C and -20°C.

All reagents can be used until the expiration date indicated on the kit.

Many freezing/thawing cycles (> 3x) must be avoided, and could lead to decrease in sensitivity.

CAUTIONS AND NOTES

Read carefully instructions before starting.

- ◇ The experiment must be performed by competent staff.
- ◇ Instruments must have been properly installed, calibrated and maintained according to the manufacturer's recommendations.
- ◇ Clinical samples are potentially infectious and must be processed under a laminar flow hood.
- ◇ The experiment must be performed according to good laboratory practices.
- ◇ Do not use this kit after expiration date.
- ◇ The kit is shipped with dry ice, and the components of the kits must arrive frozen. If one or more components are defrosted, or if the tubes have been damaged, contact Eurobio-Scientific.
- ◇ Once defrosted, spin down briefly the tubes before use.
- ◇ It is recommended to define three distinct working areas: 1) Isolation of DNA, 2) Preparation of the reaction mix and 3) Amplification / Detection of amplified products.
- ◇ Pipettes, reagents and other materials must not cross each area.
- ◇ Specific caution is required to preserve the purity of the reagents and reaction mixtures. The use of powder-free gloves is mandatory.
- ◇ Appropriate methods of preparation/extraction of DNA to produce high quality DNA and to be followed by a real-time PCR application should be used, particularly avoiding all sources of DNase contamination.
- ◇ Always use filtered tips for micropipettes.
- ◇ Use specific lab coat and powder-free gloves in each working area.
- ◇ Do not pipette with mouth and do not eat, drink or smoke in the area.
- ◇ Avoid sprays.

COLLECTION OF SAMPLES, TRANSPORT AND STORAGE

- ◇ Collect samples in sterile tubes.
- ◇ Do not use calcium alginate swabs, as they may contain substances that inhibit PCR testing.
- ◇ It is the responsibility of the user to master its own conditions of collection, transport and storage of samples, and extraction of DNA by suitable systems to produce DNA of good quality.
- ◇ The samples should be extracted immediately or stored following the recommendations in the table below (Table 3).

Table 3 :

Recommendations of maximum storage of samples before extraction	
Room temperature	2 h
4°C	16h
-20°C	Long term storage

- ◇ The user can refer to the World Health Organization or High Health Authority for storing samples.
- ◇ Transport of clinical samples must follow local regulations for this type of infectious agents

PROCEDURE

I- DNA extraction

It is the user's responsibility that the extraction system used be compatible with downstream real time PCR technology. For this kit, we recommend to use extraction methods of bacterial DNA from nasopharyngeal aspirations, and refer to manufacturer's instructions of the extraction kit.

Note: It is not necessary to add the human cellular control at the extraction step. The human cellular control on the Cy5 channel ensures that a negative result is not due to a bad DNA extraction resulting from no extraction or the presence of PCR inhibitors and to check the quality of samples collection.

Note: Acceptable specimen types include nasopharyngeal swabs (NPS) in liquid media (e.g. UTM®, eSwabs® ...).

II- REAL-TIME PCR

General comment:

The *Bordetella* positive control contains a very high concentration of DNA. Manipulations must be performed with caution to avoid contamination.

II-1/ CONTROLS OF THE REAL-TIME PCR

To control the steps of extraction and real-time PCR amplification, it is necessary to test a positive control (CP) and a negative control (PCR water: CN-H₂O). See II-2/ Real Time PCR protocol for details.

II-2/ REAL-TIME PCR PROTOCOL

Diagram of the procedure

1 - PREPARATION OF THE MASTERMIX

Number of reaction(s)	N+3*
Enzyme	(N+3)x12.5 µl
Oligomix	(N+3)x7.5 µl
Total Volume Mastermix	(N+3)x20 µl

* N+1 on T-COR 8®-IVD

2 - PREPARATION OF THE CONTROLS AND REACTIONS

Sample

20 µl Mastermix
+
5µl DNA sample

Positive Control

20 µl Mastermix
+
5µl CP

Negative Control

20 µl Mastermix
+
5µl water supplied (CN-H2O)

3 – REAL TIME PCR INSTRUMENT (except T-COR 8®-IVD - automatic program using Barcodes- see page 37-43)

Program	Temperature	Duration	Cycle(s)	
Denaturation	95°C	3 min	1	-
Amplification	95°C	10 sec	45	-
	60°C	30 sec		Acquisition de fluorescence

Detailed Procedure

- 1) Homogenize the enzyme tube and vortex Oligomix and CP and centrifuge briefly.
- 2) Prepare Mastermix as below, N= number of reaction(s). Prepare Mastermix amount for a minimum of N+3 reactions (see part 1) of scheme above).

Number of reaction(s)	N+3*
Enzyme	(N+3) x 12.5 µl
Oligomix	(N+3) x 7.5 µl
Total Volume Mastermix	(N+3) x 20 µl

* For use on T-COR 8®-IVD, we recommend a Master Mix for N + 1 reactions for one-shot use of the kit

- 3) Mix mastermix prepared in 2) and centrifuge briefly
- 4) Distribute 20 µL Master mix using a micropipette and filtered sterile tips in each tube/well of microplate for real time PCR.
- 5) Add 5 µL of extracted DNA.
- 6) In parallel, test the following controls:
 - Positive control:
 - 20µL of mastermix + 5µL of CP
 - Negative control:
 - 20 µL mastermix + 5 µL water supplied (CN-H2O)
- 7) Close immediately with an adhesive film or transparent caps to avoid all contamination.
- 8) Centrifuge briefly to collect all the PCR reaction mix at the bottom of the tubes or plate.
- 9) Program the Real time PCR instrument as follows (except T-COR 8®-IVD - automatic program using Barcodes- see page 37-43)
Approximate duration is 1h:

Program	Temperature	Duration	Cycle(s)	
Denaturation	95°C	3 min	1	-
Amplification	95°C	10 sec	45	-
	60°C	30 sec		Acquisition of fluorescence

Note 1: On Light Cycler ® 480 systems (Roche), two optical systems are available: only "System II" is compliant with use of the kit. Apply color compensation for the FAM and HEX/VIC wavelengths (465-510, 533-580).

Note 2: For the Applied Biosystems systems, select "ROX" in "PASSIVE REFERENCE".

Note 3: On Rotorgene™, please calibrate the signal by clicking on "GAIN optimization".

VALIDATION OF THE EXPERIMENT**

If a signal appears in the negative control at the end of the cycle, the threshold must be positioned just above the negative signal's maximum signal. For the assay to be valid, the Ct values for the controls must be the following (Table 4). Outside of these values, the experiment is not valid. On T-COR 8®-IVD, we recommend that the positive and negative controls be tested when a new kit is being opened.

Table 4:

Controls	NC	CP
FAM	Ct >38 or not determined	Ct ≤ 30
HEX	Ct >38 or not determined	Ct ≤ 30
Cy5	Ct not determined or determined*	Ct <35

* Note: Ct could be determined in Cy5 channel in the negative control and it does not invalidate the run. The presence of human cellular gene in negative control could be operator dependent.

** On T-COR 8®-IVD, interpretation is automatically generated with Barcodes.

DATA ANALYSIS AND INTERPRETATION

1/ A signal (< 38 Ct) detected on channel FAM -> **The result is positive: the sample contains DNA from *Bordetella pertussis*.** For competition reasons, the signal of the internal human control on channel Cy5 can be strongly decreased (even undetectable) in presence of large amounts of *Bordetella pertussis* DNA.

2/ A signal (< 38 Ct) detected on channel HEX -> **The result is positive: the sample contains DNA from *Bordetella parapertussis*.** For competition reasons, the signal of the internal human control on channel Cy5 can be strongly decreased (even undetectable) in presence of large amounts of *Bordetella parapertussis* DNA.

3/ No signal on channels FAM or HEX (or signal > 38 Ct) and, at the same time, a positive signal on channel Cy5 -> **The sample does not contain DNA from *Bordetella pertussis* nor DNA from *Bordetella parapertussis*. The sample can be considered negative.**

4/ No signal on channels FAM or HEX or Cy5 (or signal > 38 Ct) -> **No diagnosis can be given.** DNA extraction was inhibited. It is recommended to dilute the DNA or redo a DNA extraction of the sample.

PERFORMANCES ANALYSIS

➤ Detection threshold:

Bordetella Pertussis: 100 copies/μL on positive control

Bordetella Parapertussis: 100 copies/μL on positive control

➤ Reproducibility (coefficient of variation):

Reproducibility of Ct within experiment:

B.Pertussis: 1.65%; ***B.Parapertussis***:1.61%

Reproducibility of Ct between batch:

B.Pertussis: 5.586%; ***B.Parapertussis***: 5.396%

Example:

The table below indicates mean Ct values for the IS481 plasmid:

dilutions	Ct parapertussis	SD	Coefficient of variation (%)
dilution 1	25,50	0,59	2,31
dilution 2	25,76	0,60	2,32
dilution 3	26,46	0,28	1,06

The table below indicates mean Ct values for the IS1001 plasmid:

dilutions	Ct parapertussis	SD	Coefficient of variation (%)
dilution 1	25,65	0,55	2,15
dilution 2	25,88	0,26	1,00
dilution 3	28,80	0,18	0,63

➤ **Sensitivity and specificity for *Bordetella Pertussis* and *Bordetella Parapertussis*:**

Study 1: In total 87 pretested samples were used to determine the performance of the kit including: 60 pretested positive for *B.Pertussis*, 2 pretested positive for *Bordetella Holmesii*, 26 pretested positive for *B.Parapertussis*. Samples pretested positive for *B.Pertussis* were negative for *B.Parapertussis* and conversely, except one sample co infected.

<i>Bordetella Pertussis</i>		EBX-001	
		POS	NEG
Pre-tested	POS	58	2
	NEG	0	22
<i>Bordetella Parapertussis</i>		EBX-001	
		POS	NEG
Pre-tested	POS	22	1
	NEG	0	59

Specificity: 100%.

No cross reactivity was observed between *B. Pertussis* and *B.Parapertussis*.

This test does not distinguish between infection by Bordetella pertussis and Bordetella holmesii. Both germs contain the IS481 gene. The incidence of Bordetella holmesii in clinical specimens has been reported to be very low (<0.5%), and it is still unclear whether Bordetella holmesii can cause respiratory disease in otherwise in good health.

Thereby the two positive sample for B. Holmesii were positive on the Bordetella Pertussis target.

Sensitivity:

96.7% samples pretested for *Bordetella Pertussis* were detected, 91.7% with a Ct less than 38.

95.6% samples pretested for *Bordetella Parapertussis* were detected.

Study 2: Comparison between Focus *Bordetella* Simplexa and Eurobioplex PCR:

Tests were performed on 20 clinical samples from respiratory swabs previously tested by the laboratory routine test (Focus *Bordetella* Simplexa), with the EurobioPlex *Bordetella pertussis/Bordetella parapertussis* on Smartcycler II ® (Cepheid).

Reference PCR	Eurobioplex PCR			
		Positive	Negative	Total
	Positive	11	0	11
	Negative	0	9	9
	Total	11	9	20

Sensitivity and Specificity between the 2 techniques are 100%.

Study 3: Quality control report from Reference Center Pasteur Institute CQBORD-2012:

The quality controls from Reference Center are composed of 6 tubes. Tests are performed on LC480® (Roche).

CQ-2012	Ct IS481	Ct IS1001
CQ1	24,1	-
CQ2	-	33,7
CQ3	27,1	-
CQ4	33	-
CQ5	-	-
CQ6	22,7	-

The results are consistent with expected results from Reference Center.

SPECIFICITIES OF THE REAL-TIME PCR T-COR 8®-IVD INSTRUMENT

T-COR 8®-IVD is a real-time PCR instrument with 8 independent wells, which can be programmed independently, patient by patient, regarding thermic protocol, assays, and time of start of the run.

The 8 wells can also be used simultaneously with the same assay.

Note: Mix the sample by pipetting up and down, and avoiding to make bubbles

Controls

On T-COR 8®-IVD, we recommend that the positive and negative controls be tested when a new kit is being opened.

After this initial control, the cellular control/Human DNA allows to check that the PCR is working properly, without PCR inhibition and that sample extraction work well.

Patients' Tests

<i>With validation of positive and negative controls once, during first use of the kit.</i>	24 reactions
Number of possible patients tests, patient by patient, with a maximum of 12 freezing/thawing cycles (3 times per tube/4 tubes per reagent)	12 patients

Fluorophores

Four combinations of fluorophores are available on T-COR 8®-IVD.

For all EBX, such as EBX-001, FAM/HEX/Texas Red/Cy5 combination is used. If some channels are not used, do not consider this channel for results analysis.

Use of Barcodes (available on page 41 to 43)

- 1- Select Menu > New Run
- 2- Select Barcode
- 3- Scan the appropriate Barcode on the right side of the instrument:
 - For the positive control (Barcode EBX-001 Pos Ctrl),
 - For the negative control (Barcode EBX-001 Neg Ctrl),
 - For a patient sample (Barcode Assay: EBX-001)

The instrument randomly selects one of the 8 available wells. Follow the instructions

on the instrument.

4- Lift up the lid of the selected well (well x), check that the blue LED light is on, and select "Yes"

5- Place the tube in the corresponding well and select "Next".

6- If you wish to name the sample (optional), select "sample x", name it and select "Accept".

7- Select "Next"

8- To add/test another control or sample, select « Add well » and go back to point 3-.

9- Select "Start run".

Note: The name of a sample or of a run cannot be modified or added after the run is finished. It must be done before or during the run.

Presentation and Interpretation of results

Ct values and amplification curves are available in the SmartCT™ Table (Ct values on each channel), and "Ct versus PCR cycles" graphs, both being available in real-time while the run is ongoing.

An automatic analysis for positive and negative controls as well as patients' samples is available in "Interpretations", at the end of the run, in the "View" window.

For positive and negative controls, the following results are possible:

« Neg Ctrl Fail »: Not valid

« Neg Ctrl Valid »: Valid

« Pos Ctrl Fail »: Non valid

« Pos Ctrl Valid »: Valid

For patients status:

"Detected": Positive for at least one target * → green box

"Not detected": Negative → red box

* The results for the specific bacteria are specified.

CAUTION !

When the box is green, it is important to read the status for each target as some targets may be negative.

« Invalid »: Invalid result -> retest → orange box

Example of automatic interpretation of results on T-COR 8®-IVD instrument screen and associated report:

- For valid positive and negative controls, « Detected » for its corresponding targets, for an invalid positive control to retest, for four negative samples for the two targets mentioned as “Not Detected” (samples 1 to 4) and for one positive samples for only one target indicated “Detected B.pertussis and B.parapertussis: not detected” (samples 871):


❖ On the instrument screen:

		Interpretations	
1	CP	EBX-001 Pos Ctrl	
		Invalid - Pos Ctrl Fail	
2	NTC	EBX-001 Neg Ctrl	
		Detected: Neg Ctrl Valid	
3	871	EBX-001	
		Detected: B. pertussis - B. parapertussis not detected	
4	Ech 1	EBX-001	
		Not Detected - B. parapertussis, B. pertussis	
5	ECH 2	EBX-001	
		Not Detected - B. parapertussis, B. pertussis	
6	ECH 3	EBX-001	
		Not Detected - B. parapertussis, B. pertussis	
7	ECH 4	EBX-001	
		Not Detected - B. parapertussis, B. pertussis	
8	CP	EBX-001 Pos Ctrl	
		Detected: Pos Ctrl Valid	

❖ On the report:

Résumé							
	Puits	Échantillon	Dosage	FAM	HEX	TxR	Cy5
1	CP	EBX-001 Pos Ctrl					
2	NTC	EBX-001 Neg Ctrl					
3	871	EBX-001	23.8			29.1	
4	Ech 1	EBX-001				34.2	
5	ECH 2	EBX-001				28.5	
6	ECH 3	EBX-001				26.2	
7	ECH 4	EBX-001				26.6	
8	CP	EBX-001 Pos Ctrl	23.5	24.0		23.1	


Example of amplification curves:

Barcodes for EBX-001 for use on T-COR8®-IVD


EBX-001 Pos Ctrl

Positive Control CP


EBX-001 Neg Ctrl

Water = Negative Control (CN-H2O)

Assay: EBX-001

BIBLIOGRAPHY

Fry N.K., Duncan J., Wagner K., Tzivra O., Doshi N., Litt D.J., Crowcroft N., Miller E., George R.C., Harrison T.G. Role of PCR in the diagnosis of pertussis infection in infants: 5 years' experience of provision of a same-day real-time PCR service in England and Wales from 2002 to 2007. *J. of Medic. Microbiol.* 2009, 58, 1023-1029.

Walsh P., Overmeyer C., Kimmel L., Feola M., Pusavat J., Nguyen T.A., Kuan S., Emery K., Rosengreen M., Mordechai E., Adelson M.E. Prevalence of *Bordetella pertussis* and *Bordetella parapertussis* in Samples Submitted for RSV Screening. 2008, Vol IX, No 3, 135-140.

Watanabe M., Connelly B., Weiss A.A. Characterization of Serological Responses to Pertussis. *Clin. And Vacc. Imm.* 2006, Vol 13, No 3, 341-348.

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Roorda L., Buitenwerf J., Osseward J.M., Van der Zee A. A real-time PCR assay with improved specificity for detection and discrimination of all clinically relevant *Bordetella* species by the presence and distribution of three Insertion Sequence elements. *BMC Research Notes* 2011, 4:11.

Bidet P., Liguori S., De Lauzanne A., Caro V., Lorrot M., Carol A., Faye A., Guiso N., Bingen E., Bonacorsi S. Real-Time PCR Measurement of Persistence of *Bordetella pertussis* DNA in Nasopharyngeal Secretions during Antibiotic Treatment of Young Children with Pertussis. *J. of Clin. Microbiol.* 2008, Vol 46, No 11, 3636-3638.

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WASTE ELIMINATION

Be in accordance with the law on the elimination of waste of clinical infectious material.

SYMBOLS

	Reference
	Batch number
	Highest storage temperature
	Expiration Date
	Content sufficient for « N » reactions
	Keep protected from light
	Manufacturer
	CE labeled product
	In Vitro Diagnostic
	Instructions for use

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