EurobioPlex RNA QC

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BIBLIOGRAPHIE

1. Van Maarseveen NM, Wessels E, de Brouwer CS, Vossen AC and Claas EC, Diagnosis of viral gastroenteritis by simultaneous detection of adenovirus group F, astrovirus, rotavirus group A, norovirus genogroups I and II, and sapovirus in two internally controlled multiplex real-time PCR assays, J Clin virol, 2010 ; 49(3) ; 205-10.
2. Ratcliff RM, Chang G, Kok T, Sloots TP. Molecular diagnosis of medical viruses. Curr Issues Mol Biol. 2007;9(2):87-102.
3. Pierce LR, Willey JC, Crawford EL, Palsule VV, Leaman DW, Faisal M, Kim RK, Shepherd BS, Stanoszek LM, Stepien CA. A new StaRT-PCR approach to detect and quantify fish Viral Hemorrhagic Septicemia virus (VHSV): enhanced quality control with internal standards. J Virol Methods. 2013;189(1):129-42.
4. Neisters HGM, Clinical virology in real time, 2002, J Clin Virol. 30-214-217.
5. Kim K, Park J, Chung Y, Cheon D, Lee IB, Lee S, Yoon J, Cho H, Song C, Lee KH, Use of internal standard RNA molecules for the RT-PCR amplification of the faeces-borne RNA viruses, J Virol Methods, 2002, 104(2):107-15.

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7, avenue de Scandinavie
91953 Les Ulis Cedex
FRANCE


EurobioPlex RNA QC

Control of RNA extraction and amplification

Real Time RT-PCR

REF **EBX-003**

100 tests

Storage conditions: keep all reagents at -20°C until use and after first use

Instructions for use

INTENDED USE

The EurobioPlex QC Control of RNA extraction and amplification is used to control variations that may occur during the RNA extraction steps, reverse transcription and qPCR. Control RNA added before extraction or in the PCR reaction ensure that a negative result is not due to an extraction problem or due to the presence of RT-PCR inhibitors at high quantity. CI-RNA can be detected on the CY5 channel (or equivalent). This kit is not intended to be used for other RT-PCR applications.

INTRODUCTION

During the preparation of RNA from samples, problems of RNA extraction, reverse transcription or the presence of RT-PCR inhibitors can occur. Inhibitors for example are very common from samples such as blood or stool, but not only. It is therefore important to monitor these phenomena. A result can be made false negative and may be due to a bad extraction or the presence of inhibitors. It is therefore essential to add internal RNA control to each sample.

PRINCIPLE OF DETECTION

The EurobioPlex QC Control of RNA extraction and amplification aims to reverse transcribe and amplify the RNA internal control using two primers and a probe. The specific RNA internal control probe is marked by the fluorophore Cy5 and emits a specific fluorescence following its hydrolysis during elongation. Internal control RNA amplification has been validated for temperatures of hybridization of primers between 55 ° C and 60 ° C and of real time PCR comprising 3 stages (denaturation, hybridization of primers, elongation at 72 ° C) or 2 steps (denaturation, hybridization + elongation simultaneously at 60 ° C).

This control must be added before the RNA extraction.

Two results can be obtained:

- 1 / internal control RNA test is positive (see §II analysis): RNA has been properly extracted, and there is no inhibitor of RT - PCR. The result can be validated.
- 2 / the internal control RNA test is negative: RNA was not extracted, the RT did not work well, or the presence of PCR inhibitors inhibits the PCR reaction. It is then recommended to repeat the extraction or dilute the sample.

DESCRIPTION AND CONTENT OF THE KIT

The EurobioPlex QC Control of RNA extraction and amplification is ready to use.

The kit contains:

- Internal control RNA (CI-RNA), which is encapsulated, non-infectious RNA template.

- Specific primers and a fluorophore Cy5-labeled probe for CI-RNA.

The kit does not include the enzymes reverse transcriptase or Taq polymerase.

Fluorescence is emitted and measured individually by an optical system during real-time PCR. The detection of the amplified fragment of the internal control RNA is achieved by a fluorometer using the channel shown below:

Target	Fluorophore	Excitation	Emission
Internal Control RNA (CI-RNA)	Cy5	647 nm	667 nm

Fluorophores to select depending of the qPCR instrument:

- Cy5 for ABI systems, Chromo 4/CFX96, Mx systems,
- Channel Alexa647 (SmartCycler II),
- Channel 670 (LC 480),
- Channel Red (RotorGene).

Content	Label	Volume	Reconstitution
Internal Control RNA	CI ARN	4 x 140 µl	Ready to use
Primers and probe for CI-RNA amplification	PPR mix	4 x 30 µl	Ready to use

CI-ARN is present at a concentration of 5.10^7 CFU/ml, approximately 5.10^7 copies/µl. We recommend adding 5 µl of stock CI-RNA per extraction and a final concentration by PCR reaction of $1,25.10^7$ copies/reaction (5 µl of a 100 µl extraction eluate) but this setting must be determined and validated for your specific application.

Required material not provided:

- ◇ Biological Hood
- ◇ qPCR instrument
- ◇ Micro centrifuge
- ◇ Vortex
- ◇ Plates / tubes for qPCR
- ◇ Micropipettes
- ◇ Filter tips for micropipettes
- ◇ Sterile microtubes
- ◇ Gloves

STORAGE

All reagents must be stored at -20°C. Storage at + 4°C is not recommended. Once prepared ARNs should be stored at -80°C.

All reagents can be used until the expiration date indicated on the kit.
Many freezing/thawing cycles (> 3x) must be avoided, and could lead to decrease in sensitivity.

CAUTIONS AND NOTES

Read carefully instructions before starting.

- ◇ The experiment must be performed by competent staff.
- ◇ Instruments must have been properly installed, calibrated and maintained according to the manufacturer's recommendations.
- ◇ Clinical samples are potentially infectious and must be processed under a laminar flow hood.
- ◇ The experiment must be performed according to good laboratory practices.
- ◇ Do not use this kit after expiration date.
- ◇ Many freezing/thawing cycles of the reagents must be avoided, and could lead to decrease in sensitivity.
- ◇ Once thawed, spin down briefly the tubes before use using a benchtop centrifuge.
- ◇ Use of ice or cooling block is advised in case of long delay due for instance to large number of samples or high temperature.
- ◇ Define three working areas: 1) Isolation of RNA, 2) Preparation of the reaction mix 3) Amplification / Detection of amplified products.
- ◇ Pipettes, reagents and other materials must not cross each area.
- ◇ Specific caution is required to preserve the purity of the reagents and reaction mixtures. The use of gloves is mandatory. Appropriate methods of preparation of RNA for cDNA transcripts should be used particularly to avoid RNase contamination.
- ◇ Always use filtered tips for micropipettes.
- ◇ Use specific lab coat and gloves in each working area.
- ◇ Do not pipette with mouth and do not eat, drink or smoke in the area.
- ◇ Avoid sprays.

PROCEDURE

General comments:

- 1-Mix the reagents of the kit before handling (by suction/delivery).
- 2-The internal control RNA contains a high concentration of matrix RNA. Handling must be done carefully to avoid contamination.
- 3-We recommend a concentration of $1,25 \cdot 10^7$ copies per PCR reaction, but this must be determined and validated in your specific application.

I-RNA extraction from biological samples:

Add the internal RNA control to each clinical specimen before extraction. Knowing that we recommend a final amount of $1,25 \cdot 10^7$ copies/test RT - PCR, the volume added will be determined by the user according to the kit and the different volumes of

solutions in the extraction procedure (Lysis / elution) used. Many RNA extraction kits are available from several manufacturers. You can use your own extraction system or a suitable commercial system by referring to the manufacturer's instructions.

For example: for samples of viruses, with the kit InnuPREP RNA virus kit (845-KS-4800250 available at Eurobio), it is recommended for a 150 µl sample volume and a volume of elution of 100 µl, add 5 µl of RNA control prior to extraction and use 5 µl of the eluate in the qRT - PCR.

It is recommended to add the internal control in the Lysis/sample buffer mixture and not directly in the Lysis buffer.

II- RT-PCR and Results analysis

RT-PCR PROTOCOL

- Prepare the RT - PCR reagents according to the Protocol that you are using in your real-time RT-PCR application, using 1 µl of PPR mix for a volume of PCR reaction of 25 µl.
- Perform the qRT-PCR test.
- In parallel tests the following controls:
 - Positive control: 1 µl PPR mix + 1 µl CI-ARN + volume of RT-PCR Mastermix of your protocol (according to providers' recommendations) + water for a total volume of 25 µl
 - Negative control: 1 µl PPR mix + volume of RT-PCR Mastermix of your protocol (according to providers' recommendations) + water for a total volume of 25 µl
- Once the qRT - PCR real-time has been performed, analyze the RNA amplification curve for CI-RNA.

Results analysis:

Threshold must be adjusted in the exponential phase to determine the value of Ct.

The positive control should have a Ct<35.

The negative control should have a Ct>40.

1-If Ct<35, RNA was properly extracted, and the results are validated.

2-If Ct >35, RNA was not properly extracted, or RT - PCR was inhibited. It is then recommended to repeat the extraction or to dilute the sample.

BIBLIOGRAPHY

- 1-Van Maarseveen NM, Wessels E, de Brouwer CS, Vossen AC and Claas EC, Diagnosis of viral gastroenteritis by simultaneous detection of adenovirus group F, astrovirus, rotavirus group A, norovirus genogroups I and II, and sapovirus in two internally controlled multiplex real-time PCR assays, J Clin virol, 2010 ; 49(3) ; 205-10.
- 2-Ratcliff RM, Chang G, Kok T, Sloots TP. Molecular diagnosis of medical viruses. Curr Issues Mol Biol. 2007;9(2):87-102.
- 3-Pierce LR, Willey JC, Crawford EL, Palsule VV, Leaman DW, Faisal M, Kim RK, Shepherd BS, Stanoszek LM, Stepien CA. A new StaRT-PCR approach to detect and quantify fish Viral Hemorrhagic Septicemia virus (VHSV): enhanced quality control with internal standards. J Virol Methods. 2013;189(1):129-42.
- 4-Neisters HGM, Clinical virology in real time, 2002, J Clin Virol. 30-214-217.
- 5-Kim K, Park J, Chung Y, Cheon D, Lee IB, Lee S, Yoon J, Cho H, Song C, Lee KH, Use of internal standard RNA molecules for the RT-PCR amplification of the faeces-borne RNA viruses, J Virol Methods, 2002, 104(2):107-15.

WASTE ELIMINATION

Be in accordance with the law on the elimination of waste of clinical infectious material.

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