EurobioPlex

Leptospire

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96/48/24 reactions

Approved on:

- CFX96™ Real Time PCR detection system (Bio-Rad) with analysis on CFX Manager version 3.1 (Bio-Rad)
- LightCycler®480 Instrument II (Roche) with analysis on LightCycler® 480 software v1.5 (Roche)
- T-COR8™-IVD (Tetracore Inc.) with T-COR 8SmartCT™ (auto v1) software

Operating instructions

Available at www.eurobio-scientific.com

The Summary of Safety and Performance Characteristics (SPC) will be made available by the certifying Notified Body on EUDAMED once it is operational. It can also be obtained on request.

Instructions for use EurobioPlex Leptospire

Table of contents

1. GENERAL INFORMATION	27
2. PURPOSE OF THE DEVICE	29
3. SYMBOLS.....	30
4. PRINCIPLE.....	31
5. KIT COMPONENTS.....	32
6. CONSERVATION AND STORAGE	32
7. EQUIPMENT REQUIRED NOT PROVIDED:.....	32
8. REAL-TIME PCR INSTRUMENT	33
9. WARNINGS AND PRECAUTIONS.....	33
10. PROTOCOL.....	34
i. COLLECTION OF SAMPLES.....	34
ii. NUCLEIC ACID EXTRACTION.....	34
iii. PERFORMING REAL-TIME RT-PCR.....	34
iv. DETAILED PROCEDURE.....	36
11. VALIDATION OF THE EXPERIMENT.....	37
12. DATA ANALYSIS AND INTERPRETATION.....	38
13. PERFORMANCE ANALYSIS.....	39
14. SPECIFICITIES OF THE T-COR 8™-IVD REAL-TIME PCR INSTRUMENT	40
15. QUALITY CONTROL	47
16. BIBLIOGRAPHY	47
17. WASTE DISPOSAL	48
18. INCIDENT DECLARATION.....	48
19. TECHNICAL SUPPORT	48

Instructions for use EurobioPlex Leptospire

1. GENERAL INFORMATION

Leptospirosis, sometimes called "*rat fever*", is an infectious disease of bacterial origin, caused by leptospires. Leptospires are a group of bacteria often grouped under the species *Leptospira interrogans*. This disease is ubiquitous. There has recently been a sharp increase in the number of reported cases. As a zoonosis it affects animals living in close proximity to humans, such as rats and pigs. The leptospire is a human pathogen.

In most parts of the world small mammals (rats, mice, etc.) are presumed to be the natural reservoir. Their urine almost always seems to be the direct or indirect source of human infections. But as with echinococcus, the natural infection of dogs (*canine leptospirosis*) and other canids and other canids also exists, probably very underestimated, and constitutes a risk for human beings. Leptospirosis occurs in different parts of the world, in Reunion, Wallis and Futuna Islands, the West Indies, French Guiana, tropical countries, South East Asia, Cambodia, Vietnam and Laos. Leptospires have even been detected in Western France, in the south and south-west of France near pig farms in humid regions.

A wide polymorphism of non-specific symptoms, and the diversity of organs affected, make clinical diagnosis very difficult. In addition, there are acute, sub-acute and chronic forms of the disease which evolve very differently over time, depending on the individuals infected, and on the serovars (possibly new or emerging), and which can be attenuated by vaccination. There are many serovars (e.g., icterohaemorrhagiae, canicola, pomona) that do not have a homogeneous antigenic signature, making it difficult to design effective vaccines. The real-time PCR test is of interest for studying the epidemiology of leptospirosis.

Leptospirosis manifests itself in a variety of forms, which makes its diagnosis difficult because it can be confused with a strong flu (high fever and aches). It may start with spread-out pain or localized pain (e.g., head pain) which, if not diagnosed in time, leads to a loss of speech and reasoning (due to high urea which increases in the blood due to kidney failure). In some cases, it is accompanied by jaundice (about 40% of cases in the trenches during the First World War) and/or nephritis. After dialysis and antibiotics, the patient recovers in 5-6 weeks if the disease was moderate (although bacteria can still be found in the patient's urine several weeks after the symptoms have disappeared). Otherwise, leptospirosis may progress to haemorrhaging haemorrhage - sometimes lethal - or severe kidney disease. Multi-organ forms are possible in weakened patients.

According to WHO recommendations, the biological confirmation of leptospirosis is based on either the isolation of the bacteria or the identification of its nucleic acids in biological samples, or on positive serology in an indicative clinical and epidemiological context. The use of bacterial culture is very limited due to the difficulty of creating it and the time needed for the bacteria to grow.

Serology is the most commonly used test for the diagnosis of leptospirosis. There is no standardised antigen either nationally or internationally for the detection of antibodies to leptospires.

Numerous tests for the detection of leptospirosis are available; they can be grouped into two types: direct detection tests (identification of the bacterium or bacterial molecules including DNA by real-time PCR) and tests detecting the patient's humoral response (serological tests). The microscopic agglutination test (MAT) is the reference serological technique for the confirmation of leptospirosis. The serological test is used to detect IgM class antibodies either by ELISA or on strips.

Early detection of leptospirosis can be done by real-time PCR, usually from plasma prepared from the patient's blood, which gives better results than serum. The optimal time to collect blood is before 8 days of illness. There is no point in performing the real-time PCR test in blood beyond 10 days after the symptomatological onset of the disease (end of the bacteremic period). The French

Instructions for use EurobioPlex Leptospire

National Authority for Health (HAS) has indicated that the simple PCR test (not in real time) is of no interest in the specific biological diagnosis of leptospirosis. The HAS considers that the real-time PCR technique is currently the most interesting for the early specific biological diagnosis of leptospirosis. It should be reserved for the first 8 days after the onset of fever and should be performed on a blood sample taken preferably before the start of antibiotic therapy. Because of the usual low concentration of leptospires in blood, any negative real-time PCR should be investigated serologically.

Instructions for use EurobioPlex Leptospire

2. PURPOSE OF THE DEVICE

The EurobioPlex Leptospire is a real-time reverse transcription-polymerase chain reaction (RT-PCR) assay designed for the qualitative detection of the presence or absence of leptospires by amplification of leptospire nucleic acids in a nucleic acid extract. The test is indicated for making a presumptive diagnosis of leptospirosis infection in humans, or to complement a proven or undetermined diagnosis. The extract is the starting material for the EurobioPlex Leptospire kit.

This system has been approved on the following types of samples:

- Plasma
- Serum
- Urine

The EurobioPlex Leptospire test is a medical device for *in vitro* diagnosis and must be used by qualified medical laboratory personnel. This device should not be recycled after use.

This device is not intended for screening of leptospires in blood or organ banks.

Instructions for use EurobioPlex Leptospire

3. SYMBOLS

Reference

Batch number

Temperature limit

Expiry date

Sufficient content for "N" reactions

Manufacturer

Date of manufacture

CE marked product

In vitro Diagnostics

Instructions for use

Warning

Do not use if packaging is damaged

Store in a dark place

Instructions for use EurobioPlex Leptospire

4. PRINCIPLE

The EurobioPlex Leptospire is a leptospire nucleic acid amplification test, together with an encapsidated RT-PCR extraction and inhibition control, which uses real-time RT-PCR amplification. The test is performed on the sample extract by means of a single reaction in a single well/tube.

The RT-PCR Extraction and Inhibition Control (IC-RNA) is used to determine possible variations that may occur during the RNA and DNA extraction from biological samples and real-time PCR amplification steps. This ensures that a negative result cannot be due to poor DNA and RNA extraction and/or the presence of too many PCR inhibitors.

Leptospire nucleic acids are detected using a TEXAS RED labelled probe. The RT-PCR extraction and inhibition control is detected using a Cy5 labelled probe. During the elongation of the amplification product, the probes emit specific fluorescence following their hydrolysis. The measurement of fluorescence intensities in real time is related to the accumulation of specific amplification products.

Table 1: Detection of targets by fluorophores

Target	Fluorophore	Excitement	Emission
Leptospire	Texas Red	485 nm	605 nm
PCR extraction and inhibition control (IC-RNA)	Cy5	647 nm	667 nm

Equivalent channels on different PCR instruments:

- **Texas Red** Channel (ABI Systems, SmartCycler II, Chromo4/CFX96™, Mx Systems, T-COR8™-IVD), LC Red 610 (LC480®), Orange Channel (RotorGene)
- Channel **Cy5** (ABI Systems, Chromo4/CFX96™, Mx Systems, T-COR8™-IVD), Channel Alexa647 (SmartCycler II), Channel 660 (LC 480), Channel Red (RotorGene)

Note: On LC480 instrument II: Apply colour compensation for TexasRed-Cy5 wavelengths (533-610, 618-660).

Instructions for use EurobioPlex Leptospire

5. KIT COMPONENTS

The EurobioPlex Leptospire real-time RT-PCR kit is ready to use for the specific detection of leptospire.

Fluorescence is emitted and measured individually by an optical system during the PCR process. Detection of the amplified fragments is performed by a fluorometer using the channels shown in Table 1.

The kit contains the reagents and enzymes required for the amplification of leptospire nucleic acids, and the RT-PCR extraction and inhibition control (Table 2).

Table 2: Kit components

Cap Colour	Contents of the kit	EBX-005A (96 tests)	EBX-005-48 (48 tests)	EBX-005-24 (24 tests)	Reconstruction
Red	Enzymes	2 X 1100µl	1100µl	4 x 130µl	Ready to use
Transparent	Oligomix	2 X 400µl	400µl	4 x 50µl	Ready to use
Blue	Water = negative control (CN-H2O)	2 x 1 ml	1 ml	1 ml	Ready to use
Yellow	Leptospire positive control (PC)	2 x 150 µl	150 µl	4 x 20 µl	Ready to use
White	RT-PCR extraction and inhibition control (IC-RNA)	2 x 650 µl	650 µl	4 x 80 µl	Ready to use

6. CONSERVATION AND STORAGE

All reagents should be stored between -15°C and -22°C.

All reagents can be used until the expiry date indicated on the kit label.

Multiple freeze/thaw cycles (>3x) should be avoided as this may reduce the sensitivity of the analysis.

7. EQUIPMENT REQUIRED NOT PROVIDED:

- Biohazard hood
- Real-time PCR apparatus
- Centrifuge for microtubes
- Vortex
- Plates/tubes for real-time PCR reaction
- Micropipettes
- DNase-free and RNase-free filter tips for micropipettes
- Sterile microtubes
- Gloves (talc-free)
- MICROLAB NIMBUS IVD® pipetting machine (Hamilton).

Instructions for use EurobioPlex Leptospire

8. REAL-TIME PCR INSTRUMENT

The EurobioPlex Leptospire kit has been developed and approved for use with the following real-time PCR instruments:

- CFX96™ Real Time PCR detection system (Bio-Rad) with analysis on CFX Manager version 3.1 (Bio-Rad)
- LightCycler®480 Instrument II (Roche) with analysis on LightCycler® 480 software v1.5 (Roche)
- T-COR8™-IVD (Tetracore Inc.) with T-COR 8SmartCT™ (auto v1) software

9. WARNINGS AND PRECAUTIONS

Read these instructions carefully before starting the procedure.

- ◇ This testing must be carried out by medical analysis laboratory technicians.
- ◇ Ensure that instruments have been installed, calibrated and maintained in accordance with the manufacturer's recommendations.
- ◇ It is the user's responsibility to use equipment other than those approved, and in this case, performance specifications are not guaranteed.
- ◇ Clinical specimens should be considered as potentially infectious material and should be prepared in a laminar flow hood.
- ◇ This testing must be carried out according to good laboratory practice.
- ◇ Do not use this kit after the expiry date indicated on the label.
- ◇ Avoid prolonged exposure of the reagents to light, which should be limited to the technical time needed to prepare the PCR plate, otherwise performance is not guaranteed.
- ◇ The kit is shipped in dry ice, and the kit components must arrive frozen. If one or more of the components arrives thawed, or if the tubes have been damaged during transport, contact Eurobio Scientific.
- ◇ Once the reagents have thawed, centrifuge the tubes briefly before use.
- ◇ It is recommended to define three distinct working areas: 1) Isolation of RNA/DNA, 2) Preparation of reaction mixture and 3) Amplification / Detection of amplified products.
- ◇ Wear separate gowns and gloves (without talc) in each work area.
- ◇ Pipettes, reagents and other working materials must not be moved between these areas.
- ◇ Particular care must be taken to maintain the purity of reagents and reaction mixtures.
- ◇ Appropriate nucleic acid preparation/extraction methods for quality RNA/DNA production and real-time RT-PCR application should be used, especially to avoid any risk of contamination with RNases/DNases.
- ◇ Use filter tips for micropipettes, DNase-free and RNase-free.
- ◇ Do not pipette by mouth. Do not eat, drink or smoke in the laboratory.
- ◇ Do not use the device if the packaging is damaged.
- ◇ Avoid aerosols.

Instructions for use EurobioPlex Leptospire

10. PROTOCOL

I. COLLECTION OF SAMPLES

- ◇ Collect the samples in sterile tubes.
- ◇  **The use of heparin as an anticoagulant is prohibited.**
- ◇ It is the responsibility of the user to control their own sample collection, transport, storage and extraction conditions to ensure that RNA/DNA extraction by suitable systems produces quality RNA/DNA.
- ◇ It is recommended that samples are extracted immediately, or stored according to the recommendations for sample storage before extraction (Table 3).

Table 3: Storage recommendations

Recommendations for maximum storage of samples before extraction	
Ambient temperature	2 h
+4°C	1 week
-20°C (preferably -80°C)	Long-term storage

Warning	
	◇ The user can refer to the recommendations issued by the World Health Organization or the French National Authority for Health for the proper storage of samples. ◇ The extracted RNAs should be stored at -80°C ◇ The transport of clinical samples is subject to local regulations for the transport of infectious agents.

II. NUCLEIC ACID EXTRACTION

It is the responsibility of users to ensure that the nucleic acid extraction system used is compatible with real-time RT-PCR technology. For this kit, we recommend using methods for extracting nucleic acids from plasma or urine samples and refer to the instructions of the supplier of the extraction kit used.

In the Leptospire kit, the IC-RNA on the CY5 channel can be added before extraction or in the PCR reaction. It ensures that a negative result cannot be due to an extraction problem or the presence of too much RT-PCR inhibitor to validate the test. We recommend the addition of 10 µl of IC-RNA per extraction and an elution volume of 50 µl. If the IC-RNA is only added to control the RT-PCR, it is added to the reaction mixture (1 microlitre per PCR reaction). See Real-time PCR Protocol for more details.

The IC-RNA is also available from Eurobio Scientific (Ref EurobioPlex EBX-003).

III. PERFORMING REAL-TIME RT-PCR


General remark:

The positive control as well as the extraction and RT-PCR amplification control (IC-RNA) contain high concentrations of template. Handling must be done carefully to avoid contamination. To check the correct functioning of the PCR, it is necessary to test the positive control PC as well as a negative control (supplied water = CN-H₂O + IC-RNA) (see II-2/6) of the real-time RT-PCR protocol).

On T-COR8™-IVD, for proper use of the kit, it is necessary that these checks are carried out the first time a new kit is used, only regular checks can be carried out thereafter.

Instructions for use EurobioPlex Leptospire

Diagram of the procedure:

Instructions for use EurobioPlex Leptospire

IV. DETAILED PROCEDURE

- 1) Homogenise the Enzyme tube, and vortex Oligomix, PC and IC-RNA, then centrifuge.
- 2) Prepare the Mastermix as follows, where N is the number of reactions (including positive and negative controls), allow for the amount of Mastermix for at least N+3 reactions* (refer to part 1-a or 1-b of the previous diagram as appropriate).

Case (a): Extraction of sample WITHOUT addition of IC-RNA

Number of Reaction(s)	1	N+3*
Enzymes	15 µl	(N+3) x 15 µl
Oligomix	5 µl	(N+3) x 5 µl
IC-RNA	1 µl	(N+3) x 1 µl
Total volume Mastermix	21 µl**	(N+3) x 21 µl

Case (b): Extraction of the sample with the addition of IC-RNA

Number of Reaction(s)	1	N+3*
Enzymes	15 µl	(N+3) x 15 µl
Oligomix	5 µl	(N+3) x 5 µl
Total volume Mastermix	20 µl**	(N+3) x 20 µl

* For use on the T-COR8™-IVD, we recommend making a Mastermix for N+1 reactions for a one-by-one usage of the kit.

**The difference in reaction volume between case (a) and case (b) is negligible and has no effect on performance.

- 10) Homogenise the Mastermix prepared in 2) and centrifuge briefly.
- 11) Dispense 20 µL of the Mastermix** using a micropipette and filter tips into each real-time PCR microplate tube/well.
- 12) Add 5 µL of the extracted samples
- 13) In parallel, carry out the following checks (see this specific point for T-COR8™-IVD pages 14 to 17):
 - Positive control: 20 µL Mastermix + 5 µL CP positive control
 - Negative control:
 - o Case (a) Extraction of sample WITHOUT addition of IC-RNA:
20µl Mastermix + 5µl supplied water (CN-H2O)
 - o Case (b) Extraction of sample WITH addition of IC-RNA:
20µl Mastermix + 4µl supplied water (CN-H2O) + 1µl IC-RNA
- 14) Immediately seal the tubes or plate with adhesive film to avoid contamination.
- 15) Briefly centrifuge to collect the reaction mixture at the bottom of the tubes or microplate wells.
- 16) Perform the following program on the real-time PCR instrument (on the T-COR8™-IVD no manual programming is required due to the barcodes - see pages 18 to 20). The programme lasts about 1 hour and 15 minutes.

Programme	Temperature	Duration	Cycle(s)	
Reverse Transcription	50°C	10 min	1	-
Denaturation	95°C	3 min	1	-
Amplification	95°C	10 sec	45	-
	60°C	30 sec		Fluorescence acquisition

Instructions for use EurobioPlex Leptospire

Note 1: Two optical systems are available for the LightCycler® 480 (Roche): only the "System II" is compatible with the use of the kit. Apply colour compensation for Texas Red-Cy5 wavelengths (533-610,618-660).

Note 2: For Applied Biosystems devices, select "NONE" in "PASSIVE REFERENCE".


Note 3: In the Rotorgene™, calibrate the signal by clicking on "GAIN OPTIMIZATION".

Note 4: On CFX96™ (Bio-Rad), launch the run from version 1.6 or later, then analyse with version 3.1 (see § Validation of the experiment).


11. VALIDATION OF THE EXPERIMENT

Post-acquisition data analysis on a CFX96™ PCR instrument (Bio-Rad) must be performed using CFX Manager software version 3.1 (Bio-Rad). In order to upgrade to this version from a run started on an earlier version, please follow the following procedure: at the end of the run, the data file with the suffix .pcrd must be opened and processed with CFX Manager version 3.1 (Bio-Rad).

If the run was started with the CFX Manager v1.6 software for example, to open a data file with the CFX Manager v3.1 software, click on the CFX Manager v3.1 icon. The home screen appears.


- Click on File and select Open then Data File.

- Select the file you wish to analyse and click Open.

The "drift correction" option must be applied in the "Settings" tab as shown in the image below: click on the Settings tab, then on Baseline Setting and Apply Fluorescence Drift Correction.

Instructions for use EurobioPlex Leptospire

Once this is done, the analysis can begin.

For the assay to be valid, the Ct values for the controls must be as follows (Table 4); otherwise the experiment cannot be validated.

On T-COR8™-IVD, for proper use of the kit, it is necessary that these checks are carried out the first time a new kit is used only, regular checks can be carried out thereafter.

Table 4:

Positive control	
TEXAS RED	Ct < 27
Negative control	
TEXAS RED	Ct not determined
CY5	Ct ≤ 45

On T-COR8™-IVD, the interpretation is automatically generated with the barcodes

12. DATA ANALYSIS AND INTERPRETATION

For clinical samples, the following results are possible:

***Ct threshold for sample positivity:**
TEXAS RED Channel *Leptospire*: + Positive => Ct ≤ 45

PCR signal			Validity of the test/comment
TEXAS RED	CY5	Leptospire	
+	+	Positive	Valid
-	+	Negative	Valid
+	-	Positive	Possible extraction or RT-PCR inhibition that does not prevent detection of the bacteria; valid
-	-	Not interpretable	RT-PCR extraction or inhibition problem -dilute 5 x sample; if same result re-extract sample

Instructions for use EurobioPlex Leptospire

13. PERFORMANCE ANALYSIS


Detection limit/analytical sensitivity :

5 CP copies/ μ l (1 copy of leptospire = 5 genome copy equivalents (Bourhy et al. 2011)).

On the bacterial strain *Leptospire interrogans*: 100 bacteria/ml in plasma

250 bacteria/ml in serum, and 500 bacteria/ml in urine.

Example of an experiment performed on a CFX96™ real-time PCR thermocycler (Bio-Rad):

Signal reproducibility on the Leptospire channel (Texas Red)

Reproducibility	Coefficient of Variation (%) on the Positive Control	Coefficient of Variation (%) on RNA
Intra-run	1.35	2.44
Inter-runs	5.15	2.04
Inter-batch	1.83	2.45

Specificity: > 99 %

Specificity was assessed on a total of 98 extracted samples pre-tested and confirmed negative for leptospirosis. 45 of these samples were known to be positive for :

- *Neisseria meningitidis* (n=5)
- *Streptococcus pneumoniae* (n=5)
- *Haemophilus influenzae b* (n=1)
- *Pseudomonas aeruginosa* (n=4)
- *Streptococcus agalactiae* (n=1)
- *Escherichia coli* (n=4)
- *Legionella pneumophila* (n=1)
- *Bartonella* (n=5)
- *Coxiella* (n=3)
- *Staphylococcus aureus* (n=1)
- *Mycobacterium avium* (n=1)
- *Staphylococcus oralis/mitis/sanguis* (n=1)
- *Francisella tularensis* (n=1)
- *Klebsiella kingae* (n=2)

Instructions for use EurobioPlex Leptospire

- *Tropheryma whipplei* (n=1)
- *Chikungunya* (n=1)
- Dengue (n=1 Dengue 1+ n=1 Dengue 2+ n=1 Dengue 3+ n=1 Dengue 4+ n=3 of unknown serotype)
- Zikavirus (n=1)

No cross-reactivity was therefore observed with these pathogens.

Sensitivity: > 99 %

Sensitivity was assessed on a total of 50 pre-tested and confirmed leptospirosis-positive samples, of which 39 were plasma from blood/EDTA, and 11 were urine.

	Eurobio PCR		
	Positive	Negative	Total
Reference PCR			
Positive	50	0	50
Negative	0	98	98
Total	50	98	148

The overall concordance is > 99%.

The tests were performed on CFX96™.

The extraction was performed with the NucliSENS essayMAG extractor (Biomérieux) following the supplier's instructions.

14. SPECIFICITIES OF THE T-COR 8™-IVD REAL-TIME PCR INSTRUMENT

T-COR8™-IVD is a real-time PCR instrument with 8 independently programmable wells that can be used, patient by patient, independently of each other, in terms of the tests/assays performed, the thermal programs, and the timing of the PCR run.

All 8 wells can also be used simultaneously with the same "Test/Assay".

Note: Mix thoroughly after placing the sample in the tube with repeated pipetting, avoiding the formation of bubbles and ensuring that the volume of liquid is properly located at the bottom of the tube. Close each tube tightly after each deposit of control or sample to avoid contamination.

Checks

On the T-COR8™-IVD, for correct usage of the kit, it is necessary that these checks are carried out the first time a new kit is used only, regular checks can be carried out thereafter.

Instructions for use EurobioPlex Leptospire

Subsequently, the RT-PCR inhibition check (IC-RNA) ensures that within the sample tested, there is no RT-PCR inhibition, that the extraction has been successful and that the PCR is functioning correctly.

Patient tests

<i>With validation of positive and negative controls only once on first use</i>	24 reactions
Number of patient tests possible patient by patient with thawing/freezing 12 times maximum (3 times per tube/4 tubes per reagent)	12 patients

Fluorophores

There are 4 fluorophore combinations available on the T-COR8™-IVD.

For all EBX, including EBX-005, the FAM/HEX/Texas Red/Cy5 combination is the one chosen. If some channels are not used in some kits, disregard them when analysing the results.

Use of Barcodes (available on pages 18 to 20)

1- Select Menu > New Analysis/New Run

2- Select Barcode

3- Scan the desired barcode on the right side of the device:

- or for the positive control (Barcode EBX-005 Pos Ctrl),
- or for the negative control (Barcode EBX-005 Neg Ctrl),
- or for a sample (Barcode Assay: EBX-005)

The instrument automatically selects one of the 8 available wells. Follow the instructions given by the instrument.

4- Lift the valve of the well x selected by the device, and check that the blue LED light is on, then select "Yes".

5- Position the tube in the corresponding well and fold the valve back, then select "Next".

6- If you wish to name it (optional), select "sample x", and name it. Select "Accept".

7- Select "Next"

8- To add/test another control or sample, select "Add well" and start again at point 3-

9- Select "Start run"

***Note:** the name of a sample or run cannot be changed or added to after the run has ended. It must be done before or during the run.*

Instructions for use EurobioPlex Leptospire

Presentation and Interpretation of Results

Ct results and amplification curves are available, from the SmartCT™ Values Table (Ct values assigned on each channel), and Ct vs. cycle graphs on each channel, both of which can be viewed in real time on the machine.

An analysis for the negative and positive controls, as well as for the samples, is generated automatically. This is available in "Interpretations" at the end of the run in the "View" window.

For the positive and negative control, the following results are possible:

"Neg Ctrl Fail": Invalid
"Neg Ctrl Valid": Valid
"Pos Ctrl Fail": Not valid
"Pos Ctrl Valid": Valid


Determination of status for samples:

"Detected": Positive → green box
"Not detected": Negative → red box
"Invalid -> retest → orange box

Example of presentation of automatic interpretation of results on the instrument and in the report:

- Example for valid negative control (well 4) and positive control (well 6) indicated "Detected", for Invalid negative control (well 5), for a negative sample indicated "Not Detected" (sample well 7) and for a positive sample detected for the Leptospira target indicated (Detected Leptospira) (sample well 8)

❖ On the instrument display:



❖ In the report:

Instructions for use EurobioPlex Leptospire

Summary							
Well	Sample	Assay	FAM	HEX	TxR	Cy5	Call
4	CNeg	EBX-005 Neg Ctrl				36.1	Detected • Neg Ctrl Valid
5	CNeg	EBX-005 Neg Ctrl					Invalid
6	CPos	EBX-005 Pos Ctrl			21.4	34.7	Detected • Pos Ctrl
7	Negative sample	EBX-005				34.8	Not Detected
8	Positive sample	EBX-005			30.2	34.5	Detected • Leptospira
							Leptospira Neg


Example of amplification curves


Instructions for use EurobioPlex Leptospire

Barcodes on T-COR8™-IVD for EBX-005


EBX-005 Pos Ctrl
Positive control PC

EBX-005 Neg Ctrl
Water = negative control (CN-H2O)

Assay: EBX-005

Instructions for use EurobioPlex Leptospire

15. QUALITY CONTROL

In accordance with Eurobio Scientific's ISO EN 13485 certified quality management system, each batch of EurobioPlex Leptospire is tested according to predefined specifications to ensure consistent product quality.

16. BIBLIOGRAPHY

- Giry C, Roquebert B, Li-Pat-Yuen G, Gasque P, Jaffar-Bandjee M.C. Simultaneous detection of chikungunya virus, dengue virus and human pathogenic *Leptospira* genomes using a multiplex TaqMan® assay. *BMC Microbiol.* 2017;17(1):105.
- Perez J, Goarant C. Rapid *Leptospira* identification by direct sequencing of the diagnostic PCR products in New Caledonia. *BMC Microbiol.* 2010;10:325.
- Goarant C, Laumond-Barny S, Perez J, Vernel-Pauillac F, Chanteau S, Guigon A. Outbreak of leptospirosis in New Caledonia: diagnosis issues and burden of disease. *Trop Med Int Health.* 2009;14(8):926-9.
- Tubiana S, Mikulski M, Becam J, Lacassin F, Lefèvre P, Gourinat AC, Goarant C, D'Ortenzio E. Risk factors and predictors of severe leptospirosis in New Caledonia. *PLoS Negl Trop Dis.* 2013;7(1):e1991.
- van de Werve C, Perignon A, Jauréguiberry S, Bricaire F, Bourhy P, Caumes E. Travel-related leptospirosis: a series of 15 imported cases. *J Travel Med.* 2013;20(4):228-31.
- Hochedez P, Escher M, Decoussy H, Pasgrimaud L, Martinez R, Rosine J, Théodose R, Bourhy P, Picardeau M, Olive C, Ledrans M, Cabie A. Outbreak of leptospirosis among canyoning participants, Martinique, 2011. *Euro Surveill.* 2013;18(18):20472.
- Bourhy P, Herrmann Storck C, Theodose R, Olive C, Nicolas M, Hochedez P, Lamaury I, Zinini F, Brémont S, Landier A, Cassadou S, Rosine J, Picardeau M. Serovar diversity of pathogenic *Leptospira* circulating in the French West Indies. *PLoS Negl Trop Dis.* 2013;7(3):e2114.
- Bourhy P, Brémont S, Zinini F, Giry C, Picardeau M. Comparison of real-time PCR assays for detection of pathogenic *Leptospira* spp. in blood and identification of variations in target sequences. *J Clin Microbiol.* 2011;49(6):2154-60.
- Legkoby T, David DJ, Lee Robin SH. Biological diagnosis of leptospirosis, Avis sur les actes, Haute autorité de santé, service évaluation des actes professionnels, www.has-santé.fr
- Merien F, Portnoi D, Bourhy P, Charavay F, Berlioz-Arthaud A, Baranton G. A rapid and quantitative method for the detection of *Leptospira* species in human leptospirosis. *FEMS Microbiol Lett.* 2005;249(1):139-47.
- Levett PN. Sequence-based typing of *Leptospira*: epidemiology in the genomic era. *PLoS Negl Trop Dis.* 2007;1(2):e120.
- Levett PN, Morey RE, Galloway RL, Turner DE, Steigerwalt AG, Mayer LW. Detection of pathogenic leptospires by real-time quantitative PCR. *J Med Microbiol.* 2005 Jan;54(Pt 1):45-9.
- Slack A. Leptospirosis. *Aust Fam Physician.* 2010 Jul;39(7):495-8. Review.

Instructions for use EurobioPlex Leptospire

- Stoddard RA. Detection of pathogenic *Leptospira* spp. through real-time PCR (qPCR) targeting the LipL32 gene. *Methods Mol Biol.* 2013;943:257-66.
- Ren SX, Fu G, Jiang XG, Zeng R, Miao YG, Xu H, Zhang YX, Xiong H, Lu G, Lu LF, Jiang HQ, Jia J, Tu YF, Jiang JX, Gu WY, Zhang YQ, Cai Z, Sheng HH, Yin HF, Zhang Y, Zhu GF, Wan M, Huang HL, Qian Z, Wang SY, Ma W, Yao ZJ, Shen Y, Qiang BQ, Xia QC, Guo XK, Danchin A, Saint Girons I, Somerville RL, et al. Unique physiological and pathogenic features of *Leptospira interrogans* revealed by whole-genome sequencing. *Nature.* 2003;422(6934):888-93.
- Cryopréservation de tissus, cellules et liquides biologiques issus du soin, Recommandations de bonnes pratiques, www.has-santé.fr

17. WASTE DISPOSAL

Dispose of all waste in accordance with the legislation on DASRI.

18. INCIDENT DECLARATION

Any serious incident occurring in connection with the device shall be notified to EUROBIO SCIENTIFIC and to the competent authority of the Member State in which the user and/or patient is established.

19. TECHNICAL SUPPORT

For assistance with our products, please contact our technical support.


EUROBIO SCIENTIFIC's customer service can be reached by e-mail at adv@eurobio-scientific.com or by phone at +33 (0)1.69.07.94.77.

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