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---

## Validated on:

- CFX96™ Real Time PCR detection system (Bio-Rad) and analysis with the CFX Manager version 3.1 (Bio-Rad)
- LightCycler®480 Instrument II (Roche) and analysis with the LightCycler® 480 software v1.5 (Roche)
- Applied Biosystems® 7500 Real-Time PCR Systems (Applied Biosystems) and analysis with the 7500 Software v2
- T-COR 8®-IVD (Tetracore Inc.) with a T-COR 8 SmartCT™ (auto v1) software

Instructions for use

Available at [www.eurobio-scientific.com](http://www.eurobio-scientific.com)

*A Summary of Safety and Performance Characteristics (SSCP/RCSP) will be made available by the Notified Body on EUDAMED once the later is operational. It can also be obtained upon request.*

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---

## 1. General informations

Chikungunya and Dengue viruses are Arthropod-Borne viruses (Arbovirus), transmitted by arthropods (mosquitoes, ticks, sandflies, etc).

Dengue virus belongs to the *Flaviviridae* family and Chikungunya virus to the *Alphaviridae* from the *Togaviridae* family. They are transmitted by various mosquitoes species of the genus *Aedes* such as *aegypti* (mainly in tropical and sub-tropical America and in the south of the United States), *albopictus* (mainly in Asia), *furcifer*, *africanus*, *taylori* and *luteocephalus* (with a worldwide distribution, especially in the tropical regions). The first two species are identified to date as the main epidemic vectors, because of their adaptation to the human settlement areas. In France, the southern regions are more particularly concerned.

**Dengue:** There are 4 serotypes of Dengue (DENV 1-4). The Dengue genome is a single-stranded RNA. The virus replication and virus transcription occur in the cytoplasm of the host cell. The illness begins with fever, muscle and joint pain. In most of the cases, the infection do not present complications, but in 20% of the cases it can be more severe and sometimes life-threatening, especially in the case of secondary infections, with haemorrhagic shock, as an infection by one strain does not protect the host from being infected by another strain.

**Chikungunya:** The Chikungunya virus (CHIKV) is a small enveloped spherical single-stranded virus of about 12 kbp. The viraemic phase of the Chikungunya infection usually last 3 to 10 days. The major clinical signs are the sudden onset of high fever (> 40 °C), which lasts 24 to 48 hours, with arthralgia and myalgia, very painful polyarthralgia that mainly affect the joint of the extremities such as wrists, ankles or phalanges, with a pronounced joint swelling, cutaneous symptoms such as maculo-papular or erythematous rash or oedema. Infections with CHIKV can evolve into a chronic disease. Symptoms are similar to those of malaria, of dengue infection, and it is not unusual to find them both simultaneously. For a certain number of patients (10 to 70 %), clinical recurrences of the disease are observed several months after the acute phase. There is no specific treatment for chikungunya.

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---

## 2. Intended purpose

The Chikungunya/Dengue test is a real-time amplification test by reverse transcription- polymerase chain reaction (RT-PCR), designed for the qualitative detection of the presence or absence of Chikungunya and of the 4 Dengue virus serotypes, without differentiation, in an RNA nuclear acid extract.

The Chikungunya/Dengue test is indicated to assist in making a diagnosis of suspected infection in humans, or to complete a proven or undetermined serological diagnosis. The RNA extract is the starting material for the Eurobioplex Chikungunya/Dengue kit. The RNA is extracted from the patient's plasma, or serum.

This amplification system has been validated on plasma and serum samples, as well as on encapsulated RNA bacteriophage (CI-ARN).

The EurobioPlex Chikungunya/Dengue test is an in vitro diagnostic medical device. It must be used by qualified medical analysis laboratory personnel. For single use only, it should not be recycled after use.

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---

## 3. Symbols

	Reference
	Batch number
	Lower and upper limits of storage temperature
	Expiry date
	Enough content for « N » reactions
	Manufacturer
	CE marked product
	In Vitro Diagnostic
	Instruction for use
	Caution
	Date of manufacture
	Keep away from sunlight
	Do not use in case of damaged packaging

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## 4. Principle of the test

The Chikungunya/Dengue test is a ribonucleic acid (RNA) amplification test for Chikungunya and for the 4 Dengue serotypes, without any differentiation, as well as for an RNA extraction and RT-PCR inhibition control, encapsulated, which uses the real-time RT-PCR amplification. The test is performed on the RNA extracted from a sample with a single reaction, in a single well.

The primers and probes of the kit allow for the RNA detection of the 4 serotypes (screening) without typing.

The RNA extraction and RT-PCR inhibition control (CI-ARN) allows to determine possible variations that may occur during the course of RNA extraction steps from biological samples and real-time PCR amplification. It is a way to check that a negative result cannot be the consequence of a bad RNA extraction and/or the presence of PCR inhibitors in too high quantities.

The Chikungunya virus RNA is detected using a FAM labeled probe and the Dengue one using an HEX labeled probe. The extraction and PCR inhibition control is detected with a Cy5 labeled probe. They all emit a specific fluorescence.

**Table 1: Target detection according to the fluorophores**

Target	Fluorophore	Excitation	Emission
CHIKUNGUNYA	FAM	495 nm	515 nm
DENGUE	HEX	535 nm	555 nm
RNA Control (CI-ARN)	Cy5	650 nm	670 nm

Alternative channels on different PCR instruments:

- Channel FAM (Systems ABI, SmartCycler II, Systems Mx, Chromo4/CFX96, T-COR 8@-IVD), Channel 510 (LC 480), Channel Green (RotorGene)
- Channel HEX (Chromo 4/CFX96, Systems Mx, T-COR 8@-IVD), Channel VIC (Systems ABI), Channel Alexa532 (SmartCycler II), Channel 580 (LC 480), Channel Yellow (RotorGene)
- Channel Cy5 (Systems ABI, Systems Mx, Chromo4/CFX96, T-COR 8@-IVD), Channel Alexa647 (SmartCycler II), Channel 660 (LC 480), Channel Red (RotorGene)

**Note:** On the LC480 instrument II, apply a color compensation for the wavelengths FAM-HEX/VIC-Cy5 (465-510,533-580,618-660).

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## 5. Kit components

The Eurobioplex Chikungunya/Dengue real-time RT-PCR kit is ready to use for these viruses specific detection.

The fluorescence is emitted and measured individually by an optical system during the PCR. The detection of the amplified fractions is performed by a fluorimeter using the channels as mentioned in Table 1.

The kit contains the reagents and enzymes necessary for the viral RNA amplification from Chikungunya and Dengue, and for the extraction and amplification by RT-PCR of the control (Table 2).

Both EBX-009A (96 tests) and EBX-009A-48 (48 tests) formats are suitable for use with the pipetting system MICROLAB NIMBUS IVD (Hamilton) as they provide sufficient volume, including dead volume, to carry out respectively 96 or 48 tests, on that type of machine.

**Table 2 : Kit components**

Cap colour	Kit content	EBX-009A (96 tests)	EBX-009A-48 (48 tests)	EBX-009 (96 tests)	EBX-009-48 (48 tests)	EBX-009-24 (24 tests)	Reconstitution
<b>Red</b>	Enzymes	2 X 1080 µl	1080 µl	4 x 410 µl	2 x 410 µl	4 x 100 µl	Ready to use
<b>Transparent</b>	Oligomix	2 X 360 µl	360 µl	4 x 135 µl	2 x 135 µl	4 x 60 µl	Ready to use
<b>Blue</b>	Water = Negative control (CN-H2O)	2 x 1 ml	1 ml	1 ml	1 ml	1 ml	Ready to use
<b>Yellow</b>	Positive control for Chikungunya and Dengue (CP- chikDeng)	2 x 145 µl	145 µl	4 x 40 µl	2 x 40 µl	4 x 40 µl	Ready to use
<b>White</b>	RNA control (CI-ARN)	2 x 625 µl	625 µl	4 x 300 µl	2 x 300 µl	4 x 75 µl	Ready to use

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---

## 6. Storage

All reagents should be stored between entre -15°C and -22°C.

All reagents can be used until the expiration date mentioned on the kit label.

Several freeze / thaw cycles (>3x) should be avoided, as it could decrease the analysis sensitivity.

## 7. Materials required but not provided

- Biological hood
- Real-time PCR instrument
- Centrifuge suitable for microtubes
- Vortex
- Plates / tubes for real-time PCR
- Micropipettes
- DNase-free and RNase-free tips containing filters for micropipettes
- Sterile microtubes
- Disposable gloves (talc-free)
- Automatic pipetting machine MICROLAB NIMBUS IVD (Hamilton) if EBX-009A (96 tests) and EBX-009A-48 (48 tests) are used in an automated version

## 8. Real time PCR Instrument

The EurobioPlex Chikungunya/Dengue kit has been developed and validated for use with the following real-time PCR instruments:

- CFX96™ Real Time PCR detection system (Bio-Rad) and analysis with the CFX Manager version 3.1 (Bio-Rad)
- LightCycler®480 Instrument II (Roche) and analysis with the LightCycler® 480 software v1.5 (Roche)
- Applied Biosystems® 7500 Real-Time PCR Systems (Applied Biosystems) and analysis with the 7500 Software v2
- T-COR 8®-IVD (Tetracore Inc.) with the T-COR 8 SmartCT™ (auto v1) software

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## 9. General precautions

**Carefully read these instructions before starting the procedure.**

- ◇ This experimentation should be carried out by qualified medical analysis laboratory technicians.
- ◇ Check that the instruments have been installed, calibrated and maintained according to the manufacturer's instructions.
- ◇ Clinical samples should be considered as potentially infectious material and should be prepared under a laminar air flow cabinet.
- ◇ This experimentation must be performed according to good laboratory practice
- ◇ Do not use this kit after the expiration date mentioned on the label.
- ◇ The kit is shipped under dry ice and the kit components should stay frozen. If one or more of these components arrive unfrozen, or if the tubes are damaged during transportation, please contact Eurobio Scientific.
- ◇ Several freeze / thaw cycles should be avoided, as it could decrease the analysis sensitivity.
- ◇ Avoid prolonged exposure of reagents to light. It should be limited to the technical time necessary to prepare the PCR plate, otherwise performances are not guaranteed.
- ◇ Once the reagents are thawed, briefly centrifuge the tubes before use.
- ◇ It is recommended to perform the activities in three separate areas: 1) RNA isolation,  
2) Preparation of the reaction mixture and 3) Amplification/Detection of the amplified products
- ◇ Wear separate lab coats and gloves (talc-free) in each working area.
- ◇ Pipettes, reagents and other materials should not be moved around between these areas.
- ◇ Special care must be taken to maintain the purity of reagents and reaction mixtures.
- ◇ Appropriate nucleic acids preparation/extraction methods for quality RNA production and RT-PCR application should be used, in order to avoid any risk of contamination with RNases.
- ◇ Use disposable tips with filters for the micropipettes, RNase-free and DNase-free.
- ◇ Do not pipette by mouth. Do not eat, drink or smoke in the laboratory.
- ◇ Avoid aerosols.

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- ◇ This kit is not intended for single use only. When re-using the device, it is necessary to follow the recommendations concerning the authorised freeze/thaw cycles, and the precautions to prevent reagents contamination.
- ◇ The kit components should not be used separately (neither with other reagents nor with reagents from other batches).
- ◇ Thawing of the reagents should be done carefully in order not to impair the system performance (at +2°C/+8°C or room temperature).

## 10. Protocol

### 10.1 Sample collection

- ◇ Collect samples in sterile tubes.
  - For an RNA extraction from serum, collect samples on dry tubes
  - For an RNA extraction from plasma, collect samples on EDTA tubes.
  -  - Use of heparin as an anticoagulant is prohibited.
- ◇ It is up to the user to control his own conditions of collection, transport, storage and extraction of samples, so that RNA extraction by suitable systems produces quality RNA.
- ◇ It is recommended that samples be extracted immediately or stored according to samples storage recommendations before extraction (Table 3).
- ◇ Bibliographic references in Section "16. References" provide guidelines on samples and RNA stability.

**Table 3: Blood storage recommendations, before plasma or serum preparation and before extraction**

Recommendations for storage before plasma and serum preparation from blood	
<24h	
Recommendations for maximum sample storage before extraction	
Room temperature	2h
+2°C/+8°C	5 days
<-70°C (preferred over -20°C)	Long term storage (>5days and maximum 2 months at -20°C)

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<b>Caution</b>	
	<ul style="list-style-type: none"><li>◇ The user can refer to the recommendations issued by the World Health Organization or the 'Haute Autorité de Santé' for the proper storage of samples.</li><li>◇ The extracted RNA must be stored at &lt;-70°C in order to ensure their stability. After one year storage of the RNA at &lt;-70°C, an increase of the Ct obtained may be observed. It is recommended to re-extract a biological sample which has been stored for more than one year. It is recommended to limit to 3 the number of RNA thawings.</li><li>◇ Clinical samples transport should follow local regulations for infectious agents.</li></ul>

## 10.2 RNA extraction

It is the user responsibility to ensure that the nucleic acids extraction system used is compatible with the real-time RT-PCR technology. For this kit, we recommend using virus RNA extraction methods from plasma or serum samples and referring to the manufacturer instructions for the extraction kit used.

Please refer to Table 3 recommendations for the RNA storage.

In the Chikungunya/Dengue kit, the CI-ARN on channel CY5 can be added before the extraction or in the PCR reaction. It is a way to check that a negative result cannot be the consequence of an extraction problem or of the presence of RT-PCR inhibitors in too high quantities to validate the test.

We recommend adding 10 µl of CI-ARN per extraction and using a 50 µl elution volume at the end of the extraction. If the CI-ARN is only added in order to check the RT-PCR, it is added to the reaction mixture (1 microlitre per PCR reaction). See the real-time PCR Protocol for more details.

The CI-ARN is also available at Eurobio Scientific (Ref EurobioPlex EBX-003).

Performances, as described in Section 13 of this document, have been obtained with the following extraction techniques: Qiagen columns (QIAamp Ultrasens Virus kit Ref: 53704), or extractor EZ1 (Qiagen) or Maelstrom TANbeads (Taiwan Adadvanced Nanotech).

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## 10.3 Real-time PCR procedure

### General remark:

The positive control as well as the extraction and RT-PCR amplification control (CI-ARN) contain high matrix concentrations. Handling must be done with care in order to prevent contamination. In order to check the proper PCR functioning, it is necessary to test the positive control CP-ChikDeng as well as a negative control (water supplied = CN-H2O + CI-ARN)

(See II-2/6 in the real time RT-PCR protocol). On the T-COR 8™-IVD, for a proper use of the kit, it is necessary that these controls checks are performed when a new kit is first used, regular checks can be done later on.

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## Procedure summary:

### 1 - PREPARATION OF REACTION MIXTURE / MASTERMIX

(a) If samples extracted without CI-ARN (b) If samples extracted with CI-ARN

Number of reactions	N+3**	Number of reactions	N+3**
Enzymes	(N+3) x 15 µl	Enzymes	(N+3) x 15 µl
Oligomix	(N+3) x 5 µl	Oligomix	(N+3) x 5 µl
CI-ARN	(N+3) x 1 µl		
Total volume Mastermix	(N+3) x 21 µl*	Total volume Mastermix	(N+3) x 20 µl*

\*\* N+1 on T-COR 8®-IVD

### 2 - PREPARATION OF REACTIONS

#### Sample

20 µl Mastermix  
+5 µl RNA sample

#### Positive control

20 µl Mastermix + 5 µl CP-ChikDeng

#### Negative control

20 µl Mastermix  
+ 5 µl Water molecular biology (CN-H2O)

#### Sample

20 µl Mastermix+  
5 µl RNA sample

#### Positive control

20 µl Mastermix+ 5 µl CP-ChikDeng

#### Negative control

20 µl Mastermix  
+ 4 µl Molecular biology grade water (CN-H2O) +1 µl CI-ARN

### 3 - REAL-TIME PCR INSTRUMENT (except T-COR 8®-IVD - automatic programming with barcodes - see p.26 to 27)

Program	Temperature	Time	Cycle(s)	
Reverse	50 °C	10 min	1	-
Denaturation	95°C	3 min	1	-
Amplification	95°C	10 sec	45	-
	60°C	30 sec		Fluorescence acquisition

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## 10.4 Procedure in details

- 3) Homogenize the enzymes tube, vortex Oligomix, CP-ChikDeng and CI-ARN, then centrifuge.
- 4) Prepare Mastermix as indicated below: N being the reaction number, provide the amount of Mastermix for at least N+3 reactions

(Refer to part 1a or 1b from the previous table according to condition).

*Condition (a): Sample extraction WITHOUT adding CI-ARN*

Number of reactions	1	N+3**
Enzymes	15 µl	(N+3) x 15 µl
Oligomix	5 µl	(N+3) x 5 µl
CI-ARN	1 µl	(N+3) x 1 µl
Total volume Mastermix	21 µl*	(N+3) x 21 µl

\* The difference in reaction volume between condition (a) and condition (b) is negligible and has no effect on the performance.

*Condition (b): Sample extraction WITH addition of CI-ARN*

Number of reactions	1	N+3**
Enzymes	15 µl	(N+3) x 15 µl
Oligomix	5 µl	(N+3) x 5 µl
Total volume Mastermix	20 µl*	(N+3) x 20 µl

\*\* For use on the T-COR 8®-IVD, we recommend preparing a Master Mix for N+1 reactions for a single sample utilization.

- 3) Homogenize the Mastermix prepared in 2) and briefly centrifuge.
- 11) Distribute 20 µl of Mastermix\* with a micropipette and filter tips into each real time PCR tube/microplate well.

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---

12) Add 5 µl extracted RNA samples.

13) In parallel, prepare the following controls:

- Positive Control : 20 µl Mastermix + 5 µl Positive Control CP-ChikDeng.
- Negative Control :
  - Condition (a) Sample extraction WITHOUT adding CI-ARN:
    - 20 µl Mastermix + 5 µl water provided (CN-H2O)
  - Condition (b) Sample extraction WITH addition of CI-ARN
    - 20 µl Mastermix + 4µl water provided (CN-H2O) + 1 µl CI-ARN

14) Close immediately with an adhesive foil or transparent caps in order to prevent any contamination.

15) Briefly centrifuge to collect the reaction mixture at the bottom of the tubes or microplate wells.

16) Prepare the following program on the real-time PCR device (on the T-COR 8®-IVD no manual programming required due to the barcodes - see page 23).

17) Approximative duration: 1h10.

Programme	Temperature	Time	Cycle(s)	
<b>Reverse Transcription</b>	50 °C	10 min	1	-
<b>Denaturation</b>	95°C	3 min	1	-
<b>Amplification</b>	95°C	10 sec	45	-
	60°C	30 sec		Fluorescence acquisition

**Note 1:** On the LightCycler® 480, two optical systems exist. Only the “System II” is compatible with the kit use. Apply a color compensation for the wavelengths FAM-HEX/VIC-Cy5 (465-510,533-580,618-660).

**Note 2:** For the Applied Biosystems devices, select “ROX” in “PASSIVE REFERENCE”.

**Note 3:** On the Rotorgene™, calibrate the signal by clicking on “GAIN OPTIMISATION”.


**Note 4:** On the CFX96 (Bio-Rad), launch the run from version 1.6 or later, then analyze with version 3.1 (See Section “Experiment validation”).

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
## 11 Experiment validation

Post-acquisition data analysis on the CFX96 PCR device (Bio-Rad) must be done with the CFX Manager software version 3.1 (Bio-Rad). In order to switch to this version from a run launched on a previous version, please follow the following procedure: at the end of the run, open and treat the data file with the .pcrd suffix with the CFX Manager version 3.1 (Bio-Rad).

If the run has been launched with the CFX Manager Software v1.6 for example, in order to open a data file with the CFX Manager Software v3.1, click on the CFX Manager v3.1 icon. The homepage appears.

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
- Select the file you want to analyse and click on "Open".

-The option "Drift correction" must be applied in the tab "Settings" as shown in the figure below: Click on the "Settings" tab then on "Baseline Setting" and on "Apply Fluorescence Drift Correction".

Once this step is completed, the analysis can start.

**In order for the test to be valid**, Ct values for the controls should be as mentioned in Table 4. Outside of these values, the experiment cannot be validated.

On the T-COR 8<sup>TM</sup>-IVD, for a proper use of the kit, it is necessary that these controls checks are performed when a new kit is first used, regular checks can be done later on.

**Table 4:**

Positive control	
Channel FAM	Ct ≤ 30
Channel HEX	Ct ≤ 30
Negative control	
Channel FAM	Ct ≥ 40 or not defined
Channel HEX	Ct ≥ 40 or not defined
Channel CY5	Ct ≤ 45

On the T-COR 8<sup>®</sup>-IVD, interpretation is automatically generated with the barcodes.

## 12 Data analysis and interpretation

RNA extraction and RT-PCR inhibition control in the samples:

Two results may be obtained:

1/ The RNA extraction and RT-PCR inhibition control (CI-ARN) test is positive (Ct ≤45): The RNA has been properly extracted and there are no RT-PCR inhibitors. The result can be validated.

2/ The RNA extraction and RT-PCR inhibition control (CI-ARN) test is negative (Ct ≥45): either the RNA has not been extracted, or the RT-PCR did not work well, or the presence of PCR inhibitors inhibit the PCR reaction. It is recommended to repeat the extraction or to dilute the sample, unless a specific signal appears in channel FAM or HEX.

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For clinical samples, the following results are possible:

\*Ct threshold for sample: + FAM Positive => Ct positive < 40

+ HEX Positive => Ct positive < 40

\*\*Ct threshold for CI-ARN: + CY5 Positive => Ct positive ≤ 45

PCR signal			Presence of virus		Validity of the test/comment
FAM	HEX	CY5	CHIKUNGUNYA	DENGUE	
+	-	+	Yes	No	Valid
-	+	+	No	Yes	Valid
+	+	+	Yes	Yes	Valid
+	+	-	Yes	Yes	Valid; possible extraction or RT-PCR inhibition or high competition which does not prevent the virus detection
-	-	+	No	No	Valid
+	-	-	Yes	Uninterpretable	Invalid; extraction or RT-PCR inhibition problem or high competition – dilute the sample 5x; if same results, re-extract the sample
-	+	-	Uninterpretable	Yes	Invalid; extraction or RT-PCR inhibition problem or high competition - dilute the sample 5x; if same results, re-extract the sample
-	-	-	Uninterpretable	Uninterpretable	Invalid; extraction or RT-PCR inhibition problem - dilute the sample 5x; if same results, re-extract the sample

## 13 Analytical Performance

Detection limit/analytical sensitivity on CP:

Chikungunya and Dengue	CFX-96™ (Bio-Rad)	LightCycler®480 Instrument II (Roche)	Applied Biosystems® 7500 Real-Time PCR Systems (Applied Biosystems)	T-COR 8®-IVD (Tetracore)
LOD <sub>100%</sub> on CP	20 copies/μl	200 copies/μl	200 copies/μl	20 copies/μl

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
**On the CFX-96 (Bio-rad):**

**Detection limit on the Chikungunya International Standard** (1st World Health Organization International Standard for Chikungunya virus RNA for Nucleic acid amplification techniques (NAT)-based assays) (PEI code 11785/16):

385 IU/ml of serum or plasma, and 2,5 IU/μl eluate

**Linearity:** CHIKUNGUNYA and DENGUE: 20 to 2.107 copies/μl CP

## Signal variability in the Chikungunya (FAM) and Dengue (HEX) viruses channels

	Variation coefficient (CV %) intra-run	Variation coefficient (CV %) inter-runs	Variation coefficient (CV %) inter-lots
<b>Chikungunya</b>	2,31	3,97	1,28
<b>Dengue</b>	2,82	6,13	2,41

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## Clinical performances

**Trials have been carried out on pre-characterised positive or negative samples**

<b>EBX-009 Chikungunya</b>				
		POS	NEG	TOTAL
Pre-characterised chikungunya	POS	26	0	26
	NEG	2	74	76
TOTAL		28	74	102
<b>EBX-009 Dengue</b>				
		POS	NEG	TOTAL
Pre-characterised dengue	POS	71	2	73
	NEG	0	52	52
TOTAL		71	54	125

These tests have been performed on CFX96 (Bio-Rad).

For this study, samples extraction was performed using Qiagen columns (QIAamp Ultra-sens Virus kit Ref: 53704), or Extractor EZ1 (Qiagen) or Maelstrom TANbeads (Taiwan Adadvanced Nanotech).

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The sensitivity and specificity performances are:

<b>Diagnostic specificity</b>	<b>Chikungunya</b>	<b>&gt;99%</b>
<b>Diagnostic specificity</b>	<b>Chikungunya</b>	<b>97,4% [95% confidence interval (CI): 98,8% - 101%]</b>
<b>Diagnostic sensitivity</b>	<b>Dengue</b>	<b>97,2% [95% confidence interval (CI): 93,5% - 101%]</b>
<b>Diagnostic sensitivity</b>	<b>Dengue</b>	<b>&gt;99%</b>

No cross-reactivity has been observed between Chikungunya and Dengue, or with Zika virus. Furthermore, an *in silico* analysis has shown the primers and probes specificity for the two viruses, and no potential cross-reactivity for, by example, Usutu viruses, West Nile virus, Japanese encephalitis virus, Tick encephalitis virus, Yellow fever virus, Saint Louis encephalitis virus, Hepatitis A virus, or Hepatitis C virus.

## 14T-COR 8®-IVD real-time PCR instrument related specificities

The T-COR 8®-IVD is a real time PCR device with 8 independently programmable wells which can be used, patient per patient, independently from each other, in terms of performed tests, thermal programs and PCR launch time.

The 8 wells may also be used simultaneously with the same 'Dosage' (Assay).

Note: Mix well after depositing the sample in the tube, by pipetting up and down, avoiding bubbles, and check that the liquid volume is well located at the bottom of the tube.

### **Controls**

On the T-COR 8™-IVD, for a proper use of the kit, it is necessary that the controls are performed only when a new kit is first used.

Later on, the CI-ARN ensures that within the tested sample, there is no PCR inhibition and that the PCR is working properly.

### **Patients tests**

<b><i>With a positive and negative controls validation, only once when used for the first time</i></b>	<b>24 tests</b>
<b>Number of patients tests possible patient per patient with a maximum of 12 freeze/thaw cycles (3 times per tube/4 tubes per reagents)</b>	<b>12 patients</b>

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---

## **Fluorophores**

There are 4 fluorophores combinations available on the T-COR 8®-IVD.

For all EBX, including EBX-009, the FAM/HEX/Texas Red/Cy5 combination is selected. If some of the channels are not used in some kits, ignore them in the results analysis.

## **Using the barcodes (available in pages 27 to 29)**

1- Select 'Menu' > 'Nouvelle Analyse' (New Run)

2- Select 'Code-barres' (Barcode)

3- Scan on the right of the device the desired barcode:

- Either for the Positive Control (Barcode EBX-009 Pos Ctrl),
- Or for the Negative Control (Barcode EBX-009 Neg Ctrl),
- Or for a sample (Barcode Assay: EBX-009)

*The device will automatically select one out of the 8 available wells. Follow the instructions given by the device.*

4- Lift the flap of the well x selected by the device, and check that the blue LED light is on, then select 'Oui' (Yes).

5- Place the tube in the corresponding well and fold down the flap, then select 'Suivant' (Next).

6- If you wish to name it (optional), select 'échantillon x' (sample x), and name it. Select 'Accepter' (Accept).

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8- In order to add/test another control or sample, select 'Ajouter un puits' (Add well) and start again at point 3-

9- Select 'Démarrer l'analyse' (Start the run)

*Note: Naming of a sample or a run cannot be modified or added after the end of the run. It has to be done before, or during the run.*

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---

## **Presentation and Interpretation of the results**

Ct results and amplification curves are available in the Table 'Valeurs SmartCT™' (Ct values attributed for each channel), as well as the Ct graphs according to the cycles in each channels, both being viewable in real-time on the instrument.

An analysis for the negative and positive controls, and for the samples as well, is automatically generated. It is available in 'Interprétations', at the end of the run, in the window 'Vue'.

### **For the positive and the negative controls, the following results are possible:**

"Neg Ctrl Fail" : Invalid

"Neg Ctrl Valid" : Valid

"Pos Ctrl Fail" : Invalid

"Pos Ctrl Valid" : Valid

### **Sample status determination:**

'DéTECTÉ' (Detected): Positive for at least one target \* → **framed in green**

'Non DÉTECTÉ' (Not detected): Negative → **framed in red**

\* Detected or non-detected bacteria are specified.

### **CAUTION!**

**When the frame is green, it is important to check the status for each target, as some of the targets may be negative.**

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Example of an automatic interpretation presentation of the results on the instrument and in the report:

- For the valid negative and positive controls stated as “Detected” (Détecté), for a CP “Invalid” (Invalide) to be retested, for three negative samples for the two targets, stated “Not Detected” (Non détecté) (Samples 1 to 3) and for two samples detected for one of the targets stated as:
  - “Detected Dengue” (Dengue détecté)/ “Chikungunya Not detected” (Chikungunya Non détecté) for the sample Dengue and
  - “Detected Chikungunya” (Chikungunya Détecté)/ “Dengue Not detected” (Dengue Non détecté) for the sample Chikungunya

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



---

## 15 Quality control

According to the quality management system from Eurobio, which is certified ISO EN 13485, each batch of EurobioPlex Chikungunya /Dengue is tested according to defined specifications, in order to guarantee a constant quality of the products.

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## 17 Waste disposal

All the waste should be disposed of according to DASRI regulations.

## 18 Incident report

Any serious incident that has occurred in relation with the product shall be reported to EUROBIO SCIENTIFIC and the competent authority of the country in which the user and/or the patient is established.

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---

## 19 Technical assistance

For assistance with our products, please contact our technical support.

You can contact EUROBIO SCIENTIFIC Customer Service:

- by mail at: [adv@eurobio-scientific.com](mailto:adv@eurobio-scientific.com)

or

- by phone at: +33 (0)1.69.07.94.77.

7, avenue de Scandinavie

ZA de Courtabœuf

91940 Les Ulis

FRANCE