EurobioPlex HEV

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7 avenue de Scandinavie
ZA de Courtaboeuf
91940 Les Ulis
FRANCE

Instructions for use EurobioPlex HEV


EurobioPlex HEV

For **qualitative/quantitative** real time RT-PCR

REF **EBX-010-25 / EBX-010-50 / EBX-010-100**

 **25/50/100 reactions**

EBX-010 IFU - Version 8.03 of 10/01/2024

Validated on :

CFX96™ Real Time PCR detection system (Bio-Rad) with analysis on CFX Manager version 3.1 (Bio-Rad)

Instructions for use
Available on sur www.eurobio-scientific.com

The Summary of Safety and Performance Characteristics (SPSC) will be made available by the Notified Body on EUDAMED once it is operational. It can also be obtained on request.

Instructions for use EurobioPlex HEV

Table of contents

1. Introduction.....	23
2. Purpose of the system	31
3. Symbols	25
4. Principle of detection.....	26
5. Content of the kit.....	27
6. Storage.....	27
7. Materials required not provided.....	27
8. Real-time PCR instrument.....	28
9. Cautions and note.....	28
10. Procedure.....	29
11. Validation of the experiment	33
12. Data analysis and interpretation.....	35
13. Performance analysis.....	36
14. Bibliography	53
15. Quality control	53
16. Waste disposal.....	54
17. Incident report.....	41
18. Technical assistance	42

Instructions for use EurobioPlex HEV

1. Introduction

Hepatitis E Virus (HEV) belongs to the family *Hepesviridae*, genus *HepEvirus*. It is the single hepatotropic virus with an animal reservoir known to date. It is a non-enveloped single-stranded RNA virus which is approximatively 7.5 kb long. It is considered to be a major causal agent of transmitted viral hepatitis, very little by vertical transmission or by contaminated blood transfusion but more specifically by fecal-oral route through consumption of contaminated food. It is recognized as an infectious agent responsible for global public health problem.

HEV is responsible for approximately 50% of acute viral hepatitis in developing countries with low level of hygiene, in Asia, Africa and South America causing endemic-epidemic infections. In these regions the seroprevalence can reach up to 60%. In contrast, the virus exists rather as a zoonotic agent responsible for indigenous sporadic or imported infections in developed countries with high level of hygiene. In some geographical areas of these countries, HEV infection is actually endemic since the use of far more sensitive serodiagnostic techniques.

HEV infection is usually asymptomatic and benign. It affects people of all ages, but it has two still unresolved clinical features: infected people are more often men than women, and the mortality rate is particularly high (up to 30-40%) for pregnant women especially in case of fulminant form. Chronic infection can also be observed in particular in immunocompromised people, transplanted with solid organs or suffering from hematological malignancies.

On the basis of the genomic sequences of many human and animal isolates, the HEV is classified into four major genotypes divided into many subtypes. Genotypes 1 and 2 mainly infecting humans are the leading cause of HEV endemic-epidemic infections in Asia, Africa and Mexico. With relatively well conserved nucleotide sequences, they are respectively subdivided into 5 (1a-1e) and 2 (2a-2b) sub-types. Genotypes 3 and 4 are primarily detected in pigs. Genotype 3 circulates mainly in the United States but also in some European countries and Japan. Genotype 4 circulates primarily in Asia: China, Japan and Taiwan. Genotypes 3-4 are also detected in human acute viral hepatitis and are respectively subdivided into 10 (3a-3j) and 7 (4a-4g) subtypes. Their zoonotic transmission has been clearly established by comparison of genomic sequences (e.g. only 0,2 - 2% difference between human and animal genotype 3 isolates) and also by assessing the anti-HEV seroprevalence between hog producers and the general population, going respectively from 11 to 51% against 2 to 25 %.

Diagnosis strategy of HEV infection is based on its detection from stool samples or from serum or plasma by the RT-PCR technique and by the serological markers (IgM and IgG), present early in the disease, which is characterized by jaundice and asthenia. IgM appears early and may persist beyond three months while IgG which follow IgM very quickly, persist for several years after their appearance.

The quantitative RT-PCR technique is a highly sensitive and specific technic for HEV detection.

Instructions for use EurobioPlex HEV

2. Purpose of the system

Eurobioplex HEV is a real-time polymerase chain reaction (PCR) assay designed for the qualitative and quantitative detection of the presence or absence of this virus in an RNA nucleic acid extract. The test is indicated for the diagnosis of presumptive HEV infection in humans, or to supplement a proven or undetermined diagnosis. The test is not intended for screening blood or organs for HEV. The RNA extract is the starting material for the Eurobioplex HEV kit.

It is up to the user to use extraction methods adapted to the tested samples.

The kit has been tested on plasma.

The EurobioPlex HEV test is an in vitro diagnostic medical device and must be used by qualified medical laboratory personnel. The device must not be recycled after use.

Instructions for use EurobioPlex HEV

3. Symbols

Catalogue reference

Batch code

Upper and lower storage temperature limits

Expiration date

Content sufficient for « N » tests

Manufacturer

Manufacturing date

CE marked product

In Vitro Diagnostic Medical Device

Instructions for use

Caution

Do not use if package is damaged

Keep away from light

Instructions for use EurobioPlex HEV

4. Principle of detection

The Eurobioplex HEV kit is a test using reverse-transcription and real-time amplification of viral RNA of HEV. The test is performed from extracted RNA from a sample using a reaction in a single tube/well.

The RNA extraction and RT-PCR inhibition control allows to check for variations that may occur during the RNA extraction step of biological samples and real-time PCR amplification. Thus, it ensures that a negative result may be due to a bad RNA extraction, and/or due to the presence of PCR inhibitors in too large quantities.

RNA of HEV is detected using a HEX labeled probe. The endogenous control is detected using a CY5 labeled probe. Probes emit a specific fluorescence following their hydrolysis during the elongation of the amplification product. The measurement of the intensity of real-time fluorescence correlates with the accumulation of amplification products.

Fluorescence is emitted and individually recorded through optical measurements during the PCR. The detection of the amplified fragment is performed by a fluorimeter using the channels shown in the Table 1.

Table 1: Detection of targets by fluorophores

Target	Fluorophore	Excitation	Emission
HEV	HEX	535 nm	555 nm
Endogenous gene	Cy5	647 nm	667 nm

Equivalent channels on different real time PCR cyclers:

- Channel **HEX** (Chromo 4/CFX96, Mx Systems), Channel VIC (ABI Systems), Channel Alexa532 (SmartCycler II), Channel 580 (LC 480), Channel Yellow (RotorGene),
- Channel **Cy5** (CFX96/Chromo4, Systems ABI, Systems Mx), Channel Alexa647 (SmartCycler II), Channel 660 (LC 480), Channel Red (RotorGene)

Instructions for use EurobioPlex HEV

5. Content of the kit

The EurobioPlex HEV real-time RT-PCR kit is ready to use and contains the necessary reagents and enzymes for the detection and quantification of this virus (Table 2).

Table 2: Content of the kit

Cap color	Content of the kit	Volume for 25 reactions	Volume for 50 reactions	Volume for 100 reactions	Reconstitution
Red	Enzymes	420 µL	2 x 420 µL	4 x 420 µL	Ready to use
Transparent	Oligomix *	140 µL	2 x 140 µL	4 x 140 µL	Ready to use
Yellow	4 Standards HEV **: S-1 : 10 ⁶ copies/µl S-2 : 10 ⁵ copies/µl S-3 : 10 ⁴ copies/µl S-4 : 10 ³ copies/µl	50 µL of each	2 x 50 µL of each	4 x 50 µL of each	Ready to use
Orange	Positive Control CP	40 µL	80 µL	160 µL	Ready to use
Blue	Water = negative control (CN- H2O)	1000 µL	1000 µL	1000 µL	Ready to use

*Oligomix: contains primers and probes for HEV and for the endogenous control.

** The EBX-010 kit contains quantification standards in 4 different concentrations (see chapter 10.3 Quantitative RT-PCR).

6. Storage

All reagents should be stored between -15°C and -22°C.

All reagents can be used until the expiration date indicated on the kit label.

The sensitivity of the assay may decrease if the kit components undergo multiple freeze/thaw cycles. The kit can be used after the initial opening for up to 5 freeze/thaw cycles.

7. Materials required not provided

- Biological cabinet
- Real time PCR instrument
- Centrifuge for microtubes
- Vortex
- Plates / tubes for real-time PCR reaction
- Micropipettes
- DNase-free and RNase-free filter tips for micropipettes
- Sterile microtubes
- Gloves (talc-free)

Instructions for use EurobioPlex HEV

8. Real-time PCR instrument

The EurobioPlex HEV kit has been developed and validated for use with the following real-time PCR instrument:

CFX96™ Real Time PCR detection system (Bio-Rad) with analysis on CFX Manager version 3.1 (Bio-Rad)

9. Cautions and note

Read these instructions carefully before starting the procedure.

- The experiment must be performed by competent staff, trained to technical and safety techniques.
- Instruments must have been properly installed, calibrated and maintained according to the manufacturer's recommendations.
- Clinical samples are potentially infectious and must be processed under a laminar flow hood.
- It is the user's responsibility if he/she uses other non-validated equipment, and if so, the performances are not guaranteed.
- The experiment must be performed according to good laboratory practices.
- Do not use this kit after expiration date.
- The kit is shipped with dry ice, and the components of the kits must arrive frozen. If one or more components are defrosted, or if the tubes have been damaged, contact Eurobio-Scientific.
- Avoid freezing/thawing reagents, as this can lead to a reduction in test sensitivity.
- Once defrosted, spin down briefly the tubes before use.
- It is recommended to define three working areas: 1) Isolation of RNA, 2) Preparation of the reaction mix and 3) Amplification / Detection of amplified products.
- Oligomix containing oligonucleotides and especially dye-labelled oligonucleotides should be protected from light for optimal stability. Material's exposure to light is valid during the time required for the PCR setup.
- Positive controls should be opened and handled separately from biological samples to be tested and other kit components, to avoid cross-contamination.
- Use specific lab coat and gloves (powder free) in each working area.
- Pipettes, reagents and other materials must not cross each area.
- Specific caution is required to preserve the purity of the reagents and reaction mixtures.
- The internal endogenous control detects a cellular target detectable in all human samples but not in the CN-H2O negative control provided in the kit. The absence of a signal in the negative control prevents cross-contamination from occurring.
- Appropriate methods of preparation/extraction of RNA to produce high quality RNA and to be followed by a real-time RT-PCR application should be used, particularly avoiding all sources of RNase contamination.
- Always use RNase-free DNase-free filtered tips for micropipettes.
- Do not pipette with mouth and do not eat, drink or smoke in the area.
- Avoid sprays.
- The kit is not intended for single use. When reusing the device, follow the recommendations concerning authorized freeze-thaw cycles, and precautions to prevent contamination of reagents.
- Kit components must not be used separately (neither with other reagents, nor with reagents from other batches).

Instructions for use EurobioPlex HEV

- Reagents must be thawed carefully so as not to impair device performance (at +2°C/+8°C or room temperature).
- The device is not automated. It is the user's responsibility to use equipment other than that validated, in which case performance is not guaranteed.

10. Procedure

10.1 Samples collection

- Collect samples in sterile tubes.
- It is the responsibility of the user to control their own sample collection, transport, storage and extraction conditions to ensure that RNA/DNA extraction using appropriate systems produces RNA/DNA of high quality.
- It is recommended that samples are extracted immediately, or stored according to the sample storage recommendations before extraction (Table 3).
- The bibliographical references in section "6. Bibliography" provide indicative data on sample and RNA stability.

Table 3 : Storage recommendations before use

Recommendations for storage at room temperature before preparing EDTA plasma from blood	
< 24 h	
Storage recommendations at +4°C before preparing EDTA plasma from blood	
< 72 h	
Recommendations for maximum storage of samples before extraction	
Ambient temperature	2 h
+2°C/+8°C	5 days
<-70°C (preferred over -20°C)	Long term storage (>5 days and maximum 2 months at -20°C)

Caution	
	- The user may refer to the recommendations issued by the World Health Organization or the French National Authority for Health for the proper storage of samples. - Extracted RNAs should be stored at <-70°C to ensure stability. After one year's storage of RNAs at <-70°C, the Ct values obtained may increase. It is advisable to re-extract a biological sample that has been stored for more than one year. It is advisable to limit the number of times RNAs are thawed to 3. - The transport of clinical samples is subject to local regulations governing the transport of infectious agents.

Instructions for use EurobioPlex HEV

10.2 RNA/DNA extraction

It is the responsibility of the user to ensure that the nucleic acid extraction system used is compatible with real-time RT-PCR technology. We recommend using RNA/DNA extraction methods suitable to the samples used, and to refer to the instructions of the supplier of the extraction kit used.

In the EBX-010 kit, the endogenous internal control on the Cy5 channel is already present in the clinical sample. Properly detected after amplification, the endogenous internal control is used to ensure the quality of the sample and the extraction.

The performance described in this document was obtained with the following extraction technique: Maelstrom 9600 (TanBEAD) and the OptiPure Viral Auto Plate kit (ref: W665A46).

10.3 Real-time RT-PCR

General comment:

- HEV Standards and the positive control contain a very high concentration of matrix. Manipulations must be performed carefully to avoid contamination.
- To control the functioning of the RT-PCR and steps of extraction and real-time RT-PCR amplification, it is necessary to test a positive control, CP, as well as the negative control (water PCR provided = CN-H2O) (see II-3/6) of real-time RT-PCR procedure for details).

Qualitative RT-PCR

No standard range is performed. Only the CP positive control and the negative control are tested.

Quantitative RT-PCR

In addition to the negative and positive controls (CP), a standard curve must be produced using the four ready-to-use HEV standard tubes S-1 to S-4. Standards do not need to be extracted.

In order to correct for possible experimental variation, we recommend testing standards HEV in triplicates.

To generate a standard curve on a real-time PCR System, the concentrations of 10^6 , 10^5 , 10^4 and 10^3 copies/ μ l must be used and defined as standard HEV by specifying the corresponding number of copies/ μ l (table 4a).

Table 4a : Standards concentrations HEV

Standards	Concentration
S-1	10^6 copies/ μ l
S-2	10^5 copies/ μ l
S-3	10^4 copies/ μ l
S-4	10^3 copies/ μ l

Concentrations in copies/ μ l can be converted to IU (International Units), where 2.25 copies/mL correspond to 1 IU/mL, or **0.00225 copies/ μ l** correspond to 1 IU/mL (table 4b).

Instructions for use EurobioPlex HEV

Table 4b: Calculating HEV standards

HEV standard value in cp/ μ L (table 4a)	$\div$	0.00225 (cp/ μ L) =	HEV standard value in UI/mL
--	--------	-------------------------	--

Instructions for use EurobioPlex HEV

Diagram of the procedure:

1 - PREPARATION OF REACTION MIX / MASTERMIX

Number of reaction(s)	N+3
Enzymes	(N+3) x 15 µl
Oligomix	(N+3) x 5 µl
Total Volume Mastermix	(N+3) x 20 µl

2 - PREPARATION OF REACTIONS

Sample

20 µl Mastermix

+

5µl RNA sample

Positive Control(s)

Qualitative :

20 µl Mastermix

+

5µl CP

Quantitative :

20 µl Mastermix

+

5µl S-1 to S-4

Negative Control

20 µl Mastermix

+

5µl water supplied (CN-H2O)

3 - REAL TIME PCR INSTRUMENT

Program	Temperature	Duration	Cycle(s)	
Reverse Transcription	50 °C	10 min	1	-
Denaturation	95°C	3 min	1	-
Amplification	95°C	10 sec	45	-
	60°C	30 sec		Fluorescence Acquisition

Instructions for use EurobioPlex HEV

10.4 Detailed procedure

- 1) Homogenize the tube of Enzymes, and vortex Oligomix, CP, and HEV standards tubes before starting, and centrifuge.
- 2) Prepare a Mastermix as below. N is the number of PCR reactions. Plan to prepare enough reagents for at least N+3 reactions.

Number of reaction(s)	1	N+3
Enzymes	15 µl	(N+3) x 15 µl
Oligomix	5 µl	(N+3) x 5 µl
Total Volume Mastermix	20 µl	(N+3) x 20 µl

- 3) Homogenize the Mastermix prepared in 2) and centrifuge briefly.
- 4) Distribute 20 µl of the Mastermix with a micropipette and filtered tips in each tube or well of microplate for real time PCR.
- 5) Add 5 µl of extracted RNA sample.
- 6) In parallel test:
 - Positive control(s):
 - o Qualitative Test: 20µl Mastermix + 5µl CP
 - o Quantitative Test: 20µl Mastermix + 5µl S-1 to S-4 to generate HEV standard range
 - Negative control:
 - 20 µl Mastermix + 5 µl water supplied (CN-H2O)
- 7) Close immediately with an adhesive film or transparent caps to avoid all contamination.
- 8) Centrifuge briefly to collect all the PCR reaction mix at the bottom of the tubes or plate.
- 9) Program the real-time PCR instrument as follows.

Program	Temperature	Duration	Cycle(s)	
Reverse Transcription	50 °C	10 min	1	-
Denaturation	95°C	3 min	1	-
Amplification	95°C	10 sec	45	-
	60°C	30 sec		Fluorescence acquisition


Note 4: On CFX96 (Biorad), start the run using the v1.6 or later version of CFX Manager software, and analyse with v 3.1 (see § 11. Validation of the Experiment)

11 Validation of the experiment


Post-acquisition data analysis on a CFX96 PCR instrument (Bio-Rad) must be performed using version 3.1 of the CFX Manager software (Bio-Rad). In order to change to this version from a run started on an earlier version, please follow the procedure below: At the end of the run, the data file with the extension .pcrd must be opened and processed with CFX Manager (Bio-Rad) version 3.1.

Instructions for use EurobioPlex HEV

If the run was started with CFX Manager v1.6 for example, to open a data file with CFX Manager v3.1, click on the CFX Manager v3.1 icon. The home screen appears.


Click on "File" and select "Open" then "Data File".

Select the file you need to analyse and click on "Open".

The "Drift correction" option must be applied in the "Settings" tab as shown in the image below: click on the "Settings" tab, then on "Baseline Setting" and on "Apply Fluorescence Drift Correction".

Instructions for use EurobioPlex HEV

Once this is done, the analysis can start.

For the assay to be valid, the Ct values for the controls must be the following (Table 5). Otherwise, the experiment is not valid.

Table 5:

Positive control	
HEX	Ct ≤ 30
CY5	Ct ≤ 30
Negative control	
HEX	Undetermined Ct
CY5	Undetermined Ct

12 Data analysis and interpretation

A) Qualitative analysis

For clinical samples, the following results are possible:

Ct threshold for sample positivity (HEX channel): + Positive => Ct < 45

PCR Signal		Presence of HEV virus	Test validity/comment
HEX	CY5		
-	+	No	VALID
+	+	Yes	VALID
+	-	Yes	VALID Inhibition of RT-PCR or problem with extraction or strong competition that does not prevent the detection of virus.
-	-	No possible interpretation	Inhibition of RT-PCR or problem with extraction → dilute 5 x the sample; otherwise extract once again the sample

B) Quantitative analysis

Define the concentrations of the standards HEV (channel HEX) (10^6 - 10^5 - 10^4 - 10^3 copies/μl for S-1, S-2, S-3 and S-4, respectively) during the setting of the thermal cycler. The software will draw the standard curve and will calculate directly the corresponding number of copies/μl for the samples.

Instructions for use EurobioPlex HEV

Depending on the volume of corporal fluid/sample extracted, and the elution volume, the concentration of HEV RNA in the eluate can be translated in the original specimen sample as the concentration/ml of specimen.

To determine the viral load per ml original specimen, the following formula can be applied.

$$\frac{\text{Viral load per microlitre eluate} \times \text{Volume eluate (microliters)}}{\text{Volume of extracted sample (milliliters)}} = \text{Viral load per milliliter sample}$$

It can be converted to IU (International Units) knowing that 2.25 copies/ml correspond to 1IU/ml.

The results below quantification threshold cannot be quantitatively validated, but retain their quality of detectability (positivity validated by the HEX channel).

Moreover, for quantitative PCR, in addition to the valid controls of Table 5, efficiency should be between 90 and 110% on the HEX channel, and the coefficient of correlation > 0.98.

Limitations on use and interpretation:

- ❖ The EurobioPlex HEV kit is used for first-line diagnostic purposes.
- ❖ All samples should be treated as potentially infected with bacteria of interest, and local biosafety regulations should be carefully followed.
- ❖ Interpretation of results should consider the possibility of false negatives and false positives.

False negatives may be due to:

- Inadequate collection of samples, or incorrect storage,
- Inappropriate extraction conditions or use of non-validated PCR instruments,
- Experimentation not respecting all elements of these instructions for use.

False positives may be due to:

- Contamination due to mishandling of high positive samples, the positive control, or the PCR amplification products,
- Failure to follow the procedure described in these instructions for use, especially to avoid sources of contamination.
- ❖ All results should be interpreted by medical professionals in the context of the patient's clinical situation, history and symptoms.
- ❖ This assay does not exclude the presence of pathogens other than those targeted in the kit.
- ❖ A negative result of this assay does not absolutely exclude a possible infection with the virus targeted in this kit.

13. Performance analysis

Analytical sensitivity

The analytical sensitivity of the EurobioPlex EBX-010 HEV kit was determined on the positive control of the kit. A dilution range was performed using a plasmid mixture ranging from 10⁵ to 1 copy/μL.

Instructions for use EurobioPlex HEV

A dilution range was performed using an eluate of the WHO HEV International Standard (Paul Ehrlich Institute PEI code: 6329/10) ranging from 2083 IU/mL to 16.2 IU/mL. This dilution range was used to determine the detection limit of the EBX-010 kit. A probit assay was performed to determine the 95% positivity rate.

Detection limitation/analytical sensitivity on CP: 5 copies/ μ L

Limit of quantification on international WHO IS (1st world health organization international standard for Hepatitis E RNA Nucleic acid amplification (NAT) assay, Paul Ehrlich Institute PEI code : 6329/10) at 95%:

HEX : 63,9 IU/mL

Limit of quantification


A dilution range was performed using an eluate of the WHO HEV International Standard (Paul Ehrlich Institute PEI code: 6329/10) ranging from 2083 IU/ml to 16,2 IU/mL. This dilution range was used to determine the limit of quantification of the EBX-010 kit.

Limit of quantification on international WHO IS (1st world health organization international standard for Hepatitis E RNA Nucleic acid amplification (NAT) assay, Paul Ehrlich Institute PEI code : 6329/10) :

HEX : 260 IU/mL

Linearity range

The linear measurement range of the EurobioPlex HEV kit was determined by analyzing serial dilutions of a sample titrated to 3.54 log IU/mL. Each dilution was tested in duplicate from 1.13 to 3.54 log IU/mL.

Instructions for use EurobioPlex HEV

Precision

The accuracy of the kit was established by comparing the mean of the measured values of the International WHO HEV standard (Paul Ehrlich Institute PEI code: 6329/10) to its known true value. This method was also used on a clinical sample with a known viral load.

True value of HEV standard = 3.02 log IU/mL

Bias (%) = $((3.02-3.00) / 3.02) * 100 = 0.62\%$

True value of titrated sample = 4.92 log IU/ml

Bias (%) = $((4.92-4.78) / 4.92) * 100 = 2.74\%$

Accuracy

- Intra-experimental variability:

	EBX-010	
	CY 5	HEX
Mean CV intralot CP	0.77	0.70

CV: coefficient of variation

- Inter-batch variability:

	EBX-010	
	CY 5	HEX
Mean CV inter lots CP in %	3.27	2.76

CV: coefficient of variation

Analytical specificity

The validation was performed on 10 samples negative for HEV but positive for other liver viruses, to test for cross-reactivity.

Number of samples	Sample positive for	Hepatitis E	HEV status on EBX-010
3	Hepatitis C	Negative	Negative
3	Hepatitis B	Negative	Negative
2	Hepatitis D	Negative	Negative
2	Hepatitis A	Negative	Negative

Instructions for use EurobioPlex HEV

As a reinforcement to this study, an in silico analysis was used which demonstrates the specificity of the primers and probes for HEV in the EBX-010 kit.

No aspecificity or cross reaction was detected.

The kit is 100% specific.

Diagnostic sensitivity and specificity

This validation was carried out on 66 HEV positive samples and 51 HEV negative samples.

For this study, RNA extraction was performed on Maelstrom 9600 (TANBead) with kit OptiPure Viral Auto Plate (réf : W665A46) kit, and the RT-PCR on CFX96 (Bio-Rad).

		EBX-010	
		Positive HEV	Negative HEV
Pre-tested with Altona kit	Positive HEV	66	0
	Negative HEV	0	51

Overall sensitivity: > 99% (66/66)

Overall specificity : > 99% (51/51)

Overall Concordance : > 99% (117/117)

Interfering substances (including various PCR inhibitors) were tested to ensure that they would not distort the diagnostic result obtained with the kit. A panel of AccroMetrix™ Inhibition Panel REF 956400 containing 7 substances were tested at concentrations known to be the highest possible.

All samples at each concentration were HEV positive, in HEX, in the presence of these interfering substances. There is therefore no impact of these substances on the EBX-010 kit.

Cross-contamination

The cross-contamination study on 28 CP and 28 NTC revealed sensitivity, specificity and overall concordance > 99% with intra- and inter-plate coefficients of variation < 5% on CP, for all targets.

14. Bibliography

Meng X.-J., Purcell R. H., Halbur P. G., Lehman J. R., Webb D. M., Tsareva T. T., Haynes J. S., Thacker B. J., and Emerson S. U. A novel virus in swine is closely related to the human hepatitis E virus . Proc Natl Acad Sci 1997; 94: 9860-9865

Instructions for use EurobioPlex HEV

Wang Y. , Zhang H. , Ling R. , Li H. , and Harrisson T. J. The complete sequence of hepatitis E virus genotype 4 reveals an alternative strategy for translation of open reading frames 2 and 3 . *Journal of General Virology* 2000; 81: 1675-1666

Drobeniuc J. , Favorov M. O. , Shapiro C. N. , Bell B. P. , Mast E. E. , Dadu A. , Culver D. , Iarovoi P. , Robertson B H and Margolis H S. Hepatitis E virus antibody prevalence among persons who work with swine. *The Journal of Infectious Diseases* 2001; 184: 1594-1597

Meng X-J. , Wiseman B. , Elvinger F. , Guenette D. K. Toth T. E. , Engle R. E. , Emerson S. U. , and Purcell R. H. Prevalence of antibodies to hepatitis E virus in veterinarians working with swine and in normal blood donors in the United States and other countries . *Journal of clinical microbiology* 2002; 40 (1): 117-122

Orru G , Masia G , Orru L , Piras V , Coppola R C . Detection and quantitation of hepatitis E virus in human faeces by real-time quantitative PCR. *Journal of Virological Methods* 2004; 118: 77-82

Olsen B., Axelsson-Olsson D., Thelin A. and Weiland O. Unexpected high prevalence of IgG-antibodies to hepatitis E virus in Swedish pig farmers and controls . *Scand J Infect Dis* 2006; 38 (1): 55-58

Bouwknegt M. , Engel B. , Herremans M. M. P. T. , Widdowson M. A., Worm H. C. , Koopmans M. P. G. , Frankena K. , De Roda Husman A. M. , De Jong M. C. M. , and van der Poel W. H. M. Bayesian estimation of hepatitis E virus seroprevalence for populations with different exposure levels to swine in The Netherlands . *Epidemiol Infect* 2008; 136: 567-576

Kuniholm M. H. , Purcell R. H. , McQuillan G. M. , Engle R. E. , Wasley A. , and Nelson K. E. Epidemiology of hepatitis E virus in the United States: Results from the third National Health and nutrition examination survey, 1988-1994 *Ht Journal of Infectious Diseases* 2009; 200: 48-56

Ward P , Poitras E , Leblanc D , Letellier A , Brassard J , Plante D , Houde A . Comparative analysis of different Taqman real-time RT-PCR assays for detection of swine Hepatitis E virus and integration of Feline calicivirus as internal control. *Journal of Applied Microbiology* 2009; 106: 1360-1369

Zhang W , Yang S , Ren L , Shen Q , Cui L , Fan K , Huang F , Kang Y , Shan T , Wei J , Xiu H , Lou Y , Liu J , Yang Z , Zhu J , Hua X . Hepatitis E virus infection in Central China reveals no evidence of cross-species transmission between human and swine in this area. *PLoS one* 2009 ,4(12) eb156: 1-8

Meng X. J. Hepatitis E virus : Animal reservoirs and zoonotic risk. *Vet Microbiol* 2010; 140 (3-4): 256-265

Legrand-Abravanel F. , Kamar N. , Sandres-Saune K. , Garrouste C. , Dubois M. , Mansuy J-M. , Muscari F. , Sallusto F. , Rostaing L. , and Izopet J. Characteristics of autochthonous hepatitis E virus infection in solid-organ transplant recipients in France . *The Journal of Infectious Diseases* 2010; 202 (6): 835-844

Purcell R. H. and Emerson S. U. Hidden danger :The raw facts about hepatitis E virus. *The Journal of Infectious Diseases* 2010; 202 (6): 819-821

Baylis S A , Hanschmann K-M , Blümel J , and Micha Nübling C et al . Standardization of hepatitis E virus (HEV) nucleic acid amplification technique-based assays : an initial study to evaluate a panel of HEV strains and investigate Laboratory performance . *J. Clin. Microbiol.* 2011, 49(4): 1234-1239

Hakze-van der Honing R. W. , van Coillie E. , Antonis A. F. G. , van der Poel W. H.M. First isolation of hepatitis E virus genotype 4 in Europe through swine surveillance in the Netherlands and Belgium. *PloS One* 2011; 6 (8) : e22673

Instructions for use EurobioPlex HEV

Mansuy J-M., Bendall R., Legrand-Abravanel F., Sauné K., Miédouge M., Ellis V., Rech H., Destruel F., Kamar N., Dalton H. R. and Izopet J. Hepatitis E virus antibodies Emerging Infectious Diseases 2011; 17 (12): 2309-2312

Faber M. S., Wenzel J. J., Jilg W., Thamm M., Höhle M., and Stark K. Hepatitis E virus seroprevalence among adults ,Germany . Emerging Infectious Diseases 2012; 18 (10): 1654-1657

Ruggeri F. M., Di Bartolo I., Ponterio E., Angeloni G., Trevisani M., Ostanello F. Zoonotic transmission of hepatitis E virus in industrialized countries. New Microbiol 2013; 36 (4): 331-344

Smith D B , Purdy M A , Simmonds P. Genetic variability and the classification of Hepatitis E Virus. J Virol. 2013 , 87(8): 4161-4169

Dalton H.R. Hepatitis E: The 'new kid on the block' or an old friend ?. Transfus Med Hemother. 2014 , 41: 6-9

Smith D B , Simmonds P , Jameel S , Emerson S U , Harisson T J , Meng X-J , Okamoto H , Van der Poel W H M and Purdy M A. Consensus proposals for classification of the family Hepeviridae Journal of General Virology 2014, 95: 2223-2232

Izopet J , Labrique A B , Basnyat B , Dalton H R , Kmush B , Heaney C D et al; Hepatitis E virus seroprevalence in three hyperendemic areas ; Nepal, Bangladesh and southwest France . J. Clin. Virology. 2015, 70: 39-42

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15. Quality control

In accordance with Eurobio's ISO EN 13485 certified quality management system, each batch of EurobioPlex HEV is tested according to predefined specifications to ensure consistent product quality.

16. Waste disposal

Eliminate all waste in accordance with the legislation on DASRI.

17. Incident report

Serious incidents related to the system must be reported to Eurobio Scientific and to the competent authority of the Member States in which the user and/or the patient is registered.

Instructions for use EurobioPlex HEV

18. Technical assistance

For assistance with our products, please contact our technical support.

Eurobio Scientific's customer service can be contacted by e-mail at adv@eurobio-scientific.com or by phone at +33 (0)1.69.79.64.80.

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Wang Y., Zhang H., Ling R., Li H., and Harrisson T. J. The complete sequence of hepatitis E virus genotype 4 reveals an alternative strategy for translation of open reading frames 2 and 3 . Journal of General Virology 2000; 81: 1675-1666

Drobeniuc J., Favorov M. O., Shapiro C. N., Bell B. P., Mast E. E., Dadu A., Culver D., Iarovoi P., Robertson B H and Margolis H S. Hepatitis E virus antibody prevalence among persons who work with swine. The Journal of Infectious Diseases 2001; 184: 1594-1597

Meng X-J., Wiseman B., Elvinger F., Guenette D. K. Toth T. E., Engle R. E., Emerson S. U., and Purcell R. H. Prevalence of antibodies to hepatitis E virus in veterinarians working with swine and

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in normal blood donors in the United States and other countries . Journal of clinical microbiology 2002; 40 (1): 117-122

Orru G , Masia G , Orru L , Piras V , Coppola R C . Detection and quantitation of hepatitis E virus in human faeces by real-time quantitative PCR. Journal of Virological Methods 2004; 118: 77-82

Olsen B., Axelsson-Olsson D., Thelin A. and Weiland O. Unexpected high prevalence of IgG-antibodies to hepatitis E virus in Swedish pig farmers and controls . Scand J Infect Dis 2006; 38 (1): 55-58

Bouwknegt M. , Engel B. , Herremans M. M. P. T. , Widdowson M. A., Worm H. C. , Koopmans M. P. G. , Frankena K. , De Roda Husman A. M. , De Jong M. C. M. , and van der Poel W. H. M. Bayesian estimation of hepatitis E virus seroprevalence for populations with different exposure levels to swine in The Netherlands . Epidemiol Infect 2008; 136: 567-576

Kuniholm M. H. , Purcell R. H. , McQuillan G. M. , Engle R. E. , Wasley A. , and Nelson K. E. Epidemiology of hepatitis E virus in the United States: Results from the third National Health and nutrition examination survey, 1988-1994 Ht Journal of Infectious Diseases 2009; 200: 48-56

Ward P , Poitras E , Leblanc D , Letellier A , Brassard J , Plante D , Houde A . Comparative analysis of different Taqman real-time RT-PCR assays for detection of swine Hepatitis E virus and integration of Feline calicivirus as internal control. Journal of Applied Microbiology 2009; 106: 1360-1369

Zhang W , Yang S , Ren L , Shen Q , Cui L , Fan K , Huang F , Kang Y , Shan T , Wei J , Xiu H , Lou Y , Liu J , Yang Z , Zhu J , Hua X . Hepatitis E virus infection in Central China reveals no evidence of cross-species transmission between human and swine in this area. PLoS one 2009 ,4(12) eb156: 1-8

Meng X. J. Hepatitis E virus : Animal reservoirs and zoonotic risk. Vet Microbiol 2010; 140 (3-4): 256-265

Legrand-Abravanel F. , Kamar N. , Sandres-Saune K. , Garrouste C. , Dubois M. , Mansuy J-M. , Muscari F. , Sallusto F. , Rostaing L. , and Izopet J. Characteristics of autochthonous hepatitis E virus infection in solid-organ transplant recipients in France . The Journal of Infectious Diseases 2010; 202 (6): 835-844

Purcell R. H. and Emerson S. U. Hidden danger :The raw facts about hepatitis E virus. The Journal of Infectious Diseases 2010; 202 (6): 819-821

Baylis S A , Hanschmann K-M , Blümel J , and Micha Nübling C et al . Standardization of hepatitis E virus (HEV) nucleic acid amplification technique-based assays : an initial study to evaluate a panel of HEV strains and investigate Laboratory performance . J. Clin. Microbiol. 2011, 49(4): 1234-1239

Hakze-van der Honing R. W. , van Coillie E. , Antonis A. F. G. , van der Poel W. H.M. First isolation of hepatitis E virus genotype 4 in Europe through swine surveillance in the Netherlands and Belgium. PLoS One 2011; 6 (8) : e22673

Mansuy J-M. , Bendall R. , Legrand-Abravanel F. , Sauné K. , Miédouge M. , Ellis V. , Rech H., Destruel F., Kamar N., Dalton H. R. and Izopet J. Hepatitis E virus antibodies Emerging Infectious Diseases 2011; 17 (12): 2309-2312

Faber M. S., Wenzel J. J., Jilg W., Thamm M., Höhle M., and Stark K. Hepatitis E virus seroprevalence among adults ,Germany . Emerging Infectious Diseases 2012; 18 (10): 1654-1657

Ruggeri F. M., Di Bartolo I., Ponterio E. , Angeloni G. , Trevisani M. , Ostanello F. Zoonotic transmission of hepatitis E virus in industrialized countries. New Microbiol 2013; 36 (4): 331-344

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Wang Y., Zhang H., Ling R., LI H., and Harrison T. J. The complete sequence of hepatitis E virus genotype 4 reveals an alternative strategy for translation of open reading frames 2 and 3. *Journal of General Virology* 2000; 81: 1675-1666

Drobeniuc J., Favorov M. O., Shapiro C. N., Bell B. P., Mast E. E., Dadu A., Culver D., Iarovoi P., Robertson B H and Margolis H S. Hepatitis E virus antibody prevalence among persons who work with swine. *The Journal of Infectious Diseases* 2001; 184: 1594-1597

Meng X-J., Wiseman B., Elvinger F., Guenette D. K. Toth T. E., Engle R. E., Emerson S. U., and Purcell R. H. Prevalence of antibodies to hepatitis E virus in veterinarians working with swine and in normal blood donors in the United States and other countries. *Journal of clinical microbiology* 2002; 40 (1): 117-122

Orru G, Masia G, Orru L, Piras V, Coppola R C. Detection and quantitation of hepatitis E virus in human faeces by real-time quantitative PCR. *Journal of Virological Methods* 2004; 118: 77-82

Olsen B., Axelsson-Olsson D., Thelin A. and Weiland O. Unexpected high prevalence of IgG-antibodies to hepatitis E virus in Swedish pig farmers and controls. *Scand J Infect Dis* 2006; 38 (1): 55-58

Bouwknegt M., Engel B., Herremans M. M. P. T., Widdowson M. A., Worm H. C., Koopmans M. P. G., Frankena K., De Roda Husman A. M., De Jong M. C. M., and van der Poel W. H. M. Bayesian estimation of hepatitis E virus seroprevalence for populations with different exposure levels to swine in The Netherlands. *Epidemiol Infect* 2008; 136: 567-576

Kuniholm M. H., Purcell R. H., McQuillan G. M., Engle R. E., Wasley A., and Nelson K. E. Epidemiology of hepatitis E virus in the United States: Results from the third National Health and nutrition examination survey, 1988-1994 *Ht Journal of Infectious Diseases* 2009; 200: 48-56

Ward P, Poitras E, Leblanc D, Letellier A, Brassard J, Plante D, Houde A. Comparative analysis of different Taqman real-time RT-PCR assays for detection of swine Hepatitis E virus and integration of Feline calicivirus as internal control. *Journal of Applied Microbiology* 2009; 106: 1360-1369

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Zhang W, Yang S, Ren L, Shen Q, Cui I, Fan K, Huang F, Kang Y, Shan T, Wei J, Xiu H, Lou Y, Liu J, Yang Z, Zhu J, Hua X. Hepatitis E virus infection in Central China reveals no evidence of cross-species transmission between human and swine in this area. PLoS one 2009 ,4(12) eb156: 1-8

Meng X. J. Hepatitis E virus: Animal reservoirs and zoonotic risk. Vet Microbiol 2010; 140 (3-4): 256-265

Legrand-Abravanel F., Kamar N., Sandres-Saune K., Garrouste C., Dubois M., Mansuy J-M., Muscari F., Sallusto F., Rostaing L., and Izopet J. Characteristics of autochthonous hepatitis E virus infection in solid-organ transplant recipients in France. The Journal of Infectious Diseases 2010; 202 (6): 835-844

Purcell R. H. and Emerson S. U. Hidden danger: The raw facts about hepatitis E virus. The Journal of Infectious Diseases 2010; 202 (6): 819-821

Baylis S A, Hanschmann K-M, Blümel J, and Micha Nübling C et al. Standardization of hepatitis E virus (HEV) nucleic acid amplification technique-based assays: an initial study to evaluate a panel of HEV strains and investigate Laboratory performance. J. Clin. Microbiol. 2011, 49(4): 1234-1239

Hakze-van der Honing R. W., van Coillie E., Antonis A. F. G., van der Poel W. H.M. First isolation of hepatitis E virus genotype 4 in Europe through swine surveillance in the Netherlands and Belgium. PLoS One 2011; 6 (8): e22673

Mansuy J-M., Bendall R., Legrand-Abravanel F., Sauné K., Miédouge M., Ellis V., Rech H., Destruel F., Kamar N., Dalton H. R. and Izopet J. Hepatitis E virus antibodies Emerging Infectious Diseases 2011; 17 (12): 2309-2312

Faber M. S., Wenzel J. J., Jilg W., Thamm M., Höhle M., and Stark K. Hepatitis E virus seroprevalence among adults, Germany. Emerging Infectious Diseases 2012; 18 (10): 1654-1657

Ruggeri F. M., Di Bartolo I., Ponterio E., Angeloni G., Trevisani M., Ostanello F. Zoonotic transmission of hepatitis E virus in industrialized countries. New Microbiol 2013; 36 (4): 331-344

Smith D B, Purdy M A, Simmonds P. Genetic variability and the classification of Hepatitis E Virus. J Virol. 2013, 87(8): 4161-4169

Dalton H.R. Hepatitis E: The 'new kid on the block' or an old friend? Transfus Med Hemother. 2014, 41: 6-9

Smith D B, Simmonds P, Jameel S, Emerson S U, Harisson T J, Meng X-J, Okamoto H, Van der Poel W H M and Purdy M A. Consensus proposals for classification of the family Hepeviridae Journal of General Virology 2014, 95: 2223-2232

Izopet J, Labrique A B, Basnyat B, Dalton H R, Kmush B, Heaney C D et al; Hepatitis E virus seroprevalence in three hyperendemic areas; Nepal, Bangladesh and southwest France. J. Clin. Virology. 2015, 70: 39-42

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7 avenue de Scandinavie
ZA de Courtaboeuf
91940 Les Ulis
FRANCE

Description of changes:

Version	Description
8.02	Change of the address according to the K-bis update
8.03	Correction of the enzyme volume, table 10.4 of the IFU in CZ and EN.