EurobioPlex

*Bartonella spp/Bartonella henselae/Bartonella
quintana*

PCR EN TEMPS REEL

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Version 5.01 2024/08/06

Validated on:

- CFX96™ Real Time PCR detection system (Biorad) with analysis on CFX Manager v 3.1 (Biorad)
- LightCycler® 480 Instrument II (Roche) with analysis on LightCycler® 480 software v1.5 (Roche)
- Applied Biosystems® 7500 Real-Time PCR Systems (Applied Biosystems) with analysis on 7500 Software v2
- T-COR 8®-IVD (Tetracore Inc.) with T-COR 8 SmartCT™ (auto v1) software

Storage conditions:

Keep all reagents between -15°C and -22°C until use and after first use

Instructions for use

TABLE OF CONTENTS

Intended use	27
Introduction	27
Principle of detection	29
Description and content of the kit	29
Storage	30
Cautions and notes	31
Samples collection, transport and storage	31
Procedure	32
I- DNA extraction	32
II- Real-time PCR procedure	32
II-1/ Diagram of the procedure	33
II-2/ Detailed procedure	34
Validation of the experiment	35
Data analysis and Interpretation	37
Performance analysis.....	38
Specificities of the Real-Time PCR T-COR 8®-IVD instrument	41
Barcodes for EBX-013 for use on T-COR8®-IVD.....	45
Bibliography	48
Waste disposal	48
Symbols	49

INTENDED USE

The EBX-013 test uses real-time polymerase chain reaction (PCR) amplification and is designed for the qualitative detection of *Bartonella* spp. (all *Bartonella*), *Bartonella henselae* and *Bartonella quintana* from a DNA nucleic acids extract. This test is indicated to confirm a diagnosis of presumption of infection in patients or complete a proven or indeterminate diagnosis by standard techniques (serology, histology and culture). Extracted DNA is the starting material for the Eurobioplex *Bartonella* spp./ *B. henselae*/ *B. quintana* kit.

The Eurobioplex EBX-013 has been validated on the following specimen:

- Serum
- Blood
- Lymph node biopsy
- Spleen biopsy
- Cardiac biopsy
- Bone marrow extract

INTRODUCTION

Bartonella are bacteria that infect many species of mammals, including humans. Chronic infection of erythrocytes is a specific adaptation of *Bartonella* to mammals. They can be also found in the skin, the bone tissue and endothelial cells. The main species of *Bartonella* described in the literature have been isolated in the 1990s. About half of them are implicated as agents of zoonoses. The major reservoir is then domestic and wild animals living close to humans. For example in France, it is mainly cats. Human diseases caused by *Bartonella* spp. are called the Bartonellosis. The two main *Bartonella* species are *Bartonella henselae* responsible for “cat scratch disease”, and *B. Quintana* responsible for “trench fever” and endocarditis. Other less common species of *Bartonella* are likely to cause other human pathologies: *Bartonella bacilliformis*, *B. vinsonii berkhoffi*, *B. clarridgeiae*, *B. tamiae*, *B. rochalimae*, *B. elizabethae*, *B. koehlerae*, *B. graham* and *B. alsatica*. Endocarditis are the most serious complication.

Depending on the species *Bartonella* have various modes of transmission. Humans can be infected with the animal directly by scratch, bite, or simple caresses if the skin of the hands is scratched. The infection can also be transmitted indirectly through a blood-sucking arthropod vector. Each species of *Bartonella* usually has a single natural reservoir host which is generally asymptomatic due to a long-term bacteremia. The accidental hosts however develop symptoms generally more severe.

It is difficult to propose a treatment for bacteria causing so many variable clinical manifestations. Moreover, the evolution of the disease varies according to the bacteria and immune system of the host. It can be acute, recurrent or even chronic. The treatment has to be adapted to each clinical situation. There is so far no vaccine prophylaxis.

Cat scratch disease (CSD)

The pathogen agent is *Bartonella henselae*. The CSD is characterized by infectious conjunctivitis accompanied by a strong lymph node reaction. It was described in 1950 by R. Debré. The prevalence of the CSD is 6.6 cases per 100,000 population. The peak of incidence is reported during the cold period of the year, between September and April when pets live more in contact with their master at home.

The typical symptoms are lymphadenopathy, fever (38-41°C), fatigue, headaches, splenomegaly and cutaneous manifestations such as transitional maculo-papular rash. The disease vectors are fleas with a contamination by the way of a scratch or a bite of an infected cat or through the bite of an infected flea. It is necessary to adopt a number of preventive measures for exposure to animals, especially domestic pets. The prognosis is good in immunocompetent patients, with spontaneous resolution in 90% of cases. On the other hand, antibiotic therapy is often necessary for children and immunocompromised people, as well as in cases of neurological and cardiac forms.

Trench fever

The pathogen agent is *Bartonella Quintana*. It is an intracellular and ubiquitous germ transmitted by body lice causing periodic fever with headaches and severe pain of the lower limbs. Fever trenches is the clinical principal form of this infection that has been essentially described during the First World War causing then epidemics with more than a million cases estimated in Europe. It reappeared during the Second World War and finally disappeared. Currently it is a re-emerging disease which primarily concerns populations exposed to precarious hygiene and the presence of body lice. It notably affects homeless persons in the industrialized countries with a high prevalence (16 to 30%).

Diagnosis

Diagnosis of these infections relies on serology, histology, blood culture and the research of bacterial DNA from different samples. *Bartonella* are small gram-negative bacilli, aerobic, polymorphic, non-acido-alcohol-resistant, often immobile, from 1 to 1,2µm long and 0,5µm in diameter. They are slow-growing, as difficult to cultivate to isolate. However, these bacteria can be grown on media with fresh blood, humid atmosphere enriched in CO₂, incubated at 35°C, but growth of bacteria on these environments remains particularly tedious and costly. Diagnosis is based first on the clinical description and then on serology (indirect immunofluorescence or immunoassay enzyme with titration of IgG or IgM antibodies). The use of the techniques of molecular biology (PCR and sequencing) allows a rapid diagnosis from tissues (lymph node, skin, heart valve,...). It also enables the diagnosis of species. The diagnosis by qPCR brings sensitivity, reproducibility and early detection of the presence of the bacteria.

PRINCIPLE OF DETECTION

The Eurobioplex EBX-013 is a test using real-time amplification of bacterial DNA of *Bartonella spp.*, *Bartonella henselae* and *Bartonella Quintana*, as well as a DNA extraction and PCR inhibition control. The test is performed from the DNA extracted from the sample in a single well/tube.

The DNA extraction and PCR inhibition control allows to checks for variations that may occur during the DNA extraction step of biological samples and real-time PCR amplification. Thus, it ensures that a negative result may not be due to a bad DNA extraction, and/or the presence of PCR inhibitors in too large quantities.

This kit is composed of one quadriplex: *Bartonella spp.* (Bspp.), *Bartonella henselae* (Bh) and *Bartonella quintana* (Bq) + control of DNA extraction and PCR inhibition (CI-PCR).

DNA of Bspp., Bh, and Bq is respectively detected using FAM, HEX, and TEXAS RED labeled probes. DNA extraction and PCR inhibition control is detected using a CY5 labeled probe. All probes emit a specific fluorescence following their hydrolysis during the elongation of the amplification product. The measurement of the intensities of real-time fluorescence correlates with the accumulation of specific amplification products.

DESCRIPTION AND CONTENT OF THE KIT

The real-time PCR EBX-013 kit is ready to use for the specific detection of *Bartonella spp. / henselae / quintana*.

Fluorescence is emitted and individually recorded through optical measurements during the PCR. The detection of the amplified fragment is performed by a fluorimeter using the channels shown in Table 1 below.

The kit contains reagents and enzyme for the amplification of DNA of these bacteria, and the DNA extraction and PCR inhibition control (see Table 2).

Table 1 :

Target	Fluorophore	Excitation	Emission
<i>Bartonella spp</i>	FAM	495 nm	515 nm
<i>Bartonella henselae</i>	HEX	535 nm	555 nm
<i>Bartonella quintana</i>	Texas Red	585 nm	605 nm
PCR inhibition control (CI-PCR)	Cy5	647 nm	667 nm

Equivalent channels on different real time PCR cyclers:

- Channel **FAM** (ABI Systems, SmartCycler II, Chromo 4/CFX96, Systems Mx, T-COR8®-IVD), Channel 510 (LC 480), Channel Green (RotorGene)
- Channel **HEX** (Chromo 4/CFX96, Mx Systems, T-COR8®-IVD), Channel VIC (ABI Systems), Channel Alexa532 (SmartCycler II), Channel 580 (LC 480), Channel Yellow (RotorGene),

- Channel **Texas Red** (Chromo 4/CFX96, Mx Systems, T-COR8®-IVD), LC Red 610 (LC480), Channel Orange (RotorGene)
 - Channel **Cy5** (ABI Systems, Chromo 4/CFX96, Mx Systems, T-COR8®-IVD), Channel 660 (LC 480), Channel Alexa647 (SmartCyclerII), Channel Red (RotorGene)
 Note: On LC480 instrument II, apply color compensation for the following channels: FAM-HEX/VIC – Texas Red-Cy5 (465-510, 533-580, 533-610, 618-660).

Table 2 :

Cap color	Components of the kit	24 reactions	48 reactions	96 reactions	Reconstitution
Red	ENZYME	4 x 100 µl	2 x 350 µl	4 x 350 µl	Ready to use
Transparent	OLIGOMIX	4 x 60 µl	2 x 210 µl	4 x 210 µl	Ready to use
Yellow	Positive Control CP-Barto	4 x 50 µl	2 x 50 µl	4 x 50 µl	Ready to use
White	PCR inhibition control (CI-PCR)	4 x 280 µl	2 x 280 µl	4 x 280 µl	Ready to use
Blue	Water = negative control (CN-H ₂ O)	1 ml	1 ml	1 ml	Ready to use

Oligomix : contains the primers and probes of the quadriplex (Bspp+Bh+Bq+CI-PCR)

Required material not provided:

- ◇ Biological Hood
- ◇ Real-time PCR instrument
- ◇ Micro centrifuge
- ◇ Vortex
- ◇ Plates / tubes for real-time PCR
- ◇ Micropipettes
- ◇ DNase-free RNase-free filter tips for micropipettes
- ◇ Sterile microtubes
- ◇ Gloves (powder free)

STORAGE

All reagents must be stored between -15 and -22°C.

All reagents can be used until the expiration date indicated on the kit.

Many freezing / defrosting cycles (> 3x) must be avoided, and could lead to decrease in sensitivity.

CAUTIONS AND NOTES

Read carefully instructions before starting.

- ◇ The experiment must be performed by competent staff.
- ◇ Instruments must have been properly installed, calibrated and maintained according to the manufacturer's recommendations.
- ◇ Clinical samples are potentially infectious and must be processed under a laminar flow hood.
- ◇ The experiment must be performed according to good laboratory practices.
- ◇ Do not use this kit after expiration date.
- ◇ The kit is shipped with dry ice, and the components of the kits must arrive frozen. If one or more components are defrosted, or if the tubes have been damaged, contact Eurobio Scientific.
- ◇ Once defrosted, spin down briefly the tubes before use.
- ◇ It is recommended to define three working areas: 1) Isolation of DNA, 2) Preparation of the reaction mix and 3) Amplification / Detection of amplified products.
- ◇ Use specific lab coat and gloves (powder free) in each working area.
- ◇ Pipettes, reagents and other materials must not cross each area.
- ◇ Specific caution is required to preserve the purity of the reagents and reaction mixtures.
- ◇ Appropriate methods of preparation/extraction of DNA to produce high quality DNA and to be followed by a real-time PCR application should be used, particularly avoiding all sources of DNase contamination.
- ◇ Always use RNase-free DNase-free filtered tips for micropipettes.
- ◇ Do not pipette with mouth and do not eat, drink or smoke in the area.
- ◇ Avoid sprays.

SAMPLES COLLECTION, TRANSPORT AND STORAGE

- ◇ Collect samples in sterile tubes. **The use of heparin as an anticoagulant is prohibited.**
- ◇ It is the responsibility of the user to master its own conditions of collection, transport and storage of samples, and extraction of DNA by suitable systems to produce DNA of good quality.
- ◇ The samples should be extracted immediately or stored following the recommendations in the table below (Table 3).

Table 3 :

Recommandations of maximum storage of samples before extraction	
Room temperature	2h
4°C	16h
-20°C	Long term storage

- ◇ The user can refer to the World Health Organization or High Health Authority for storing samples.

- ◇ Transport of clinical samples must obey local regulations for this type of infectious agents.

PROCEDURE

I- DNA Extraction

It is the user's responsibility that the extraction system used be compatible with downstream real time PCR technology. For this kit, we recommend to use extraction methods of bacterial DNA adapted to the type of sample tested (e.g. serum, blood, biopsy or bone marrow extract), and refer to manufacturer's instructions of the extraction kit.

In the EBX-013 kit, CI-PCR on the CY5 channel can be added before extraction or in the PCR reaction. It ensures that a negative result is not due to an extraction problem or due to the presence of PCR inhibitors at high quantity.

We recommend the addition of 10 µl of CI-PCR per extraction. If the CI-PCR is added to control the PCR, CI-PCR is added to the reaction mix (1µl per PCR reaction). See real-time PCR protocol for details.

CI-PCR is available from Eurobio (Ref EurobioPlex EBX-002).


II- Real-time PCR procedure

General comment:

Positive control and the DNA extraction and PCR inhibition control (CI-PCR) contain a very high concentration of matrix. Manipulations must be performed carefully to avoid contamination. To control the functioning of the PCR and steps of extraction and real-time PCR amplification, it is necessary to test the positive control CP-Barto as well as the negative control (water supplied = CN-H₂O + CI-PCR) (see II-2/6) of real-time PCR procedure).

On T-COR 8®-IVD, we recommend that the positive and negative controls be tested when a new kit is being opened.

II-1/ Diagram of the procedure

II-2/ Detailed procedure

- 1) Homogenize the tube of Enzyme, and vortex Oligomix, CP-Barto, and CI-PCR tubes before starting, and centrifuge.
- 2) Prepare Mastermix as below; prepare the reaction mix as follows by multiplying by the number of samples N to test (including positive and negative controls). Prepare enough reagents for at least N+3 reactions^α. (Refer to part 1-(a) or 1-(b) of the previous diagram according to the condition).

Case (a): for a DNA sample extracted without CI-PCR

Number of reaction(s)	1	N+3 ^α
Enzyme	12.5 µl	(N+3) x 12.5 µl
Oligomix	7.5 µl	(N+3) x 7.5 µl
CI-PCR	1 µl	(N+3) x 1 µl
Total volume mastemix	21 µl*	(N+3) x 21 µl

* The volume difference between condition (a) or (b) has no effect on performance.

^α N+1 on T-COR 8®-IVD

Case (b): for a DNA sample extracted WITH CI-PCR

Number of reaction(s)	1	N+3 ^α
Enzyme	12.5 µl	(N+3) x 12.5 µl
Oligomix	7.5 µl	(N+3) x 7.5 µl
Total volume mastemix	20 µl*	(N+3) x 20 µl

* The volume difference between condition (a) or (b) has no effect on performance.

^α N+1 on T-COR 8®-IVD

- 3) Homogenize mastermix prepared in 2) and centrifuge briefly
- 4) Distribute 20 µL Mastermix* using a micropipette and filtered tips in each tube or well of microplate for real time PCR.
- 5) Add 5 µL of extracted DNA sample.
- 6) In parallel, test the following controls (see this specific point for T-COR 8®-IVD on pages 41-47):
 - Positive control:
 - 20µL of mastermix + 5µL of CP-Barto
 - Negative control:
 - Case (a): for a DNA sample extracted without CI-PCR
 - 20 µL mastermix + 5 µL supplied water (CN-H2O)
 - Case (b): for a DNA sample extracted with CI-PCR
 - 20 µL mastermix + 4 µL water supplied (CN-H2O) + 1 µL CI-PCR
- 7) Close immediately with an adhesive film or transparent caps to avoid contamination.
- 8) Centrifuge briefly to collect all the PCR reaction mix at the bottom of the tubes or plate.
- 9) Program the real-time PCR instrument as follows (on T-COR 8®-IVD, no programming is needed thanks to the use of Barcodes- see pages 41-47) :

Program	Temperature	Duration	Cycle(s)	
Activation of Taq polymerase	95°C	3 min	1	-
Amplification	95°C	10 sec	45	-
	60°C	30 sec		Fluorescence acquisition

Note 1: On LightCycler ® 480 systems (Roche), two optical systems are available: only "System II" is compliant with use of the kit. Apply color compensation for the following channels: FAM-HEX/VIC-Texas Red-Cy5 (465-510,533-580,533-610,618-660).

Note 2: For the Applied Biosystems systems, select "NONE" in "PASSIVE REFERENCE".


Note 3: On Rotorgene ™, please calibrate the signal by clicking on "GAIN optimization".

Note 4: On CFX96 (Biorad), start the run using the v1.6 or later version of CFX Manager software, and analyse with v 3.1 (see § Validation of the Experiment)


VALIDATION OF THE EXPERIMENT

The analysis of data post-acquisition on a CFX96 PCR instrument (Biorad) must be done with version 3.1 of CFX Manager Software (Biorad). In order to use this v3.1 version from a run started with an older version, follow the procedure below: at the end of the run, the data file with .pcrd suffix must be open and treated with version 3.1 of CFX Manager (Biorad).

If the run was done with CFX Manager v1.6 for instance, to open the data file with CFX Manager v3.1, click on CFX Manager v3.1 icon. The screen below appears.


- Click on File and select Open, then Data File

- Select the file you want to analyze and click on Open.

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Once this is done, the analysis can start.

For the assay to be valid, the results for the controls must be the following (Table 5). Otherwise, the experiment is not valid. On T-COR 8®-IVD, we recommend that the positive and negative controls be tested when a new kit is being opened.

Table 4 :

Positive control	
FAM	Ct ≤ 32
HEX	Ct ≤ 32
Texas Red	Ct ≤ 32
Negative control	
FAM	Ct not determined
HEX	Ct not determined
Texas Red	Ct not determined
CY5	Ct ≤ 35

DATA ANALYSIS AND INTERPRETATION

DNA extraction and PCR inhibition control in biological samples:

Two results can be obtained:

1/ the DNA extraction and PCR inhibition control (CI-PCR; channel Cy5) test is positive: DNA has been properly extracted, and there is no inhibitor of PCR. The result can be validated.

2/ the DNA extraction and PCR inhibition control (CI-PCR; channel Cy5) test is negative: either DNA was not extracted, or the presence of PCR inhibitors inhibits the PCR reaction. It is then recommended to repeat the extraction or dilute the sample, unless a specific signal appears in the targeted channels (FAM, HEX, Texas Red).

For clinical samples, the following results are possible:

*** Ct cut off for positivity of samples:**

Channel FAM *Bartonella spp*: + Positive => Ct positive (≤ 45)

Channel HEX *Bartonella henselae*: + Positive => Ct positive (≤ 45)


Channel TEXAS RED *Bartonella quintana*: + Positive => Ct positive (≤ 45)

(Ct cut off for positivity of samples for T-COR 8®-IVD is: + Positive => Ct positive (≤ 40) for the 3 channels).

PCR Signal		Presence of bacteria	Test validity/comment
FAM or HEX or Texas Red	CY5	Bartonella spp Bartonella henselae Bartonella quintana	
+	+	Yes	VALID
-	+	No	VALID
+	-	Yes	VALID Possible inhibition of PCR or problem with extraction that doesn't prevent the detection of the bacteria.
-	-	No possible interpretation	Inhibition of PCR or problem with extraction → dilute 5 x the sample; otherwise extract once again the sample

PERFORMANCE ANALYSIS

Example of experiment performed on real-time PCR thermocycleur CFX96 (Biorad):

Analytical sensitivity: CP-Barto:

12,5 copies/μl for *B.spp.*,

5 copies/μl for *B.henselae*,

5 copies/μl for *B.quintana*

Linearity for quantification: CP-Barto: from the limit of detection to 10^9 copies/μl

Coefficient of correlation and efficiency

	CHANNEL	Efficiency	Coefficient of correlation R 2	Slope
Barto spp	FAM	95.5 %	0.999	-3.436
Barto henselae	HEX	98.5 %	0.995	-3.358
Barto quintana	TEXAS RED	95.3 %	0.998	-3.440

Variability of the signal

Studies of reproducibility of the kit EurobioPlex Bartonella EBX-013 have been achieved on plasmid BspPHQ (positive control of Bartonella spp., Bartonella henselae and Bartonella quintana) on CFX96 (Biorad). The tables below indicate the average of the coefficient of variation (CV) depending on the concentration of the plasmid.

Reproducibility within an experiment (intra-experiment):

copies/μl	CV % CHANNEL FAM Bartonella spp	CV % CHANNEL HEX Bartonella henselae	CV % CHANNEL TEXAS RED Bartonella quintana
10 8	0.50	0,05	0,35
10 7	0.35	0,48	0,44
10 6	0.48	0,22	0,19
10 5	0.30	0,09	0,12
10 4	0.50	0,29	0,30
10 3	0.27	0,28	0,28
100	0.35	0,50	0,33
50	1.56	2,10	1,63
25	2.00	2,16	2,32
12.5	2.80	2,88	2,52
10			1,60
5	/	/	2,91
AVERAGE CV %	0,91	0,91	1,08

Reproducibility between experiments (inter-experiments):

copies/μl	CV % CHANNEL FAM Bartonella spp	CV % CHANNEL HEX Bartonella henselae	CV % CHANNEL TEXAS RED Bartonella quintana
10 8	8,91	3,96	2,80
10 7	6,36	1,99	4,04
10 6	3,77	1,32	4,62
10 5	3,11	1,30	3,49
10 4	2,58	1,09	3,27
10 3	0,19	0,48	4,86
10 2	1,38	2,08	3,89
10	1,24	2,40	4,00
5			4,26
AVERAGE CV %	3,44	1,83	3,91

Validation on clinical samples

The tests were performed on 109 DNA extracts from: 3 sera, 6 blood samples, 90 lymph node biopsies, 1 spleen biopsy, 1 extract from bone marrow, 7 cardiac biopsies, 1 sample of unknown origin). These samples were previously tested by the routine laboratory technique from a French National Reference Center (CNR) and compared with the kit EurobioPlex Bartonella EBX-013 on CFX96 (Biorad).

The specificity and sensitivity of the EBX-013 kit have been evaluated by analysis of :

- 60 Bartonella spp POSITIVES
- 49 Bartonella spp NEGATIVES
- 45 Bartonella henselae POSITIVES/quintana NEGATIVES
- 64 Bartonella henselae NEGATIVES, 15 of them are the 15 quintana positives ; no cross-reactivity
- 15 Bartonella quintana POSITIVES/henselae NEGATIVES
- 94 Bartonella quintana NEGATIVES, 45 of them are the 45 henselae positives

		EBX-013 Quadriplex Bartonella EUROBIO								
		B spp			B henselae			B quintana		
		POS	NEG	TOTAL	POS	NEG	TOTAL	POS	NEG	TOTAL
Pretest	POS	60	0	60	41	4	45	15	0	15
	NEG	2	47	49	0	64	64	0	94	94
	TOTAL	62	47	109	41	68	109	15	94	109

Sensitivity

Spp : > 99 % (n=60)

Henselae : 41/45= 91,1 %

Quintana : > 99 % (n=15)

Specificity

Spp : 47/49= 95,9%

Henselae : > 99 % (n=64)

Quintana : > 99 % (n=94)

The tests were performed on CFX96.

Bacterial DNA extraction was done using the kits detailed below, following manufacturer's instructions:

- EZ1 Blood DNA kit (Qiagen) for blood and serum samples,
- EZ1 Tissue DNA kit (Qiagen) for tissue samples (detailed above).

SPECIFICITIES OF THE REAL-TIME PCR T-COR 8®-IVD INSTRUMENT

T-COR 8®-IVD is a real-time PCR instrument with 8 independent wells, which can be programmed independently, patient by patient, regarding thermic protocol, assays, and time of start of the run.

The 8 wells can also be used simultaneously with the same assay.

Note: Vortex the tubes with the provided T-core vortex before putting the tubes in the T-COR 8®-IVD, and always check that there is no bubble, and that the liquid is all located at the bottom of the tube.

Controls

On T-COR 8®-IVD, we recommend that the positive and negative controls be tested when a new kit is being opened.

After this initial control, the internal control allows to check that the PCR is working properly, and that there is no PCR inhibition.

Fluorophores

Four combinations of fluorophores are available on T-COR 8®-IVD.

For all EBX, such as EBX-013, FAM/HEX/Texas Red/Cy5 combination is used. If some channels are not used, do not consider this channel for results analysis.

Use of Barcodes (available on pages 45-47)

1- Select Menu > New Run

2- Select Barcode

3- Scan the appropriate Barcode on the right side of the instrument:

- For the positive control (Barcode EBX-013 Pos Ctrl),
- For the negative control (Barcode EBX-013 Neg Ctrl),
- For a patient sample (Barcode EBX-013)

The instrument randomly selects one of the 8 available wells. Follow the instructions on the instrument.

4- Lift up the lid of the selected well (well x), check that the blue LED light is on, and select "Yes"

5- Place the tube in the corresponding well and select "Next".

6- If you wish to name the sample (optional), select "sample x", name it and select "Accept".

7- Select "Next"

8- To add/test another control or sample, select « Add well » and go back to point 3-.

9- Select "Start run".

Note: *The name of a sample or of a run cannot be modified or added after the run is finished. It must be done before or during the run.*

Presentation and Interpretation of results

Ct values and amplification curves are available in the SmartCT™ Table (Ct values on each channel), and “Ct versus PCR cycles” graphs, both being available in real-time while the run is ongoing.

An automatic analysis for positive and negative controls as well as patients' samples is available in “Interpretations”, at the end of the run, in the “View” window.

For positive and negative controls, the following results are possible:

« Neg Ctrl Fail »: Not valid

« Neg Ctrl Valid »: Valid

« Pos Ctrl Fail »: Non valid

« Pos Ctrl Valid »: Valid

For patients status:

"Detected": Positive for at least one target * → green box

"Not detected": Negative → red box

* The results for the specific species are specified.

CAUTION !

When the box is green, it is important to read the status for each target as some targets may be negative.

« Invalid »: Invalid result -> retest

Example of automatic interpretation of results on T-COR 8®-IVD instrument screen and associated report:

➤ For valid positive and negative controls, for the limit of detection (LOD) and two positive samples (samples 3 and 4) « Detected » for their corresponding targets, and for two negative samples indicated “Not Detected” (samples 1 and 2):


❖ On the instrument screen:

Interpretations		
1	LOD Repl1	EBX-013
	Detected: B. quint, B. hens, B. spp	
2	LOD Repl2	EBX-013
	Detected: B. quint, B. hens, B. spp	
3	ech1	EBX-013
	Not Detected - B. quint, B. hens, B. spp	
4	ech2	EBX-013
	Not Detected - B. quint, B. hens, B. spp	
5	ech3	EBX-013
	Detected: B. hens, B. spp - B. quint not detected	
6	ech4	EBX-013
	Detected: B. quint, B. spp - B. hens not detected	
7	CN	EBX-013 Neg Ctrl
	Détecté(e-s): B. hens, B. spp - B. quint not detected	
8	CP	EBX-013 Pos Ctrl
	Non détecté - B. quint B. hens B. spp	

❖ On the report:

Well	Sample	Assay	FAM	HEX	TxR	Cy5	Call	Note
1	LOD Rep11	EBX-013	30.7	31.1	30.5	25.7	Detected • B. quint • B. hens • B. spp	
2	LOD Rep12	EBX-013	32.0	31.9	31.2	25.6	Detected • B. quint • B. hens • B. spp	
3	ech1	EBX-013				25.9	Not Detected	B. quint, B. hens, B. spp
4	ech2	EBX-013				25.9	Not Detected	B. quint, B. hens, B. spp
5	ech3	EBX-013	28.3	30.0		25.3	Detected • B. hens • B. spp	B. quint not detected
6	ech4	EBX-013	22.5		22.9	25.1	Detected • B. quint • B. spp	B. hens not detected
7	CN	EBX-013 Neg Ctrl				26.0	Detected • Neg Ctrl Valid	
8	CP	EBX-013 Pos Ctrl	24.3	24.0	24.0	25.1	Detected • Pos Ctrl Valid	

Example of amplification curves:

Barcodes for EBX-013 for use on T-COR8®-IVD

EBX-013 Pos Ctrl
Positive Control CP-Barto

EBX-013 Neg Ctrl

Water = Negative Control (CN-H2O)

EBX-013

BIBLIOGRAPHY

Angelakis E, Raoult D. Pathogenicity and treatment of Bartonella infections. Int J Antimicrob Agents. 2014 Jul;44(1):16-25.

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WASTE DISPOSAL

Be in accordance with the law on the elimination of waste of clinical infectious material.

SYMBOLS

	Reference
	Batch number
	Limits of storage temperature
	Expiration Date
	Content sufficient for « N » reactions
	Keep protected from light
	Manufacturer
	Instructions for use
	CE labeled product
	In vitro diagnostic

eurobio
SCIENTIFIC
 7, avenue de Scandinavie
 ZA de Courtaboeuf
 91940 Les Ulis
 FRANCE