EurobioPlex

Coxiella burnetii

PCR EN TEMPS REEL

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
Version 5.00 - 2020/07/07

Validated on:

- CFX96™ Real Time PCR detection system (Bio-Rad) with analysis on CFX Manager v 3.1 (Bio-Rad)
- LightCycler®480 Instrument II (Roche) with analysis on LightCycler® 480 software v1.5 (Roche)
- Applied Biosystems® 7500 Real-Time PCR Systems (Applied Biosystems) with analysis on 7500 Software v2
- T-COR 8™-IVD (Tetracore Inc.) with T-COR 8 Per-Well SmartCT™ software

Storage conditions:

Keep all reagents between -15°C and -22°C until use and after first use

Instructions for use

TABLE OF CONTENTS

Intended use	25
Introduction	25
Principle of detection	26
Description and content of the kit	26
Storage	28
Cautions and notes	28
Samples collection, transport and storage	28
Procedure	29
I- DNA extraction	29
II- Real-time PCR procedure	29
II-1/ Diagram of the procedure.....	30
II-2/ Detailed procedure	31
Validation of the experiment	32
Data analysis and Interpretation	34
Performance analysis.....	35
Specificities of the Real-Time PCR T-COR 8™-IVD instrument	36
Barcodes for EBX-017 for use on T-COR 8™-IVD	40
Bibliography	43
Waste disposal	43
Symbols	44

INTENDED USE

The EBX-017 test uses real time polymerase chain reaction (PCR) amplification and is designed for the qualitative detection of *uburnetii* from a DNA nucleic acids extract. This kit also contains a control of inhibition of PCR. This test is indicated to confirm a diagnosis of presumption of infection in patients or complete a proven or indeterminate diagnosis by standard techniques (diagnosis by serology, by culture, immunohistochemistry or gene amplification by PCR). Extracted DNA is the starting material for the EurobioPlex *Coxiella burnetii*.

The EurobioPlex EBX-017 has been validated on the following specimen:

- Serum
- Blood
- Cardiac biopsy
- Pharyngeal biopsy
- Bone biopsy
- Puncture fluid (pus)
- Lymph node biopsy
- Negative culture supernatants
- Biopsy of aneurysm

INTRODUCTION

Q fever (Q = query) was so named to highlight the uncertainties about its etiology and epidemiology. It is a zoonotic disease caused by an intracellular gram-negative bacterium: coccobacilli *Coxiella burnetii*. First described in 1935 in employees of a slaughterhouse in Brisbane, Australia, this disease is nevertheless ubiquitous but mostly found in rural areas. Indeed, the bacterium infects most domestic mammals (cattle, sheep, goats, horses, dogs, cats, rabbits) and wild mammals (deer, foxes, rodents). Infected animals transmit the bacteria by the placenta, urine, feces and milk. Inter-humans transmission remains anecdotal, most often by the way of obstetrical manipulations of infected women. Contamination in humans occurs mainly through inhalation of infected fine particles from herds with duration of incubation of one to three weeks. Individuals of all ages and both sexes can be affected. The prevalence is not known accurately and varies greatly from one region to another depending on exposure to ruminants. Nevertheless, regional data such as those of the National Reference Center allow an estimation of the incidence of acute Q fever of 50 per 100,000 inhabitants per year in the South of France. In Europe, the highest incidence is during spring and early summer.

The primo-infection can be asymptomatic (in 60% of cases) or symptomatic with the acute form of Q fever (40% of cases). The clinical presentation of acute Q fever is a flu syndrome, atypical pneumonia or hepatitis. It is primarily described by high fever (91%), sometimes prolonged, major headache (51%), myalgia (37%), arthralgia (27%) and cough (34%). In healthy subjects, spontaneous evolution is a complete healing. In some people, on the other hand, *Coxiella burnetii* is able to multiply despite the response triggered by the primary infection. These subjects at risk are pregnant women, patients with valvulopathy, with a valvular prosthesis, and immuno-compromised patients. When the immune system is unable to control the infection, chronic Q fever can develop. Endocarditis is the main manifestation of chronic Q fever (found in 80% of cases). Endocarditis is the severe factor of the infection with a

mortality of 25 to 60% of the cases in the absence of treatment. Q fever is 5% of all cases of endocarditis in France. The curative treatment consists in an antibiotic therapy. Doxycycline is administered for a period of two weeks to several months or even years in the case of endocarditis (associated with another antibiotic). The effectiveness of the treatment is evaluated by the follow up of antibody levels. Serological surveillance must be of at least five years.

The diagnosis of this disease is often tricky because these bacteria are difficult to cultivate. The diagnostic approach most often requires the use of several complementary methods depending on the immune status of the patient and the symptomatology. Specific diagnosis relies both on the indirect diagnosis (serology), direct identification of the bacterium on histological sections by Immunohistochemistry using monoclonal antibody, and gene amplification by PCR from blood or serum or different other samples such as heart valve, placenta, aneurysm, cerebrospinal fluid. Indirect immunofluorescence reaction is to date the technique of reference for the serodiagnosis of Q fever, with quantification of IgA, IgG and IgM for both phases (antigens of phase II in the acute form and phase I in the chronic form). **Direct PCR amplification is the most sensitive and most specific technique for the diagnosis of infection with *Coxiella burnetii*.**

PRINCIPLE OF DETECTION

The EurobioPlex EBX-017 is a test using real-time amplification of bacterial DNA of *Coxiella burnetii*, as well as a DNA extraction and PCR inhibition control. The test is performed from extracted DNA in a single well/tube.

The DNA extraction and PCR inhibition control allows to check for variations that may occur during the DNA extraction step of biological samples and real-time PCR amplification. Thus, it ensures that a negative result may not be due to a bad DNA extraction, and/or the presence of PCR inhibitors.

This kit is composed of one duplex: *Coxiella burnetii* (Coxb) + control of DNA extraction and PCR inhibition (CI-PCR).

DNA of *Coxiella burnetii* is detected using HEX labeled probe. DNA extraction and PCR inhibition control is detected using a CY5 labeled probe. All probes emit a specific fluorescence following their hydrolysis during the elongation of the amplification product. The measurement of the intensities of real-time fluorescence correlates with the accumulation of specific amplification products.

DESCRIPTION AND CONTENT OF THE KIT

The real-time PCR EBX-017 kit is ready to use for the specific detection of *Coxiella burnetii*.

Fluorescence is emitted and individually recorded through optical measurements during the PCR. The detection of the amplified fragment is performed by a fluorimeter using the channels shown in Table 1 below.

The kit contains reagents and enzyme for the amplification of DNA of this bacterium, and the DNA extraction and PCR inhibition control (see Table 2).

Table 1:

Target	Fluorophore	Excitation	Emission
<i>Coxiella burnetii</i>	HEX	535 nm	555 nm
PCR inhibition control (CI-PCR)	Cy5	647 nm	667 nm

Equivalent channels on different real time PCR cyclers:

- Channel **HEX** (Chromo 4/CFX96™, Mx Systems, T-COR 8™-IVD), Channel VIC (ABI Systems), Channel Alexa532 (SmartCycler II), Channel 580 (LC 480), Channel Yellow (RotorGene),

- Channel **Cy5** (CFX96™ /Chromo4, Systems ABI, Systems Mx, T-COR 8™-IVD), Channel Alexa647 (SmartCycler II), Channel 660 (LC 480), Channel Red (RotorGene)

Table 2:

Cap color	Components of the kit	EBX-017-24 (24 tests)	EBX-017-48 (48 tests)	EBX-017 (96 tests)	Reconstitution
Red	ENZYME	4 x 100 µl	2 x 350 µl	4 x 350 µl	Ready to use
Transparent	OLIGOMIX	4 x 60 µl	2 x 210µl	4 x 210 µl	Ready to use
Yellow	Positive Control (CP-Coxb)	4 x 25 µl de CP	2 x 25 µl de CP	4 x 25 µl de CP	Ready to use
White	PCR inhibition control (CI-PCR)	4 x 280 µl	2 x 280 µl	4 x 280 µl	Ready to use
Blue	Water = negative control (CN-H2O)	1 ml	1 ml	1 ml	Ready to use

Oligomix: contains the primers and probes of the duplex (Coxb + CI-PCR)

Required material not provided:

- ◇ Biological Hood
- ◇ Real-time PCR instrument
- ◇ Micro centrifuge
- ◇ Vortex
- ◇ Plates / tubes for real-time PCR
- ◇ Micropipettes
- ◇ DNase-free RNase-free filter tips for micropipettes
- ◇ Sterile microtubes
- ◇ Gloves (powder free)

STORAGE

All reagents must be stored between -15 and -22°C.

All reagents can be used until the expiration date indicated on the label of the kit.

Many freezing / defrosting cycles (> 3x) must be avoided, and could lead to decrease in sensitivity.

CAUTIONS AND NOTES

Read carefully instructions before starting

- ◇ The experiment must be performed by competent staff.
- ◇ Instruments must have been properly installed, calibrated and maintained according to the manufacturer's recommendations.
- ◇ If a non validated equipment is used, it is the user's responsibility, and if so, the performances are not guaranteed.
- ◇ Clinical samples are potentially infectious and must be processed under a laminar flow hood.
- ◇ The experiment must be performed according to good laboratory practices.
- ◇ Do not use this kit after expiration date, specified on the label.
- ◇ The kit is shipped with dry ice, and the components of the kits must arrive frozen. If one or more components are defrosted, or of the tubes have been damaged, contact Eurobio Scientific.
- ◇ Once defrosted, spin down briefly the tubes before use.
- ◇ It is recommended to define three working areas: 1) Isolation of DNA, 2) Preparation of the reaction mix and 3) Amplification / Detection of amplified products.
- ◇ Use specific lab coat and gloves (powder free) in each working area.
- ◇ Pipettes, reagents and other materials must not cross each area.
- ◇ Specific caution is required to preserve the purity of the reagents and reaction mixtures.
- ◇ Appropriate methods of preparation/extraction of DNA to produce high quality DNA and to be followed by a real-time PCR application should be used, particularly avoiding all sources of DNase contamination.
- ◇ Always use RNase-free DNase-free filtered tips for micropipettes.
- ◇ Do not pipette with mouth and do not eat, drink or smoke in the area.
- ◇ Avoid sprays.

SAMPLES COLLECTION, TRANSPORT AND STORAGE

- ◇ Collect samples in sterile tubes. **The use of heparin as an anticoagulant is prohibited.**
- ◇ It is the responsibility of the user to master its own conditions of collection, transport and storage of samples, and extraction of DNA by suitable systems to produce DNA of good quality.
- ◇ The samples should be extracted immediately or stored following the recommendations in the table below (Table 3).

Table 3 :

Recommendations of maximum storage of samples before extraction	
Room temperature	2 h
+4°C	16h
-20°C	Long term storage

- ◇ The user can refer to the World Health Organization or High Health Authority for storing samples.
- ◇ Transport of clinical samples must follow local regulations for this type of infectious agents.

PROCEDURE

I- DNA extraction

It is the user's responsibility that the extraction system used be compatible with downstream real time PCR technology. For this kit, we recommend to use extraction methods of DNA from plasma or serum for instance, and refer to manufacturer's instructions of the extraction kit.

In the EBX-017 kit, CI-PCR on the CY5 channel can be added before extraction or in the PCR reaction. It ensures that a negative result is not due to an extraction problem or due to the presence of PCR inhibitors at high quantity.

We recommend the addition of 10 µl of CI-PCR per extraction and elution in 50 µl. If the CI-PCR is added to control the PCR, CI-PCR is added to the reaction mix (1 µl per PCR reaction). See real-time PCR protocol for details.

CI-PCR is available from Eurobio Scientific (Ref EurobioPlex EBX-002).


II- Real-time PCR Procedure

General comment:

Positive control and the DNA extraction and PCR inhibition control (CI-PCR) contain a very high concentration of matrix. Manipulations must be performed carefully to avoid contamination. To control the functioning of the PCR and steps of extraction and real-time PCR amplification, it is necessary to test a positive control, CP-Coxb, as well as the negative control (water supplied = CN-H₂O + CI-PCR) (see II-2/6) of real-time PCR procedure).

On T-COR 8TM-IVD, we recommend that the positive and negative controls be tested when a new kit is being opened.

II-1/ Diagram of the procedure

II-2/ Detailed Procedure

- 1) Homogenize the tube of Enzyme, and vortex Oligomix, CP-Coxb and CI-PCR tubes, and centrifuge briefly.
- 2) Prepare a Mastermix as below. N is the number of PCR reactions. Plan to prepare enough reagents for at least N+3 reactions** (Refer to part 1-(a) or 1-(b) of the previous diagram according to the condition).

Case (a): For DNA samples extracted without CI-PCR

Number of reactions	1	N+3**
Enzyme	12.5 µl	(N+3) x 12.5 µl
Oligomix	7.5 µl	(N+3) x 7.5 µl
CI-PCR	1 µl	(N+3) x 1 µl
Total Volume Mastermix	21 µl*	(N+3) x 21 µl

Case (b): For DNA samples extracted with CI-PCR

Number of reactions	1	N+3**
Enzyme	12.5 µl	(N+3) x 12.5 µl
Oligomix	7.5 µl	(N+3) x 7.5 µl
Total Volume Mastermix	20 µl*	(N+3) x 20 µl

* The volume difference between case (a) or (b) has no effect on performance.

** For use on T-COR 8™-IVD, we recommend a Mastermix for N + 1 reactions for patient by patient use of the kit

- 3) Homogenize the Mastermix prepared in 2) and centrifuge briefly.
- 4) Distribute 20 µl of Mastermix* using a micropipette and filtered tips in each tube or well of microplate for real time PCR.
- 5) Add 5 µl of sample (DNA extracted or non-extracted serum).
- 6) In parallel test the following controls (see this specific point for T-COR 8™-IVD on pages 36-39):

- Positive control:
 - 20µl Mastermix + 5µl CP-Coxb

- Negative control:

Case (a): For DNA samples extracted without CI-PCR

- 20 µl Mastermix + 5 µl water supplied (CN-H2O)

Case (b): For DNA samples extracted with CI-PCR

- 20 µl Mastermix + 4 µl water supplied (CN-H2O) + 1 µl CI-PCR

- 7) Close immediately the tubes, or the plate with an adhesive film, to avoid all contamination.
- 8) Centrifuge briefly to collect all the PCR reaction mix at the bottom of the tubes or plate.
- 9) Program the real-time PCR instrument as follows (on T-COR 8™-IVD, no programming is needed thanks to the use of Barcodes- see pages 40-42):

Program	Temperature	Duration	Cycle(s)	
Denaturation	95°C	3 min	1	-
Amplification	95°C	10 sec	45	-
	60°C	30 sec		Fluorescence Acquisition

Note 1: On LightCycler ® 480 systems (Roche), two optical systems are available: only "System II" is compliant with use of the kit.

Note 2: For the Applied Biosystems systems, select "ROX" in "PASSIVE REFERENCE".


Note 3: On Rotorgene™, please calibrate the signal by clicking on "GAIN optimization".

Note 4: On CFX96™ (Bio-Rad), start the run using the v1.6 or later version of CFX Manager software, and analyse with v 3.1 (see § Validation of the Experiment)


VALIDATION OF THE EXPERIMENT

The analysis of data post-acquisition on a CFX96™ PCR instrument (Bio-Rad) must be done with version 3.1 of CFX Manager Software (Bio-Rad). In order to use this v3.1 version from a run started with an older version, follow the procedure below: at the end of the run, the data file with .pcrd suffix must be open and treated with version 3.1 of CFX Manager (Bio-Rad).

If the run was done with CFX Manager v1.6 for instance, click on CFX Manager v3.1 icon to open the data file with CFX Manager v3.1. The screen below appears.


- Click on “File” and select “Open”, then “Data File”:

- Select the file you want to analyze and click on “Open”.

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Once this is done, the analysis can start.

For the assay to be valid, the results for the controls must be the following (Table 4). Otherwise, the experiment is invalid. On T-COR 8™-IVD, we recommend that the positive and negative controls be tested when a new kit is being opened.

Table 4:

Positive control	
HEX	Ct ≤ 28
Negative control	
HEX	Undetermined
CY5	Ct ≤ 30

On T-COR 8™-IVD, results are automatically interpreted thanks to the Barcodes.

DATA ANALYSIS AND INTERPRETATION

DNA extraction and PCR inhibition control in samples:

Two results can be obtained:

1/ the DNA extraction and PCR inhibition control (CI-PCR, channel Cy5) test is positive. The result can be validated.

2/ the DNA extraction and PCR inhibition control (CI-PCR, channel Cy5) test is negative: either DNA was not extracted, the PCR did not work well, or the presence of PCR inhibitors inhibits the PCR reaction. It is then recommended to repeat the extraction or dilute the sample, unless a specific signal appears in the targeted channel (HEX).

For clinical samples, the following results are possible:


***Ct cut off for positivity of samples: HEX channel:
+ Positive => Ct ≤ 45**

PCR Signal		Presence of <i>Coxiella burnetii</i>	Test validity/comment
HEX	CY5		
+*	+	Yes	VALID
-	+	No	VALID
+*	-	Yes	VALID Possible inhibition of PCR or problem with extraction that does not prevent the detection of the <i>Coxiella burnetii</i> .
-	-	No possible interpretation	Inhibition of PCR or problem with extraction → dilute 5 x the sample; otherwise extract once again the sample

On T-COR 8TM-IVD, results are automatically interpreted thanks to the Barcodes.

PERFORMANCE ANALYSIS

Example of experiment performed on real-time PCR thermocycleur CFX96™ (Bio-Rad):

Analytical sensitivity: CP-Coxb: 10 copies/μl

Linearity for quantification: CP-Coxb: from 10 copies/μl to 10⁹ copies/μl

Variability of the signal

Studies of reproducibility of the kit EurobioPlex *Coxiella burnetii* EBX-017 have been achieved on plasmid Coxb (positive control of *Coxiella burnetii*) on CFX96™ (Bio-Rad). The table below indicates the average of the coefficient of variation (CV) depending on the concentration of the plasmid.

Coefficient of variation (%)	<i>Coxiella burnetii</i> (CHANNEL HEX)	CI-PCR (CHANNEL CY5)
REPRODUCIBILITY WITHIN-EXPERIMENT	0.90 %	0.45 %
REPRODUCIBILITY BETWEEN-EXPERIMENTS	3.08 %	3.92 %

Validation on clinical samples

The tests were performed on 103 DNA extracts from: 13 sera, 14 blood samples, 7 cardiac biopsies, 2 biopsies of the pharynx, 30 bone biopsies, 6 puncture fluids (pus), 12 lymph node biopsies, 14 negative culture supernatants, and 5 biopsies of aneurysm. These samples were previously tested by the routine laboratory technique from a French National Reference Center, and results were compared with the kit EurobioPlex *Coxiella burnetii* EBX-017 on CFX96™ (Bio-Rad).

The specificity and sensitivity of the EBX-017 kit have been evaluated on the following samples:

- 48 *Coxiella burnetii* POSITIVES
- 55 *Coxiella burnetii* NEGATIVES

		EBX-017 <i>Coxiella burnetii</i>		
		POSITIVE	NEGATIVE	TOTAL
Pretest	POSITIVE	48	0	48
	NEGATIVE	0	55	55
	TOTAL	48	55	103

Specificity : > 99 %

Moreover, no cross-reactivity has been observed with the bacteria detailed below:

- *Mycoplasma* (n=1)=0%
- *Bartonella spp* (n=2)=0%
- *Bartonella henselae* (n=2)=0%
- *Klebsiella kingae* (n=2)=0%
- *Mycobacterium tuberculosis* (n=1)=0%
- *Tropheryma whipplei* (n=1)=0%
- *Staphylococcus aureus* (n=2)=0%
- *Mycobacterium avium* (n=1)=0%
- *Staphylococcus oralis/mitis/sanguis* (n=1)=0%
- *Francisella tularensis* (n=1)=0%

Sensitivity : > 99 %

Concordance is > 99 %.

The tests were performed on CFX96™ (Bio-Rad).

Bacterial DNA extraction was done using the kits detailed below, following manufacturer's instructions:

- EZ1 Blood DNA kit (Qiagen) for blood and serum samples,
- EZ1 Tissue DNA kit (Qiagen) for tissue samples (detailed above).

SPECIFICITIES OF THE REAL-TIME PCR T-COR 8™-IVD INSTRUMENT

T-COR 8™-IVD is a real-time PCR instrument with 8 independent wells, which can be programmed independently, patient by patient, regarding thermic protocol, assays, and time of start of the run.

The 8 wells can also be used simultaneously with the same assay.

Note: Mix the sample by pipetting up and down, and always check that there is no

bubble, and that the liquid is all located at the bottom of the tube.

Controls

On T-COR 8™-IVD, we recommend that the positive and negative controls be tested when a new kit is being opened.

After this initial control, the internal control allows to check that the PCR is working properly, and that there is no inhibition.

Fluorophores

For all EBX, such as EBX-017, FAM/HEX/Texas Red/Cy5 combination is used. Do not consider FAM and Texas Red for results analysis.

Patients' Tests

With testing of negative and positive controls only when first opening the kit and with thawing 3 times maximum each tube.	24 reactions format
Maximum total number of patient tests	22 patients
Maximum total number of patient tests, one patient at a time in each run	12 patients

Use of Barcodes (available on pages 40-42)

1- Select Menu > New Run

2- Select Barcode

3- Scan the appropriate Barcode on the right side of the instrument:

- For the positive control (Barcode EBX-017 Pos Ctrl),
- For the negative control (Barcode EBX-017 Neg Ctrl),
- For a patient sample (Barcode Assay: EBX-017)

The instrument randomly selects one of the 8 available wells. Follow the instructions on the instrument.

4- Lift up the lid of the selected well (well x), check that the blue LED light is on, and select "Yes"

5- Place the tube in the corresponding well and select "Next".

6- If you wish to name the sample (optional), select "sample x", name it and select "Accept".

7- Select "Next"

8- To add/test another control or sample, select « Add well » and go back to point 3-.

9- Select "Start run".

Note: The name of a sample or of a run cannot be modified or added after the run is finished. It must be done before or during the run.

Presentation and Interpretation of results

Ct values and amplification curves are available in the SmartCT™ Table (Ct values on each channel), and “Ct versus PCR cycles” graphs, both being available in real-time.

An automatic analysis for positive and negative controls as well as patients' samples is available in “Interpretations”, at the end of the run, in the “View” window.

For positive and negative controls, the following results are possible:

« Neg Ctrl Fail »: Invalid

« Neg Ctrl Valid »: Valid

« Pos Ctrl Fail »: Invalid

« Pos Ctrl Valid »: Valid

For patients status:

« Detected »: Positive → green box

« Not detected »: Negative → red box

« Invalid »: Invalid result -> retest → yellow box

Example of automatic interpretation of results on T-COR 8™-IVD instrument screen and associated report:

- Example: tests of valid positive and negative controls (wells 7 and 8 respectively), invalid negative control (well 1), 1 negative sample (sample 1-well 2) indicated “Not Detected”, and 1 positive sample (sample 2-well 3) indicated “Detected”:

❖ On the instrument screen:

Interpretations	
1	Neg Ctrl EBX-017 Neg Ctrl Invalid - Neg Ctrl Fail
2	Sample 1 EBX-017 Not Detected - Coxb Neg
3	Sample 2 EBX-017 Detected: Coxiella burnetii
7	Pos Ctrl EBX-017 Pos Ctrl Detected: Pos Ctrl
8	Neg Ctrl EBX-017 Neg Ctrl Detected: Neg Ctrl



❖ On the report:

Summary								
Well	Sample	Assay	FAM	HEX	TxR	Cy5	Call	Note
1	Neg Ctrl	EBX-017 Neg Ctrl					Invalid	Neg Ctrl Fail
2	Sample 1	EBX-017				25.4	Not Detected	Coxb Neg
3	Sample 2	EBX-017		26.6		25.3	Detected • <i>Coxiella burnetii</i>	
7	Pos Ctrl	EBX-017 Pos Ctrl		22.3		25.0	Detected • Pos Ctrl	
8	Neg Ctrl	EBX-017 Neg Ctrl				25.3	Detected • Neg Ctrl	

Example of amplification curves:

Graphs

Legend	
Well 1	
Well 2	
Well 3	
Well 7	
Well 8	


Barcodes for EBX-017 for use on T-COR 8™-IVD


EBX-017 Pos Ctrl

EBX-017 Neg Ctrl

Assay: EBX-017

BIBLIOGRAPHY

De Bruin A. Detection of *Coxiella burnetii* by qPCR : a comparison of available assays. RIVM Letter report 330071002/2011. National Institute for Public Health and the Environment. Ministry of Health, Welfare and Sport.

Raoult D. Q fever : still a mysterious disease. QJM 2002;95:491-492

Delsing CE, Warris A, Bleeker-Rovers CP. Q fever: still more queries than answers. Adv Exp Med Biol. 2011;719:133-43. Review.

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WASTE DISPOSAL

Be in accordance with the law on the elimination of waste of clinical infectious material.

SYMBOLS

	Reference
	Batch number
	Limits of storage temperature
	Expiration Date
	Content sufficient for « N » reactions
	Keep protected from light
	Manufacturer
	CE labeled product
	In vitro Diagnostic
	Instructions for use

eurobio
SCIENTIFIC

7, avenue de Scandinavie
ZA de Courtaboeuf
91940 Les Ulis
FRANCE