EurobioPlex
Dengue 1234
RT-PCR EN TEMPS REEL

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Callahan et al. Development and evaluation of serotype and group-specific fluorogenic reverse transcriptase PCR (TaqMan) assays for Dengue virus, J Clin. Microbiol, 2001; 41:19-4124.

Leparc-Goffart I, Baragatti M, Temmam S, Tuiskunen A, Moureau G, Charrel R, de Lamballerie X. Development and validation of real-time one-step reverse transcription-PCR for the detection and typing of dengue viruses. J Clin Virol. 2009;45(1):61-6.

Cecilia D, Kakade M, Alagarasu K, Patil J, Salunke A, Parashar D, Shah PS. Development of a multiplex real-time RT-PCR assay for simultaneous detection of dengue and chikungunya viruses. Arch Virol. 2015;160(1):323-7.

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24/48/96 reactions

Validated on:

- CFX96™ Real Time PCR detection system (Biorad) with analysis on CFX Manager v 3.1 (Biorad)
- LightCycler®480 Instrument II (Roche) with analysis on LightCycler® 480 software v1.5 (Roche)
- Applied Biosystems® 7500 Real-Time PCR Systems (Applied Biosystems) with analysis on 7500 Software v2

Storage conditions:

Keep all reagents between -15°C and -22°C until use and after first use

Instructions for use

TABLE OF CONTENTS

Intended Use	18
Introduction	18
Principle of detection	19
Description and content of the kit	19
Storage	21
Cautions and notes	21
Samples collection, transport and storage	22
Procedure	22
I- RNA extraction	22
II- Real-time RT-PCR procedure	22
II-1/ Protocol of real-time PCR	23
II-2/ Detailed procedure	24
Validation of the experiment	25
Data analysis and Interpretation	27
Performance analysis.....	28
Bibliography	29
Waste disposal	29
Symbols	30

INTENDED USE

The EBX-018 test uses reverse transcription-real time polymerase chain reaction (RT-PCR) amplification and is designed for the qualitative detection of DENGUE virus and the identification of serotypes 1, 2, 3 and 4 from a RNA nucleic acids extract. This test is indicated to confirm a diagnosis of presumption of infection in patients or to complete a proven or indetermined serological diagnosis. Extracted RNA is the starting material for the Eurobioplex Dengue 1234 kit.

The Eurobioplex Dengue 1234 has been validated on the following specimen

- Plasma

This amplification system has been validated on biological samples (plasma), and on a RNA bacteriophage (CI-ARN).

INTRODUCTION

Dengue virus is an Arthropod-Borne virus (Arboviruses), transmitted by arthropods (mosquitoes, ticks, sandflies, etc.). The Dengue virus belongs to the *Flaviviridae* virus family. They are transmitted by various species of mosquitoes of the *Aedes* genus as *aegypti* (mainly in tropical and subtropical America, and the South of the United States) and *albopictus* (mainly in Asia), *furcifer*, *africanus*, *taylori* and *luteophaeus* with a worldwide distribution, especially in the tropics. The first two species are identified to this day as the main vectors of epidemic, because of their adaptation to human habitat areas. In France, the southern regions are more particularly affected.

There are 4 serotypes of Dengue (DENV 1-4). The genome of Dengue is a RNA single-strand. The viral envelope has an E protein responsible for binding to a receptor and which constitutes the main antigen recognized by the immune system. Virus replication and transcription occurs in the cytoplasm of the host cell. The disease begins with fever and muscle and joint pains. In general case pathological complications don't occur but in 20% of cases infection can be more severe and sometimes take a life-threatening, particularly in the case of secondary infections with hemorrhagic shock. Keep in mind that infection by a strain does not protect the host from infection with another strain. No vaccine or curative treatment currently exists. Only a treatment of symptoms is provided. Dengue hemorrhagic fever has been recognized in 50's in the Philippines and Thailand. It is estimated that 100 million people are infected each year, and this number continues to be continuously growing. In several Asian countries, 4 serotypes coexist.

Diagnosis of Dengue strategy of the Haute Autorité de santé (HAS) is based on the detection by RT-PCR from blood samples or serum between the beginning and the 7th day and IgM and IgG serology from the 5th day. However, there are cross-reactions between the members of the family of *flaviviruses*, which makes difficult the analysis of serological data in the acute phase.

The HAS report of January 2013 indicates that the diagnosis of the disease is of great importance for the clinical management, surveillance and epidemic control measures. Biological diagnosis depends when sampling is performed to the date of onset of clinical signs. For OMS, viral isolation is the reference test. Detection of the virus is by viral isolation by cell culture. Detection of the virus can be made between the 1st and

the 7th day of the illness. The detection of the virus genome is done by RT-PCR that allows direct diagnosis in the early phase of the disease. Serological diagnosis is based on detection of IgM identified on average from the 5th day after the onset of clinical signs. IgG are positive between the 7th and 10th day and reached a maximum to J15.

PRINCIPLE OF DETECTION

The Eurobioplex Dengue 1234 kit is a test using reverse-transcription and real-time amplification of viral RNA of Dengue, as well as an encapsidated control of RNA extraction and RT-PCR inhibition. The test is performed from extracted RNA from a sample using two reactions in two wells/tubes.

The RNA extraction and RT-PCR inhibition control allows to check for variations that may occur during the RNA extraction step of biological samples and real-time PCR amplification. Thus, it ensures that a negative result may not be due to a bad RNA extraction, and/or due to the presence of PCR inhibitors in too large quantities.

This kit contains two different oligomixes which permit to detect the different serotypes listed below:

Oligomix 1: Dengue 1 and Dengue 4 (D14) + control of RNA extraction and RT-PCR inhibition (CI-ARN);

Oligomix 2: Dengue 2 and Dengue 3 (D23) + control of RNA extraction and RT-PCR inhibition (CI-ARN);

RNA of D1 and D4 is respectively detected using FAM and HEX labeled probes. RNA of D2 and D3 is also respectively detected using FAM and HEX labeled probes. For the two multiplexes, RNA extraction and RT-PCR inhibition control is detected using a CY5 labeled probe. All probes emit a specific fluorescence following their hydrolysis during the elongation of the amplification product. The measurement of the intensities of real-time fluorescence correlates with the accumulation of specific amplification products.

DESCRIPTION AND CONTENT OF THE KIT

The real-time RT-PCR Dengue 1234 kit is ready to use for the specific detection of this virus and the typing of its four serotypes.

Fluorescence is emitted and individually recorded through optical measurements during the PCR. The detection of the amplified fragment is performed by a fluorimeter using the channels shown in Table 1 below.

The kit contains reagents and enzymes for the amplification of each serotype of Dengue viral RNA, and the RNA extraction and RT-PCR inhibition control (Table 2).

Table 1: Serotypes of the Dengue virus detected by each oligomix

Target	Oligomix 1	Oligomix 2	Fluorophore	Excitation	Emission
Dengue 1	X	-	FAM	495 nm	515 nm
Dengue 4	X	-	HEX	535 nm	555 nm
Dengue 2	-	X	FAM	495 nm	515 nm
Dengue 3	-	X	HEX	535 nm	555 nm
RNA control (CI-ARN)	X	X	Cy5	647 nm	667 nm

Equivalent channels on different real time PCR cyclers:

- Channel **FAM** (Systems ABI, SmartCycler II, Systems Mx, CFX96/Chromo4), Channel 510 (LC 480), Channel Green (RotorGene)
- Channel **HEX** (Chromo 4/CFX96, Mx Systems), Channel VIC (ABI Systems), Channel Alexa532 (SmartCycler II), Channel 580 (LC 480), Channel Yellow (RotorGene),
- Channel **Cy5** (CFX96/Chromo4, Systems ABI, Systems Mx), Channel Alexa647 (SmartCycler II), Channel 660 (LC 480), Channel Red (RotorGene)

Note: On LC480 instrument II, apply color compensation for the following channels: FAM-HEX/VIC-Cy5 (465-510,533-580,618-660).

Table 2 :

Cap color	Components of the kit	96 reactions	48 reactions	24 reactions	Reconstitution
Red	Enzymes	4 x 410 µl	2 x 410 µl	410 µl	Ready to use
Transparent	Oligomix 1	4 x 70 µl	2 x 70 µl	70 µl	Ready to use
Green	Oligomix 2	4 x 70 µl	2 x 70 µl	70 µl	Ready to use
Blue	Water = negative control (CN-H2O)	1 ml	1 ml	1 ml	Ready to use
Yellow	Dengue 1 and 4 Positive Control (CP-D14)	4 x 50 µl	2 x 50 µl	50 µl	Ready to use
Orange	Dengue 2 and 3 Positive Control (CP-D23)	4 x 50 µl	2 x 50 µl	50 µl	Ready to use
White	RNA control (CI-ARN)	4 x 300 µl	2 x 300 µl	300 µl	Ready to use

Oligomix 1 : contains the primers and probes of the triplex (D1 + D4 + CI-ARN)

Oligomix 2: contains the primers and probes of the triplex (D2 + D3 + CI-ARN)

Required material not provided:

- ◇ Biological Hood
- ◇ Real-time PCR instrument
- ◇ Micro centrifuge
- ◇ Vortex
- ◇ Plates / tubes for real-time PCR
- ◇ Micropipettes
- ◇ DNase-free RNase-free filter tips for micropipettes
- ◇ Sterile microtubes
- ◇ Gloves (powder free)

STORAGE

All reagents must be stored between -15 and -22°C.

All reagents can be used until the expiration date indicated on the kit.

Many freezing / defrosting cycles (> 3x) must be avoided, and could lead to decrease in sensitivity.

CAUTIONS AND NOTES

Read carefully instructions before starting

- ◇ The experiment must be performed by competent staff.
- ◇ Instruments must have been properly installed, calibrated and maintained according to the manufacturer's recommendations.
- ◇ Clinical samples are potentially infectious and must be processed under a laminar flow hood.
- ◇ The experiment must be performed according to good laboratory practices.
- ◇ Do not use this kit after expiration date.
- ◇ The kit is shipped with dry ice, and the components of the kits must arrive frozen. If one or more components are defrosted, or if the tubes have been damaged, contact Eurobio-Ingen.
- ◇ Once defrosted, spin down briefly the tubes before use.
- ◇ Use of ice or cooling block is advised in case of long delay due for instance to large number of samples or high temperature.
- ◇ It is recommended to define three working areas: 1) Isolation of RNA, 2) Preparation of the reaction mix and 3) Amplification / Detection of amplified products.
- ◇ Use specific lab coat and gloves (powder free) in each working area.
- ◇ Pipettes, reagents and other materials must not cross each area.
- ◇ Specific caution is required to preserve the purity of the reagents and reaction mixtures.
- ◇ Appropriate methods of preparation/extraction of RNA to produce high quality RNA and to be followed by a real-time RT-PCR application should be used, particularly avoiding all sources of RNase contamination.
- ◇ Always use RNase-free DNase-free filtered tips for micropipettes.
- ◇ Do not pipette with mouth and do not eat, drink or smoke in the area.
- ◇ Avoid sprays.

SAMPLES COLLECTION, TRANSPORT AND STORAGE

- ◇ Collect samples in sterile tubes. **The use of heparin as an anticoagulant is prohibited.**
- ◇ It is the responsibility of the user to master its own conditions of collection, transport and storage of samples, and extraction of RNA by suitable systems to produce RNA of good quality.
- ◇ The samples should be extracted immediately or stored following the recommendations in the table below (Table 3).

Table 3 :

Recommendations of maximum storage of samples before extraction	
Room temperature	2 h
4°C	16h
-20°C (preferentially -80°C)	Long term storage

- ◇ The user can refer to the World Health Organization or High Health Authority for storing samples.
- ◇ Extracted RNAs have to be stored at -80°C.
- ◇ Transport of clinical samples must follow local regulations for this type of infectious agents.

PROCEDURE

I- RNA extraction

It is the user's responsibility that the extraction system used is compatible with downstream real time RT-PCR technology. For this kit, we recommend to use extraction methods of viral RNA for plasma, and refer to manufacturer's instructions of the extraction kit.

In the Dengue 1234 kit, CI-ARN on the CY5 channel can be added before extraction or in the PCR reaction. It ensures that a negative result is not due to an extraction problem or due to the presence of RT-PCR inhibitors at high quantity.

We recommend the addition of 10 µl of CI-ARN per extraction. If the CI-ARN is added to control the RT-PCR, CI-ARN is added to the reaction mix (1 µl per PCR reaction). See real-time RT-PCR protocol for details.

CI-ARN is available from Eurobio (Ref Eurobioplex EBX-003).


II- Real-time RT-PCR Procedure

General comment:

The positive control and the RNA extraction and RT-PCR inhibition control (CI-ARN) contain a very high concentration of matrix. Manipulations must be performed carefully to avoid contamination. To control the functioning of the PCR and steps of extraction and real-time RT-PCR amplification, it is necessary to test the CP-D14 and CP-D23 positive controls, as well as the negative controls (water supplied = CN-H₂O + CI-ARN), with the corresponding oligomix. See II-2/6) of real-time RT-PCR protocol for details.

II-1/ Protocol of real-time PCR

Diagram of the procedure

II-2/ Detailed Procedure

- 1) Homogenize the tube of Enzymes, and vortex Oligomixes 1 and 2, CP-D14 and CP-D23, and CI-ARN tubes before starting, and centrifuge.
- 2) Prepare the Mastermixes 1 and 2 as below. N is the number of PCR reactions. Plan to prepare enough reagents for at least N+3 reactions. (Refer to part 1-(a) or 1-(b) of the previous diagram according to the condition).

Case (a): For samples extracted without CI-ARN

Number of reaction(s)	1	N+3
Enzymes	15 µl	(N+3) x 15 µl
Oligomix 1 or 2	5 µl	(N+3) x 5 µl
CI-ARN	1 µl	(N+3) x 1 µl
Total Volume of Mastermix 1 or 2	21 µl*	(N+3) x 21 µl

Case (b): For samples extracted with CI-ARN

Number of reaction(s)	1	N+3
Enzymes	15 µl	(N+3) x 15 µl
Oligomix 1 or 2	5 µl	(N+3) x 5 µl
Total Volume of Mastermix 1 or 2	20 µl*	(N+3) x 20 µl

* The volume difference between condition (a) or (b) has no effect on performance.

- 3) Homogenize the Mastermixes prepared in 2) and centrifuge briefly.
- 4) Distribute 20 µL of the Mastermix* 1 or 2 with a micropipette and filtered tips in each tube or well of a microplate for real time PCR.
- 5) Add 5 µL of extracted RNA sample
- 6) In parallel test the following controls:

- Positive controls:

20 µl Mastermix 1 + 5 µl CP-D14

20 µl Mastermix 2 + 5 µl CP-D23

- Negative controls:

Case (a): For samples extracted without CI-ARN

- 20 µl Mastermix 1 + 5 µl water supplied (CN-H2O)
- 20 µl Mastermix 2 + 5 µl water supplied (CN-H2O)

Case (b): For samples extracted with CI-ARN

- 20 µl Mastermix 1 + 1 µl CI-ARN + 4 µl water supplied (CN-H2O)
- 20 µl Mastermix 2 + 1 µl CI-ARN + 4 µl water supplied (CN-H2O)

- 7) Close immediately with an adhesive film or transparent caps to avoid all contamination.
- 8) Centrifuge briefly to collect all the PCR reaction mix at the bottom of the tubes or plate.
- 9) Program the real-time PCR instrument as follows.

Program	Temperature	Duration	Cycle(s)	
Reverse Transcription	50 °C	10 min	1	-
Denaturation	95°C	3 min	1	-
Amplification	95°C	10 sec	45	-
	60°C	30 sec		Fluorescence Acquisition

Note 1: On LightCycler ® 480 systems (Roche), two optical systems are available: only "System II" is compliant with use of the kit. Apply color compensation for the following channels: FAM-HEX/VIC-Cy5 (465-510,533-580,618-660).

Note 2: For the Applied Biosystems systems, select "ROX" in "PASSIVE REFERENCE".


Note 3: On Rotorgene ™, please calibrate the signal by clicking on "GAIN optimization".

Note 4: On CFX96 (Biorad), start the run using the v1.6 or later version of CFX Manager software, and analyse with v 3.1 (see § Validation of the Experiment)


VALIDATION OF THE EXPERIMENT

The analysis of data post-acquisition on a CFX96 PCR instrument (Biorad) must be done with version 3.1 of CFX Manager software (Biorad). In order to use this v3.1 version from a run started with an older version, follow the procedure below: at the end of the run, the data file with .pcrd suffix must be open and treated with version 3.1 of CFX Manager (Biorad).

If the run was done with CFX Manager v1.6 for instance, to open the data file with CFX Manager v3.1, click on CFX Manager v3.1 icon. The screen below appears.


- Click on File and select Open, then Data File

- Select the file you want to analyze and click on Open.

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Once this is done, the analysis can start.

For the assay to be valid, the Ct values for the controls must be the following (Table 4). Otherwise, the experiment is not valid.

Table 4:

Positive Controls		
Channel	CP-D14	CP-D23
FAM	Ct ≤ 30	Ct ≤ 30
HEX	Ct ≤ 30	Ct ≤ 30
Negative Controls		
Channel	Negative Control Mastermix 1	Negative Control Mastermix 2
FAM	Ct ≥ 40 or undetermined	Ct ≥ 40 or undetermined
HEX	Ct ≥ 40 or undetermined	Ct ≥ 40 or undetermined
CY5	Ct ≤ 45	Ct ≤ 45

CP-D14: Positive control for Dengue, serotypes 1 and 4

CP-D23: Positive control for Dengue, serotypes 2 and 3

DATA ANALYSIS AND INTERPRETATION

RNA extraction and RT-PCR inhibition control in samples:

Two results can be obtained:

1/ the RNA extraction and RT-PCR inhibition control (CI-ARN; channel Cy5) test is positive: RNA has been properly extracted, and there is no inhibitor of RT-PCR. The result can be validated.

2/ the RNA extraction and RT-PCR inhibition control (CI-ARN; channel Cy5) test is negative: either RNA was not extracted, the reverse transcription did not work well, or the presence of PCR inhibitors inhibits the PCR reaction. It is then recommended to repeat the extraction or dilute the sample, unless a specific signal appears in the targeted channel (FAM, HEX).

For clinical samples, the following results are possible: example for serotypes 1 and 4 in the table below. The same principle applies to serotypes 2 (FAM) and 3 (HEX).

***Ct cut off for positivity of samples:**


On all channels (FAM or HEX) : + Positive => Positive Ct and < 40

PCR Signal			Presence of virus		Test validity/comment
FAM	HEX	CY5	DENGUE 1	DENGUE 4	
+	-	+	Yes	No	VALID
-	+	+	No	Yes	VALID
+	+	+	Yes	Yes	VALID
-	-	+	No	No	VALID
+	+	-	Yes	Yes	VALID Inhibition of RT-PCR or problem with extraction or strong competition that does not prevent the detection of virus.
+	-	-	Yes	No possible interpretation	Inhibition of RT-PCR or problem with extraction → dilute 5 x the sample; otherwise extract once again the sample
-	+	-	No possible interpretation	Yes	Inhibition of RT-PCR or problem with extraction → dilute 5 x the sample; otherwise extract once again the sample
-	-	-	No possible interpretation	No possible interpretation	Inhibition of RT-PCR or problem with extraction → dilute 5 x the sample; otherwise extract once again the sample


PERFORMANCE ANALYSIS

Example of experiment performed on real-time PCR thermocycleur CFX96 (Biorad):

Dengue serotype 1 (FAM) and serotype 4 (HEX)

Dengue serotype 2 (FAM) and serotype 3 (HEX)

Analytical sensitivity: each DENGUE serotype: 10 copies/ μ l CP

Linearity of quantification: each DENGUE serotype: 10 to 10⁹ copies/ μ l CP

Variability of the signal

Serotype	Coefficient of Variation Within experiment (%)	Coefficient of Variation Between experiments (%)
Dengue 1	2.02	3.32
Dengue 2	2.32	1.65
Dengue 3	0.90	5.19
Dengue 4	1.84	1.31

Sensitivity and Specificity

12 Dengue samples, characterized for their serotype, were analyzed and the following results were obtained.

- Sensitivity: > 99%
- Specificity: > 99%
- Concordance: > 99%

No cross-reactivity between serotype or with Chikungunya virus has been detected.

The extraction of the samples was performed on Qiagen columns: QIAamp Ultrasens Virus kit (Ref: 53704).

The tests were performed on CFX96 (Biorad).

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Callahan et al. Development and evaluation of serotype and group-specific fluorogenic reverse transcriptase PCR (TaqMan) assays for Dengue virus, J Clin. Microbiol, 2001;41:19-4124.

Leparc-Goffart I, Baragatti M, Temmam S, Tuiskunen A, Moureau G, Charrel R, de Lamballerie X. Development and validation of real-time one-step reverse transcription-PCR for the detection and typing of dengue viruses. J Clin Virol. 2009;45(1):61-6.

Cecilia D, Kakade M, Alagarasu K, Patil J, Salunke A, Parashar D, Shah PS. Development of a multiplex real-time RT-PCR assay for simultaneous detection of dengue and chikungunya viruses. Arch Virol. 2015;160(1):323-7.

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WASTE DISPOSAL

Be in accordance with the law on the elimination of waste of clinical infectious material.

SYMBOLS

	Reference
	Batch number
	Limits of storage temperature
	Expiration Date
	Content sufficient for « N » reactions
	Keep protected from light
	Manufacturer
	CE labeled product
	In Vitro Diagnostic
	Instructions for use

eurobio ingen

7, avenue de Scandinavie
91953 Les Ulis Cedex
FRANCE