EurobioPlex

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


EBX-020 olig2

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BIBLIOGRAPHIE

Van de Beek et al (2006). Community-acquired bacterial meningitis in adults in the Netherlands, 2006–14: a prospective cohort study. In *The Lancet Infectious Diseases*. March 2016 16(3):339-347.

Cohen-Bacrie S, Ninove L, Nougairé`de A, Charrel R, Richet H, et al. (2011) Revolutionizing Clinical Microbiology Laboratory Organization in Hospitals with In Situ Point-of-Care. PLoS ONE 6(7): e22403.

Hotomi M, Togawa A, Kono M, Sugita G, Sugita R, et al. (2013) An Application of Outer Membrane Protein P6-Specific Enzyme-Linked Immunosorbent Assay for Detection of Haemophilus influenzae in Middle Ear Fluids and Nasopharyngeal Secretions. PLoS ONE 8(8): e71774.

Brouwer MC, Tunkel AR, van de Beek D. (2010) Epidemiology, diagnosis, and antimicrobial treatment of acute bacterial meningitis. Clin Microbiol Rev. 2010 Jul;23(3):467-92.

Corless C.E, Guiver M, Borrow R, Edwards-Jones V, Fox A.J and Kaczmarek E.B. (2001). Simultaneous Detection of Neisseria meningitidis, Haemophilus influenzae, and Streptococcus pneumoniae in Suspected Cases of Meningitis and Septicemia Using Real-Time PCR. J. Clin. Microbiol. 2001, 39(4):1553.

<http://www.sante.ujf-grenoble.fr/SANTE/corpus/disciplines/pedia/malinfped/96/leconimprim.pdf>. BOST-BRU C. et PLANTAZ D.

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
ZA de Courtaboeuf
7, avenue de Scandinavie
91940 Les Ulis Cedex
FRANCE

EurobioPlex

*Community-acquired
bacterial meningitis*
REAL-TIME PCR

For **qualitative** real-time PCR

REF **EBX-020-36**
EBX-020-09

9 reactions
36 reactions

Version 5.00 of 2019/10/03

Validated on:

- CFX96™ Real Time PCR detection system (Biorad) with analysis on CFX Manager v 3.1 (Biorad)
- LightCycler®480 Instrument II (Roche) with analysis on LightCycler® 480 software v1.5 (Roche)
- Applied Biosystems® 7500 Real-Time PCR Systems (Applied Biosystems) with analysis on 7500 Software v2
- T-COR 8®-IVD (Tetracore Inc.) with T-COR 8 SmartCT™ (auto v1) software

Storage conditions:

Keep all reagents between -15°C and -22°C until use and after first use

Instructions for use

TABLE OF CONTENTS

Intended use	33
Introduction	33
Principle of detection	34
Description and content of the kit	35
Storage	36
Cautions and notes	36
Samples collection, transport and storage	37
Procedure	38
I- DNA extraction	38
II- Real-time PCR procedure	38
II-1/ Diagram of the procedure	39
II-2/ Detailed procedure	40
Validation of the experiment	41
Data analysis and Interpretation	43
Performance analysis.....	44
Specificities of the Real-Time PCR T-COR 8®-IVD	48
Barcodes for EBX-020 for use on T-COR8®-IVD.....	52
Bibliography	58
Waste disposal	58
Symbols	59

INTENDED USE

The EBX-020 test uses real-time polymerase chain reaction (PCR) amplification, and is designed for the qualitative detection from cerebrospinal fluid (CSF) of the three most frequently involved bacteria in the etiology of Community-acquired Bacterial Meningitis (MBcom): *Neisseria meningitidis* (meningococcus), *Streptococcus pneumoniae* (pneumococcus), and *Haemophilus influenzae*. **Real-time PCR can be carried out directly on CSF without extraction, which** can answer the need of diagnostic urgency in this pathology, but **also on extracted DNA**. This test is indicated to confirm a diagnosis of presumption of infection in patients or complete a proven or indeterminate diagnosis by other techniques (clinical diagnosis with semiological examination, bacteriological study and culture of CSF).

The Eurobioplex EBX-020 has been validated on the following specimen:

- Cerebrospinal fluid (CSF)

INTRODUCTION

Meningitis is an inflammation of meninges surrounding the central nervous system tissue, and is due to many bacterial or viral pathogens.

Bacterial meningitis affects 22 cases per million people, about 1400 cases/year in France. Two peaks of incidence are observed: one in less than one year old children and the other affects adults over age 60. Bacterial meningitis is rare (20% of cases) but is extremely serious compared to viral meningitis. Maternal-fetal transmission is most often the source of infection in the newborn. As for bacterial meningitis in adults, a substantial part would be due to nosocomial infections. Bacterial meningitis is thus often subdivided in three subcategories: community-acquired bacterial meningitis, maternal-fetal bacterial meningitis, and nosocomial bacterial meningitis.

Community-acquired bacterial meningitis are worldwide responsible for 170 000 annual deaths (WHO), and more than a million people would be affected. The epidemiological characteristics vary according to climatic zones and countries. Meningococcal outbreaks are mainly found in sub-Saharan Africa at the band of the Sahel called "meningitis belt". Endemic areas are found in some countries of South America (Brazil, Ecuador), Africa (Western Sahara, Morocco, Egypt), and Asia (Saudi Arabia, India and China). All countries are concerned for endemo-sporadic areas. Mortality is particularly high in developing countries. 20% of patients with community-acquired bacterial meningitis die and 30% suffer from sequelae. Patients' life is at stake if targeted antibiotic treatment is delayed. Whether bacterial meningitis is proven or suspected, it remains an absolute diagnostic and therapeutic emergency.

There are four main sources of surveillance in France: 1) the EPIBAC/InVS network which comprises approximately 300 hospital laboratories, 2) the national observatory for childhood bacterial meningitis: GPIP-Activ which includes approximately 250 pediatric services and 150 microbiology laboratories, 3) the regional pneumococcal observatories, 4) the national reference centers.

There is no obvious diagnostic assessment by the clinician between viral and bacterial meningitis because the symptoms are common: fever, headache, stiffness in the neck, vomiting, limited elevation of the lower limbs, and involuntary flexion of lower limbs as soon as the neck is forced to bend. The ORL entry ways are the most common. Only the cyto-bacteriological examination of the cerebrospinal fluid (CSF) enables to make

the diagnosis. Culture of CSF is rapidly and systematically done. A fraction of the sampling of CSF is separately kept in anticipation of microbiological research by polymerase chain reaction (PCR). Count of leukocytes, GRAM stain, and detection of soluble antigens can orient the diagnosis but the sensitivity of these techniques is strongly altered if the patient is already treated by antibiotic therapy.

The main bacteria involved in community-acquired bacterial meningitis (in adults) are in order of frequency: *Streptococcus pneumoniae* (Pneumococcus), *Neisseria meningitidis* (meningococcus, mandatory reporting), and *Haemophilus influenzae*. The pneumococcal and meningococcal bacteria are responsible for more than 80% of cases. The rapid and targeted diagnosis by polymerase chain reaction in real-time from the DNA of bacterial pathogens of community-acquired bacterial meningitis is the technique with the best sensitivity and specificity.

PRINCIPLE OF DETECTION

The Eurobioplex EBX-020 is a test using real-time PCR amplification of bacterial DNA of *Neisseria meningitidis*, *Streptococcus pneumoniae* and *Haemophilus influenzae*, as well as a DNA extraction and PCR inhibition control. The test is performed directly from non-extracted samples of CSF or from extracted DNA using two reactions in two wells/tubes.

The DNA extraction and PCR inhibition control allows to check for variations that may occur during the DNA extraction step of biological samples and real-time PCR amplification. Thus, it ensures that a negative result may not be due to a lack of availability of DNA in CSF, and/or a bad DNA extraction, and/or the presence of PCR inhibitors in too large quantities.

This kit contains two different oligomixes which permit to detect the different pathogens listed below:

Oligomix 1 : *Neisseria meningitidis* (NM), 2 genes + control of DNA extraction and PCR inhibition (CI-PCR),

Oligomix 2 : *Streptococcus pneumoniae* (SP) and *Haemophilus influenzae* (HI) + control of DNA extraction and PCR inhibition (CI-PCR).

For the first one, two genes of NM are respectively detected using FAM and HEX labeled probes. For the second triplex, DNA of SP and HI is also detected using respectively FAM and HEX labeled probes. For the two triplexes, DNA extraction and PCR inhibition control is detected using a CY5 labeled probe. All probes emit a specific fluorescence following their hydrolysis during the elongation of the amplification product. The measurement of the intensities of real-time fluorescence correlates with the accumulation of specific amplification products.

DESCRIPTION AND CONTENT OF THE KIT

The real-time PCR EBX-020 kit is ready to use for the specific detection of *Neisseria meningitidis*, *Streptococcus pneumoniae*, and *Haemophilus influenzae*.

Fluorescence is emitted and individually recorded through optical measurements during the PCR. The detection of the amplified fragment is performed by a fluorimeter using the channels shown in Table 1 below.

The kit contains reagents and enzyme for the amplification of DNA of these bacteria, and the DNA extraction and PCR inhibition control (see Table 2).

Table 1: Pathogens detected by each oligomix

Target	Oligomix 1	Oligomix 2	Fluorophore
<i>Neisseria meningitidis</i> – Gene 1	X	-	FAM
<i>Neisseria meningitidis</i> – Gene 2	X	-	HEX
<i>Streptococcus pneumoniae</i>	-	X	FAM
<i>Haemophilus influenzae</i>	-	X	HEX
PCR inhibition control (CI-PCR)	X	X	Cy5

Wavelengths of excitation/emission for FAM (495/515 nm), HEX (535/555 nm), CY5 (647/667 nm).

Equivalent channels on different real time PCR cyclers:

- Channel **FAM** (ABI Systems, SmartCycler II, Chromo 4/CFX96, Systems Mx, T-COR8®-IVD), Channel 510 (LC 480), Channel Green (RotorGene)
- Channel **HEX** (Chromo 4/CFX96, Mx Systems, T-COR8®-IVD), Channel VIC (ABI Systems), Channel Alexa532 (SmartCycler II), Channel 580 (LC 480), Channel Yellow (RotorGene),
- Channel **Cy5** (ABI Systems, Chromo 4/CFX96, Mx Systems, T-COR8®-IVD), Channel 660 (LC 480), Channel Alexa647 (SmartCyclerII), Channel Red (RotorGene)

Note: On LC480 instrument II, apply color compensation for the following channels: FAM-HEX/VIC -Cy5 (465-510, 533-580-618-660).

Table 2:

Cap color	Components of the kit	9 reactions	36 reactions	Reconstitution
Red	Enzyme	4 x 100 µl	16 x 100 µl	Ready to use
Transparent	Oligomix 1	2 x 60 µl	8 x 60 µl	Ready to use
Green	Oligomix 2	2 x 60 µl	8 x 60 µl	Ready to use
Blue	Water = negative control (CN-H ₂ O)	2 x 500 µl	8 x 500 µl	Ready to use
Yellow	Positive Control CP-MBcom1	2 x 15 µl	8 x 15 µl	Ready to use
Orange	Positive Control CP-MBcom2	2 x 15 µl	8 x 15 µl	Ready to use
White	PCR inhibition control (CI-PCR)	2 x 50 µl	8 x 50 µl	Ready to use

Oligomix 1: contains the primers and probes of the first triplex (2 genes of NM, and CI-PCR).
Oligomix 2: contains the primers and probes of the second triplex (SP, HI, and CI-PCR).
CP-MBcom1: positive control of the community-acquired bacterial meningitis for 2 genes of NM, triplex1
CP-MBcom2: positive control of the community-acquired bacterial meningitis for SP et HI, triplex2

Table 3: Number of patients tests according to the EBX format (9 or 36 reactions) on all instruments except T-COR8®-IVD (see page 48)

	9 Reactions	36 Reactions
Number of possible patients tests patient-by-patient with maximum 3 defrost/refreezing cycles	6 patients	24 patients
Maximum number of patients tests in 1 run without refreezing cycles	7 patients	34 patients

Required material not provided:

- ◇ Biological Hood
- ◇ Real-time PCR instrument
- ◇ Micro centrifuge
- ◇ Vortex
- ◇ Plates / tubes for real-time PCR
- ◇ Micropipettes
- ◇ DNase-free RNase-free filter tips for micropipettes
- ◇ Sterile microtubes
- ◇ Gloves (powder free)

STORAGE

All reagents must be stored between -15 and -22°C.
All reagents can be used until the expiration date indicated on the kit.
Many freezing / defrosting cycles (> 3x) must be avoided, and could lead to decrease in sensitivity.

CAUTIONS AND NOTES

Read carefully instructions before starting.

- ◇ The experiment must be performed by competent staff.
- ◇ Instruments must have been properly installed, calibrated and maintained according to the manufacturer's recommendations.
- ◇ Clinical samples are potentially infectious and must be processed under a laminar flow hood.
- ◇ The experiment must be performed according to good laboratory practices.
- ◇ Do not use this kit after expiration date.
- ◇ The kit is shipped with dry ice, and the components of the kits must arrive frozen. If one or more components are defrosted, or of the tubes have been damaged, contact Eurobio Scientific.
- ◇ Once defrosted, spin down briefly the tubes before use.

- ◇ It is recommended to define three working areas: 1) Isolation of DNA, 2) Preparation of the reaction mix and 3) Amplification / Detection of amplified products.
- ◇ Use specific lab coat and gloves (powder free) in each working area.
- ◇ Pipettes, reagents and other materials must not cross each area.
- ◇ Specific caution is required to preserve the purity of the reagents and reaction mixtures.
- ◇ Appropriate methods of preparation/extraction of DNA to produce high quality DNA and to be followed by a real-time PCR application should be used, particularly avoiding all sources of DNase contamination.
- ◇ Always use RNase-free DNase-free filtered tips for micropipettes.
- ◇ Do not pipette with mouth and do not eat, drink or smoke in the area.
- ◇ Avoid sprays.

SAMPLES COLLECTION, TRANSPORT AND STORAGE

- ◇ Collect samples in sterile tubes.
- ◇ It is the responsibility of the user to master its own conditions of collection, transport and storage of samples, and if an extraction of DNA is performed then use suitable systems to produce DNA of good quality.
- ◇ Lumbar puncture is performed after a rigorous asepsis of the area and a local anesthesia. CSF is collected by puncture at the lumbar vertebrae level of the L3-L4, L4-L5 or L5-S1 interspace. Once collected, the CSF is immediately brought to the laboratory for analysis.
- ◇ A CSF showing the presence of blood can also give a negative result. In this case, it may be wise to perform the test after extraction of DNA and not on direct LCR. This can also be done in a second attempt, if the result is negative on direct LCR.
- ◇ It is recommended that samples be stored according to the recommendations of storage of samples (Table 4).

Table 4 :

Recommendations of maximum storage of samples before extraction	
Room temperature	2 h
4°C	16h
-20°C	Long term storage

- ◇ The user can refer to the World Health Organization or High Health Authority for storing samples.
- ◇ Transport of clinical samples must obey local regulations for this type of infectious agents.

PROCEDURE

I- DNA Extraction

This part is optional because this kit can be used directly on CSF (without extraction of bacterial DNA).

It is the user's responsibility that the extraction system used be compatible with downstream real time PCR technology. For this kit, we recommend to use extraction methods of bacterial DNA from CSF samples, and refer to manufacturer's instructions of the extraction kit.

In the EBX-020 kit, CI-PCR on the CY5 channel can be added into the mastermix in case of direct CSF testing, or before extraction if one is performed. It ensures that a negative result is not due to the presence of PCR inhibitors at high quantity, or potentially due to a problem with extraction if it has been performed.

If an extraction is performed, we recommend the addition of 10 µl of CI-PCR to the biological sample for a final volume of elution of 50 µl after extraction.

If the CI-PCR is added to control the real-time PCR or in the case of direct CSF testing, CI-PCR is added to the reaction mix (1 µl per PCR reaction).

See real-time PCR protocol for details.

CI-PCR is also available separately from Eurobio (Ref EurobioPlex EBX-002).


II- Real-time PCR procedure

General comment:

Positive controls and the DNA extraction and PCR inhibition control (CI-PCR) contain a very high concentration of matrix. Manipulations must be performed carefully to avoid contamination. To control the functioning of the PCR and steps of extraction and real-time PCR amplification, it is necessary to test the positive controls CP-MBcom1 and CP-MBcom2 with their corresponding triplex 1 and 2 respectively, as well as the negative controls (water supplied = CN-H2O + CI-PCR) (see II-2/6) of real-time PCR procedure).

On T-COR 8®-IVD, we recommend that the positive and negative controls be tested when a new kit is being opened.

II-1/ Diagram of the procedure

II-2/ Detailed procedure: Proceed in the same way for the detection of NM with the oligomix 1, and for the detection of SP and HI with the oligomix 2.

- 1) Homogenize the tube of Enzyme, and vortex Oligomix 1, Oligomix 2, CP-MBcom1, CP-MBcom2, and CI-PCR tubes before starting, and centrifuge.
- 2) Prepare Mastermixes 1 and 2 as below; prepare the reaction mixes 1 and 2 as follows by multiplying by the number of samples N to test (including positive and negative controls). On average, prepare enough reagents depending on the size of the sampling for N+2 reactions (up to 3 patients tested) or N+3 reactions (beyond 3 patients tested).
(Refer to part 1-(a) or 1-(b) of the previous diagram according to the condition).

Case (a): for CSF direct testing, or for DNA extracted without CI-PCR

Number of reaction(s)	1	N+3	Example for 3 patients N = 7
Enzyme	12.5 µl	(N+3) x 12.5 µl	87.5 µl
Oligomix 1 or 2	7.5 µl	(N+3) x 7.5 µl	52.5 µl
CI-PCR	1 µl	(N+3) x 1 µl	7 µl
Total volume mastermix 1 or 2	21 µl*	(N+3) x 21 µl	147 µl

* The volume difference between condition (a) or (b) has no effect on performance.

Case (b): for DNA extracted with CI-PCR

Number of reaction(s)	1	N+3	Example for 3 patients N = 7
Enzyme	12.5 µl	(N+3) x 12.5 µl	87.5 µl
Oligomix 1 or 2	7.5 µl	(N+3) x 7.5 µl	52.5 µl
Total volume mastermix 1 or 2	20 µl*	(N+3) x 20 µl	140 µl

- 3) Homogenize the mastermixes prepared in 2) and centrifuge briefly
- 4) Distribute 20 µL Mastermix* 1 or 2 using a micropipette and filtered tips in each tube or well of microplate for real time PCR.
- 5) Add 5 µL of CSF or extracted DNA sample.
- 6) In parallel, test the following controls (see this specific point for T-COR 8®-IVD on page 48-57)

- Positive controls:

- 20µL of mastermix 1 + 5µL of CP-MBcom1
- 20µL of mastermix 2 + 5µL of CP-MBcom2

- Negative controls:

○ Case (a): CSF direct testing or for DNA extracted without CI-PCR

- 20 µL mastermix 1 + 5 µL supplied water (CN-H2O)
- 20 µL mastermix 2 + 5 µL supplied water (CN-H2O)

○ Case (b) : for DNA extracted with CI-PCR

- 20 µL mastermix 1 + 4 µL water supplied (CN-H2O) + 1 µL CI-PCR
- 20 µL mastermix 2 + 4 µL water supplied (CN-H2O)+ 1 µL CI-PCR

- 7) Close immediately with an adhesive film or transparent caps to avoid contamination.
- 8) Centrifuge briefly to collect all the PCR reaction mix at the bottom of the tubes or plate.
- 9) Program the real-time PCR instrument as follows (on T-COR 8®-IVD, no programming is needed thanks to the use of Barcodes- see page 48-57):

Program	Temperature	Duration	Cycle(s)	
Heating	95°C	10 min	1	-
Denaturation	95°C	3 min	1	-
Amplification	95°C	10 sec	45	-
	60°C	30 sec		Fluorescence acquisition

Note 1: On LightCycler ® 480 systems (Roche), two optical systems are available: only "System II" is compliant with use of the kit. Apply color compensation for FAM-HEX/VIC- Cy5 (465-510, 533-580, 618-660).

Note 2: For the Applied Biosystems systems, select "ROX" in "PASSIVE REFERENCE".


Note 3: On Rotorgene™, please calibrate the signal by clicking on "GAIN optimization".

Note 4: On CFX96 (Biorad), start the run using the v1.6 or later version of CFX Manager software, and analyse with v 3.1 (see § Validation of the Experiment)


VALIDATION OF THE EXPERIMENT

The analysis of data post-acquisition on a CFX96 PCR instrument (Biorad) must be done with version 3.1 of CFX Manager Software (Biorad). In order to use this v3.1 version from a run started with an older version, follow the procedure below: at the end of the run, the data file with .pcrd suffix must be open and treated with version 3.1 of CFX Manager (Biorad).

If the run was done with CFX Manager v1.6 for instance, to open the data file with CFX Manager v3.1, click on CFX Manager v3.1 icon. The screen below appears.


- Click on File and select Open, then Data File

- Select the file you want to analyze and click on Open.

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Once this is done, the analysis can start.

For the assay to be valid, the results for the controls must be the following (Table 5). Otherwise, the experiment is not valid. On T-COR 8®-IVD, we recommend that the positive and negatives controls be tested when a new kit is being opened.

Table 5:

Positive controls		
Channel	CP-MBcom1	CP-MBcom2
FAM	Ct ≤ 22	Ct ≤ 21
HEX	Ct ≤ 22	Ct ≤ 21
Negative controls		
Channel	Negative control Mastermix 1	Negative control Mastermix 2
FAM	Ct not determinated	Ct not determinated
HEX	Ct not determinated	Ct not determinated
CY5	Ct ≤ 35	Ct ≤ 35

CP-MBcom1 and 2: positive controls of the community-acquired bacterial meningitis (Triplex1 and Triplex2)

DATA ANALYSIS AND INTERPRETATION

DNA extraction and PCR inhibition control in samples:

Two results can be obtained:

1/ the DNA extraction and PCR inhibition control (CI-PCR; channel Cy5) test is positive: the result can be validated.

2/ the DNA extraction and PCR inhibition control (CI-PCR; channel Cy5) test is negative: either bacterial DNA is not reachable within the CSF, or DNA was not extracted (in case of extraction), the PCR did not work well, or the presence of PCR inhibitors inhibits the PCR reaction. A CSF showing the presence of blood can also give a negative result.

It is then recommended to perform an extraction (if CSF was initially tested directly), or to repeat the extraction, or dilute the sample, unless a specific signal appears in all other channels (FAM and HEX).

OLIGOMIX 1:

For clinical samples tested with Oligomix 1, the following results are possible:

*** Ct cut off for samples:**

Channel FAM-NM gene 1 and HEX-NM gene 2: + Positive => Positive Ct (≤ 45)

PCR Signal			Presence of NM	Test validity/comment
FAM	HEX	CY5		
+	+	+	Yes	valid
+	-	+	Yes	valid
-	+	+	Yes	valid
-	-	+	No	valid
+	-	-	Yes	Valid: possible problem with PCR inhibition, or with extraction if performed, which does not prevent the detection of the bacterium.
-	+	-	Yes	Valid: possible problem with PCR inhibition, or with extraction if performed, which does not prevent the detection of the bacterium.
-	-	-	Not interpretable	Not valid: possible problem with PCR inhibition, or with extraction if performed - dilute first 5 x the sample and repeat PCR; if necessary redo an extraction.

NM: *Neisseria meningitidis*

OLIGOMIX 2:

For clinical samples tested with Oligomix 2, the following results are possible:

* Ct cut off for samples:

Channel FAM-SP and HEX-HI: + Positive => Positive Ct (≤ 45)


PCR Signal			Presence of SP	Presence of HI	Test validity/comment
FAM	HEX	CY5			
+	+	+	Yes	Yes	valid
+	-	+	Yes	No	valid
-	+	+	No	Yes	valid
-	-	+	No	No	valid
+	-	-	Yes	Not interpretable	Valid for SP: possible problem with PCR inhibition, or with extraction if performed, which does not prevent the detection of SP. BUT that could prevent the detection of HI -dilute 5 x the sample; not valid for HI.
-	+	-	Not interpretable	Yes	Valid for HI: possible problem with PCR inhibition, or with extraction if performed, which does not prevent the detection of HI. BUT that could prevent the detection of SP -dilute 5 x the sample; not valid for SP.
-	-	-	Not interpretable	Not interpretable	Not valid: possible problem with PCR inhibition, or with extraction if performed - dilute first 5 x the sample and repeat PCR; if necessary redo an extraction.

SP: *Streptococcus pneumoniae*; HI : *Haemophilus influenzae*

PERFORMANCE ANALYSIS

Example of experiment performed on real-time PCR thermocycler CFX96 (Biorad):


- With Oligomix 1:

Analytical sensitivity: CP-MBcom1: 10 copies/μl

Linearity for quantification: CP-MBcom1 from 10 copies/μl to 10⁹ copies/μl

- With Oligomix 2:

— FAM : *Streptococcus pneumoniae*

— HEX : *Haemophilus influenzae*

Analytical sensitivity: CP-MBcom2: 10 copies/μl

Linearity for quantification: CP-MBcom2 from 10 copies/μl to 10⁹ copies/μl

Coefficient of correlation and efficiency

EBX-020 Triplex 1	EBX-020 Triplex 2
Bacteria and abbreviations: <i>Neisseria meningitidis</i> (NM) 2 genes	Bacteria and abbreviations: <i>Streptococcus pneumoniae</i> (SP) and <i>Haemophilus influenzae</i> (HI)
Efficiency FAM: NM gene 1: 98,2 % HEX: NM gene 2: 101,9 % Coefficient of correlation FAM: NM gene 1: 0.996 HEX: NM gene 2: 0.995	Efficiency FAM: SP: 100,3 % HEX: HI: 96,3 % Coefficient of correlation FAM: SP: 0.993 HEX: HI: 0.994

Variability of the signal

EBX-020		INTRA-experiment	INTER-experiments
	Bacterial targets	CV %	CV %
Triplex 1	<i>Neisseria meningitidis</i> Gene 1 (FAM)	0.98	1.87
	<i>Neisseria meningitidis</i> Gene 2 (HEX)	0.61	2.01
Triplex 2	<i>Streptococcus pneumoniae</i> (FAM)	0.62	2.65
	<i>Haemophilus influenzae</i> (HEX)	0.54	3.52

Validation on clinical samples

Specificity and sensitivity were analysed from samples of CSF (DNA not extracted) and bacterial strains (HI positive strains diluted in HI-negative CSF) previously tested by the reference laboratory.

TRIPLEX 1: EBX-020 Com1

1/ Analytical specificity

Neisseria meningitidis gene 1: 95.5 %*

Neisseria meningitidis gene 2: 95.5 %*

2/ Analytical sensitivity

Neisseria meningitidis gene 1: > 98 %

Neisseria meningitidis gene 2: > 98 %

3/ Concordance

Neisseria meningitidis gene 1: 96 % (n = 50)

Neisseria meningitidis gene 2: 96 % (n = 50)

TRIPLEX 2: EBX-020 Com2

1/ Analytical specificity

Streptococcus pneumoniae: 95.1 %*

Haemophilus influenzae: > 98 %

2/ Analytical sensitivity

Streptococcus pneumoniae: > 98 %

Haemophilus influenzae: > 98 %

3/ Concordance

Streptococcus pneumoniae: 95.7 % (n = 46)

Haemophilus influenzae: > 98 % (n = 14)

* Note that samples not detected by the reference method but detected by the Eurobio method are at the limit of detection. It's likely that the EBX-020 Eurobioplex is more sensitive.

Study of potentially interfering microorganisms: none of the bacteria listed below was detected by the EBX-020 kit (n=1 sample of each bacterium).

No cross-reactivity of the Triplex 1 with:

- *Streptococcus pneumoniae*
- *Haemophilus influenzae*
- *Staphylococcus aureus*
- *Pseudomonas aeruginosa*
- *Listeria monocytogenes*
- *Escherichia coli*
- *Streptococcus agalactiae*
- *Bartonella spp*
- *Coxiella burnetii*

- *Mycobacterium avium*
- *Streptococcus oralis mitis sanguis*
- *Francisella tularensis*
- *Mycoplasma*
- *Klebsiella kingae*
- *Mycobacterium tuberculosis*
- *Tropheryma whipplei*

No cross-reactivity of the Triplex 2 with:

- *Neisseria meningitidis*
- *Staphylococcus aureus*
- *Listeria monocytogenes*
- *Escherichia coli*
- *Streptococcus agalactiae*
- *Bartonella spp*
- *Coxiella burnetii*
- *Mycobacterium avium*
- *Mycoplasma*
- *Klebsiella kingae*
- *Mycobacterium tuberculosis*
- *Tropheryma whipplei*

The tests were performed on CFX96.

SPECIFICITIES OF THE REAL-TIME PCR T-COR 8®-IVD INSTRUMENT

T-COR 8®-IVD is a real-time PCR instrument with 8 independent wells, which can be programmed independently, patient by patient, regarding thermic protocol, assays, and time of start of the run.

The 8 wells can also be used simultaneously with the same assay.

Note: Vortex the tubes with the provided T-core vortex before putting the tubes in the T-COR 8®-IVD, and always check that there is no bubble, and that the liquid is all located at the bottom of the tube.

Controls

On T-COR 8®-IVD, we recommend that the positive and negative controls be tested when a new kit is being opened.

After this initial control, the internal control allows to check that the PCR is working properly, and that there is no PCR inhibition.

Patients' Tests

Patient's tests require 2 wells, one for each oligomix tested.

<i>With validation of positive and negative controls once, during first use of the kit.</i>	9 reactions	36 reactions
Number of possible patients tests, patient by patient, with a maximum of 3 freezing/thawing cycles	7 patients	34 patients

Fluorophores

Four combinations of fluorophores are available on T-COR 8®-IVD.

For all EBX, such as EBX-020, FAM/HEX/Texas Red/Cy5 combination is used. If some channels are not used (Texas Red for EBX-020), do not consider this channel for results analysis.

Use of Barcodes (available on page 52 to 57)

1- Select Menu > New Run

2- Select Barcode

3- Scan the appropriate Barcode on the right side of the instrument:

- For positive control 1 or 2 (Barcode EBX-020 CP1 or EBX-020 CP2),

- For negative control 1 or 2 (Barcode EBX-020 CN1 or EBX-020 CN2),
- For a patient sample tested with oligomix 1 (Barcode EBX-020 olig1) or oligomix 2 (Barcode EBX-020 olig2)

The instrument randomly selects one of the 8 available wells. Follow the instructions on the instrument.

4- Lift up the lid of the selected well (well x), check that the blue LED light is on, and select “Yes”

5- Place the tube in the corresponding well and select “Next”.

6- If you wish to name the sample (optional), select “sample x”, name it and select “Accept”.

7- Select “Next”

8- To add/test another control or sample, select « Add well » and go back to point 3-.

9- Select “Start run”.

***Note:** The name of a sample or of a run cannot be modified or added after the run is finished. It must be done before or during the run.*

Presentation and Interpretation of results

Ct values and amplification curves are available in the SmartCT™ Table (Ct values on each channel), and “Ct versus PCR cycles” graphs, both being available in real-time while the run is ongoing.

An automatic analysis for positive and negative controls as well as patients’ samples is available in “Interpretations”, at the end of the run, in the “View” window.

For positive and negative controls, the following results are possible:

« Neg Ctrl 1 (for oligomix 1) or Neg Ctrl 2 (for oligomix 2) Fail »: Not valid

« Neg Ctrl 1 (for oligomix 1) or Neg Ctrl 2 (for oligomix 2) Valid »: Valid

« Pos Ctrl 1 (for oligomix 1) or Pos Ctrl 2 (for oligomix 2) Fail »: Non valid

« Pos Ctrl 1 (for oligomix 1) or Pos Ctrl 2 (for oligomix 2) Valid »: Valid

For patients status:

"Detected": Positive for at least one target * → green box

"Not detected": Negative → red box

* For each oligomix, the results for the specific bacteria or genes are specified.

CAUTION !

When the box is green, it is important to read the status for each target as some targets may be negative.

« Invalid »: Invalid result -> retest

Example of automatic interpretation of results on T-COR 8®-IVD instrument screen and associated report:

- For valid positive and negative controls, for the limit of detection (LOD) indicated as « Detected » for its corresponding targets, and for a negative sample for both oligomix, indicated “Not Detected” for its corresponding targets:




❖ On the instrument screen:

Interpretations		
1	CN2	EBX-020 CN2
Detected: Neg Ctrl 2 Valid		
2	LOD OLIGO1	EBX-020 olig1
Detected: N. men gene 2, N. men gene 1		
3	Sample 1	EBX-020 olig1
Not Detected - N. men		
4	LOD OLIGO2	EBX-020 olig2
Detected: H. inf, S. pneumo		
5	Sample 1	EBX-020 olig2
Not Detected - S. pneumo, H. inf		
6	CP1	EBX-020 CP1
Detected: Pos Ctrl 1 Valid		
7	CP2	EBX-020 CP2
Detected: Pos Ctrl 2 Valid		
8	CN1	EBX-020 CN1
Detected: Neg Ctrl 1 Valid		

❖ On the report:

Summary							
Well	Sample	Assay	FAM	HEX	TxR	Cy5	Call
1	CN2	EBX-020 CN2				23.9	Detected • Neg Ctrl 2 Valid
2	LOD OLIGO1	EBX-020 olig1	31.9	33.3		23.4	Detected • N. men gene 2 • N. men gene 1
3	Sample 1	EBX-020 olig1				23.2	Not Detected
4	LOD OLIGO2	EBX-020 olig2	31.2	30.9		23.6	Detected • H. inf • S. pneumo
5	Sample 1	EBX-020 olig2				23.7	Not Detected
6	CP1	EBX-020 CP1	17.2	19.0		22.1	Detected • Pos Ctrl 1 Valid
7	CP2	EBX-020 CP2	15.8	15.3		23.4	Detected • Pos Ctrl 2 Valid
8	CN1	EBX-020 CN1				23.4	Detected • Neg Ctrl 1 Valid


Example of amplification curves:

Barcodes for EBX-020 for use on T-COR8®-IVD


EBX-020 CN1

**Water = Negative control (CN-H2O)
(to test with oligomix 1)**


EBX-020 CN2

**Water = Negative control (CN-H₂O)
(to test with oligomix 2)**


EBX-020 CP1

**Positive control CP-MBcom1
(to test with oligomix 1)**


EBX-020 CP2

**Positive control CP-MBcom2
(to test with oligomix 2)**


EBX-020 olig1

(to test with oligomix 1)

EBX-020 olig2

(to test with oligomix 2)

BIBLIOGRAPHY

Van de Beek et al (2006). Community-acquired bacterial meningitis in adults in the Netherlands, 2006–14: a prospective cohort study. In *The Lancet Infectious Diseases*. March 2016 16(3):339-347.

Cohen-Bacrie S, Ninove L, Nougairé de A, Charrel R, Richet H, et al. (2011) Revolutionizing Clinical Microbiology Laboratory Organization in Hospitals with In Situ Point-of-Care. *PLoS ONE* 6(7): e22403.

Hotomi M, Togawa A, Kono M, Sugita G, Sugita R, et al. (2013) An Application of Outer Membrane Protein P6-Specific Enzyme-Linked Immunosorbent Assay for Detection of *Haemophilus influenzae* in Middle Ear Fluids and Nasopharyngeal Secretions. *PLoS ONE* 8(8): e71774.

Brouwer MC, Tunkel AR, van de Beek D. (2010) Epidemiology, diagnosis, and antimicrobial treatment of acute bacterial meningitis. *Clin Microbiol Rev*. 2010 Jul;23(3):467-92.

Corless C.E, Guiver M, Borrow R, Edwards-Jones V, Fox A.J and Kaczmarek E.B. (2001). Simultaneous Detection of *Neisseria meningitidis*, *Haemophilus influenzae*, and *Streptococcus pneumoniae* in Suspected Cases of Meningitis and Septicemia Using Real-Time PCR. *J. Clin. Microbiol*. 2001, 39(4):1553.

<http://www.sante.ujf-grenoble.fr/SANTE/corpus/disciplines/pedia/malinfped/96/leconimprim.pdf>. BOST-BRU C. et PLANTAZ D.

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WASTE DISPOSAL

Be in accordance with the law on the elimination of waste of clinical infectious material.

SYMBOLS


	Reference
	Batch number
	Limits of storage temperature
	Expiration Date
	Content sufficient for « N » reactions
	Keep protected from light
	Manufacturer
	Instructions for use
	CE labeled product
	In vitro diagnostic

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