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

Version 4.00 of 2019/10/03

**Validated on:**

- CFX96™ Real Time PCR detection system (Biorad) with analysis on CFX Manager v 3.1 (Biorad)
- LightCycler®480 Instrument II (Roche) with analysis on LightCycler® 480 software v1.5 (Roche)
- Applied Biosystems® 7500 Real-Time PCR Systems (Applied Biosystems) with analysis on 7500 Software v2

**Storage conditions:**

Keep all reagents between -15°C and -22°C until use and after first use

Instructions for use

|                          |
|--------------------------|
| <b>TABLE OF CONTENTS</b> |
|--------------------------|

|                                                        |           |
|--------------------------------------------------------|-----------|
| <b>Intended use .....</b>                              | <b>20</b> |
| <b>Introduction .....</b>                              | <b>20</b> |
| <b>Principle of detection .....</b>                    | <b>21</b> |
| <b>Description and content of the kit .....</b>        | <b>22</b> |
| <b>Storage .....</b>                                   | <b>23</b> |
| <b>Cautions and notes .....</b>                        | <b>23</b> |
| <b>Samples collection, transport and storage .....</b> | <b>24</b> |
| <b>Procedure .....</b>                                 | <b>25</b> |
| <b>I- DNA extraction .....</b>                         | <b>25</b> |
| <b>II- Real-time PCR procedure .....</b>               | <b>25</b> |
| II-1/ Diagram of the procedure .....                   | 26        |
| II-2/ Detailed procedure .....                         | 27        |
| <b>Validation of the experiment .....</b>              | <b>28</b> |
| <b>Data analysis and Interpretation .....</b>          | <b>29</b> |
| <b>Performance analysis.....</b>                       | <b>30</b> |
| <b>Bibliography .....</b>                              | <b>32</b> |
| <b>Waste disposal .....</b>                            | <b>32</b> |
| <b>Symbols .....</b>                                   | <b>33</b> |

## INTENDED USE

The EBX-021 test uses real-time polymerase chain reaction (PCR) amplification and is designed for the qualitative detection from cerebrospinal fluid (CSF) of the two most frequently involved bacteria in the etiology of nosocomial Bacterial Meningitis (MBnoso): *Staphylococcus aureus* and *Pseudomonas aeruginosa*. **Real-time PCR can be carried out directly on CSF without extraction, which** can answer the need of diagnostic urgency in this pathology, but **also on extracted DNA**. This test is indicated to confirm a diagnosis of presumption of infection in patients or complete a proven or indeterminate diagnosis by other techniques (clinical diagnosis with semiological examination, bacteriological study and culture of CSF).

The Eurobioplex EBX-021 has been validated on the following specimen:

- Cerebrospinal fluid (CSF)

## INTRODUCTION

Meningitis is an inflammation of meninges surrounding the central nervous system tissue, and is due to many bacterial or viral pathogens.

Bacterial meningitis affects 22 cases per million people, about 1400 cases/year in France. There are four main sources of surveillance in France: 1) the EPIBAC/InVS network which comprises approximately 300 hospital laboratories, 2) the national observatory for childhood bacterial meningitis: GPIP-Activ which includes approximately 250 pediatric services and 150 microbiology laboratories, 3) the regional pneumococcal observatories, 4) the national reference centers.

Two peaks of incidence are observed: one in less than one year old children and the other affects adults over age 60. Bacterial meningitis is rare (20% of cases) but is extremely serious compared to viral meningitis. Maternal-fetal transmission is most often the source of infection in the newborn. As for bacterial meningitis in adults, a substantial part would be due to nosocomial infections. Bacterial meningitis is thus often subdivided in three subcategories: community-acquired bacterial meningitis, maternal-fetal bacterial meningitis, and nosocomial bacterial meningitis.

Meningitis is known as nosocomial if it occurs at the Hospital (more than 48 hours after the beginning of hospitalization) or if it follows a gesture that leads to contamination (lumbar puncture, neurosurgical intervention). The relative frequency of nosocomial meningitis compared with those of Community-acquired meningitis origin increased since the 1960s. It is about 50% currently. Factors that predispose to nosocomial bacterial meningitis are mainly recent neurosurgical interventions (CSF leakage, wound infection), the presence of a neurosurgical device, but also the immunosuppression related to diabetes or undernutrition. Postoperative nosocomial meningitis represents half of the infections after neurosurgery. The local inflammatory reaction induced by a foreign body promotes infections by reducing the defense systems of the host. These infections are serious with variable morbidity, neurological sequelae in up to 50% cases, on average 15% mortality, up to 38% in patients over age 60. The main cause of mortality remains infection (primitive and secondary infections) and increases in case of multidrug resistant germs. Patients' life is at stake if targeted antibiotic treatment is delayed and if the infectious material is not removed. Whether bacterial meningitis is proven or suspected, it remains an absolute diagnostic and therapeutic emergency.

### Diagnosis:

There is no obvious diagnostic assessment by the clinician between viral and bacterial meningitis because the symptoms are common: fever, headache, stiffness in the neck, vomiting, limited elevation of the lower limbs, and involuntary flexion of lower limbs as soon as the neck is forced to bend. Moreover, the clinical symptoms of nosocomial bacterial meningitis are not often clear, especially in a post-operative or post-traumatic context. Only the cyto-bacteriological examination of the cerebrospinal fluid (CSF) enables to make the diagnosis. Culture of CSF is rapidly and systematically done. A fraction of the sampling of CSF is separately kept in anticipation of microbiological research by polymerase chain reaction (PCR). Count of leukocytes, GRAM stain, and detection of soluble antigens can orient the diagnosis but the sensitivity of these techniques is strongly altered if the patient is already treated by antibiotic therapy.

The bacteria found in the case of nosocomial meningitis are different from those found in community-acquired meningitis. Gram-negative bacilli and *Staphylococcus* are most often encountered. These are mainly *Pseudomonas aeruginosa* (40% of cases) and *Staphylococcus aureus* (10% of cases). The rapid and targeted diagnosis by polymerase chain reaction in real-time from the DNA of bacterial pathogens is the technique with the best sensitivity and specificity.

## **PRINCIPLE OF DETECTION**

The Eurobioplex EBX-021 is a test using real-time PCR amplification of bacterial DNA of *Staphylococcus aureus* (SA) and *Pseudomonas aeruginosa* (PA), as well as a DNA extraction and PCR inhibition control. The test is performed directly from non-extracted samples of CSF or from extracted DNA or in a single well/tube.

The DNA extraction and PCR inhibition control allows to checks for variations that may occur during the DNA extraction step of biological samples and real-time PCR amplification. Thus, it ensures that a negative result may not be due to a lack of availability of DNA in CSF, and/or a bad DNA extraction, and/or the presence of PCR inhibitors in too large quantities.

This kit is composed of one triplex: *Staphylococcus aureus* (SA), *Pseudomonas aeruginosa* (PA) + control of DNA extraction and PCR inhibition (CI-PCR). DNA of SA and PA is respectively detected using FAM and HEX labeled probes. DNA extraction and PCR inhibition control is detected using a CY5 labeled probe. All probes emit a specific fluorescence following their hydrolysis during the elongation of the amplification product. The measurement of the intensities of real-time fluorescence correlates with the accumulation of specific amplification products.

## DESCRIPTION AND CONTENT OF THE KIT

The real-time PCR EBX-021 kit is ready to use for the specific detection of *Staphylococcus aureus* and *Pseudomonas aeruginosa*.

Fluorescence is emitted and individually recorded through optical measurements during the PCR. The detection of the amplified fragment is performed by a fluorimeter using the channels shown in Table 1 below.

The kit contains reagents and enzyme for the amplification of DNA of these bacteria, and the DNA extraction and PCR inhibition control (see Table 2).

**Table 1 :**

| Target                          | Fluorophore | Excitation | Emission |
|---------------------------------|-------------|------------|----------|
| <i>Staphylococcus aureus</i>    | FAM         | 495 nm     | 515 nm   |
| <i>Pseudomonas aeruginosa</i>   | HEX         | 535 nm     | 555 nm   |
| PCR inhibition control (CI-PCR) | Cy5         | 647 nm     | 667 nm   |

Equivalent channels on different real time PCR cyclers:

- Channel **FAM** (ABI Systems, SmartCycler II, Chromo 4/CFX96, Systems Mx), Channel 510 (LC 480), Channel Green (RotorGene)
- Channel **HEX** (Chromo 4/CFX96, Mx Systems), Channel VIC (ABI Systems), Channel Alexa532 (SmartCycler II), Channel 580 (LC 480), Channel Yellow (RotorGene),
- Channel **Cy5** (ABI Systems, Chromo 4/CFX96, Mx Systems), Channel 660 (LC 480), Channel Alexa647 (SmartCyclerII), Channel Red (RotorGene)

Note: On LC480 instrument II, apply color compensation for the following channels: FAM-HEX/VIC -Cy5 (465-510, 533-580,618-660).

**Table 2 :**

| Cap color          | Components of the kit             | 9 reactions | 36 reactions | Reconstitution |
|--------------------|-----------------------------------|-------------|--------------|----------------|
| <b>Red</b>         | ENZYME                            | 2 x 100 µl  | 8 x 100 µl   | Ready to use   |
| <b>Transparent</b> | OLIGOMIX                          | 2 x 60 µl   | 8 x 60 µl    | Ready to use   |
| <b>Yellow</b>      | Positive Control CP-MBnoso        | 2 x 15 µl   | 8 x 15 µl    | Ready to use   |
| <b>White</b>       | PCR inhibition control (CI-PCR)   | 2 x 50 µl   | 8 x 50 µl    | Ready to use   |
| <b>Blue</b>        | Water = negative control (CN-H2O) | 2 x 500 µl  | 8 x 500 µl   | Ready to use   |

Oligomix : contains the primers and probes of the triplex (SA+PA+CI-PCR)

**Table 3: Number of patients tests according to the EBX format (9 or 36 reactions)**

|                                                                                                        | 9 Reactions | 36 Reactions |
|--------------------------------------------------------------------------------------------------------|-------------|--------------|
| * Number of possible patients tests patient-by-patient with <b>maximum 3 defrost/refreezing cycles</b> | 6 patients  | 24 patients  |
| * Maximum number of patients tests in 1 run without refreezing cycles                                  | 7 patients  | 34 patients  |

**Required material not provided:**

- ◇ Biological Hood
- ◇ Real-time PCR instrument
- ◇ Micro centrifuge
- ◇ Vortex
- ◇ Plates / tubes for real-time PCR
- ◇ Micropipettes
- ◇ DNase-free RNase-free filter tips for micropipettes
- ◇ Sterile microtubes
- ◇ Gloves (powder free)

**STORAGE**

All reagents must be stored between -15 and -22°C.

All reagents can be used until the expiration date indicated on the kit.

Many freezing / defrosting cycles (> 3x) must be avoided, and could lead to decrease in sensitivity.

**CAUTIONS AND NOTES**

**Read carefully instructions before starting.**

- ◇ The experiment must be performed by competent staff.
- ◇ Instruments must have been properly installed, calibrated and maintained according to the manufacturer's recommendations.
- ◇ Clinical samples are potentially infectious and must be processed under a laminar flow hood.
- ◇ The experiment must be performed according to good laboratory practices.
- ◇ Do not use this kit after expiration date.
- ◇ The kit is shipped with dry ice, and the components of the kits must arrive frozen. If one or more components are defrosted, or if the tubes have been damaged, contact Eurobio-Ingen.
- ◇ Once defrosted, spin down briefly the tubes before use.
- ◇ Use of ice or cooling block is advised in case of long delay due for instance to large number of samples or high temperature.
- ◇ It is recommended to define three working areas: 1) Isolation of DNA, 2) Preparation of the reaction mix and 3) Amplification / Detection of amplified products.

- ◇ Use specific lab coat and gloves (powder free) in each working area.
- ◇ Pipettes, reagents and other materials must not cross each area.
- ◇ Specific caution is required to preserve the purity of the reagents and reaction mixtures.
- ◇ Appropriate methods of preparation/extraction of DNA to produce high quality DNA and to be followed by a real-time PCR application should be used, particularly avoiding all sources of DNase contamination.
- ◇ Always use RNase-free DNase-free filtered tips for micropipettes.
- ◇ Do not pipette with mouth and do not eat, drink or smoke in the area.
- ◇ Avoid sprays.

## SAMPLES COLLECTION, TRANSPORT AND STORAGE

- ◇ Collect samples in sterile tubes.
- ◇ It is the responsibility of the user to master its own conditions of collection, transport and storage of samples, and if an extraction of DNA is performed then use suitable systems to produce DNA of good quality.
- ◇ Lumbar puncture is performed after a rigorous asepsis of the area and a local anesthesia. CSF is collected by puncture at the lumbar vertebrae level of the L3-L4, L4-L5 or L5-S1 interspace. Once collected, the CSF is immediately brought to the laboratory at room temperature for analysis.
- ◇ A CSF showing the presence of blood can also give a negative result. In this case, it may be wise to perform the test after extraction of DNA and not on direct LCR. This can also be done in a second attempt, if the result is negative on direct LCR.
- ◇ It is recommended that samples be stored according to the recommendations of storage of samples (Table 4).

**Table 4 :**

| Recommendations of maximum storage of samples before extraction |                   |
|-----------------------------------------------------------------|-------------------|
| Room temperature                                                | 2 h               |
| 4°C                                                             | 16h               |
| -20°C (preferentially -80°C)                                    | Long term storage |

- ◇ The user can refer to the World Health Organization or High Health Authority for storing samples.
- ◇ Transport of clinical samples must obey local regulations for this type of infectious agents.

## PROCEDURE

### **I- DNA extraction**

***This part is optional because this kit can be used directly on CSF (without extraction of bacterial DNA).***

*It is the user's responsibility that the extraction system used be compatible with downstream real time PCR technology. For this kit, we recommend to use extraction methods of bacterial DNA from CSF samples, and refer to manufacturer's instructions of the extraction kit.*

In the EBX-021 kit, CI-PCR on the CY5 channel can be added into the mastermix in case of direct CSF testing, or before extraction, if one is performed. It ensures that a negative result is not due to the presence of PCR inhibitors at high quantity, or potentially due to a problem with extraction if it has been performed.

If an extraction is performed, we recommend the addition of 10 µl of CI-PCR to the biological sample for a final volume of elution of 50 µl after extraction.

If the CI-PCR is added to control the real-time PCR or in the case of direct CSF testing, CI-PCR is added to the reaction mix (1 µl per PCR reaction).

See real-time PCR protocol for details.


CI-PCR is also available separately from Eurobio (Ref EurobioPlex EBX-002).

### **II- Real-time PCR procedure**

#### **General comment:**

Positive control and the DNA extraction and PCR inhibition control (CI-PCR) contain a very high concentration of matrix. Manipulations must be performed carefully to avoid contamination. To control the functioning of the PCR and steps of extraction and real-time PCR amplification, it is necessary to test the positive control CP-MBnos0 as well as the negative controls (water supplied = CN-H<sub>2</sub>O + CI-PCR) (see II-2/6) of real-time PCR procedure).

## II-1/ Diagram of the procedure

## II-2/ Detailed procedure

- 1) Homogenize the tube of Enzyme, and vortex Oligomix, CP-MBnoso, and CI-PCR tubes before starting, and centrifuge.
- 2) Prepare Mastermix as below; prepare the reaction mix as follows by multiplying by the number of samples N to test (including positive and negative controls). On average, prepare enough reagents depending on the size of the sampling for N+2 reactions (up to 3 patients tested) or N+3 reactions (beyond 3 patients tested).  
(Refer to part 1-(a) or 1-(b) of the previous diagram according to the condition).

### Case (a): for CSF direct testing, or for DNA extracted without CI-PCR

| Number of reaction(s) | 1       | N+3             | Example for 3 patients N = 7 |
|-----------------------|---------|-----------------|------------------------------|
| Enzyme                | 12.5 µl | (N+3) x 12.5 µl | 87.5 µl                      |
| Oligomix              | 7.5 µl  | (N+3) x 7.5 µl  | 52.5 µl                      |
| CI-PCR                | 1 µl    | (N+3) x 1 µl    | 7 µl                         |
| Total volume mastemix | 21 µl*  | (N+3) x 21 µl   | 147 µl                       |

\* The volume difference between condition (a) or (b) has no effect on performance.

### Case (b): for DNA extracted with CI-PCR

| Number of reaction(s) | 1       | N+3             | Example for 3 patients N = 7 |
|-----------------------|---------|-----------------|------------------------------|
| Enzyme                | 12.5 µl | (N+3) x 12.5 µl | 87.5 µl                      |
| Oligomix              | 7.5 µl  | (N+3) x 7.5 µl  | 52.5 µl                      |
| Total volume mastemix | 20 µl*  | (N+3) x 20 µl   | 140 µl                       |

- 3) Homogenize mastermix prepared in 2) and centrifuge briefly
- 4) Distribute 20 µL Mastermix\* using a micropipette and filtered tips in each tube or well of microplate for real time PCR.
- 5) Add 5 µL of CSF or extracted DNA sample.
- 6) In parallel, test the following controls:
  - Positive control:
    - 20µL of mastermix + 5µL of CP-MBnoso
  - Negative control:
    - Case (a): for CSF direct testing or for DNA extracted without CI-PCR
      - 20 µL mastermix + 5 µL supplied water (CN-H2O)
    - Case (b): for DNA extracted with CI-PCR
      - 20 µL mastermix + 4 µL water supplied (CN-H2O) + 1 µL CI-PCR
- 7) Close immediately with an adhesive film or transparent caps to avoid contamination.
- 8) Centrifuge briefly to collect all the PCR reaction mix at the bottom of the tubes or plate.
- 9) Program the real-time PCR instrument as follows:

| Program       | Temperature | Duration | Cycle(s) |                          |
|---------------|-------------|----------|----------|--------------------------|
| Heating       | 95°C        | 10 min   | 1        | -                        |
| Denaturation  | 95°C        | 3 min    | 1        | -                        |
| Amplification | 95°C        | 10 sec   | 45       | -                        |
|               | 60°C        | 30 sec   |          | Fluorescence acquisition |

Note 1: On LightCycler ® 480 systems (Roche), two optical systems are available: only "System II" is compliant with use of the kit. Apply color compensation for the following channels: FAM-HEX/VIC-Cy5 (465-510,533-580,618-660).

Note 2: For the Applied Biosystems systems, select "ROX" in "PASSIVE REFERENCE".


Note 3: On Rotorgene ™, please calibrate the signal by clicking on "GAIN optimization".

Note 4: On CFX96 (Biorad), start the run using the v1.6 or later version of CFX Manager software, and analyse with v 3.1 (see § Validation of the Experiment)


## VALIDATION OF THE EXPERIMENT

The analysis of data post-acquisition on a CFX96 PCR instrument (Biorad) must be done with version 3.1 of CFX Manager Software (Biorad). In order to use this v3.1 version from a run started with an older version, follow the procedure below: at the end of the run, the data file with .pcrd suffix must be open and treated with version 3.1 of CFX Manager (Biorad).

If the run was done with CFX Manager v1.6 for instance, to open the data file with CFX Manager v3.1, click on CFX Manager v3.1 icon. The screen below appears.


- Click on File and select Open, then Data File

- Select the file you want to analyze and click on Open.

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Once this is done, the analysis can start.

**For the assay to be valid**, the results for the controls must be the following (Table 5). Otherwise, the experiment is not valid.

**Table 5 :**

| Positive control |                             |
|------------------|-----------------------------|
| FAM              | Ct ≤ 21                     |
| HEX              | Ct ≤ 21                     |
| Negative control |                             |
| FAM              | Ct not determinated         |
| HEX              | Ct ≥ 38 or not determinated |
| CY5              | Ct ≤ 35                     |

## DATA ANALYSIS AND INTERPRETATION

DNA extraction and PCR inhibition control in samples:

Two results can be obtained:

1/ the DNA extraction and PCR inhibition control (CI-PCR; channel Cy5) test is positive: the result can be validated.

2/ the DNA extraction and PCR inhibition control (CI-PCR; channel Cy5) test is negative: either bacterial DNA is not reachable within the CSF, or DNA was not extracted (in case of extraction), the PCR did not work well, or the presence of PCR inhibitors inhibits the PCR reaction. A CSF showing the presence of blood can also give a negative result.

It is then recommended to perform an extraction (direct CSF), or to repeat the extraction (extracted DNA), or dilute the sample.

For clinical samples, the following results are possible:

\* Ct cut off for samples:

Channel FAM *Staphylococcus aureus* : + Positive => Positive Ct ( $\leq 45$ )


Channel HEX *Pseudomonas aeruginosa* : + Positive => Positive Ct and  $< 38$

| PCR Signal |     |     | Presence of SA    | Presence of PA    | Test validity/comment                                                                                                                                                                                                                                                            |
|------------|-----|-----|-------------------|-------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| FAM        | HEX | CY5 |                   |                   |                                                                                                                                                                                                                                                                                  |
| +          | +   | +   | Yes               | Yes               | valid                                                                                                                                                                                                                                                                            |
| +          | -   | +   | Yes               | No                | valid                                                                                                                                                                                                                                                                            |
| -          | +   | +   | No                | Yes               | valid                                                                                                                                                                                                                                                                            |
| -          | -   | +   | No                | No                | valid                                                                                                                                                                                                                                                                            |
| +          | -   | -   | Yes               | Not interpretable | <b>Valid for SA:</b> possible PCR inhibition, or problem with extraction if one was performed, that does not prevent the detection of SA <b>BUT</b> that could prevent the detection of PA - <b>not valid for PA.</b> Dilute 5 x the sample; if same result, re-extract the CSF  |
| -          | +   | -   | Not interpretable | Yes               | <b>Valid for PA:</b> possible PCR inhibition, or problem with extraction if one was performed, that does not prevent the detection of PA. <b>BUT</b> that could prevent the detection of SA - <b>not valid for SA.</b> Dilute 5 x the sample; if same result, re-extract the CSF |
| -          | -   | -   | Not interpretable | Not interpretable | <b>Not valid:</b> possible PCR inhibition, or problem with extraction if one was performed. Dilute 5 x the sample; if same result, re-extract the CSF                                                                                                                            |

SA : *Staphylococcus aureus*, PA : *Pseudomonas aeruginosa*

## PERFORMANCE ANALYSIS

Example of experiment performed on real-time PCR thermocycleur CFX96 (Biorad):

— **HEX : *Pseudomonas aeruginosa***

**Analytical sensitivity:** CP-MBno: 10 copies/μl

**Linearity for quantification:** CP-MBno: from 10 copies/μl to 10<sup>9</sup> copies/μl

**Coefficient of correlation and efficiency**

**Efficiency**

FAM: SA: 96.0 %

HEX: PA: 98.9 %

**Coefficient of correlation**

FAM: SA: 0.998

HEX: PA: 0.998

*Staphylococcus aureus* (SA) and *Pseudomonas aeruginosa* (PA)

**Variability of the signal**

| EBX-021                             | INTRA-experiment | INTER-experiments |
|-------------------------------------|------------------|-------------------|
| Bacterial targets                   | CV %             | CV %              |
| <i>Staphylococcus aureus</i> (FAM)  | 1.31             | 0.64              |
| <i>Pseudomonas aeruginosa</i> (HEX) | 1.89             | 4.84              |

CV : Coefficient of Variation

**Validation on clinical samples**

Specificity and sensitivity were analysed from samples of CSF (DNA not extracted) previously tested by the reference laboratory.

**1/ Analytical specificity**

*Staphylococcus aureus* : > 98%

*Pseudomonas aeruginosa* : > 98%

**2/ Analytical sensitivity**

*Staphylococcus aureus* : > 98%

*Pseudomonas aeruginosa* : > 98%

**3/ Concordance**

*Staphylococcus aureus* : > 98% (n = 48)

*Pseudomonas aeruginosa* : > 98% (n = 49)

**Study of potentially interfering microorganisms: none of the bacteria listed below was detected by the EBX-021 kit (n=1 sample of each bacterium):**

- *Neisseria meningitidis*
- *Streptococcus pneumoniae*
- *Haemophilus influenza*
- *Streptococcus agalactiae*

- *Listeria monocytogenes*
- *Escherichia coli*
- *Bartonella spp*
- *Coxiella burnetii*
- *Mycobacterium avium*
- *Klebsiella kingae*

The tests were performed on CFX96.

## BIBLIOGRAPHY

Van de Beek et al (2006). Community-acquired bacterial meningitis in adults in the Netherlands, 2006–14: a prospective cohort study. In *The Lancet Infectious Diseases*. March 2016 16(3):339-347.

Cohen-Bacrie S, Ninove L, Nougairé A, Charrel R, Richet H, et al. (2011) Revolutionizing Clinical Microbiology Laboratory Organization in Hospitals with In Situ Point-of-Care. *PLoS ONE* 6(7): e22403.

Hotomi M, Togawa A, Kono M, Sugita G, Sugita R, et al. (2013) An Application of Outer Membrane Protein P6-Specific Enzyme-Linked Immunosorbent Assay for Detection of *Haemophilus influenzae* in Middle Ear Fluids and Nasopharyngeal Secretions. *PLoS ONE* 8(8): e71774.

Brouwer MC, Tunkel AR, van de Beek D. (2010) Epidemiology, diagnosis, and antimicrobial treatment of acute bacterial meningitis. *Clin Microbiol Rev*. 2010 Jul;23(3):467-92.

Corless C.E, Guiver M, Borrow R, Edwards-Jones V, Fox A.J and Kaczmarek E.B. (2001). Simultaneous Detection of *Neisseria meningitidis*, *Haemophilus influenzae*, and *Streptococcus pneumoniae* in Suspected Cases of Meningitis and Septicemia Using Real-Time PCR. *J. Clin. Microbiol*. 2001, 39(4):1553.

Charvet A, Garcin F, Albanèse J, Martin C (2009). Méningites nosocomiales. *Antibiotiques*, Volume 11, Issue 1, February 2009, Pages 18-28

<http://www.sante.ujf-grenoble.fr/SANTE/corpus/disciplines/pedia/malinfped/96/leconimprim.pdf>. BOST-BRU C. et PLANTAZ D.

[http://opac.invs.sante.fr/doc\\_num.php?explnum\\_id=8208](http://opac.invs.sante.fr/doc_num.php?explnum_id=8208)  
Epidémiologie des méningites aiguës en France. Parent I. Paris. Institut de Veille Sanitaire. 2012.

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## WASTE DISPOSAL

Be in accordance with the law on the elimination of waste of clinical infectious material.

## SYMBOLS

|                                                                                     |                                        |
|-------------------------------------------------------------------------------------|----------------------------------------|
|    | Reference                              |
|    | Batch number                           |
|    | Limits of storage temperature          |
|    | Expiration Date                        |
|    | Content sufficient for « N » reactions |
|    | Keep protected from light              |
|   | Manufacturer                           |
|  | Instructions for use                   |
| <b>CE</b>                                                                           | CE labeled product                     |
| <b>IVD</b>                                                                          | In vitro diagnostic                    |

eurobio

SCIENTIFIC

7, avenue de Scandinavie  
ZA de Courtaboeuf  
91940 Les Ulis Cedex  
FRANCE