EurobioPlex

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
7, avenue de Scandinavie
ZA de Courtaboeuf
91940 Les Ulis
FRANCE

EurobioPlex

Neonatal bacterial meningitis REAL-TIME PCR

For **qualitative** real-time PCR

REF EBX-022-36
EBX-022-09

9 reactions
36 reactions

Version 6.00 of 2019/10/03

Validated on:

- CFX96™ Real Time PCR detection system (Biorad) with analysis on CFX Manager v 3.1 (Biorad)
- LightCycler®480 Instrument II (Roche) with analysis on LightCycler® 480 software v1.5 (Roche)
- Applied Biosystems® 7500 Real-Time PCR Systems (Applied Biosystems) with analysis on 7500 Software v2
- T-COR 8®-IVD (Tetracore Inc.) with T-COR 8 SmartCT™ (auto v1) software

Storage conditions:

Keep all reagents between -15°C and -22°C until use and after first use

Instructions for use

TABLE OF CONTENTS

Intended use	29
Introduction	29
Principle of detection	30
Description and content of the kit	30
Storage	32
Cautions and notes	32
Samples collection, transport and storage	33
Procedure	33
I- DNA extraction	33
II- Real-time PCR procedure	34
II-1/ Diagram of the procedure	35
II-2/ Detailed procedure	36
Validation of the experiment	37
Data analysis and Interpretation	39
Performance analysis.....	40
Specificities of the Real-Time PCR T-COR 8®-IVD instrument	43
Barcodes for EBX-022 for use on T-COR8®-IVD.....	47
Bibliography	50
Waste disposal	50
Symbols	51

INTENDED USE

The EBX-022 test uses real-time polymerase chain reaction (PCR) amplification and is designed for the qualitative detection from cerebrospinal fluid (CSF) of the three most frequently involved bacteria in the etiology of neonatal/maternal-fetal Bacterial Meningitis (MBmf): *Streptococcus agalactiae*, *Listeria monocytogenes* and *Escherichia coli*. **Real-time PCR can be carried out directly on non-extracted samples of CSF**, which can answer the need of diagnostic urgency in this pathology, but also on **extracted DNA**. This test is indicated to confirm a diagnosis of presumption of infection in patients or complete a proven or indeterminate diagnosis by other techniques (clinical diagnosis with semiological examination, bacteriological study and culture of CSF).

The Eurobioplex EBX-022 has been validated on the following specimen:

- Cerebrospinal fluid (CSF)

INTRODUCTION

Meningitis is an inflammation of meninges surrounding the central nervous system tissue, and is due to many bacterial or viral pathogens.

Bacterial meningitis affects 22 cases per million people, about 1400 cases/year in France. There are four main sources of surveillance in France: 1) the EPIBAC/InVS network which comprises approximately 300 hospital laboratories, 2) the national observatory for childhood bacterial meningitis: GPIP-Activ which includes approximately 250 pediatric services and 150 microbiology laboratories, 3) the regional pneumococcal observatories, 4) the national reference centers.

Two peaks of incidence are observed: one in less than one year old children and the other affects adults over age 60. Bacterial meningitis is rare (20% of cases) but is extremely serious compared to viral meningitis. Maternal-fetal transmission is most often the source of infection in the newborn. As for bacterial meningitis in adults, a substantial part would be due to nosocomial infections. Bacterial meningitis is thus often subdivided in three subcategories: community-acquired bacterial meningitis, maternal-fetal bacterial meningitis, and nosocomial bacterial meningitis.

Mother to child infections are relatively rare (incidence = 0.5%) but are very serious (10% of neonatal mortality). Premature babies, low birth weight children (less than 2 kg), and babies whose immune system is deficient or immature are the most vulnerable to meningitis. It requires immediate treatment in an intensive care unit until antibiotics act. Children's life is at stake if targeted antibiotic treatment is delayed. Whether bacterial meningitis is proven or suspected, it remains an absolute diagnostic and therapeutic emergency.

Diagnosis :

There is no obvious diagnostic assessment by the clinician between viral and bacterial meningitis because the symptoms are common. Maternal-fetal bacterial meningitis can occur in different forms. The order of appearance of the symptoms varies among babies or young children and some may never appear. The following symptoms are most frequently noted: screams or unusual whimpers, rapid breathing, restlessness, irritability when the child is being manipulated, vomiting, refusal to eat, pale or mottled skin, hypotonic child, drowsiness, cold hands and feet, protruding fontanels, spots or rash. Several symptoms of meningitis appear only when the disease is already well advanced and most of them are also observed in several diseases such as flu for

example. This makes the diagnosis even more difficult. Only the cyto-bacteriological examination of the cerebrospinal fluid (CSF) enables to make the diagnosis. Culture of CSF is rapidly and systematically done. A fraction of the sampling of CSF is separately kept in anticipation of microbiological research by polymerase chain reaction (PCR). Count of leukocytes, GRAM stain, and detection of soluble antigens can orient the diagnosis but the sensitivity of these techniques is strongly altered if the patient is already treated by antibiotic therapy.

Neonatal meningitis is most often due to bacteria: *Streptococcus* Group B (*Streptococcus agalactiae*), *Escherichia coli* or *Listeria monocytogenes*. The rapid and targeted diagnosis by polymerase chain reaction in real-time from the DNA of bacterial pathogens is the technique with the best sensitivity and specificity.

PRINCIPLE OF DETECTION

The Eurobioplex EBX-022 is a test using real-time amplification of bacterial DNA of *Streptococcus agalactiae*, *Listeria monocytogenes*, *Escherichia coli*, as well as a DNA extraction and PCR inhibition control. The test is performed directly from non-extracted samples of CSF, or from extracted DNA or in a single well/tube.

The DNA extraction and PCR inhibition control allows to checks for variations that may occur during the DNA extraction step of biological samples and real-time PCR amplification. Thus, it ensures that a negative result may not be due to a lack of availability of DNA in CSF, and/or a bad DNA extraction, and/or the presence of PCR inhibitors in too large quantities.

This kit is composed of one quadriplex: *Streptococcus agalactiae* (Sag), *Listeria monocytogenes* (LM), *Escherichia coli* (EC) + control of DNA extraction and PCR inhibition (CI-PCR). DNA of Sag, LM and EC is respectively detected using FAM, HEX and Texas Red labeled probes. DNA extraction and PCR inhibition control is detected using a CY5 labeled probe. All probes emit a specific fluorescence following their hydrolysis during the elongation of the amplification product. The measurement of the intensities of real-time fluorescence correlates with the accumulation of specific amplification products.

DESCRIPTION AND CONTENT OF THE KIT

The real-time PCR EBX-022 kit is ready to use for the specific detection of *Streptococcus agalactiae*, *Listeria monocytogenes*, and *Escherichia coli*.

Fluorescence is emitted and individually recorded through optical measurements during the PCR. The detection of the amplified fragment is performed by a fluorimeter using the channels shown in Table 1 below.

The kit contains reagents and enzyme for the amplification of DNA of these bacteria, and the DNA extraction and PCR inhibition control (see Table 2).

Table 1 :

Target	Fluorophore	Excitation	Emission
--------	-------------	------------	----------

<i>Streptococcus agalactiae</i>	FAM	495 nm	515 nm
<i>Listeria monocytogenes</i>	HEX	535 nm	555 nm
<i>Escherichia coli</i>	Texas Red	585 nm	605 nm
PCR inhibition control (CI-PCR)	Cy5	647 nm	667 nm

Equivalent channels on different real time PCR cyclers:

- Channel **FAM** (ABI Systems, SmartCycler II, Chromo 4/CFX96, Mx Systems, T-COR8®-IVD), Channel 510 (LC 480), Channel Green (RotorGene)

- Channel **HEX** (Chromo 4/CFX96, Mx Systems, T-COR8®-IVD), Channel VIC (ABI Systems), Channel Alexa532 (SmartCycler II), Channel 580 (LC 480), Channel Yellow (RotorGene),

- Channel **Texas Red** (Chromo 4/CFX96, Mx Systems, T-COR8®-IVD), LC Red 610 (LC480), Channel Orange (RotorGene)

- Channel **Cy5** (ABI Systems, Chromo 4/CFX96, Mx Systems, T-COR8®-IVD), Channel 660 (LC 480), Channel Alexa647 (SmartCyclerII), Channel Red (RotorGene)

Note: On LC480 instrument II, apply color compensation for the following channels: FAM-HEX/VIC – Texas Red-Cy5 (465-510, 533-580, 533-610, 618-660).

Table 2 :

Cap color	Components of the kit	9 reactions	36 reactions	Reconstitution
Red	ENZYME	2 x 100 µl	8 x 100 µl	Ready to use
Transparent	OLIGOMIX	2 x 60 µl	8 x 60 µl	Ready to use
Yellow	Positive Control CP-MBmf	2 x 15 µl	8 x 15 µl	Ready to use
White	PCR inhibition control (CI-PCR)	2 x 50 µl	8 x 50 µl	Ready to use
Blue	Water = negative control (CN-H2O)	2 x 500 µl	8 x 500 µl	Ready to use

Oligomix : contains the primers and probes of the quadriplex (SAg+LM+EC+CI-PCR)

Table 3: Number of patients tests according to the EBX format (9 or 36 reactions) on all instruments except T-COR8®-IVD (see page 43)

	9 Reactions	36 Reactions
* Number of possible patients tests patient-by-patient with maximum 3 defrost/refreezing cycles	6 patients	24 patients
* Maximum number of patients tests in 1 run without refreezing cycles	7 patients	34 patients

Required material not provided:

- ◇ Biological Hood
- ◇ Real-time PCR instrument

- ◇ Micro centrifuge
- ◇ Vortex
- ◇ Plates / tubes for real-time PCR
- ◇ Micropipettes
- ◇ DNase-free RNase-free filter tips for micropipettes
- ◇ Sterile microtubes
- ◇ Gloves (powder free)

STORAGE

All reagents must be stored between -15 and -22°C.

All reagents can be used until the expiration date indicated on the kit.

Many freezing / defrosting cycles (> 3x) must be avoided, and could lead to decrease in sensitivity.

CAUTIONS AND NOTES

Read carefully instructions before starting.

- ◇ The experiment must be performed by competent staff.
- ◇ Instruments must have been properly installed, calibrated and maintained according to the manufacturer's recommendations.
- ◇ Clinical samples are potentially infectious and must be processed under a laminar flow hood.
- ◇ The experiment must be performed according to good laboratory practices.
- ◇ Do not use this kit after expiration date.
- ◇ The kit is shipped with dry ice, and the components of the kits must arrive frozen. If one or more components are defrosted, or of the tubes have been damaged, contact Eurobio Scientific.
- ◇ Once defrosted, spin down briefly the tubes before use.
- ◇ It is recommended to define three working areas: 1) Isolation of DNA, 2) Preparation of the reaction mix and 3) Amplification / Detection of amplified products.
- ◇ Use specific lab coat and gloves (powder free) in each working area.
- ◇ Pipettes, reagents and other materials must not cross each area.
- ◇ Specific caution is required to preserve the purity of the reagents and reaction mixtures.
- ◇ Appropriate methods of preparation/extraction of DNA to produce high quality DNA and to be followed by a real-time PCR application should be used, particularly avoiding all sources of DNase contamination.
- ◇ Always use RNase-free DNase-free filtered tips for micropipettes.
- ◇ Do not pipette with mouth and do not eat, drink or smoke in the area.
- ◇ Avoid sprays.

SAMPLES COLLECTION, TRANSPORT AND STORAGE

- ◇ Collect samples in sterile tubes.

- ◇ It is the responsibility of the user to master its own conditions of collection, transport and storage of samples, and if an extraction of DNA is performed then use suitable systems to produce DNA of good quality.
- ◇ Lumbar puncture is performed after a rigorous asepsis of the area and a local anesthesia. CSF is collected by puncture at the lumbar vertebrae level of the L3-L4, L4-L5 or L5-S1 interspace. Once collected, the CSF is immediately brought to the laboratory for analysis.
- ◇ A CSF showing the presence of blood can give a negative result. In this case, it may be wise to perform the test after extraction of DNA and not on direct LCR. This can also be done in a second attempt, if the result is negative on direct LCR.
- ◇ It is recommended that samples be stored according to the recommendations of storage of samples (Table 4).

Table 4 :

Recommendations of maximum storage of samples before extraction	
Room temperature	2 h
4°C	16h
-20°C (preferentially -80°C)	Long term storage

- ◇ The user can refer to the World Health Organization or High Health Authority for storing samples.
- ◇ Transport of clinical samples must obey local regulations for this type of infectious agents.

PROCEDURE

I- DNA extraction

This part is optional because this kit can be used directly on CSF (without extraction of bacterial DNA).

It is the user's responsibility that the extraction system used be compatible with downstream real time PCR technology. For this kit, we recommend to use extraction methods of bacterial DNA from CSF samples, and refer to manufacturer's instructions of the extraction kit.

In the EBX-022 kit, CI-PCR on the CY5 channel can be added into the mastermix in case of direct CSF testing or before extraction, if one is performed. It ensures that a negative result is not due to the presence of PCR inhibitors at high quantity, or potentially due to a problem with extraction if it has been performed.

If an extraction is performed, we recommend the addition of 10 µl of CI-PCR to the biological sample for a final volume of elution of 50 µl after extraction.

If the CI-PCR is added to control the real-time PCR or in the case of direct CSF testing, CI-PCR is added to the reaction mix (1 µl per PCR reaction).

See real-time PCR protocol for details.

CI-PCR is also available separately from Eurobio (Ref EurobioPlex EBX-002).


II- Real-time PCR procedure

General comment:

Positive control and the DNA extraction and PCR inhibition control (CI-PCR) contain a very high concentration of matrix. Manipulations must be performed carefully to avoid contamination. To control the functioning of the PCR and steps of extraction and real-time PCR amplification, it is necessary to test the positive control CP-MBmf as well as the negative controls (water supplied = CN-H₂O + CI-PCR) (see II-2/6) of real-time PCR procedure).

On T-COR 8®-IVD, we recommend that the positive and negative controls be tested when a new kit is being opened.

II-1/ Diagram of the procedure

II-2/ Detailed procedure

- 1) Homogenize the tube of Enzyme, and vortex Oligomix, CP-MBmf, and CI-PCR tubes before starting, and centrifuge.
- 2) Prepare Mastermix as below; prepare the reaction mix as follows by multiplying by the number of samples N to test (including positive and negative controls). On average, prepare enough reagents depending on the size of the sampling for N+2 reactions (up to 3 patients tested) or N+3 reactions (beyond 3 patients tested).
(Refer to part 1-(a) or 1-(b) of the previous diagram according to the condition).

Case (a): for CSF direct testing, or for DNA extracted without CI-PCR

Number of reaction(s)	1	N+3	Example for 3 patients N = 7
Enzyme	12.5 µl	(N+3) x 12.5 µl	87.5 µl
Oligomix	7.5 µl	(N+3) x 7.5 µl	52.5 µl
CI-PCR	1 µl	(N+3) x 1 µl	7 µl
Total volume mastemix	21 µl*	(N+3) x 21 µl	147 µl

* The volume difference between condition (a) or (b) has no effect on performance.

Case (b): for DNA extracted with CI-PCR

Number of reaction(s)	1	N+3	Example for 3 patients N = 7
Enzyme	12.5 µl	(N+3) x 12.5 µl	87.5 µl
Oligomix	7.5 µl	(N+3) x 7.5 µl	52.5 µl
Total volume mastemix	20 µl*	(N+3) x 20 µl	140 µl

- 3) Homogenize mastermix prepared in 2) and centrifuge briefly
- 4) Distribute 20 µL Mastermix* using a micropipette and filtered tips in each tube or well of microplate for real time PCR.
- 5) Add 5 µL of CSF or extracted DNA sample.
- 6) In parallel, test the following controls (see this specific point for T-COR 8®-IVD on page 43-49):
 - Positive control:
 - 20µL of mastermix + 5µL of CP-MBmf
 - Negative control:
 - Case (a): for CSF direct testing or for DNA extracted without CI-PCR
 - 20 µL mastermix + 5 µL supplied water (CN-H2O)
 - Case (b): for DNA extracted with CI-PCR
 - 20 µL mastermix + 4 µL water supplied (CN-H2O) + 1 µL CI-PCR
- 7) Close immediately with an adhesive film or transparent caps to avoid contamination.
- 8) Centrifuge briefly to collect all the PCR reaction mix at the bottom of the tubes or plate.

- 9) Program the real-time PCR instrument as follows (on T-COR 8®-IVD, no programming is needed thanks to the use of Barcodes- see page 43-49):

Program	Temperature	Duration	Cycle(s)	
Heating	95°C	10 min	1	-
Denaturation	95°C	3 min	1	-
Amplification	95°C	10 sec	45	-
	60°C	30 sec		Fluorescence acquisition

Note 1: On LightCycler ® 480 systems (Roche), two optical systems are available: only "System II" is compliant with use of the kit. Apply color compensation for the following channels: FAM-HEX/VIC-Texas Red-Cy5 (465-510,533-580,533-610,618-660).

Note 2: For the Applied Biosystems systems, select "NONE" in "PASSIVE REFERENCE".


Note 3: On Rotorgene ™, please calibrate the signal by clicking on "GAIN optimization".

Note 4: On CFX96 (Biorad), start the run using the v1.6 or later version of CFX Manager software, and analyse with v 3.1 (see § Validation of the Experiment)

VALIDATION OF THE EXPERIMENT

The analysis of data post-acquisition on a CFX96 PCR instrument (Biorad) must be done with version 3.1 of CFX Manager Software (Biorad). In order to use this v3.1 version from a run started with an older version, follow the procedure below: at the end of the run, the data file with .pcrd suffix must be open and treated with version 3.1 of CFX Manager (Biorad).

If the run was done with CFX Manager v1.6 for instance, to open the data file with CFX Manager v3.1, click on CFX Manager v3.1 icon. The screen below appears.


- Click on File and select Open, then Data File

- Select the file you want to analyze and click on Open.

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Once this is done, the analysis can start.

For the assay to be valid, the results for the controls must be the following (Table 5). Otherwise, the experiment is not valid. On T-COR 8®-IVD, we recommend that the positive and negative controls be tested when a new kit is being opened.

Table 5:

Positive Control	
FAM	Ct ≤ 28
HEX	Ct ≤ 28
TEXAS RED	Ct ≤ 28
Negative Control	
FAM	Ct not determined
HEX	Ct not determined
TEXAS RED	Ct ≥ 38 or not determined*
CY5	Ct ≤ 35

* On T-COR 8®-IVD, Ct ≥ 34 or not determined. Interpretation is automatically generated with Barcodes.

DATA ANALYSIS AND INTERPRETATION

DNA extraction and PCR inhibition control in samples:

Two results can be obtained:

1/ the DNA extraction and PCR inhibition control (CI-PCR; channel Cy5) test is positive: the result can be validated.

2/ the DNA extraction and PCR inhibition control (CI-PCR; channel Cy5) test is negative: either bacterial DNA is not reachable within the CSF, or DNA was not extracted (in case of extraction), the PCR did not work well, or the presence of PCR inhibitors inhibits the PCR reaction. A CSF showing the presence of blood can also give a negative result.

It is then recommended to perform an extraction (if CSF was initially tested directly), or to repeat the extraction, or dilute the sample, unless a signal is detected on the target channel (FAM, HEX, Texas Red).

For clinical samples, the following results are possible:

*** Ct cut off for positivity of samples:**

Channel FAM *Streptococcus agalactiae*: + Positive => Ct positive (≤ 45)

Channel HEX *Listeria monocytogenes*: + Positive => Ct positive (≤ 45)

Channel TEXAS RED *Escherichia coli*: + Positive => Ct positive and < 38

(Specificity of T-COR 8®-IVD, Canal TEXAS RED *Escherichia coli*: + Positive => Ct positive and < 34)

- PCR signal in the Cy5 channel (Ct ≤ 45)

PCR Signal				Presence of SAg	Presence of LM	Presence of EC	Test validity / comment
FAM	HEX	Texas Red	CY5				
+	+	+	+	Yes	Yes	Yes	VALID
+	+	-	+	Yes	Yes	No	
+	-	+	+	Yes	No	Yes	
-	+	+	+	No	Yes	Yes	
+	-	-	+	Yes	No	No	
-	+	-	+	No	Yes	No	
-	-	+	+	No	No	Yes	
-	-	-	+	No	No	No	

- No PCR signal in the Cy5 channel (underdetermined Ct) but Positive for one or more pathogens (Ct < Ct cut off of positivity):

PCR Signal				Presence of SAg	Presence of LM	Presence of EC	Test validity / comment
FAM	HEX	Texas Red	CY5				
+	-	-	-	Yes	Not interpretable	Not interpretable	Possible PCR inhibition or problem with extraction if performed, that does not prevent the detection of one pathogen (valid)
-	+	-	-	Not interpretable	Yes	Not interpretable	
-	-	+	-	Not interpretable	Not interpretable	Yes	

							for this pathogen) BUT that could prevent the detection of the other(s) -dilute 5 x the sample and repeat PCR; if same result redo an extraction
--	--	--	--	--	--	--	--


- No PCR signal in the Cy5 channel (undetermined Ct), and Ct signal > Ct cut off for positivity for the pathogens:

PCR Signal				Presence of SAg	Presence of LM	Presence of EC	Test validity / comment
FAM	HEX	Texas Red	CY5				
-	-	-	-	Not interpretable	Not interpretable	Not interpretable	Possible PCR inhibition or problem with extraction if performed - dilute first 5 x and repeat PCR; if same result redo an extraction

SAg: *Streptococcus agalactiae*, LM: *Listeria monocytogenes*, EC: *Escherichia coli*

PERFORMANCE ANALYSIS

Example of experiment performed on real-time PCR thermocycleur CFX96 (Biorad):

Analytical sensitivity: CP-MBmf: 10 copies/μl

Linearity for quantification: CP-MBmf: from 10 copies/μl to 10⁹ copies/μl

Coefficient of correlation and efficiency

Efficiency

FAM: SAg: 95.4 %

HEX: LM: 100.9 %

Texas Red: EC: 99.1 %

Coefficient of correlation

FAM: SAg: 0.993

HEX: LM: 0.994

Texas Red: EC: 0.997

Streptococcus agalactiae (SAg), *Listeria monocytogenes* (LM), *Escherichia coli* (EC)

Variability of the signal

EBX-022	INTRA-experiment	INTER-experiments
Bacterial targets	CV %	CV %
<i>Streptococcus agalactiae</i> (FAM)	0.69	2.26
<i>Listeria monocytogenes</i> (HEX)	0.57	1.56
<i>Escherichia coli</i> (TEXAS RED)	0.64	3.21

CV : Coefficient of Variation

Validation on clinical samples

Specificity and sensitivity were analyzed from samples of CSF (DNA not extracted) previously tested by the reference laboratory.

		EBX-022								
		SAg			LM			EC		
		POS	NEG	TOTAL	POS	NEG	TOTAL	POS	NEG	TOTAL
Pretest	POS	5	0	5	1	0	1	7	0	7
	NEG	1	30	31	0	42	42	1	35	36
	TOTAL	6	30	36	1	42	43	8	35	43

Bacterial targets EBX-022	Specificity	Sensitivity	Concordance
<i>Streptococcus agalactiae</i> (FAM)	96.8 % *	> 98 %	97.2 % (n = 36)
<i>Listeria monocytogenes</i> (HEX)	> 98 %	> 98 %	> 98 % (n = 43)
<i>Escherichia coli</i>	97.2 % *	> 98 %	97.7 % (n = 43)

* Note that samples not detected by the reference method but detected by the Eurobio method are at the limit of detection. It is likely that the EBX-022 Eurobioplex is more sensitive.

Study of potentially interfering microorganisms: Non interference

No cross-reaction for:

- *Streptococcus agalactiae* on *Listeria monocytogenes* (n=5) = 0%
- *Streptococcus agalactiae* on *Escherichia coli* (n=5) = 0 %
- *Listeria monocytogenes* on *Streptococcus agalactiae* (n=1) = 0%
- *Listeria monocytogenes* on *Escherichia coli* (n=1) = 0 %
- *Escherichia coli* on *Streptococcus agalactiae* (n=7) = 0%
- *Escherichia coli* on *Listeria monocytogenes* (n=7) = 0%

None of the bacteria listed below was detected by the EBX-022 kit (n=1 sample of each bacterium):

- *Neisseria meningitidis*
- *Streptococcus pneumoniae*
- *Haemophilus influenza*
- *Staphylococcus aureus*
- *Pseudomonas aeruginosa*
- *Bartonella spp*
- *Coxiella burnetii*
- *Mycobacterium avium*
- *Klebsiella kingae*

The tests were performed on CFX96.

SPECIFICITIES OF THE REAL-TIME PCR T-COR 8®-IVD INSTRUMENT

T-COR 8®-IVD is a real-time PCR instrument with 8 independent wells, which can be programmed independently, patient by patient, regarding thermic protocol, assays, and time of start of the run.

The 8 wells can also be used simultaneously with the same assay.

Note: Vortex the tubes with the provided T-core vortex before putting the tubes in the T-COR 8®-IVD, and always check that there is no bubble, and that the liquid is all located at the bottom of the tube.

Controls

On T-COR 8®-IVD, we recommend that the positive and negative controls be tested when a new kit is being opened.

After this initial control, the internal control allows to check that the PCR is working properly, and that there is no PCR inhibition.

Patients' Tests

<i>With validation of positive and negative controls once, during first use of the kit.</i>	9 reactions	36 reactions
Number of possible patients tests, patient by patient, with a maximum of 3 freezing/thawing cycles	7 patients	34 patients

Fluorophores

Four combinations of fluorophores are available on T-COR 8®-IVD.

For all EBX, such as EBX-022, FAM/HEX/Texas Red/Cy5 combination is used. If some channels are not used, do not consider this channel for results analysis.

Use of Barcodes (available on page 47 to 49)

- 1- Select Menu > New Run
- 2- Select Barcode
- 3- Scan the appropriate Barcode on the right side of the instrument:
 - For the positive control (Barcode EBX-022 Pos Ctrl),
 - For the negative control (Barcode EBX-022 Neg Ctrl),
 - For a patient sample (Barcode Assay: EBX-022)

The instrument randomly selects one of the 8 available wells. Follow the instructions on the instrument.

- 4- Lift up the lid of the selected well (well x), check that the blue LED light is on, and

select “Yes”

5- Place the tube in the corresponding well and select “Next”.

6- If you wish to name the sample (optional), select “sample x”, name it and select “Accept”.

7- Select “Next”

8- To add/test another control or sample, select « Add well » and go back to point 3-.

9- Select “Start run”.

***Note:** The name of a sample or of a run cannot be modified or added after the run is finished. It must be done before or during the run.*

Presentation and Interpretation of results

Ct values and amplification curves are available in the SmartCT™ Table (Ct values on each channel), and “Ct versus PCR cycles” graphs, both being available in real-time while the run is ongoing.

An automatic analysis for positive and negative controls as well as patients’ samples is available in “Interpretations”, at the end of the run, in the “View” window.

For positive and negative controls, the following results are possible:

« Neg Ctrl Fail »: Not valid

« Neg Ctrl Valid »: Valid

« Pos Ctrl Fail »: Non valid

« Pos Ctrl Valid »: Valid

For patients status:

"Detected": Positive for at least one target * → green box

"Not detected": Negative → red box

* The results for the specific bacteria are specified.

CAUTION !

When the box is green, it is important to read the status for each target as some targets may be negative.

« Invalid »: Invalid result -> retest

Example of automatic interpretation of results on T-COR 8®-IVD instrument screen and associated report:

- For valid positive and negative controls, for the limit of detection (LOD) « Detected » for its corresponding targets, and for two negative samples indicated “Not Detected” (samples 1 and 2):


❖ On the instrument screen:

Interpretations	
1	Sample 1 EBX-022 Not Detected - S. ag, L. mono, E. coli
2	Sample 2 EBX-022 Not Detected - S. ag, L. mono, E. coli
5	EBX-022 POS Ctrl EBX-022 Pos Ctrl Detected: Pos Ctrl Valid
6	EBX-022 NEG ctrl EBX-022 Neg Ctrl Detected: Neg Ctrl Valid
7	LOD 10c/microl repl1 EBX-022 Detected: E. coli, L. mono, S. ag
8	LOD 10c/microl repl2 EBX-022 Detected: E. coli, L. mono, S. ag

❖ On the report:


Summary							
Well	Sample	Assay	FAM	HEX	TxR	Cy5	Call
1	Sample 1	EBX-022				21.5	Not Detected
2	Sample 2	EBX-022				21.5	Not Detected
5	EBX-022 POS Ctrl	EBX-022 Pos Ctrl	21.1	20.6	21.1	21.1	Detected • Pos Ctrl Valid
6	EBX-022 NEG ctrl	EBX-022 Neg Ctrl				21.4	Detected • Neg Ctrl Valid
7	LOD 10c/microl repl1	EBX-022	31.4	31.6	31.5	21.5	Detected • E. coli • L. mono • S. ag
8	LOD 10c/microl repl2	EBX-022	31.5	31.6	31.4	21.4	Detected • E. coli • L. mono • S. ag

Example of amplification curves:


Barcodes for EBX-022 for use on T-COR8®-IVD

EBX-022 Pos Ctrl
Positive Control CP-MBmf


EBX-022 Neg Ctrl

Water = Negative Control (CN-H2O)

Assay: EBX-022

BIBLIOGRAPHY

Van de Beek et al (2006). Community-acquired bacterial meningitis in adults in the Netherlands, 2006–14: a prospective cohort study. In *The Lancet Infectious Diseases*. March 2016 16(3):339-347.

Cohen-Bacrie S, Ninove L, Nougairé A, Charrel R, Richet H, et al. (2011) Revolutionizing Clinical Microbiology Laboratory Organization in Hospitals with In Situ Point-of-Care. *PLoS ONE* 6(7): e22403.

Hotomi M, Togawa A, Kono M, Sugita G, Sugita R, et al. (2013) An Application of Outer Membrane Protein P6-Specific Enzyme-Linked Immunosorbent Assay for Detection of *Haemophilus influenzae* in Middle Ear Fluids and Nasopharyngeal Secretions. *PLoS ONE* 8(8): e71774.

Brouwer MC, Tunkel AR, van de Beek D. (2010) Epidemiology, diagnosis, and antimicrobial treatment of acute bacterial meningitis. *Clin Microbiol Rev*. 2010 Jul;23(3):467-92.

Corless C.E, Guiver M, Borrow R, Edwards-Jones V, Fox A.J and Kaczmarek E.B. (2001). Simultaneous Detection of *Neisseria meningitidis*, *Haemophilus influenzae*, and *Streptococcus pneumoniae* in Suspected Cases of Meningitis and Septicemia Using Real-Time PCR. *J. Clin. Microbiol.* 2001, 39(4):1553.

<http://www.sante.ujf-grenoble.fr/SANTE/corpus/disciplines/pedia/malinfped/96/leconimprim.pdf>. BOST-BRU C. et PLANTAZ D.

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WASTE DISPOSAL

Be in accordance with the law on the elimination of waste of clinical infectious material.

SYMBOLS

	Reference
	Batch number
	Limits of storage temperature
	Expiration Date
	Content sufficient for « N » reactions
	Keep protected from light
	Manufacturer
	Instructions for use
	CE labeled product
	In vitro diagnostic

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