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
T-COR 8®-IVD (Tetracore Inc.) with T-COR 8 Per-well SmartCT™ (auto v1) software

Storage conditions:

Keep all reagents between -15°C and -22°C until use and after the first use

Instructions for use

Available on www.eurobio-scientific.com

Fiche technique EurobioPlex Streptocoque Groupe B
Technical Sheet EurobioPlex Group B Streptococcus
Ficha técnica de EurobioPlex estreptococo grupo B

TABLE OF CONTENTS

FRENCH VERSION	2
1 – Intended Purpose	23
2 – General Information	23
3 – Principle of the test	24
4 – Kit components.....	24
5 – Conservation and Storage	25
6 – Precautions.....	26
7 - Procedure	26
I- Sample collection	26
II- Real-time PCR program	27
Controls	27
Procedure	28
Detailed Procedure	29
III- Barcodes use on T-COR 8® - IVD	29
IV- Results presentation	30
8 – Data Analysis and Interpretation.....	30
9 - Performance Analysis.....	32
10 – Quality Control	35
11 – Waste Disposal	35
12 – Incident Report.....	35
13 - Symbols	36
SPANISH VERSION	38

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1 – INTENDED PURPOSE

The EBX-039 assay is a real-time polymerase chain reaction (PCR) assay for *in vitro* use, designed for the qualitative detection of Group B Streptococcus (GBS), *Streptococcus agalactiae* (*S. agalactiae*) in a vaginal swab. The real-time PCR **can be performed directly on vaginal swabs, without extraction**, which allows for a rapid answer to the diagnostic urgency of this pathology. The test is indicated to make a diagnosis of presumed infection in women or new-borns, or to complete a diagnosis that has been proven or undetermined by other techniques.

The EurobioPlex EBX-039 has been validated on the following specimen:

- Vaginal swab

The EurobioPlex EBX-039 test is an *in vitro* medical device. It must be used by qualified medical analysis technicians. This device should not be recycled after use.

2 – GENERAL INFORMATION

Group B Streptococcus (GBS), *S. agalactiae*, belongs to the commensal flora and is responsible for serious invasive infections (bacteraemia, meningitis) in new-borns. Mother to child transmission by a vertical route is most of the time the cause of the infection in new-borns, has a high potential for severity (10% of early neonatal mortality, especially in premature or low birth-weight infants), and constitute an absolute diagnostic and therapeutic emergency. The vital prognosis is rapidly engaged if the targeted antibiotic is delayed.

The main reservoirs for GBS are the gastrointestinal and urinary tracts, from which the colonization of the female genital tract occurs. Vaginal carriage varies from 5 to 40% during pregnancy. Mother to child transmission is essentially aerodigestive, by inhalation or ingestion of contaminated amniotic fluid during delivery. The rate of transmission from the colonized mother to the new-born averages 50 %, and in women colonized by GBS, approximately 1 to 2 % of new-borns will develop an infection in the absence of an antibiotic prophylaxis at the time of delivery.

There are two types of neonatal GBS infections: early infections, which occurs in 80 % of cases during the first 24 hours of life, and late infections occurring between the first week and the third month of life.

- Early infections commonly present with respiratory distress and bacteraemia (89 % of the cases). Meningitis remains a less common clinical presentation in early forms (10-20% of the cases).
- Late neonatal infections are associated with bacteraemia and meningitis in most cases.

Several symptoms of meningitis only appear when the disease is well advanced and many of them can also be seen in several illnesses, complicating the diagnosis.

Rapid and targeted diagnosis by real-time gene amplification of *S. agalactiae*, directly on vaginal swab without extraction, is the technic offering the best sensitivity and specificity. This diagnosis allows an early detection of the presence of the bacteria and the initiation of an early antibiotic prophylaxis.

Fiche technique EurobioPlex Streptocoque Groupe B
Technical Sheet EurobioPlex Group B Streptococcus
Ficha técnica de EurobioPlex estreptococo grupo B

3 – PRINCIPLE OF THE TEST

The EurobioPlex EBX-039 is a bacterial DNA nucleic acid amplification test of *S. agalactiae* as well as of a PCR Inhibition Control, which uses real-time PCR amplification. The test is performed directly on a vaginal swab in a single well/tube.

A PCR Inhibition Control (CI-PCR) is used, in order to check that a negative result is not due to an excess of PCR inhibitors.

S. agalactiae and CI-PCR DNAs are detected using a probe labelled with FAM and CY5 respectively. Following the hydrolysis of their respective probe during the extension of the amplification product in real-time, the measurement of specific fluorescence emission is related to the accumulation of specific amplification products.

4 – KIT COMPONENTS

The real-time PCR EBX-039 kit is ready for use for the specific detection of *S. agalactiae* and a PCR Extraction/Inhibition Control (CI-PCR).

The fluorescence is emitted and measured individually by an optical system during PCR. The detection of amplified fragments is performed by a fluorometer using the channels listed in Table 1 below.

The kit contains the necessary reagents and enzyme for the amplification of the DNA of these bacteria and of the PCR Extraction/Inhibition Control (see Table 2).

Table 1:

Target	Fluorophore	Excitation	Emission
<i>Group B Streptococcus/S. agalactiae</i>	FAM	495 nm	515 nm
PCR Inhibition Control (CI-PCR)	Cy5	647 nm	667 nm

Equivalent channels on different PCR instruments:

- Channel **FAM** (Systems ABI, SmartCycler II, Chromo 4/CFX96, Systems Mx, T-COR 8[®]-IVD), Channel 510 (LC 480), Channel Green (RotorGene)
- Channel **Cy5** (Systems ABI, Chromo 4/CFX96, Systems Mx, T-COR 8[®]-IVD), Channel Alexa647 (SmartCycler II), Channel 660 (LC 480), Channel Red (RotorGene)

Fiche technique EurobioPlex Streptocoque Groupe B
Technical Sheet EurobioPlex Group B Streptococcus
Ficha técnica de EurobioPlex estreptococo grupo B

Table 2:

Cap color	Kit content	25 tests	50 tests	100 tests	Reconstitution
Red	Enzyme	5 x 90 µl	10 x 90 µl	20 x 90 µl	Ready to use
Transparent	Oligomix	3 x 75 µl	6 x 75 µl	12 x 75 µl	Ready to use
Yellow	Positive Control	3 x 25 µl	6 x 25 µl	12 x 25 µl	Ready to use
White	PCR Inhibition Control (CI-PCR)	3 x 20 µl	6 x 20 µl	12 x 20 µl	Ready to use
Blue	Negative Control (CN-H ₂ O)	500 µl	500 µl	500 µl	Ready to use

Number of patients tests within the recommended number of thawing/freezing cycles:

	25 tests	50 tests	100 tests
One by one	15	30	60
By 2 or more	23	48	98

Materials required but not provided:

- ◇ Biological safety cabinet
- ◇ Real-time PCR instrument T-COR 8®-IVD
- ◇ Sterile microtubes
- ◇ Real-time PCR microtubes suitable for T-COR 8®-IVD
- ◇ Vortex
- ◇ Centrifuge suitable for microtubes
- ◇ Micropipettes
- ◇ Disposable tips DNase-free and RNase-free containing filters for micropipettes
- ◇ Gloves (without talc)

5 – CONSERVATION AND STORAGE

All reagents should be stored between -15°C and -22°C.

All reagents may be used until the expiration date mentioned on the kit label

 Multi freeze/thaw cycles should be avoided for the enzyme (>3x) or for the other reagents (>5x), as the test sensitivity could be reduced.

6 – PRECAUTIONS

Carefully read these instructions before starting the procedure.

- ◇ This experimentation should be carried out by qualified medical analysis technicians. Check that the instruments have been installed, calibrated and maintained according to the manufacturer's instructions.
- ◇ It is the user's responsibility to use another system than those validated. In this case, performances are not guaranteed.
- ◇ Clinical samples should be considered as potentially infectious material and should be prepared under a laminar air flow cabinet.
- ◇ This experimentation must be performed according to good laboratory practice
- ◇ Do not use this kit after the expiration date mentioned on the label.
- ◇ Avoid a long exposure of the reagents to light. It should be limited to the technical time necessary to prepare the PCR plate, otherwise performances are not guaranteed.
- ◇ The kit is shipped under dry ice and the kit components should stay frozen. If one or more of these components arrive unfrozen, or if the tubes are damaged during transportation, please contact Eurobio Scientific.
- ◇ Once the reagents are thawed, briefly centrifuge the tubes before use.
- ◇ It is recommended to perform the activities in three separate areas: 1) RNA/DNA isolation, 2) Preparation of the reaction mixture and 3) Amplification/Detection of the amplified products.
- ◇ Wear separate lab coats and gloves (without talc) in each working area
- ◇ Pipettes, reagents and other material should not be moved around between these areas.
- ◇ Special care must be taken to maintain the purity of reagents and reaction mixtures.
- ◇ Appropriate nucleic acids preparation/extraction methods for quality RNA/DNA production and real-time PCR application should be used, in order to avoid any risk of contamination with RNases/DNases.
- ◇ Use disposable tips with filters for the micropipettes, DNase-free and RNase-free.
- ◇ Do not pipette by mouth. Do not eat, drink or smoke in the laboratory.
- ◇ Do not use the device if the packaging is damaged
- ◇ Avoid aerosols.

7 - PROCEDURE

I- Sample collection

- **Vaginal specimen collection should be performed according to published guidelines for the collection of clinical specimens for group B streptococcal culture.**
- It is up to the user to control his own conditions of collection, transportation and conservation of the samples, in order to preserve the nucleic acids quality.
- Samples should be kept as described in the table below (Table 4).

Table 4:

Maximum recommended storage for the samples	
Room temperature	2 h

Fiche technique EurobioPlex Streptocoque Groupe B
Technical Sheet EurobioPlex Group B Streptococcus
Ficha técnica de EurobioPlex estreptococo grupo B

4°C	16 h
-20°C	Long term storage

The user may refer to the recommendations issued by the World Health Organization or the “Haute Autorité de Santé” for a proper storage of the samples.

The transport of clinical samples should follow local regulations for the transport of infectious agents.

II- Real-time PCR program

T-COR 8[®]-IVD is a real-time PCR system with 8 wells independently programmable that can be used, patient by patient, independently of each other, in terms of tests performed, thermal programs, and timing of the PCR start.

The 8 wells may as well be used simultaneously in the same assay.

Controls

On T-COR 8[®]-IVD, we recommend testing the Positive Control and the Negative Control at least when first opening a new kit.

Thereafter, the Internal Control ensures that in the tested sample there is no PCR inhibition and that the PCR works properly.

There is a high matrix concentration in the Positive Control and the Extraction and Inhibition PCR Control (CI-PCR). Handle them with care in order to avoid any contamination.

Fiche technique EurobioPlex Streptocoque Groupe B
Technical Sheet EurobioPlex Group B Streptococcus
Ficha técnica de EurobioPlex estreptococo grupo B

Procedure

1 – REACTION MIXTURE PREPARATION / MASTERMIX

Number of test (s)	N+1*
Enzyme	(N+1) x 12.5 µl
Oligomix	(N+1) x 7.5 µl
Cl-PCR	(N+1) x 1 µl
Total Mastermix volume	(N+1) x 21 µl

*For 1 or 2 tubes, pipet the exact quantity of reagents in each tube

2 – REACTION PREPARATION

Sample

20 µl Mastermix
+
5µl sample

Positive Control*

20 µl Mastermix
+
5µl CP

Negative Control*

20 µl Mastermix
+
5µl Molecular Biology Water (CN-H₂O)

* On T-COR 8®-IVD, we recommend testing the Positive Control and the Negative Control at least when first opening a new kit.

Mix well after adding the sample in the tube, with back and forth pipetting avoiding foaming and check that the liquid volume is at the bottom of the tube. To prevent contamination, carefully close each tube after each pipetting of a control or a sample.

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Detailed Procedure

- 4) Homogenize the Enzyme tube and vortex Oligomix, CP and CI-PCR, then centrifuge.
- 5) For 1 or 2 tubes, pipet the exact quantity of reagents in each tube. Proceed directly to point III – Barcodes use.
- 6) From 3 tubes on, prepare the Mastermix as described in the table below (N = number of tests). Prepare the necessary Mastermix volume for at least N+1 tests.

Number of test (s)	1	N+1	Example for 1 patient + CP + C Neg N = 3
Enzyme	12.5 µl	(N+1) x 12.5 µl	50 µl
Oligomix	7.5 µl	(N+1) x 7.5 µl	20 µl
CI-PCR	1 µl	(N+1) x 1 µl	4 µl
Total Mastermix volume	21 µl	(N+1) x 21 µl	74 µl

- 7) Homogenize the prepared Mastermix and briefly centrifuge.
- 8) Using a micropipette and tips with filters, distribute 20 µl Mastermix in each tube for real-time PCR.
- 9) Add 5 µl non extracted sample.
- 10) In parallel, make the following controls, at least when first opening a new kit:
 - Positive Control: 20µl Mastermix + 5 µl Positive Control CP
 - Negative Control: 20µl Mastermix + 5µl water supplied (CN-H₂O)

Mix well after adding the sample in the tube, with back and forth pipetting avoiding foaming and check that the liquid volume is at the bottom of the tube. To prevent contamination, carefully close each tube after each pipetting of a control or a sample.

III- Barcodes use on T-COR 8[®] - IVD

1- Select Menu > New Run

2- Select Barcode

3- Scan the desired barcode on the right of the system:

- Either for the Positive Control (Barcode EBX-039 Pos Ctrl),
- Or for the Negative Control (Barcode EBX-039 Neg Ctrl),
- Or for a sample (Barcode Assay: EBX-039)

The system will automatically select one of the 8 available wells. Follow the instructions given by the instrument.

4- Lift the flap from the well x selected by the system, check that the blue LED light is on, then select « Yes ».

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5- Position the tube in the corresponding well and lower the flap, then select « Next ».

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7- Select « Next »

8- To add or test another control or sample, select « Add well » and start again at point 3.

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***Note:** Naming a sample or a run cannot be modified or added after the end of the run. It must be done before or during the run.*

IV- Results presentation

The C_t results and the amplification curves are available, in the SmartCT™ Table (C_t values assigned for each channel), as well as C_t graphs according to the cycles in each channel, which can both be viewed in real time on the instrument.

An analysis for the Negative and Positive Controls, as well as for the samples, is generated automatically. It is available in « Interpretations », at the end of the run, in the « View » window.

8 – DATA ANALYSIS AND INTERPRETATION

Interpretation will be done automatically on T-COR 8®-IVD

For Positive and Negative Controls, the following results may occur:

« Neg Ctrl Fail »: Not valid

« Neg Ctrl Valid »: Valid

« Pos Ctrl Fail »: Not valid

« Pos Ctrl Valid »: Valid

Having been tested at least when the kit is opened for the first time, they must be valid in order to use the kit.

Status determination for samples:

Internal Control CI-PCR automatically validates the result of each sample.

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

Presentation example of an automatic result interpretation:

Summary							
Well	Sample	Assay	FAM	HEX	TxR	Cy5	Call
2	Ctrl POS	EBX-039 Pos Ctrl	21.6			25.0	Detected • Pos Ctrl
3	Ctrl NEG	EBX-039 Neg Ctrl				24.4	Detected • Neg Ctrl
4	Ctrl NEG	EBX-039 Neg Ctrl					Invalid
5	POS sample	EBX-039	29.0			24.9	Detected • S. Agalactiae
6	POS sample	EBX-039	28.9				Detected • S. Agalactiae
7	NEG sample	EBX-039					Invalid
8	NEG sample	EBX-039				25.1	Not Detected
							Note
							Pos Ctrl Valid
							Neg Ctrl Valid
							Neg Ctrl Fail
							S. Ag Positive
							S. Ag Positive
							Retest
							S. Ag Negative

Example of amplification curves

Fiche technique EurobioPlex Streptocoque Groupe B
Technical Sheet EurobioPlex Group B Streptococcus
Ficha técnica de EurobioPlex estreptococo grupo B

9 - PERFORMANCE ANALYSIS

- **Analytical sensitivity of Positive Control:** 10 copies/microliters
- **Signal variability:** Coefficient of Variation (%)

EBX-039	INTRA-lot	INTER-lots
<i>Group B Streptococcus/S. agalactiae</i> (FAM)	0,71	1,52
<i>Internal Control (CI-PCR)</i>	0,66	1,77

- **Validation on clinical samples**

The specificity and sensitivity have been analysed on pretested vaginal samples (non-extracted).

		EBX-039		
		Group B Streptococcus/<i>S. agalactiae</i>		
		POS	NEG	TOTAL
Pre-test	POS	15	0	15
	NEG	0	18	18
	TOTAL	15	18	33

EBX-039	Specificity	Sensitivity	Agreement
<i>Streptococcus agalactiae</i> (FAM)	>99 %	>99 %	>99 %

Cross-reaction study:

Each of the following bacteria has been tested and show no cross-reactivity:

- *Listeria monocytogenes*
- *Escherichia coli*
- *Neisseria meningitidis*
- *Streptococcus pneumoniae*
- *Haemophilus influenza*
- *Staphylococcus aureus*
- *Pseudomonas aeruginosa*
- *Bartonella spp*
- *Coxiella burnetii*
- *Mycobacterium avium*
- *Klebsiella kingae*

Fiche technique EurobioPlex Streptocoque Groupe B
Technical Sheet EurobioPlex Group B Streptococcus
Ficha técnica de EurobioPlex estreptococo grupo B

Barcodes on T-COR8®-IVD for EBX-039


EBX-039 Pos Ctrl
Positive Control CP

EBX-039 Neg Ctrl
Water = Negative Control (CN-H2O)

Assay: EBX-039

10 – QUALITY CONTROL

In accordance with Eurobio Scientific quality management system, certified ISO EN 13485, each batch of EurobioPlex Streptocoque Groupe B is tested according to predefined specifications, in order to guarantee a constant quality of the products.

11 – WASTE DISPOSAL

Waste elimination should be done following PIMW rules.

12 – INCIDENT REPORT

Any serious incident related to this device should be reported to EUROBIO and to the competent authority of the member State in which the user and/or the patient is established.

13 - SYMBOLS

	Reference
	Lot number
	Limits of storage temperature
	Expiration date
	Enough for « N » tests
	Keep away from sunlight
	Manufacturer
	CE marked product
	In Vitro Diagnostic
	Instructions for use
	Caution
	Not to be used if the packaging is damaged

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