EurobioPlex SARS-CoV-2 Multiplex

RT-PCR EN TEMPS REEL

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

EBX-041

14. Bibliographie

Chan JF, Kok KH, Zhu Z, Chu H et al. Genomic characterization of the 2019 novel human-pathogenic coronavirus isolated from a patient with atypical pneumonia after visiting Wu-han. *Emerg microbes infect.* 2020 Dec; 9(1):221-236.

Cohen J, Normile D. New SRAS-like virus in China triggers alarm, *Science*. 2020 Jan 17;367(6475):234-235.

Drosten C, Muth D, Corman VM, Hussain R et al. An observational, laboratory-based study of outbreaks of middle east respiratory syndrome coronavirus in Jeddah and Riyadh, kingdom of Saudi Arabia, 2014. *Clin infect dis.* 2015 Feb 1; 60(3):369-77.

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
Reference EBX-041 Version 14.00 of January **27/01/2023**

Validated on:

- CFX96™ Real Time PCR detection system (Bio-Rad) with analysis on CFX Manager version 3.1 (Bio-Rad)
- CFX Opus 96 (Bio-Rad) with analysis on CFX Maestro version 2.2 (Bio-Rad)
- LightCycler®480 with analysis on LightCycler® 480 software v1.5 (Roche)
- T-COR 8®-IVD (Tetradcore Inc.) with T-COR 8 Per-Well SmartCT™ (auto v1) software
- QuantStudio™ 5 & QuantStudio™ 7 Pro (Applied Biosystems™ by ThermoFisher Scientific) avec QuantStudio 6/7 Pro Real-Time PCR Systems Software (Applied Biosystems™ by ThermoFisher Scientific)

Storage conditions:

Keep all reagents between -15°C and -22°C until use and after first use

Instructions for use

Available on www.eurobio-scientific.com

Fiche technique EurobioPlex SARS-CoV-2 Multiplex

Table of contents

1. Introduction.....	32
2. Purpose of the system	31
3. Symbols	34
4. Principle of detection.....	35
5. Content of the kit.....	35
6. Storage	36
7. Materials required not provided.....	37
8. Real-time PCR instrument	37
9. Cautions and note.....	37
10. Procedure.....	39
11. Validation of the experiment	43
12. Data analysis and interpretation.....	45
13. Performance analysis.....	47
14. Bibliography	53
15. Quality control	53
16. Waste disposal.....	54
17. Incident report.....	57
18. Technical assistance	57

Fiche technique EurobioPlex SARS-CoV-2 Multiplex

1. Introduction

SARS-CoV-2/Covid-19 virus appeared in China end of 2019 in the city of Wuhan. It belongs to the Coronaviridae family and the sub-genera Sarbecovirus.

Its genome consists of 29903 ribonucleic acid bases (RNA). SARS-CoV-2 is the 7th coronavirus infecting humans to be identified, since human coronavirus (HCoV) 229E, NL63, OC43, HKU1, SARS-CoV (coronavirus responsible for severe acute respiratory syndrome) and MERS-CoV (coronavirus inducing Middle-East respiratory syndrome).

Eighteen genomes have been isolated and reported, amongst which BetaCoV / Wuhan / IVDC-HB-01/2019, BetaCoV / Wuhan / IVDC-HB-04/2020, BetaCoV / Wuhan / IVDC-HB-05/2019, BetaCoV / Wuhan / WIV04 / 2019 et BetaCoV / Wuhan / IPBCAMS-WH-01/2019. These sequences of SARS-CoV-2 have similarities with the ones of betacoronaviruses found in bats. SARS-CoV-2 is genetically distinct from other human coronaviruses such as the ones related to SARS and to MERS.

Symptoms can come as a cold, fever, cough, difficulty breathing, pneumonia, to severe respiratory syndrome, which can be fatal. The level of mortality is not precise at the beginning of the epidemic, around 2 % end of January. This is much less compared to 10 % for SARS-CoV or 30% for MERS-CoV. SARS-CoV-2 is highly contagious with more than 90 000 cases worldwide beginning of March 2020.

EurobioPlex SARS-CoV-2 was designed to detect all known SARS-CoV-2 published sequences by alignment in silico with sequences of other coronaviruses.

A decision algorithm based on the detection of 3 targets (RdRp Gene target 1 and RdRp target 2, identical to the ones recommended by the World Health Organization (<https://www.who.int/emergencies/diseases/novel-coronavirus-2019/technical-guidance/laboratory-guidance>: published under Real-time RT-PCR assays for the detection of SARS-CoV-2 Institut Pasteur, Paris (2 March 2020)), and N Gene) can be used to determine SARS-CoV-2 patient status (see Data analysis and interpretation and limitations). The diagnostic must always be made by a physician and using the clinical context, history and symptoms of the patient.

Extracted RNA is the starting material for the EurobioPlex SARS-CoV-2 kit.

It is the user's responsibility to choose extraction methods relevant to the type of samples tested.

Fiche technique EurobioPlex SARS-CoV-2 Multiplex

2. Purpose of the system

EurobioPlex SARS-CoV-2 is a test based on real-time reverse-transcription and amplification designed for *in vitro* use for qualitative determination of absence or presence of SARS-CoV-2 in a RNA extract. This test is indicated to help diagnose the occurrence of this infection in humans, or complement a proven or indeterminate diagnosis.

The EurobioPlex SARS-CoV-2 test is an *In Vitro* diagnostic medical device, it must be used by qualified medical biology analysis laboratory personnel. It should not be recycled after use.

The EurobioPlex SARS-CoV-2 has been validated on the following specimen:

- Nasopharyngeal aspirate
- Nasopharyngeal swab
- Bronchoalveolar fluid
- Sputum
- Nasal swab
- Salivary samples

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Fiche technique EurobioPlex SARS-CoV-2 Multiplex

4. Principle of detection

The EBX-041 is a test using reverse-transcription and real-time amplification of viral RNA of SARS-CoV-2 based on amplification of 3 genes (RdRp Gene Target 1, RdRp Gene Target 2 and N Gene target 3). If at least one RdRp target is positive, a SARS-CoV-2 positive diagnostic can be called (see Data analysis and Interpretation).

The kit contains one oligomix to detect the 3 targets, as well as an encapsulated control of RNA extraction and RT-PCR inhibition. The test is performed from extracted RNA from a sample, using one RT-PCR reaction in a single distinct well/tube.

This test ensures sensitivity, and specificity with regard to other coronaviruses known, such as HCoV, SARS-CoV, betacoronavirus associated to bats (BtCoV) and MERS-CoV. The RNA extraction and RT-PCR inhibition control allows to check for variations that may occur during the RNA extraction step of biological samples and real-time PCR amplification. Thus, it ensures that a negative result may be due to a bad RNA extraction, and/or due to the presence of PCR inhibitors in too large quantities.

RNA of SARS-CoV-2 is detected with specific probes of each target, labeled respectively with FAM (target 1/RdRp gene), HEX (target 2/RdRp gene) and Texas red (target 3/N gene). The endogenous control is detected using a CY5 labeled probe. Probes emit a specific fluorescence following their hydrolysis during the elongation of the amplification product. The measurement of the intensity of real-time fluorescence correlates with the accumulation of amplification products.

5. Content of the kit

The RT-PCR real-time EurobioPlex SARS-CoV-2 kit is ready to use and contains all reagents and enzymes for the detection of this virus (Table 1).

Fluorescence is emitted and individually recorded through optical measurements during the PCR. The detection of the amplified fragment is performed by a fluorimeter using the channels shown in the Table 2.

Fiche technique EurobioPlex SARS-CoV-2 Multiplex

Table 1: Content of the kit

Cap color	Content of the kit	25 reactions	50 reactions	100 reactions	192 reactions	200 reactions	600 reactions	Reconstitution
Red	Enzymes	2x100 µL	400 µL	900 µl	2 x 900µl	2 x 900µl	3 x 1300 µl	Ready to use
Transparent	Oligomix	2x80 µL	320 µL	720 µl	1440µL	1440µL	3 x 1050 µl	Ready to use
Yellow	Positive control CP	2x20 µL	80 µL	160 µl	320 µl	320 µl	320 µl	Ready to use
Blue	Water = negative control (CN-H2O)	1500 µL	1500 µL	1500µl	2x1500 µl	2x1500 µl	5 x 1500 µl	Ready to use

Oligomix: contains primers and probes for the 3 targets and for the endogenous control

Table 2: Detection of targets by fluorophores

Targets	Fluorophore	Excitation	Emission
Target 1 / RdRp Gene	FAM	495 nm	515 nm
Target 2 / RdRp Gene	HEX	535 nm	555 nm
Target 3 / N Gene	Texas Red	585 nm	605 nm
Endogenous internal control	Cy5	650 nm	670 nm

Equivalent channels on different real time PCR cyclers:

- Channel **FAM** (Systems ABI, QS5, QS7, SmartCycler II, Systems Mx, CFX96™/CFX Opus 96/Chromo4, T-COR8™-IVD, Dtprime, LC 480), Channel Green (Rotor-Gene Q)
- Channel **HEX** (Chromo 4/CFX96/CFX Opus 96, Systèmes Mx, T-COR 8®-IVD, Dtprime, LC 480), Canal VIC (Systèmes ABI, QS5, QS7.), Canal Alexa532 (SmartCycler II), Canal Yellow (Rotor-Gene Q),
- Channel **Texas Red** (ABI 7500, Chromo 4/CFX96/CFX Opus 96, Systèmes Mx, T-COR 8®-IVD), LC Red 610 (LC480), Channel Orange (Rotor-Gene Q), Channel Rox (QS5, QS7, DTprime)
- Channel **Cy5** (CFX96™/CFX Opus 96/Chromo4, Systems ABI, QS5, QS7, Systems Mx, T-COR8™-IVD, Dtprime, LC 480), Channel Alexa647 (SmartCycler II), Channel Red (Rotor-Gene Q)

Note1: On LC480, apply color compensation for the following channels: FAM/HEX-VIC/Texas Red/Cy5 (FAM : 483-533 / HEX-VIC : 523-568 /Texas Red : 558-610 / Cy5 : 615-670).

Note2: During the analysis using LC480 Software release 1.5.0 apply a noise band of 2,1 for the FAM channel and 4,5 for the HEX channel.

Fiche technique EurobioPlex SARS-CoV-2 Multiplex

6. Storage

All reagents should be stored between -15°C and -22°C.

All reagents can be used until the expiration date indicated on the kit label.

The sensitivity of the assay may decrease if the kit components undergo multiple freeze/thaw cycles. The kit can be used after the initial opening for up to 5 freeze/thaw cycles.

7. Materials required not provided

- ◇ Biological cabinet
- ◇ Real time PCR instrument
- ◇ Centrifuge for microtubes
- ◇ Vortex
- ◇ Plates / tubes for real-time PCR reaction
- ◇ Micropipettes
- ◇ DNase-free and RNase-free filter tips for micropipettes
- ◇ Sterile microtubes
- ◇ Gloves (talc-free)

8. Real-time PCR instrument

The EurobioPlex SARS-CoV-2 Fast-SVD kit has been developed and validated for use with the following real-time PCR instruments:

- CFX96™ Real Time PCR detection system (Bio-Rad) with analysis on CFX Manager version 3.1 (Bio-Rad)
- CFX Opus 96 (Bio-Rad) with analysis on CFX Maestro version 2.2 (Bio-Rad)
- LightCycler® 480 with analysis on LightCycler® 480 software v1.5 (Roche)
- T-COR 8®-IVD (Tetradcore Inc.) with T-COR 8 Per-Well SmartCT™ (auto v1) software
- QuantStudio™ 5 & QuantStudio™ 7 Pro (Applied Biosystems™ by ThermoFisher Scientific) avec QuantStudio 6/7 Pro Real-Time PCR Systems Software (Applied Biosystems™ by ThermoFisher Scientific)

9. Cautions and note

Read these instructions carefully before starting the procedure.

- ◇ This experimentation should be performed by medical laboratory technicians.
- ◇ The local and national biosafety regulations in place for the detection of SARS-CoV-2 must be followed strictly at all times, especially in laboratories and with laboratory equipment in agreement.
- ◇ Ensure that instruments have been installed, calibrated and maintained in accordance with the manufacturer's recommendations.
- ◇ It is the responsibility of the user to use equipment other than that validated, and in such cases performance is not guaranteed.
- ◇ Clinical samples should be considered as potentially infectious material and should be prepared in a laminar flow cabinet.
- ◇ This experiment should be performed according to good laboratory practice.
- ◇ Do not use this kit after the stated expiration date.
- ◇ The kit is shipped under dry ice, and the kit components should arrive frozen. If one or more of the components arrive thawed, or if the tubes have been damaged in transit, contact Eurobio Scientific.
- ◇ Avoid freeze/thaw cycles of reagents as this may lead to a decrease in assay sensitivity.
- ◇ Once the reagents have been thawed, centrifuge the tubes briefly before use.
- ◇ The use of ice or an ice pack is recommended in case of long delays due to, for example, a large number of samples to be processed or high temperatures.
- ◇ It is recommended to define three separate work areas: 1) RNA isolation, 2) Preparation of the reaction mixture, and 3) Amplification/Detection of the amplified products.
- ◇ The fluorescent probes in the oligomix are light sensitive. Any prolonged exposure of the oligomix to light should be limited to the technical time required to prepare the PCR plate.
- ◇ It is recommended that the positive control be opened and handled separately from the biological samples to be tested and the other kit components to avoid cross-contamination.
- ◇ Distinct lab coat and gloves (talc-free) should be worn in each work area.

Fiche technique EurobioPlex SARS-CoV-2 Multiplex

- ◇ Pipettes, reagents and other working materials should not be moved between these areas.
- ◇ Special care should be taken to maintain the purity of reagents and reaction mixtures.
- ◇ The internal endogenous control detects a detectable cell target in all samples of human origin but not in the CN-H2O negative control provided in the kit. The absence of a signal in the negative control prevents cross-contamination from occurring.
- ◇ Appropriate RNA preparation/extraction methods for quality RNA production and RT-PCR application should be used, especially to avoid contamination with RNAs.
- ◇ Use filter tips for micropipettes, RNase-free and DNase-free.
- ◇ Do not pipette with the mouth. Do not eat, drink or smoke in the laboratory.
- ◇ Avoid aérosols.

10. Procedure

a. Samples collection

- ◇ Collect samples in sterile tubes.
- ◇ It is the responsibility of the user to control their own sample collection, transport, storage and extraction conditions to ensure that RNA extraction using appropriate systems produces RNA of high quality.
- ! ◇ It is recommended that samples are extracted immediately, or stored according to the sample storage recommendations before extraction (Table 3).

Table 3: Storage recommendations before use

Recommendations for maximum storage of samples before extraction	
Ambient temperature	2 h
4°C	72 h
-20°C (preferably -80°C)	Long-term storage

Caution	
	◇ The user may refer to the recommendations issued by the World Health Organization or the French National Authority for Health for the proper storage of samples. ◇ Extracted RNA should be stored at -80°C. ◇ The transport of clinical samples is subject to local regulations for the transport of infectious agents.

Fiche technique EurobioPlex SARS-CoV-2 Multiplex

b. RNA extraction

It is the responsibility of the user to ensure that the nucleic acid extraction system used is compatible with real-time RT-PCR technology. We recommend using virus RNA extraction methods suitable for respiratory specimens, and to refer to the instructions of the supplier of the extraction kit used.

For salivary samples, we recommend the use of a Dithiothreitol -based fluidifier (DTT) such as the product marketed by Eurobio (Digest-EUR - reference DD0DIG00-AI). In the case of the use of the latter, we recommend to use 10 µl of the reference DD0DIG00-AI for 200 µl of saliva). Without fluidizer, performance is not guaranteed.

In the EBX-041 kit, the endogenous internal control on the Cy5 channel is already present in the clinical sample. Properly detected after amplification, the endogenous internal control is used to ensure the quality of the sample and the extraction. The performance described in this document was obtained with the following extraction technique: Maelstrom 9600 (TanBEAD) and the OptiPure Viral Auto Plate kit (ref: W665A46).

c. Real-time RT-PCR

General note:

- The positive control contains high concentrations of matrix. Handling must be done carefully to avoid contamination.
- To confirm that the RT-PCR is working correctly, it is necessary to test the positive control, as well as the negative control (water supplied = CN-H2O) (see II-2/6 of the real-time RT-PCR protocol).

Fiche technique EurobioPlex SARS-CoV-2 Multiplex

Diagram of the procedure:

1 - PREPARATION OF MASTERMIX

Number of reactions	N+3
Enzymes	(N+3) x 6,25 µl
Water	(N+3) x 8,75 µl
Oligomix	(N+3) x 5 µl
Total Volume Mastermix	(N+3) x 20 µl

2 - PREPARATION OF REACTIONS

Sample

20 µl Mastermix
+
5µl RNA sample

Positive Control

20 µl Mastermix
+
5µl CP

Negative Control

20 µl Mastermix
+
5µl Molecular biology water (CN-H₂O)

3 - REAL TIME PCR INSTRUMENT

Program	Temperature	Duration	Cycle(s)	
Reverse Transcription	45°C	5 min	1	-
Denaturation	98°C	20 sec	1	-
Amplification	98°C	3 sec	40	-
	58°C	10 sec		Acquisition of fluorescence

Fiche technique EurobioPlex SARS-CoV-2 Multiplex

d. Detailed procedure

- 1) Ensure that the reagents are completely thawed. Homogenise the enzyme tubes, and vortex the Oligomix and CP for about 15 seconds, then centrifuge briefly.
- 2) Prepare the Mastermix as below. N is the number of reactions (including positive and negative controls. Plan to prepare enough reagents for at least N+3 reactions.

Number of reactions	N+3
Enzymes	(N+3) x 6,25 μ l
Water	(N+3) x 8,75 μ l
Oligomix	(N+3) x 5 μ l
Total Volume Mastermix	(N+3) x 20 μ l

(≤ 10):
is sufficient.

**For smaller series
preparing for N+2*

- 3) Homogenise the Mastermix prepared in 2) and centrifuge briefly
- 4) Dispense 20 μ L of Mastermix into separate microplate tubes/wells using a micropipette and filter tips.
- 5) Add 5 μ L of extracted RNA sample.
- 6) In parallel perform the following controls:
 - Positive Control:
 - 20 μ L of Mastermix + 5 μ L of CP.
 - Negative Control:
 - 20 μ L of Mastermix + 5 μ L supplied water (CN-H2O).
- 7) Immediately close the tubes, or plate with an adhesive film to avoid contamination.
- 8) Centrifuge briefly to collect the reaction mixture at the bottom of the tubes or microplate wells.
- 9) Run the following program on the real-time PCR instrument.

Program	Temperature	Duration	Cycle(s)	
Reverse Transcription	45°C	5 min	1	-
Denaturation	98°C	20 sec	1	-
Amplification	98°C	3 sec	40	-
	58°C	10 sec		Acquisition of fluorescence

Note1: On CFX96™ (Bio-Rad), start the run using the v1.6 or later version of CFX Manager Software, and analyze with v 3.1 (see § Validation of the Experiment)

Note 2: On LightCycler ® 480 (Roche), apply color compensation for the following channels: FAM/HEX-VIC/Texas Red/Cy5 (FAM : 483-533 / HEX-VIC : 523-568 /Texas Red : 558-610 / Cy5 : 615-670).


Fiche technique EurobioPlex SARS-CoV-2 Multiplex

Note 3 : During the analysis using LC480 Software release 1.5.0 apply a noise band of 2,1 for the FAM channel and 4,5 for the HEX channel.

11. Validation of the experiment


Post-acquisition data analysis on a CFX96 PCR instrument (Bio-Rad) must be performed using version 3.1 of the CFX Manager software (Bio-Rad). In order to change to this version from a run started on an earlier version, please follow the procedure below: At the end of the run, the data file with the extension .pcrd must be opened and processed with CFX Manager (Bio-Rad) version 3.1.

If the run was started with CFX Manager v1.6 for example, to open a data file with CFX Manager v3.1, click on the CFX Manager v3.1 icon. The home screen appears.


Click on "File" and select "Open" then "Data File".

Fiche technique EurobioPlex SARS-CoV-2 Multiplex


Select the file you need to analyse and click on "Open".

The "Drift correction" option must be applied in the "Settings" tab as shown in the image below: click on the "Settings" tab, then on "Baseline Setting" and on "Apply Fluorescence Drift Correction".


To optimise the run analysis, select the "Log Scale" checkbox for each channel analysed. Then place the threshold bar above the background noise corresponding to the middle of the exponential phase. The 5 channels of interest have high RFUs and differ according to the channel considered, so the "Log Scale" option allows for a better reading and analysis of the run.

Fiche technique EurobioPlex SARS-CoV-2 Multiplex

Once this step has been completed, the analysis can proceed for the 5 channels considered **FAM, HEX, Texas Red and CY5**.

To validate the assay, the Ct values for the controls must be the following (Table 4). Outside of these values, the experiment cannot be validated.

Table 4: Run validation

Positive Control	
FAM	Ct ≤ 32
HEX	Ct ≤ 32
Texas Red	Ct ≤ 32
Cy5	Ct ≤ 32
Negative Control	
FAM	Undetermined Ct
HEX	Undetermined Ct
Texas Red	Ct ≥ 38 or Undetermined
Cy5	Undetermined Ct

12. Data analysis and interpretation

RNA extraction and RT-PCR inhibition control in samples:

The proper functioning of RT-PCR reaction can be evaluated on the Cy5 channel measuring the endogenous control.

In some cases, it is recommended to repeat the extraction or to dilute the sample 5 times, because the result cannot be interpreted (NI) (See column « validity of the test » on Table 5. All cases that can be encountered are described in Table 5.

For clinical samples, and determination of presence or absence of SARS-CoV-2:

Cut-off Ct values for positivity:
Target 1 and Target 2 RdRp Gene and Target 3 N gene: Ct < 40

Fiche technique EurobioPlex SARS-CoV-2 Multiplex

Ct cut off for positivity of samples for T-COR 8®-IVD is: + Positive => Ct positive (≤ 45) for the 3 channels, but automatic interpretation is provided with Barcodes.

The following results are possible:

Table 5: Detection of SARS-CoV-2

Target 1 RdRp	Target 2 RdRp	Target 3 N	Endogenous control	Validity of the test	Presence of SARS-CoV-2 or no possible interpretation (NI)
FAM	HEX	Texas Red	CY5		
+	+/-	+/-	+ / -	Yes	YES
+/-	+	+/-	+ / -	Yes	YES
-	-	-	+	Yes	NO
-	-	+	+	Yes	Undetermined
-	-	+/-	-	Limited	NI

NI: no possible interpretation because of RT-PCR inhibition or failed extraction: no conclusion can be given. **It is then recommended to proceed to a new sampling and/or repeat the extraction and/or to dilute 5 times the sample.**

Limitations of use and interpretation:

- ❖ All samples must be treated as potentially infected by SARS-CoV-2, and biosafety local regulations must be thoroughly followed.
- ❖ Interpretation of results must take into account the possibility of false negatives and false positives.
 - False negative can be due to:
 - Inappropriate collection of samples, or bad storage
 - Samples outside the viremic phase
 - Incorrect extraction methods or use of non-validated PCR instruments
 - Manipulations that do not rigorously follow all the indications of this manual.
 - False positive can be due to:
 - A contamination related to wrong manipulation of highly positive samples, or from the positive control, or PCR products
 - Procedures that do not rigorously follow precautions to avoid contamination described in this manual
- ❖ Results must be interpreted by medical professionals in the clinical context of the patient, its history and symptoms.
- ❖ This test does not exclude the presence of other pathogens than the SARS-CoV-2.
- ❖ A negative result for this test does not absolutely exclude a possible infection with SARS-CoV-2.

Fiche technique EurobioPlex SARS-CoV-2 Multiplex

13. Performance analysis

Limit of detection/analytical sensitivity

A dilution range was performed from a plasmid mixture ranging from 10^5 to 1 copy/ μ L. This dilution range was used to determine the detection limit (LoD) of the kit.

Detection limitation/analytical sensitivity:

- **Positive control:** 10 copies/ μ L CP.
- **A SARS-CoV-2 synthetic RNA** (Twist Bioscience, USA) was detected with a sensitivity of 12.5 copies ARN/ μ L.
- **The First WHO International Standard for SARS-CoV-2 RNA** was detected with a sensitivity of 7.5 copies/ μ L
- **Sensitivity study on 106 known positive samples, diluted in known samples**
Negatives at 3 x LOD: 100% sensitivity (detection of 106 samples / 106)

- **Interfering substances:** 7 substances from the AcroMetrixTM Inhibition Panel REF 956400 kit, Oxymetazoline and Zanamivir. A high concentration, which may be contained in some drugs, was added to the positive sample before extraction: No impact of the tested interfering substances was observed on the result obtained.
- **The limit of detection on salivary sample**, with a confidence of 95% calculated from the data reported below is under **1940 c/mL**:

Copies/ml SARS-CoV-2	Number of saliva tested	Number of positive saliva	% detection
800	20	16	80
1200	13	12	91,7
3600	20	20	100
12000	20	20	100

Signal variability on FAM, HEX, Texas Red and Cy5 channels

- Between-lots Ct variability (3 lots) of EBX-041 is the following:

Fiche technique EurobioPlex SARS-CoV-2 Multiplex

Coefficient of variation %	RdRp Gene Target 1	RdRp Gene Target 2	N Gene Target 3	Endogenous control gene
CFX96	1.21	0.37	1.67	4.44
LC480	1.24	0.91	2.18	1.41
T-COR8-IVD	1.80	1.66	1.33	2.53
QS5 / QS7 pro	3.90	4.41	3.31	6.50

Diagnostic specificity

This validation included :

- 508 SARS-CoV-2 negative samples, pretested and characterized by CE approved tests.
- 50 negative samples for known coronaviruses, some of which are positive for other respiratory viruses, in order to evaluate potential cross-reactions.

Number of samples	Positive for	Coronavirus	Status SARS-CoV-2 EBX-047
1	Adenovirus	Negative	Negative
5	Parainfluenzae	Negative	Negative
1	Rhinovirus	Negative	Negative
1	Bocavirus	Negative	Negative
2	Legionella pneumoniae	Negative	Negative
2	Metapneumovirus	Negative	Negative
3	Mycoplasma pneumoniae	Negative	Negative
3	Mycoplasma tuberculosis	Negative	Negative
2	Bordetella pertussis	Negative	Negative
1	Parechovirus Type 1	Negative	Negative
1	Parechovirus Type 2	Negative	Negative
1	Parechovirus Type 3	Negative	Negative
1	Parechovirus Type 4	Negative	Negative
1	Parechovirus Type 5	Negative	Negative
1	Parechovirus Type 6	Negative	Negative
3	Enterovirus	Negative	Negative
1	Streptococcus pyogenes	Negative	Negative
3	Streptococcus pneumoniae	Negative	Negative
2	Staph.epiderm/oralis/mitis/sanguis	Negative	Negative
3	Candida albicans	Negative	Negative
2	Pseudomonas aeruginosa	Negative	Negative

Fiche technique EurobioPlex SARS-CoV-2 Multiplex

3	CoV NL63	Positive	Negative
2	CoV OC43	Positive	Negative
1	CoV 229E	Positive	Negative
1	CoV KU1	Positive	Negative
1	HPV	Negative	Negative

In silico analysis was used and demonstrated the specificity of the primers and probes for SARS-CoV-2. For example, on the sequences of MERS and SARS viruses, no significant homology of SARS-CoV-2 primers and probes that could amplify these viruses was found.

No aspecificity or cross reaction was detected.

The kit is 100% specific.

Diagnostic sensibility

The validation of the performance was performed on:

558 negative samples	
508 SARS-CoV-2 negative samples	Characterised with Allplex 2019 n-CoV (Seegene) kit, CE marked
50 positive respiratory samples	RSV-A/B, Parainfluenzae, Rhinovirus, Metapneumovirus, CoV NL63, CoV OC43, CoV 229E, CoV KU1

270 positive samples	
270 positive samples	Characterized by EurobioPlex EBX-041 CE-IVD kit

For this study, sample extraction was performed on Maelstrom 9600 (TANBead) with the OptiPure Viral Auto Plate kit (ref: W665A46), and RT-PCR on CFX96 (Bio-Rad). And tests were performed on CFX96™ (Bio-Rad).

Fiche technique EurobioPlex SARS-CoV-2 Multiplex

		EBX-041	
		SARS-CoV-2 positive	SARS-CoV-2 negative
Pre-tested EBX-041 and Seegene kit	SARS-CoV-2 positive	270	0
	SARS-CoV-2 negative	0	558

General sensitivity: > 99% (270/270)
General specificity: > 99% (558/558)
General concordance: > 99% (828/828)

Diagnostic specificity and sensitivity on salivary samples/nasopharyngeal swabs taken concurrently on the same patients:

Evaluation was carried out on 103 salivary samples and 103 nasopharyngeal samples, taken concurrently on the same patients.

Out of the 103 nasopharyngeal samples, 37 were positive and 66 negative for SARS-CoV-2 using EBX-041.

RNA extraction was performed on Genolution system (Nextractor®). Saliva was pre-treated with Digest-EUR (Eurobio. Réf. : DD0DIG00-AI) using 12.5 µL in a mix of saliva and PBS (250µL/250µL).

➤ **Diagnostic specificity on salivary samples/nasopharyngeal swabs taken concurrently on the same patients:**

Evaluation was carried out on the concurrent 66 salivary samples from the same 66 patients that had a SARS-CoV-2 negative nasopharyngeal sample.

No specific signal was detected in any of the 66 saliva, all displaying a valid control of extraction and RT-PCR amplification.

➤ **Diagnostic sensitivity on salivary samples/nasopharyngeal swabs taken concurrently on the same patients:**

Evaluation was carried out on the concurrent 37 salivary samples from the same 37 patients that had a positive SARS-CoV-2 nasopharyngeal sample.

Fiche technique EurobioPlex SARS-CoV-2 Multiplex

Overall performances on concurrent NP and saliva samples are the following:

Performances on salivary /nasopharyngeal (NP) concurrent samples on the same patients		Salivary samples	
		Saliva SARS-CoV-2 positive	Saliva SARS-CoV-2 negative
Nasopharyngeal samples	NP SARS-CoV-2 positive	34	3
	NP SARS-CoV-2 negative	0	66

Tests were performed on CFX96™ (Bio-Rad).

Sensitivity on saliva (from patients with a NP positive result)/positive concordance:

91.9 % (34/37) (IC₉₅ : 83.1 , 100)

Specificity on saliva (from patients with a NP negative result)/negative concordance:

> 99% (66/66) (IC₉₅ : 95 - 100)

Overall concordance between results on concurrent salivary and nasopharyngeal samples: 97.1% (100/103) (IC₉₅ : 93.8 , 100)

The positive concordance of EBX-041 test between nasopharyngeal and salivary samples is consistent with the HAS recommendations, that require a minimum clinical sensitivity of the salivary RT-PCR test of at least 80% (considering the lower confidence interval limit at 95% of the estimated sensitivity).

Fiche technique EurobioPlex SARS-CoV-2 Multiplex

SPECIFICITIES OF THE REAL-TIME PCR T-COR 8®-IVD INSTRUMENT

T-COR 8®-IVD is a real-time PCR instrument with 8 independent wells, which can be programmed independently, patient by patient, regarding thermic protocol, assays, and time of start of the run.

The 8 wells can also be used simultaneously with the same assay.

Note: Mix well after adding the sample in the tube with pipetting back and forth, avoiding the formation of bubbles and make sure that the volume of liquid is located at the bottom of the tube. Close each tube tightly after each deposit of the control or sample to avoid contamination.

Controls

On T-COR 8®-IVD, the positive and negative controls must be tested at least when a new kit is being opened.

After this initial control, the internal control allows to check that the PCR is working properly, and that there is no PCR inhibition.

Fluorophores

Four combinations of fluorophores are available on T-COR 8®-IVD.

For all EBX, such as EBX-041, FAM/HEX/Texas Red/Cy5 combination is used. If some channels are not used, do not consider this channel for results analysis.

Use of Barcodes (available on pages 45-47)

1- Select Menu > New Run

2- Select Barcode

3- Scan the appropriate Barcode on the right side of the instrument:

- For the positive control (Barcode EBX-041 Pos Ctrl),
- For the negative control (Barcode EBX-041 Neg Ctrl),
- For a patient sample (Barcode EBX-041)

The instrument randomly selects one of the 8 available wells. Follow the instructions on the instrument.

4- Lift up the lid of the selected well (well x), check that the blue LED light is on, and select "Yes"

5- Place the tube in the corresponding well and select "Next".

6- If you wish to name the sample (optional), select "sample x", name it and select "Accept".

7- Select "Next"

8- To add/test another control or sample, select « Add well » and go back to point 3-.

Fiche technique EurobioPlex SARS-CoV-2 Multiplex

9- Select "Start run".

Note: The name of a sample or of a run cannot be modified or added after the run is finished. It must be done before or during the run.

Presentation and Interpretation of results

Ct values and amplification curves are available in the SmartCT™ Table (Ct values on each channel), and "Ct versus PCR cycles" graphs, both being available in real-time while the run is ongoing.

An automatic analysis for positive and negative controls as well as patients' samples is available in "Interpretations", at the end of the run, in the "View" window.

For positive and negative controls, the following results are possible:

« Neg Ctrl Fail »: Not valid

« Neg Ctrl Valid »: Valid

« Pos Ctrl Fail »: Non valid

« Pos Ctrl Valid »: Valid

For patients status:

"Detected": Positive for at least one target → green box

"Not detected": Negative → red box




« Invalid »: Invalid result → yellow box: retest

Example of visual automatic interpretation of results on T-COR 8®-IVD

Summary								
Well	Sample	Assay	FAM	HEX	TxR	Cy5	Call	Note
1	ECH POS	EBX-041	30.6	30.9	29.4	28.9	Detected • Sars-Cov-2	Sars-Cov-2 positive
2	ECH NEG	EBX-041				16.3	Not Detected	Sars-Cov-2 negative
3	ECH NEG	EBX-041					Invalid	Retest
5	Ctrl POS	EBX-041 Pos Ctrl	35.7	36.8	35.7		Invalid	Pos Ctrl Fail
6	Ctrl POS	EBX-041 Pos Ctrl	24.5	24.7	23.5	23.0	Detected • Pos Ctrl	Pos Ctrl Valid
7	Ctrl NEG	EBX-041 Neg Ctrl	37.0			32.5	Invalid	Neg Ctrl Fail
8	Ctrl NEG	EBX-041 Neg Ctrl					Detected • Neg Ctrl	Neg Ctrl valid

Exemple of amplification curves

Fiche technique EurobioPlex SARS-CoV-2 Multiplex

Fiche technique EurobioPlex SARS-CoV-2 Multiplex

Barcodes for EBX-041 for use on T-COR8®-IVD

EBX-041 Pos Ctrl
Positive Control CP

EBX-041 Neg Ctrl
Water = Negative Control (CN-H2O)

EBX-041

15. Bibliography

Chan JF, Kok KH, Zhu Z, Chu H et al. Genomic characterization of the 2019 novel human-pathogenic coronavirus isolated from a patient with atypical pneumonia after visiting Wu-han. Emerg microbes infect. 2020 Dec; 9(1):221-236.

Cohen J, Normile D. New SRAS-like virus in China triggers alarm, Science. 2020 Jan 17;367(6475):234-235.

Drosten C, Muth D, Corman VM, Hussain R et al. An observational, laboratory-based study of outbreaks of middle east respiratory syndrome coronavirus in Jeddah and Riyadh, kingdom of Saudi Arabia, 2014. Clin infect dis. 2015 Feb 1; 60(3):369-77.

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16. Quality control

In accordance with Eurobio's ISO EN 13485 certified quality management system, each batch of EurobioPlex SARS-CoV-2 is tested according to predefined specifications to ensure consistent product quality.

17. Waste disposal

Eliminate all waste in accordance with the legislation on DASRI.

18. Incident report

Serious incidents related to the system must be reported to Eurobio Scientific and to the competent authority of the Member States in which the user and/or the patient is registered.

19. Technical assistance

For assistance with our products, please contact our technical support.

Eurobio Scientific's customer service can be contacted by e-mail at adv@eurobio-scientific.com or by phone at +33 (0)1.69.79.18.26.

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