Instructions for use EurobioPlex FluCoSyn

EurobioPlex FluCoSyn


REAL-TIME RT-PCR

For **qualitative** real-time RT-PCR

REF

EBX-042

EBX-042-25 ; EBX-042-50 ; EBX-042-100 ; EBX-042-200 ; EBX-042-600

25/50/100/200/600 reactions

Reference EBX-042 Version 8.03 – 14/02/2025

Validated on:

- CFX96™ Real Time PCR detection system (Bio-Rad) with analysis on CFX Manager version 3.1 (Bio-Rad)
- LightCycler® 480 (Roche) with analysis on LightCycler® 480 software v1.5 (Roche)
- T-COR 8®-IVD (Tetracore Inc.) with T-COR 8 Per-Well SmartCT™ (auto v1) software
- CFX Opus 96 (Bio-Rad) with analysis on CFX Maestro version 2.2 (Bio-Rad)

Storage conditions:

Keep all reagents between -15°C and -22°C until use and after first use

Instructions for use

Available on request at www.info@eurobio-scientific.com

Instructions for use EurobioPlex FluCoSyn

Table of content

English.....	1
Français	27

1. Introduction	3
2. Purpose of the system	4
3. Symbols.....	5
4. Principle of detection	6
5. Content of the kit.....	7
6. Storage.....	7
7. Materials required not provided	7
8. Real-time PCR instrument	8
9. Cautions and note.....	8
10. Procedure	9
11. Validation of the experiment	13
12. Data analysis and interpretation	14
13. Performance analysis	15
14. Bibliography.....	25
15. Quality control.....	25
16. Waste disposal.....	25
17. Incident report.....	25
18. Technical assistance	26

Instructions for use EurobioPlex FluCoSyn

1. Introduction

The influenza A virus (Influenza A, Flu A) and the influenza B virus (Flu B) belong to the *Orthomyxoviridae* family. Respiratory Syncytial Virus A and B (RSV A and B) are members of the *Paramyxoviridae* family. SARS-CoV-2 virus belongs to the *Coronaviridae* family and subgenus *Sarbecovirus*. The latter emerged in China at the end of 2019, in the city of Wuhan, and caused an unprecedented worldwide pandemic in the first half of 2020. All three are RNA viruses, with human-to-human transmission mainly infecting the lungs and respiratory tract (cough and sneeze droplets from infected individuals).

The genome of the SARS-CoV-2 virus consists of 29903 ribonucleic acid (RNA) bases. The SARS-CoV-2 is the seventh coronavirus identified as infecting humans, after human coronaviruses (HCoV) 229E, NL63, OC43, HKU1, SARS-CoV (coronavirus causing severe acute respiratory syndrome) and MERS-CoV (coronavirus inducing Middle Eastern respiratory syndrome). The SARS-CoV-2 sequences present similarities with those of betacoronaviruses found in bats. The SARS-CoV-2 virus is genetically distinct from other human coronaviruses such as those related to SARS and MERS.

These viruses are a major public health problem worldwide (morbidity rate, and health and social costs). They are seasonal, mainly autumn/winter. Flu induces epidemic, with possible serious consequences especially in the elderly. RSV is the main cause of bronchiolitis in young children. SARS-CoV-2 led to a pandemic. Clinical signs usually disappear after a few days, but the symptoms are variable, commonly covering an influenza-like illness such as fever, headache, cough, runny nose, pharyngitis, myalgia, asthenia, breathing difficulty, but in the most severe cases, complications occur (pneumonia, severe respiratory syndrome, and dehydration) and can be fatal. Antivirals are currently available for influenza and RSV, but not for SARS-CoV-2, whose treatment and prophylaxis are being studied worldwide. To date, only the annual influenza vaccine has a Marketing Authorization.

In July 2020, the actual case-fatality rate associated with SARS-CoV-2 in France is between 0.5% and 1% according to the WHO, which is well above the seasonal influenza rate of 0.1%, but without any measurement compared to that associated with SARS-CoV-1 or MERS, which are 10% and 30% respectively. SARS-CoV-2 is highly contagious with more than 16 million cases worldwide at the end of July 2020. According to WHO reports, the pandemic is very active in the USA, the Middle East, and Central and South Asia, and is on the rise in countries where influenza usually circulates at this time of year (South Africa, Australia).

Respiratory viral infections are generally diagnosed clinically on the basis of symptoms and local epidemic situation. The 3 types of viruses detected by EurobioPlex FluCoSyn co-circulate. Despite the fact that the pathogenic agents lead to specific manifestations, each can lead to various symptoms of the respiratory tract, with different degree of severity, depending on the host. Moreover, the therapeutic approach is different depending on the virus. It is therefore crucial to be able to differentiate which type of virus is at the origin of the symptoms, within an epidemic and pandemic context. This is achieved with EurobioPlex EBX-042.

Instructions for use EurobioPlex FluCoSyn

2. Purpose of the system

The EurobioPlex FluCoSyn is a real-time reverse transcription polymerase chain reaction (*RT-PCR*) *amplification test in an RNA nucleic acid extract designed for in vitro* use for the qualitative determination of the presence or absence of the following respiratory viruses:

- SARS-CoV-2, an emerging virus, which appeared in China in late 2019 in the city of Wuha
- Influenza A virus (Flu A) and Influenza B virus (Flu B)
- Respiratory Syncytial Virus A and B (RSV A and B)

EurobioPlex FluCoSyn has been designed to detect Flu and RSV without differentiation of types A and B, and all known SARS-CoV-2 sequences known to date, and by *in silico* alignment, particularly allowing to exclude other benign coronaviruses or SARS-CoV-1 or MERS coronaviruses. The determination of SARS-CoV-2 status is based on the detection of 2 targets coupled to the same fluorophore: one in the RdRp gene and the other target in the N gene.

The test is indicated to help diagnose presumptive infection with these three viruses associated with respiratory syndrome in humans, or to complement a proven or undetermined diagnosis, specifically distinguishing the SARS-CoV-2 coronavirus from RSV or Influenza respiratory viruses. The diagnosis should always be made by medical personnel and in the clinical, historical and symptomatic context of the patient.

The EurobioPlex FluCoSyn test is an *In Vitro* diagnostic medical device, it must be used by qualified medical biology analysis laboratory personnel. The device should not be recycled after use.

The RNA extract from any respiratory specimen is the starting material for the EurobioPlex FluCoSyn. It is the responsibility of the user to use extraction methods appropriate to the specimens tested.

The kit has been validated on the following types of respiratory specimen:

- Nasopharyngeal swab ;
- Nasopharyngeal aspirate ;
- Bronchial aspiration ;
- Bronchoalveolar fluid ;
- Sputum.

Instructions for use EurobioPlex FluCoSyn

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Manufacturer

CE marked product

In Vitro Diagnostic Medical Device

Instructions for use

Store away from light

Do not use if packaging is damaged

Caution

Instructions for use EurobioPlex FluCoSyn

4. Principle of detection

EurobioPlex FluCoSyn is a test using reverse-transcription and real-time amplification of viral RNA of Influenza (A and B), RSV (A and B) and SARS-CoV-2 (one target in RdRp gene and the other in N gene), based on a single Multiplex RT-PCR test, detecting all 3 viruses in a signal well.

The kit contains one oligomix to detect the 3 targets, as well as an endogenous control. The test is performed from extracted RNA from a sample, using one RT-PCR reaction in a single distinct well/tube.

The endogenous control allows to check for variations that may occur during the RNA extraction step of biological samples and real-time RT-PCR amplification. Thus, it ensures that a negative result may be due to a bad RNA extraction, and/or due to the presence of PCR inhibitors in too large quantities.

Influenza A and B are detected with FAM-labeled probes, RSV A and B are detected with an HEX labeled probe, and SARS-CoV-2 with Texas Red labeled probes.

Endogenous control is detected using a Cy5 labeled probe. Probes emit a specific fluorescence following their hydrolysis during the elongation of the amplification product. The measurement of the intensity of real-time fluorescence correlates with the accumulation of amplification products.

Table 1: Detection of targets by fluorophores

Targets	Fluorophore	Excitation	Emission
Target 1 / Flu	FAM	495 nm	515 nm
Target 2 / RSV	HEX	535 nm	555 nm
Target 3 / RdRP1/RdRP2 gene	Texas Red	585 nm	605 nm
Endogenous internal control	Cy5	650 nm	670 nm

Equivalent channels on different real time PCR cyclers:

- Channel **FAM** (Systems ABI, SmartCycler II, Systems Mx, CFX96™/Chromo4, T-COR8-IVD, LC480), Channel Green (RotorGene), CFX® Opus 96 (Bio-Rad).
- Channel **HEX** (Chromo 4/CFX96, Systèmes Mx, T-COR 8-IVD, LC480), Canal VIC (Systèmes ABI), Canal Alexa532 (SmartCycler II), Canal Yellow (RotorGene), CFX® Opus 96 (Bio-Rad).
- Channel **Texas Red** (Systems ABI, Chromo 4/CFX96, Systèmes Mx, T-COR 8-IVD), LC Red 610 (LC480), Canal Orange (RotorGene), CFX® Opus 96 (Bio-Rad).
- Channel **Cy5** (CFX96™/Chromo4, Systems ABI, Systems Mx, T-COR8-IVD, LC 480), Channel Alexa647 (SmartCycler II), Channel Red (RotorGene), CFX® Opus 96 (Bio-Rad).

Instructions for use EurobioPlex FluCoSyn

Note 1: On LC480, apply color compensation for the following channels: FAM/HEX-VIC/Texas Red/Cy5 (FAM : 483-533 / HEX-VIC : 523-568 /Texas Red : 558-610 / Cy5 : 615-670).

Note 2: During the analysis using LC480 Software release 1.5.0, do not leave Noiseband Auto on Texas red and HEX channels, but apply the NoiseBan Fluor manually at 0,43 for Texas Red channel and 5,5 for HEX channel (the operator need to adjust the baseline).

5. Content of the kit

The real-time RT-PCR kit EurobioPlex FluCoSyn is ready to use and contains the reagents and enzymes necessary for the detection of these viruses (Table 2).

Fluorescence is emitted and measured individually by an optical system during the PCR process. Detection of the amplified fragments is performed by a fluorometer using the channels listed in Table 1.

Table 2: Content of the kit

Cap color	Content of the kit	25 reactions	50 reactions	100 reactions	200 reactions	600 reactions	Reconstitution
Red	Enzymes	5x25 µL	10x25 µL	525 µl	800 µl	2 x 1200 µl	Ready to use
Transparent	Oligomix	5x20 µL	10x20 µL	420 µl	700 µL	2 x 950 µl	Ready to use
Yellow	Positive Control CP	5x15 µL	10x15 µL	160 µl	320 µl	320 µl	Ready to use
Blue	Water = negative control (CN-H2O)	1500 µL	1500 µL	1500µl	1500 µl	2 x 1500 µl	Ready to use

Oligomix: contains primers and probes for the 3 targets and for the endogenous control.

6. Storage

All reagents should be stored between -15°C and -22°C.

All reagents can be used until the expiration date indicated on the kit label.

The sensitivity of the assay may decrease if the kit components undergo multiple freeze/thaw cycles. The kit can be used after the initial opening for up to 5 freeze/thaw cycles.

7. Materials required not provided

- ◇ Biological cabinet
- ◇ Real time PCR instrument
- ◇ Centrifuge for microtubes
- ◇ Vortex
- ◇ Plates / tubes for real-time PCR reaction
- ◇ Micropipettes

Instructions for use EurobioPlex FluCoSyn

- ◇ DNase-free and RNase-free filter tips for micropipettes
- ◇ Sterile microtubes
- ◇ Gloves (talc-free)

8. Real-time PCR instrument

The EurobioPlex FluCoSyn kit has been developed and validated for use with the following real-time PCR instruments:

- CFX96™ Real Time PCR detection system (Bio-Rad) with analysis on CFX Manager version 3.1 (Bio-Rad)
- LightCycler®480 with analysis on LightCycler® 480 software v1.5 (Roche)
- T-COR 8®-IVD (Tetacore Inc.) with T-COR 8 Per-Well SmartCTTM (auto v1) software
- CFX Opus 96 (Bio-Rad) with analysis on CFX Maestro version 2.2 (Bio-Rad)

9. Cautions and note

Read these instructions carefully before starting the procedure.

- ◇ This experimentation should be performed by medical laboratory technicians.
- ◇ The local and national biosafety regulations in place for the detection of SARS-CoV-2 must be followed strictly at all times, especially in laboratories and with laboratory equipment in agreement.
- ◇ Ensure that instruments have been installed, calibrated and maintained in accordance with the manufacturer's recommendations.
- ◇ It is the responsibility of the user to use equipment other than that validated, and in such cases performance is not guaranteed.
- ◇ Clinical samples should be considered as potentially infectious material and should be prepared in a laminar flow cabinet.
- ◇ This experiment should be performed according to good laboratory practice.
- ◇ Do not use this kit after the stated expiration date.
- ◇ The kit is shipped under dry ice, and the kit components should arrive frozen. If one or more of the components arrive thawed, or if the tubes have been damaged in transit, contact Eurobio Scientific.
- ◇ Avoid freeze/thaw cycles of reagents as this may lead to a decrease in assay sensitivity.
- ◇ Once the reagents have been thawed, centrifuge the tubes briefly before use.
- ◇ The use of ice or an ice pack is recommended in case of long delays due to, for example, a large number of samples to be processed or high temperatures.
- ◇ It is recommended to define three separate work areas: 1) RNA isolation, 2) Preparation of the reaction mixture, and 3) Amplification/Detection of the amplified products.

Instructions for use EurobioPlex FluCoSyn

- ◇ The fluorescent probes in the oligomix are light sensitive. Any prolonged exposure of the oligomix to light should be limited to the technical time required to prepare the PCR plate.
- ◇ It is recommended that the positive control be opened and handled separately from the biological samples to be tested and the other kit components to avoid cross-contamination.
- ◇ Distinct lab coat and gloves (talco-free) should be worn in each work area.
- ◇ Pipettes, reagents and other working materials should not be moved between these areas.
- ◇ Special care should be taken to maintain the purity of reagents and reaction mixtures.
- ◇ The internal endogenous control detects a detectable cell target in all samples of human origin but not in the CN-H2O negative control provided in the kit. The absence of a signal in the negative control prevents cross-contamination from occurring.
- ◇ Appropriate RNA preparation/extraction methods for quality RNA production and RT-PCR application should be used, especially to avoid contamination with RNAs.
- ◇ Use filter tips for micropipettes, RNase-free and DNase-free.
- ◇ Do not pipette with the mouth. Do not eat, drink or smoke in the laboratory.
- ◇ Avoid aerosols.

10. Procedure

a. Samples collection

- ◇ Collect samples in sterile tubes.
- ◇ It is the responsibility of the user to control their own sample collection, transport, storage and extraction conditions to ensure that RNA extraction using appropriate systems produces RNA of high quality.
- ! It is recommended that samples are extracted immediately, or stored according to the sample storage recommendations before extraction (Table 3).

Table 3: Storage recommendations before use

Recommendations for maximum storage of samples before extraction	
Ambient temperature	2 h
4°C	72 h
-20°C (preferably -80°C)	Long-term storage

Caution	
	◇ The user may refer to the recommendations issued by the World Health Organization or the French National Authority for Health for the proper storage of samples.

Instructions for use EurobioPlex FluCoSyn

	◇ Extracted RNA should be stored at -80°C. ◇ The transport of clinical samples is subject to local regulations for the transport of infectious agents.
--	--

b. RNA extraction

It is the responsibility of the user to ensure that the nucleic acid extraction system used is compatible with real-time RT-PCR technology. We recommend using virus RNA extraction methods suitable for respiratory specimens, and to refer to the instructions of the supplier of the extraction kit used.

In the EBX-042 kit, the endogenous internal control on the Cy5 channel is already present in the clinical sample. Properly detected after amplification, the endogenous internal control is used to ensure the quality of the sample and the extraction. The performance described in this document was obtained with the following extraction technique: Maelstrom 9600 (TanBEAD) and the OptiPure Viral Auto Plate kit (ref: W665A46).

c. Real-time RT-PCR

General note:

- The positive control contains high concentrations of matrix. Handling must be done carefully to avoid contamination.
- To confirm that the RT-PCR is working correctly, it is necessary to test the positive control, as well as the negative control (water supplied = CN-H₂O) (see II-2/6 of the real-time RT-PCR protocol).

Instructions for use EurobioPlex FluCoSyn

Diagram of the procedure:

1 - PREPARATION OF MASTERMIX

Number of reactions	N+3
Enzymes	(N+3) x 3.75 µl
Water	(N+3) x 3.25 µl
Oligomix	(N+3) x 3 µl
Total Volume Mastermix	(N+3) x 10 µl

* N+1 on T-COR 8⁺-IVD

2 - PREPARATION OF REACTIONS

Sample

10 µl Mastermix
+
5µl RNA sample

Positive Control

10 µl Mastermix
+
5µl CP

Negative Control

10 µl Mastermix
+
5µl Molecular biology water (CN-H2O)

3 – REAL TIME PCR INSTRUMENT

Program	Temperature	Duration	Cycle(s)	
Reverse Transcription	45°C	5 min	1	-
Denaturation	98°C	20 sec	1	-
Amplification	98°C	3 sec	40	-
	58°C	10 sec		Acquisition of fluorescence

Instructions for use EurobioPlex FluCoSyn

d. Detailed procedure

- 1) Ensure that the reagents are completely thawed. Homogenise the enzyme tubes, and vortex the Oligomix and CP for about 15 seconds, then centrifuge briefly.
- 2) Prepare the Mastermix as below. N is the number of reactions (including positive and negative controls. Plan to prepare enough reagents for at least N+3 reactions (*N+1 on T-COR8).

Number of reactions	N+3*
Enzymes	(N+3) x 3.75 µl
Water	(N+3) x 3.25 µl
Oligomix	(N+3) x 3 µl
Total Volume Mastermix	(N+3) x 10 µl

- 1) Homogenise the Mastermix prepared in 2) and centrifuge briefly
- 2) Dispense 10 µL of Mastermix into separate microplate tubes/wells using a micropipette and filter tips.
- 3) Add 5 µL of extracted RNA sample.
- 4) In parallel perform the following controls:
 - Positive Control:
 - 10 µL of Mastermix + 5 µL of CP.
 - Negative Control:
 - 10 µL of Mastermix + 5 µL supplied water (CN-H2O).
- 5) Immediately close the tubes, or plate with an adhesive film to avoid contamination.
- 6) Centrifuge briefly to collect the reaction mixture at the bottom of the tubes or microplate wells.
- 9) Run the following program on the real-time PCR instrument.

Program	Temperature	Duration	Cycle(s)	
Reverse Transcription	45°C	5 min	1	-
Denaturation	98°C	20 sec	1	-
Amplification	98°C	3 sec	40	-
	58°C	10 sec		Acquisition of fluorescence

Note1: On CFX96™ (Bio-Rad), start the run using the v1.6 or later version of CFX Manager Software, and analyze with v 3.1 (see § Validation of the Experiment). **Please ensure that white PCR plates are used for appropriate reading in all channels.**

Note 2: On LightCycler® 480 (Roche), apply color compensation for the following channels: FAM/HEX-VIC/Texas Red/Cy5 (FAM : 483-533 / HEX-VIC : 523-568 /Texas Red : 558-610 / Cy5 : 615-670).

Note 3 : During the analysis using LC480 Software release 1.5.0 apply a noise band of 2,1 for the FAM channel and 4,5 for the HEX channel.


Instructions for use EurobioPlex FluCoSyn

11. Validation of the experiment


Post-acquisition data analysis on a CFX96 PCR instrument (Bio-Rad) must be performed using version 3.1 of the CFX Manager software (Bio-Rad). In order to change to this version from a run started on an earlier version, please follow the procedure below: At the end of the run, the data file with the extension .pcrd must be opened and processed with CFX Manager (Bio-Rad) version 3.1.

If the run was started with CFX Manager v1.6 for example, to open a data file with CFX Manager v3.1, click on the CFX Manager v3.1 icon. The home screen appears.

If a CFX Opus 96 PCR instrument is used, the analysis should be performed on the CFX Maestro software version 2.2. The remainder of the analysis is identical to that to be performed with a CFX96, as shown below.


Click on "File" and select "Open" then "Data File".

Instructions for use EurobioPlex FluCoSyn

Select the file you need to analyse and click on "Open".

The "Drift correction" option must be applied in the "Settings" tab as shown in the image below: click on the "Settings" tab, then on "Baseline Setting" and on "Apply Fluorescence Drift Correction".

Once this is done, the analysis can start.

The results for the controls must be the following (Table 4); otherwise, the experiment is not valid.

On T-COR 8®-IVD, these controls must be performed at least during the first use of a new box of kit.

Table 4: Validation of the run

Positive Control	
FAM	Ct ≤ 30
HEX	Ct ≤ 30
Texas Red	Ct ≤ 30
Cy5	Ct ≤ 30
Negative Control	
FAM	Undetermined Ct
HEX	Undetermined Ct
Texas Red	Undetermined Ct
Cy5	Undetermined Ct

12. Data analysis and interpretation

For clinical samples, the following results are possible:

***SARS-CoV-2 (Texas Red), Flu (FAM) ou RSV (HEX): Ct < 40**

*Ct cut off for positivity of samples for T-COR 8®-IVD is: + Positive => Ct positive (≤ 45) for the 3 channels, but automatic interpretation is provided with Barcodes.

Instructions for use EurobioPlex FluCoSyn

Table 5:

PCR Signal Flu and/or RSV and/or SARS-CoV-2	CI-PCR Signal	Presence of virus	Test validity/comment
FAM and/or HEX and/or Texas Red	Cy5		
+	+	Yes	VALID
-	+	No	VALID
+	-	Yes	VALID Possible RT-PCR inhibition or problem with extraction that does not prevent the detection of viruses.
-	-	No possible interpretation	RT-PCR inhibition or problem with extraction → dilute first 5 x the sample and test again; if necessary take a new sample and/or redo an extraction

Limitations on use and interpretation:

- ❖ All samples should be treated as potentially infected with SARS-CoV-2, and local biosafety regulations should be carefully followed.
- ❖ Interpretation of results should consider the possibility of false negatives and false positives.
False negatives may be due to:
 - Inadequate collection of samples, or incorrect storage,
 - Samples outside the viremia phase,
 - Inappropriate extraction conditions or use of non-validated PCR instruments,
 - Experimentation not respecting all elements of these instructions for use.
 False positives may be due to:
 - Contamination due to mishandling of high positive samples, the positive control, or the PCR amplification products,
 - Failure to follow the procedure described in these instructions for use, especially to avoid sources of contamination.
- ❖ All results should be interpreted by medical professionals in the context of the patient's clinical situation, history and symptoms.

13. Performance analysis

A dilution range was performed using a plasmid mixture ranging from 10^5 to 1 copy/ μ L. This dilution range was used to determine the detection limit (at 100% detection) of the EBX-042 kit.

In addition, a detection limit was determined on the WHO 1st International Standard for SARS-CoV-2 RNA (NIBSC code 21/146) and on a SARS-CoV-2 synthetic RNA of Twist Biosciences.

Instructions for use EurobioPlex FluCoSyn

Detection limitation/analytical sensitivity on CP:

SARS-CoV-2 : 7.5 copies/ μ L

Influenza/Flu : 5 copies/ μ L

RSV : 5 copies/ μ L

CI Cy5 : 7.5 copies/ μ L

Detection limitation/analytical sensitivity on WHO 1st International Standard for SARS-CoV-2 RNA:

SARS-CoV-2: 12.5 copies/ μ L

Detection limitation/analytical sensitivity on a SARS-CoV-2 synthetic RNA (Twist Biosciences, USA) :

SARS-CoV-2: 6.25 copies/ μ L

Sensitivity study on 106 known positive samples diluted in known samples Negatives at 3x LOD: 100% sensitivity (detection of 106 samples / 106)

- **Interfering substances:** 7 substances from the AcroMetrixTM Inhibition Panel REF 956400 kit, Oxymetazoline and Zanamivir. A high concentration, which may be contained in some drugs, was added to the positive sample before extraction: No impact of the tested interfering substances was observed on the result obtained.

Signal variability on FAM, HEX, Texas Red and Cy5 channels

- Intra-experimental variability:

Average intra-batch CV (%)	Endogenous ctrl Gene	Influenza	RSV	SARS-CoV-2
CFX96	1,62	1,19	1,11	2,39
LC480	1,37	1,09	1,97	1,00
T-COR-8-IVD	2,97	1,94	2,37	1,96
CFX [®] Opus 96	0.39	0.33	0.76	0.64

CV: coefficient of variation

- Inter-batch variability:

Average inter-batch CV (%)	Endogenous ctrl Gene	Influenza	RSV	SARS-CoV-2
CFX96	2,75	4,06	1,18	4,85
LC480	0,45	1,80	1,95	1,13
T-COR-8-IVD	0,73	1,15	3,20	0,93
CFX [®] Opus 96	5.75	1.98	1.55	0.92

Instructions for use EurobioPlex FluCoSyn

CV: coefficient of variation

Diagnostic specificity

This validation included:

- 508 SARS-CoV-2 negative samples, pretested and characterized by CE approved tests.
- 40 samples negative for all known coronaviruses, including respiratory panel samples characterized for other viruses to test for cross-reactions.
- 3 samples positive for CoV NL63, 2 samples positive for CoV OC43, 1 sample positive for CoV 229E and 1 sample positive for CoV HKU1.

Number of sample	Positive for	Coronavirus	Status SARS-CoV-2 EBX-042
1	Adenovirus	Negative	Negative
5	ParaInfluenza	Negative	Negative
1	Rhinovirus	Negative	Negative
1	Bocavirus	Negative	Negative
2	Legionella pneumoniae	Negative	Negative
3	Mycoplasma pneumoniae	Negative	Negative
3	Mycoplasma tuberculosis	Negative	Negative
2	Bordetella pertussis	Negative	Negative
1	Parechovirus Type 1	Negative	Negative
1	Parechovirus Type 2	Negative	Negative
1	Parechovirus Type 3	Negative	Negative
1	Parechovirus Type 4	Negative	Negative
1	Parechovirus Type 5	Negative	Negative
1	Parechovirus Type 6	Negative	Negative
3	Enterovirus	Negative	Negative
1	Streptococcus pyogenes	Negative	Negative
3	Streptococcus pneumoniae	Negative	Negative
2	Staph.epiderm/oralis/mi-	Negative	Negative
3	Candida albicans	Negative	Negative
2	Pseudomonas aeruginosa	Negative	Negative
3	CoV NL63	Positive	Negative
2	CoV OC43	Positive	Negative

Instructions for use EurobioPlex FluCoSyn

1	CoV 229E	Positive	Negative
1	CoV HKU1	Positive	Negative
1	Metapneumovirus	Positive	Negative
1	Human papillomavirus	Positive	Negative

In silico analysis was used and demonstrated the specificity of the primers and probes for SARS-CoV-2. For example, on the sequences of MERS and SARS viruses, no significant homology of SARS-CoV-2 primers and probes that could amplify these viruses was found.

No aspecificity or cross reaction was detected.

The kit is 100% specific.

Diagnostic sensibility and specificity

The validation of the performance was performed on :

555 negative samples	
508 SARS-CoV-2 negative samples	Characterised with EBX-041 kit, CE marked
47 positive respiratory samples	RSV-A/B, Parainfluenzae, Rhinovirus, Metapneumovirus, CoV NL63, CoV OC43, CoV 229E, CoV KU1

335 positive samples	
270 SARS-CoV-2 positive samples	Characterized by EurobioPlex EBX-041 CE-IVD kit
29 VRS positive samples	Characterized with the Seegene Allplex SARS-Cov-2 Flu A / Flu B RSV Assay
37 Influenza positive samples	characterized with the Seegene Allplex SARS-Cov-2 Flu A / Flu B RSV Assay

For this study, sample extraction was performed on Maelstrom 9600 (TANBead) with the OptiPure Viral Auto Plate kit (ref: W665A46), and RT-PCR on CFX96 (Bio-Rad). And tests were performed on CFX96™ (Bio-Rad).

Overall performances are:

		EBX-042		
		Flu		
		POSITIVE	NEGATIVE	TOTAL
Pretested	POSITIVE	34	0	34
	NEGATIVE	3	806	809

Instructions for use EurobioPlex FluCoSyn

		RSV		
		POSITIVE	NEGATIVE	TOTAL
Pretested	POSITIVE	29	0	29
	NEGATIVE	0	814	814
		SARS-CoV-2		
		POSITIVE	NEGATIVE	TOTAL
Pretested	POSITIVE	269	0	269
	NEGATIVE	0	574	574

The tests were performed on CFX96™ (Bio-Rad).

	Sensitivity %	Specificity %	Concordance %
FLU	92	> 99	> 99
RSV	> 99	> 99	> 99
SARS-Cov-2	> 99	> 99	> 99

SPECIFICITIES OF THE REAL-TIME PCR T-COR 8®-IVD INSTRUMENT

T-COR 8®-IVD is a real-time PCR instrument with 8 independent wells, which can be programmed independently, patient by patient, regarding thermic protocol, assays, and time of start of the run.

The 8 wells can also be used simultaneously with the same assay.

Note: Mix well after adding the sample in the tube with pipetting back and forth, avoiding the formation of bubbles and make sure that the volume of liquid is located at the bottom of the tube. Close each tube tightly after each deposit of the control or sample to avoid contamination.

Controls

On T-COR 8®-IVD, the positive and negative controls must be tested at least when a new kit is being opened.

After this initial control, the internal control allows to check that the PCR is working properly, and that there is no PCR inhibition.

Fluorophores

Four combinations of fluorophores are available on T-COR 8®-IVD.

For all EBX, such as EBX-042, FAM/HEX/Texas Red/Cy5 combination is used. If some channels are not used, do not consider this channel for results analysis.

Instructions for use EurobioPlex FluCoSyn

Use of Barcodes (available on pages 45-47)

1- Select Menu > New Run

2- Select Barcode

3- Scan the appropriate Barcode on the right side of the instrument:

- For the positive control (Barcode EBX-042 Pos Ctrl),
- For the negative control (Barcode EBX-042 Neg Ctrl),
- For a patient sample (Barcode EBX-042)

The instrument randomly selects one of the 8 available wells. Follow the instructions on the instrument.

4- Lift up the lid of the selected well (well x), check that the blue LED light is on, and select "Yes"

5- Place the tube in the corresponding well and select "Next".

6- If you wish to name the sample (optional), select "sample x", name it and select "Accept".

7- Select "Next"

8- To add/test another control or sample, select « Add well » and go back to point 3-.

9- Select "Start run".

***Note:** The name of a sample or of a run cannot be modified or added after the run is finished. It must be done before or during the run.*

Presentation and Interpretation of results

Ct values and amplification curves are available in the SmartCT™ Table (Ct values on each channel), and "Ct versus PCR cycles" graphs, both being available in real-time while the run is ongoing.

An automatic analysis for positive and negative controls as well as patients' samples is available in "Interpretations", at the end of the run, in the "View" window.

For positive and negative controls, the following results are possible:

« Neg Ctrl Fail »: Not valid

« Neg Ctrl Valid »: Valid

« Pos Ctrl Fail »: Non valid

« Pos Ctrl Valid »: Valid

For patients status:

"Detected": Positive for at least one target → green box

"Not detected": Negative → red box



« Invalid »: Invalid result → yellow box: retest

Instructions for use EurobioPlex FluCoSyn



Example of automatic interpretation of results on T-COR 8®-IVD

Summary							
Well	Sample	Assay	FAM	HEX	TxR	Cy5	Call
1	Sample 1	EBX-042	22.6			16.5	Detected • Influa
2	Sample 2	EBX-042		26.2			Detected • RSV
3	Sample 3	EBX-042			24.9	19.0	Detected • SARS-CoV-2
4	Sample 4	EBX-042				15.3	Not Detected
5	Sample 5	EBX-042					Invalid
6	Ctrl pos	EBX-042 Pos Ctrl	22.5	23.2	24.3	22.3	Detected • Pos Ctrl
7	Ctrl pos	EBX-042 Pos Ctrl	24.9	28.6			Invalid
8	Ctrl neg	EBX-042 Neg Ctrl					Detected • Neg Ctrl
							Note
							RSV, SARS-CoV-2 negative
							Retest for Influa, SARS-CoV-2
							Influa, RSV negative
							Influa, RSV, SARS-CoV-2 negative
							Retest
							Pos Ctrl Valid
							Pos Ctrl Fail
							Neg Ctrl Valid

Exemple of amplification curves


Instructions for use EurobioPlex FluCoSyn


Instructions for use EurobioPlex FluCoSyn

Barcodes for EBX-042 for use on T-COR8®-IVD

EBX-042 Pos Ctrl
Positive control CP


EBX-042 Neg Ctrl
Water = Negative control (CN-H2O)

Instructions for use EurobioPlex FluCoSyn

EBX-042

Instructions for use EurobioPlex FluCoSyn

14. Bibliography

Xiaojing Wu, Ying Cai, Xu Huang, Xin Yu, Li Zhao, Fan Wang, Quanguo Li, Sichao Gu, Teng Xu, Yongjun Li, Binghuai Lu, Qingyuan Zhan, Co-infection with SARS-CoV-2 and Influenza A Virus in Patient with Pneumonia, China, *Emerg Infect Dis.* 2020 Jun;26(6):1324-1326.

Yan Li, Jiangshan Wang, Chunting Wang, Qiwen Yang, Yingchun Xu, Jun Xu, Yi Li Xuezhong Yu, Huadong Zhu, Jihai Liu Characteristics of respiratory virus infection during the outbreak of 2019 novel coronavirus in Beijing. *Int J Infect Dis.* 2020 Jul;96:266-269.

Chan JF, Kok KH, Zhu Z, Chu H *et al.* Genomic characterization of the 2019 novel human-pathogenic coronavirus isolated from a patient with atypical pneumonia after visiting Wuhan. *Emerg microbes infect.* 2020 Dec; 9(1):221-236.

Cohen J, Normile D. New SRAS-like virus in China triggers alarm, *Science.* 2020 Jan 17;367(6475):234-235.

Drosten C, Muth D, Corman VM, Hussain R *et al.* An observational, laboratory-based study of outbreaks of middle east respiratory syndrome coronavirus in Jeddah and Riyadh, kingdom of Saudi Arabia, 2014. *Clin infect dis.* 2015 Feb 1; 60(3):369-77.

Cryopreservation of tissues, cells and biological fluids from care, Recommendations for good practices, www.has-santé.fr.

15. Quality control

In accordance with Eurobio's ISO EN 13485 certified quality management system, each batch of EurobioPlex FluCoSyn is tested according to predefined specifications to ensure consistent product quality.

16. Waste disposal

Eliminate all waste in accordance with the legislation on DASRI.

17. Incident report

Serious incidents related to the system must be reported to Eurobio Scientific.

Instructions for use EurobioPlex FluCoSyn

18. Technical assistance

For assistance with our products, please contact our technical support.

Eurobio Scientific's customer service can be contacted by e-mail at adv@eurobio-scientific.com or by phone at +33 (0)1.69.79.64.80.

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