EurobioPlex Flu A/Flu B/SARS-CoV-2

RT-PCR EN TEMPS REEL

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
Version 2.00 of 2023/02/17

Validated on:

- CFX96™ Real Time PCR detection system (Bio-Rad) with analysis on CFX Manager version 3.1 (Bio-Rad)
- LightCycler®480 with analysis on LightCycler® 480 software v1.5 (Roche)
- T-COR 8®-IVD (Tetacore Inc.) with T-COR 8 Per-Well SmartCT™ (auto v1) software
- MIC® (Bio Molecular Systems) with analysis on MIC software version 2.12.6

Storage conditions:

Keep all reagents between -15°C and -22°C until use and after first use

Instructions for use
Available on www.eurobio-scientific.com

EurobioPlex Flu A/Flu B/SARS-CoV-2

Table of contents

1. Introduction.....	32
2. Purpose of the system.....	31
3. Symbols	34
4. Principle of detection	35
5. Content of the kit	36
6. Storage.....	36
7. Materials required not provided	37
8. Real-time PCR instrument	37
9. Cautions and note	37
10. Procedure	38
11. Validation of the experiment	43
12. Data analysis and interpretation	45
13. Performance analysis	46
14. Bibliography.....	53
15. Quality control.....	53
16. Waste disposal.....	54
17. Incident report.....	57
18. Technical assistance	57

1. Introduction

Influenza A and B viruses belong to the Orthomyxoviridae family. The SARS-Cov-2 virus belongs to the family Coronaviridae and the subgenus Sarbecovirus. The latter appeared in China at the end of 2019, in the city of Wuhan, is the cause of an unprecedented global pandemic in the first half of 2020. All three are RNA viruses, of human-to-human transmission infecting mainly the lungs and respiratory tract (cough droplets and sneezes of infected subjects).

The SARS-CoV-2 genome consists of 29903 ribonucleic acid (RNA) bases. SARS-CoV-2 is the seventh identified coronavirus infecting humans, after human coronaviruses (HCoV) 229E, NL63, OC43, HKU1, SARS-CoV (severe acute respiratory syndrome coronavirus) and MERS-CoV (Middle East respiratory syndrome coronavirus). The sequences of SARS-CoV-2 show similarities to those of the betacoronaviruses found in bats. SARS-CoV-2 is genetically distinct from other human coronaviruses such as those linked to SARS and MERS.

These 3 viruses are a major public health problem worldwide (morbidity rate, and health and social costs). They occur seasonally, essentially in the fall and winter. Influenza induces epidemics, which can have serious consequences, particularly in the elderly. SARS-CoV-2 is at the origin of a pandemic. Clinical signs usually disappear after a few days of evolution, but the clinical picture is variable, commonly covering an influenza-like syndrome with fever, headache, cough, runny nose, pharyngitis, myalgia, asthenia, respiratory difficulty, but in the most severe cases, complications occur (pneumonia, severe respiratory syndrome, dehydration) and can be fatal. Antivirals are currently available for influenza but not for SARS-Cov-2, whose treatment and prophylaxis are under study worldwide. To date, only the annual influenza vaccine has a Marketing Authorization.

In July 2020, the case fatality rate associated with SARS-CoV-2 in France is between 0.5% and 1% according to the WHO, which is well above that of seasonal influenza (0.1%), but in no way comparable to that associated with SARS-CoV-1 or MERS, which are 10% and 30% respectively. SARS-CoV-2 is highly contagious with more than 16 million cases worldwide at the end of July 2020. According to WHO reports, the SARS-CoV-2 pandemic is very active in the Americas, the Middle East, and Central and South Asia, and is on the rise in countries where influenza usually circulates at this time of year (South Africa, Australia).

Respiratory viral infections are usually diagnosed clinically based on symptoms and local epidemiology. All 3 types of viruses detected by EBX-043 circulate at the same time. Although the specific pathogens cause characteristic clinical manifestations, each can cause numerous viral respiratory tract syndromes, with varying degrees of severity depending on the host. In addition, the therapeutic management is different depending on the viruses responsible for the infection. It is therefore essential to be able to differentiate which type of virus is responsible for the observed symptomatology, in the context of epidemics and pandemics, which is what EBX-043 allows.

EurobioPlex Flu A/Flu B/SARS-CoV-2

2. Purpose of the system

The EurobioPlex Flu A/Flu B/SARS-CoV-2 is a real-time reverse transcription-polymerase chain reaction (RT-PCR) assay in RNA nucleic acid extract for in vitro use, designed for the detection of the following respiratory viruses:

- Influenza A virus (Influenza A or Flu A),
- Influenza B virus (Influenza B or Flu B),
- SARS-CoV-2, an emerging virus that appeared in China in late 2019 in the city of Wuhan.

The EurobioPlex Flu A/Flu B/SARS-CoV-2 has been designed to detect influenza A and B viruses without differentiating between the A H1N1 and H3N2 virus subtypes, as well as all SARS-CoV-2 sequences known to date, by in silico alignment, notably allowing the exclusion of other benign coronaviruses and SARS-CoV-1 or MERS. The determination of the SARS-CoV-2 status is based on the detection of 2 targets coupled to the same fluorophore: one in the IP2 gene (RdRp2) and the other in the IP4 gene (RdRp).

The test is indicated to make a presumptive diagnosis of infection by these three viruses or to complete a proven or undetermined diagnosis, by specifically distinguishing the SARS-CoV-2 coronavirus from the two influenza viruses. The diagnosis should always be made by medical personnel and in the clinical, historical and symptomatic context of the patient.

The RNA extract of any respiratory specimen is the starting material for the EurobioPlex EBX-043 kit. It is the responsibility of the user to use extraction methods appropriate for the specimens tested.

The EurobioPlex EBX-043 test is an in vitro diagnostic medical device and must be used by qualified medical laboratory personnel.

The kit is validated on Nasopharyngeal swab.

3. Symbols

	Reference
	Batch number
	Limits of storage temperature
	Expiration date
	Content sufficient for « N » reactions
	Manufacturer
	CE marked product
	In Vitro Diagnostic Medical Device
	Instructions for use
	Store away from light
	Do not use if packaging is damaged
	Caution

EurobioPlex Flu A/Flu B/SARS-CoV-2

4. Principle of detection

The EurobioPlex Flu A/Flu B/SARS-CoV-2 is a test for the amplification of ribonucleic acid (RNA) of viruses responsible for influenza A, influenza B and SARS-CoV-2 based on 1 multiplex RT-PCR assay (1 target gene for influenza A virus, 1 target gene for influenza B virus and two target genes for SARS-CoV-2: genes RdRp and RdRp2) in a single well.

The kit contains 1 oligomix to detect the 3 targets, as well as an endogenous control.

The presence of a human housekeeping gene serving as a control for the quality of the sample, for RNA extraction and RT-PCR inhibition allows to determine possible variations that may occur during the steps of RNA extraction from biologic samples and amplification by real-time RT-PCR. It ensures that a negative result cannot be due to poor RNA extraction and/or the presence of too many RT-PCR inhibitors.

The test is performed on RNA extracted from samples using a 1-well/tube reaction. Influenza A virus RNA is detected using a FAM-labeled probe. Influenza B virus RNA is detected using a HEX-labeled probe. Finally, the RNA of SARS-CoV-2 is detected using specific probes for each gene (RdRp and RdRp2) labeled in TEXAS RED. Endogenous control is detected using a CY5-labeled probe. During the elongation of the amplification product, the probes emit a specific fluorescence following their hydrolysis. The measurement of fluorescence intensity in real time is relative to the accumulation of amplification products.

Table 1: Detection of targets by fluorophores

Targets	Fluorophore	Excitation	Emission
Target 1 / Flu A	FAM	495 nm	515 nm
Target 2 / Flu B	HEX	535 nm	555 nm
Target 3 / RdRP1/RdRP2 genes	Texas Red	585 nm	605 nm
Endogenous internal control	Cy5	650 nm	670 nm

Equivalent channels on different real time PCR cyclers:

- Channel **FAM** (Systems ABI, SmartCycler II, Systems Mx, CFX96™/Chromo4, T-COR8 IVD, LC480), Channel Green (RotorGene), Channel Green (MIC).

EurobioPlex Flu A/Flu B/SARS-CoV-2

- Channel **HEX** (Chromo 4/CFX96, Systèmes Mx, T-COR 8-IVD, LC480), Canal VIC (Systèmes ABI), Canal Alexa532 (SmartCycler II), Canal Yellow (RotorGene), Yellow Channel (MIC).
- Channel **Texas Red** (Systems ABI, Chromo 4/CFX96, Systèmes Mx, T-COR 8-IVD), LC Red 610 (LC480), Canal Orange (RotorGene), Orange Channel (MIC).
- Channel **Cy5** (CFX96™/Chromo4, Systems ABI, Systems Mx, T-COR8-IVD, LC 480), Channel Alexa647 (SmartCycler II), Channel Red (RotorGene), Red Channel (MIC).

Note 1: On LC480, apply color compensation for the following channels: FAM/HEX-VIC/Texas Red/Cy5 (FAM : 483-533 / HEX-VIC : 523-568 /Texas Red : 558-610 / Cy5 : 615-670).

Note 2: During the analysis using LC480 Software release 1.5.0, do not leave Noiseband Auto on Texas red and HEX channels, but apply the NoiseBan Fluor manually at 0,43 for Texas Red channel and 5,5 for HEX channel (the operator need to adjust the baseline).

5. Content of the kit

The real-time RT-PCR kit EBX-043 is ready to use and contains the reagents and enzymes necessary for the detection of these viruses (Table 2).

Fluorescence is emitted and measured individually by an optical system during the PCR process. Detection of the amplified fragments is performed by a fluorimeter using the channels listed in Table 1.

Table 2: Content of the kit

Cap color	Content of the kit	25 reactions	50 reactions	100 reactions	200 reactions	600 reactions	Reconstitution
Red	Enzymes	5x25 µL	10x25 µL	525 µl	800 µl	2 x 1200 µl	Ready to use
Transparent	Oligomix	5x20 µL	10x20 µL	420 µl	700 µL	2 x 950 µL	Ready to use
Yellow	Positive Control CP	5x15 µL	10x15 µL	160 µl	320 µl	320 µl	Ready to use
Blue	Water = negative control (CN-H ₂ O)	1500 µL	1500 µL	1500µl	1500 µl	2 x 1500 µl	Ready to use

Oligomix: contains primers and probes for the 3 targets and for the endogenous control.

6. Storage

All reagents should be stored between -15°C and -22°C.

All reagents can be used until the expiration date indicated on the kit label.

EurobioPlex Flu A/Flu B/SARS-CoV-2

The sensitivity of the assay may decrease if the kit components undergo multiple freeze/thaw cycles. The kit can be used after the initial opening for up to 5 freeze/thaw cycles.

7. Materials required not provided

- ◇ Biological cabinet
- ◇ Real time PCR instrument
- ◇ Centrifuge for microtubes
- ◇ Vortex
- ◇ Plates / tubes for real-time PCR reaction
- ◇ Micropipettes
- ◇ DNase-free and RNase-free filter tips for micropipettes
- ◇ Sterile microtubes
- ◇ Gloves (talc-free)

8. Real-time PCR instrument

The EurobioPlex EBX-043 kit has been developed and validated for use with the following real-time PCR instruments:

- CFX96™ Real Time PCR detection system (Bio-Rad) with analysis on CFX Manager version 3.1 (Bio-Rad)
- LightCycler® 480 with analysis on LightCycler® 480 software v1.5 (Roche)
- T-COR 8®-IVD (Tetradcore Inc.) with T-COR 8 Per-Well SmartCTM (auto v1) software
- MIC® (Bio Molecular Systems) with analysis on MIC software version 2.12.6

9. Cautions and note

Read these instructions carefully before starting the procedure.

- ◇ This experimentation should be performed by medical laboratory technicians.
- ◇ The local and national biosafety regulations in place for the detection of SARS-CoV-2 must be followed strictly at all times, especially in laboratories and with laboratory equipment in agreement.
- ◇ Ensure that instruments have been installed, calibrated and maintained in accordance with the manufacturer's recommendations.
- ◇ It is the responsibility of the user to use equipment other than that validated, and in such cases performance is not guaranteed.

EurobioPlex Flu A/Flu B/SARS-CoV-2

- ◇ Clinical samples should be considered as potentially infectious material and should be prepared in a laminar flow cabinet.
- ◇ This experiment should be performed according to good laboratory practice.
- ◇ Do not use this kit after the stated expiration date.
- ◇ The kit is shipped under dry ice, and the kit components should arrive frozen. If one or more of the components arrive thawed, or if the tubes have been damaged in transit, contact Eurobio Scientific.
- ◇ Avoid freeze/thaw cycles of reagents as this may lead to a decrease in assay sensitivity.
- ◇ Once the reagents have been thawed, centrifuge the tubes briefly before use.
- ◇ The use of ice or an ice pack is recommended in case of long delays due to, for example, a large number of samples to be processed or high temperatures.
- ◇ It is recommended to define three separate work areas: 1) RNA isolation, 2) Preparation of the reaction mixture, and 3) Amplification/Detection of the amplified products.
- ◇ The fluorescent probes in the oligomix are light sensitive. Any prolonged exposure of the oligomix to light should be limited to the technical time required to prepare the PCR plate.
- ◇ It is recommended that the positive control be opened and handled separately from the biological samples to be tested and the other kit components to avoid cross-contamination.
- ◇ Distinct lab coat and gloves (talc-free) should be worn in each work area.
- ◇ Pipettes, reagents and other working materials should not be moved between these areas.
- ◇ Special care should be taken to maintain the purity of reagents and reaction mixtures.
- ◇ The internal endogenous control detects a detectable cell target in all samples of human origin but not in the CN-H2O negative control provided in the kit. The absence of a signal in the negative control prevents cross-contamination from occurring.
- ◇ Appropriate RNA preparation/extraction methods for quality RNA production and RT-PCR application should be used, especially to avoid contamination with RNAs.
- ◇ Use filter tips for micropipettes, RNase-free and DNase-free.
- ◇ Do not pipette with the mouth. Do not eat, drink or smoke in the laboratory.
- ◇ Avoid aerosols.

10.Procedure

10.1 Samples collection

- ◇ Collect samples in sterile tubes.

EurobioPlex Flu A/Flu B/SARS-CoV-2

- ◇ It is the responsibility of the user to control their own sample collection, transport, storage and extraction conditions to ensure that RNA extraction using appropriate systems produces RNA of high quality.
- ◇ It is recommended that samples be extracted immediately, or stored according to sample storage recommendations prior to extraction (Table 3).
- ◇ The bibliographic references in section "14. Bibliography" provide indicative data on the stability of samples and RNAs.

Table 3: Storage recommendations before use

Recommendations for maximum storage of samples before extraction	
Ambient temperature	2 hrs
+2°C/+8°C	5 days
<-70°C (preferred to -20°C)	Long term storage (>5 days and maximum 2 months at -20°C)

Caution	
	◇ The user may refer to the recommendations issued by the World Health Organization or the French National Authority for Health for the proper storage of samples : https://www.has-sante.fr/upload/docs/application/pdf/2010-02/cryopreservation_-_recommandations.pdf ◇ Extracted RNAs should be stored at <-70°C. After one year of storage of RNAs at <-70°C, the Ct obtained may increase. It is advisable to re-extract a biology sample that has been stored for more than one year. It is advisable to limit the number of times RNAs are thawed to 3. ◇ Transportation of clinical samples is subject to local regulations for the transportation of infectious agents.

10.2 RNA extraction

It is the responsibility of the user to ensure that the nucleic acid extraction system used is compatible with real-time RT-PCR technology. We recommend using virus RNA extraction methods suitable for respiratory specimens, and to refer to the instructions of the supplier of the extraction kit used.

In the EBX-043 kit, the endogenous internal control on the Cy5 channel is already present in the clinical sample. Properly detected after amplification, the endogenous internal control is

EurobioPlex Flu A/Flu B/SARS-CoV-2

used to ensure the quality of the sample and the extraction. The performance described in this document was obtained with the following extraction technique: Maelstrom 9600 (TANBead), Microlab Starlet IVD Hamilton, the OptiPure Viral Auto Plate kit (ref: W665A46) and the Starmag 96X4 kit.

10.3 Real-time RT-PCR

General note:

- The positive control contains high concentrations of matrix. Handling must be done carefully to avoid contamination.
- To confirm that the RT-PCR is working correctly, it is necessary to test the positive control, as well as the negative control (water supplied = CN-H2O) (see II-2/6 of the real-time RT-PCR protocol).

On T-COR 8®-IVD, these controls must be performed at least the first time a new kit is used.

EurobioPlex Flu A/Flu B/SARS-CoV-2

Diagram of the procedure:

1 - PREPARATION OF MASTERMIX

Number of reactions	N+3*
Enzymes	(N+3) x 3.75 µl
Water	(N+3) x 3.25 µl
Oligomix	(N+3) x 3 µl
Total Volume Mastermix	(N+3) x 10 µl

* N+1 on T-COR 8[®]-IVD

2 - PREPARATION OF REACTIONS

Sample

10 µl Mastermix
+
5µl RNA sample

Positive Control

10 µl Mastermix
+
5µl CP

Negative Control

10 µl Mastermix
+
5µl Molecular biology water (CN-H₂O)

3 – REAL TIME PCR INSTRUMENT

Program	Temperature	Duration	Cycle(s)	
Reverse Transcription	45°C	5 min	1	-
Denaturation	98°C	20 sec	1	-
Amplification	98°C	3 sec	40	-

EurobioPlex Flu A/Flu B/SARS-CoV-2

	58°C	10 sec		Acquisition of fluorescence
--	------	--------	--	-----------------------------

10.4 Detailed procedure

- 1) Ensure that the reagents are completely thawed. Homogenise the enzyme tubes, and vortex the Oligomix and CP for about 15 seconds, then centrifuge briefly.
- 2) Prepare the Mastermix as below. N is the number of reactions (including positive and negative controls. Plan to prepare enough reagents for at least N+3 reactions (*N+1 on T-COR8).

Number of reactions	N+3*
Enzymes	(N+3) x 3.75 µl
Water	(N+3) x 3.25 µl
Oligomix	(N+3) x 3 µl
Total Volume Mastermix	(N+3) x 10 µl

- 3) Homogenise the Mastermix prepared in 2) and centrifuge briefly
- 4) Dispense 10 µL of Mastermix into separate microplate tubes/wells using a micropipette and filter tips.
- 5) Add 5 µL of extracted RNA sample.
- 6) In parallel perform the following controls:
 - Positive Control:
 - 10 µL of Mastermix + 5 µL of CP.
 - Negative Control:
 - 10 µL of Mastermix + 5 µL supplied water (CN-H2O).
- 7) Immediately close the tubes, or plate with an adhesive film to avoid contamination.
- 8) Centrifuge briefly to collect the reaction mixture at the bottom of the tubes or microplate wells.
- 9) Run the following program on the real-time PCR instrument.

Program	Temperature	Duration	Cycle(s)	
Reverse Transcription	45°C	5 min	1	-
Denaturation	98°C	20 sec	1	-
Amplification	98°C	3 sec	40	-
	58°C	10 sec		Acquisition of fluorescence

Note1: On CFX96™ (Bio-Rad), start the run using the v1.6 or later version of CFX Manager Software, and analyze with v 3.1 (see § Validation of the Experiment). **Please ensure that white PCR plates are used for appropriate reading in all channels.**

EurobioPlex Flu A/Flu B/SARS-CoV-2


Note 2: On LightCycler® 480 (Roche), apply color compensation for the following channels: FAM/HEX-VIC/Texas Red/Cy5 (FAM : 483-533 / HEX-VIC : 523-568 /Texas Red : 558-610 / Cy5 : 615-670).

Note 3: When analyzing on LC480 Software release 1.5.0, do not leave the Noiseband (Auto) for Texas Red and HEX channels but apply the Fluor Noiseband manually at 0.43 for Texas Red and 5.5 for HEX.

11 Validation of the experiment


Post-acquisition data analysis on a CFX96 PCR instrument (Bio-Rad) must be performed using version 3.1 of the CFX Manager software (Bio-Rad). In order to change to this version from a run started on an earlier version, please follow the procedure below: At the end of the run, the data file with the extension .pcrd must be opened and processed with CFX Manager (Bio-Rad) version 3.1.

If the run was started with CFX Manager v1.6 for example, to open a data file with CFX Manager v3.1, click on the CFX Manager v3.1 icon. The home screen appears.


Click on "File" and select "Open" then "Data File".

EurobioPlex Flu A/Flu B/SARS-CoV-2

Select the file you need to analyse and click on "Open".

The "Drift correction" option must be applied in the "Settings" tab as shown in the image below: click on the "Settings" tab, then on "Baseline Setting" and on "Apply Fluorescence Drift Correction".

Once this is done, the analysis can start.

The results for the controls must be the following (Table 4); otherwise, the experiment is not valid. On T-COR 8®-IVD, these controls must be performed at least during the first use of a new box of kit.

Table 4: Validation of the run

Positive Control	
FAM	Ct ≤ 30
HEX	Ct ≤ 30
Texas Red	Ct ≤ 30
Cy5	Ct ≤ 30
Negative Control	
FAM	Undetermined Ct
HEX	Undetermined Ct
Texas Red	Undetermined Ct

EurobioPlex Flu A/Flu B/SARS-CoV-2

Cy5	Undetermined Ct
-----	-----------------

12 Data analysis and interpretation

For clinical samples, the following results are possible:

 *SARS-CoV-2 (Texas Red), Flu A (FAM) or Flu B (HEX): Ct < 40

*Ct cut off for positivity of samples for T-COR 8®-IVD is: + Positive => Ct positive (≤ 45) for the 3 channels, but automatic interpretation is provided with Barcodes.

Table 5:

PCR Signal Flu A and/or Flu B and/or SARS-CoV-2	Endogenous control	Presence of virus	Test validity/comment
FAM and/or HEX and/or Texas Red	Cy5	Flu A and/or Flu B and/or SARS-CoV-2	
+	+	Yes	VALID
-	+	No	VALID
+	-	Yes	VALID Possible RT-PCR inhibition or problem with extraction that does not prevent the detection of viruses.
-	-	No possible interpretation	RT-PCR inhibition or problem with extraction → dilute first 5 x the sample and test again; if necessary redo an extraction

Limitations on use and interpretation:

- ❖ All samples should be treated as potentially infected with SARS-CoV-2, and local biosafety regulations should be carefully followed.
- ❖ Interpretation of results should consider the possibility of false negatives and false positives.

False negatives may be due to:

- Inadequate collection of samples, or incorrect storage,
- Samples outside the viremia phase,
- Inappropriate extraction conditions or use of non-validated PCR instruments,
- Experimentation not respecting all elements of these instructions for use.

False positives may be due to:

EurobioPlex Flu A/Flu B/SARS-CoV-2

- Contamination due to mishandling of high positive samples, the positive control, or the PCR amplification products,
 - Failure to follow the procedure described in these instructions for use, especially to avoid sources of contamination.
- ❖ All results should be interpreted by medical professionals in the context of the patient's clinical situation, history and symptoms.

13. Performance analysis

A dilution range was performed using a plasmid mixture ranging from 10^5 to 1 copy/ μ L. This dilution range was used to determine the detection limit (at 100% detection) of the EBX-043 kit.

In addition, a detection limit was determined on the WHO 1st International Standard for SARS-CoV-2 RNA (NIBSC code 21/146) and on a SARS-CoV-2 synthetic RNA of Twist Biosciences.

Finally, the detection limit was verified on a plasmid that mimics the sequence of an H3N2 subtype that emerged in Denmark in March 2022 and whose mutations alter the efficiency of the detection system usually used on the M gene.

Jørgensen RL, Lerche CJ, Pedersen MS, Kirkby NS, Botnen AB, Trebbien R, Nilsson-Møller S, Pinholt M, Nielsen ACY, Westh H, Lisby JG, Schneider UV. Emergence of circulating influenza A H3N2 viruses with genetic drift in the matrix gene: be alert of false-negative test results. APMIS. 2022 Oct;130(10):612-617. doi: 10.1111/apm.13262. Epub 2022 Aug 10. PMID: 35836366; PMCID: PMC9544743.

Our primer and probe system used for influenza A detection targets an amplicon that is close to but offset from the design described in the Jørgensen paper.

Detection limitation/analytical sensitivity on CP:

SARS-CoV-2 : 10 copies/ μ L

Flu A : 5 copies/ μ L

Flu B : 50 copies/ μ L

CI Cy5 : 50 copies/ μ L

Detection limitation/analytical sensitivity on WHO 1st International Standard for SARS-CoV-2 RNA (cut-off 95%):

SARS-CoV-2: 5.31 copies/ μ L

EurobioPlex Flu A/Flu B/SARS-CoV-2

Detection limitation/analytical sensitivity on a SARS-CoV-2 synthetic RNA (Twist Biosciences, USA):

SARS-CoV-2: 12.5 copies/ μ L

Detection limitation/analytical sensitivity on plasmid mimicking the Danish variant (influenza AH3N2):

SARS-CoV-2 : 10 copies/ μ L

Flu A : 10 copies/ μ L

Flu B : 50 copies/ μ L

CI Cy5 : 50 copies/ μ L

Sensitivity study on 107 known positive samples diluted in known samples Negatives at 3x LOD: 100% sensitivity (detection of 107 samples / 107)

- **Interfering substances:** 7 substances from the AcroMetrix™ Inhibition Panel REF 956400 kit, Oxymetazoline and Zanamivir. A high concentration, which may be contained in some drugs, was added to the positive sample before extraction: No impact of the tested interfering substances was observed on the result obtained.

Signal variability on FAM, HEX, Texas Red and Cy5 channels

- Intra-experimental variability:

Average intra-batch CV (%)	Flu A	Flu B	SARS-CoV-2	Endogenous ctrl Gene
CFX96™	0,26	0,67	0,52	1,11
T-COR-8® IVD	1,67	1,40	1,18	1,58
LC480®	0,25	0,56	0,50	0,25
MIC®	3,20	0,39	1,26	2,12

CV: coefficient of variation

- Inter-batch variability:

Average inter-batch CV (%)	Flu A	Flu B	SARS-CoV-2	Endogenous ctrl Gene
CFX96™	0,68	0,81	1,14	3,71
T-COR-8® IVD	3,44	3,17	1,09	2,80
LC480®	1,42	1,06	1,03	0,94

EurobioPlex Flu A/Flu B/SARS-CoV-2

MIC®	2,66	4,77	5,78	6,98
------	------	------	------	------

CV: coefficient of variation

Diagnostic specificity

This validation included:

- 505 SARS-CoV-2 negative samples, pretested and characterized by CE approved tests.
- 32 samples negative for all known coronaviruses, including respiratory panel samples characterized for other viruses to test for cross-reactions.
- 1 sample positive for CoV NL63.

Number of samples	Positive sample for	Coronavirus	SARS-CoV-2 status with EBX-043
3	Adenovirus	Negative	Negative
2	Parainfluenzae	Negative	Negative
7	Rhinovirus	Negative	Negative
1	Legionella pneumoniae	Negative	Negative
1	Mycoplasma pneumoniae	Negative	Negative
3	Mycoplasma tuberculosis	Negative	Negative
1	ParaInfluenza Type 4	Negative	Negative
3	Enterovirus	Negative	Negative
3	Streptococcus pneumoniae	Negative	Negative
1	Staph.epiderm/oralis/mitis/sanguis	Negative	Negative
2	Candida albicans	Negative	Negative
2	Pseudomonas aeruginosa	Negative	Negative
1	CoV NL63	Negative	Negative
3	VRS	Negative	Negative

In silico analysis was used and demonstrated the specificity of the primers and probes for SARS-CoV-2. For example, on the sequences of MERS and SARS viruses, no significant homology of SARS-CoV-2 primers and probes that could amplify these viruses was found.

No aspecificity or cross reaction was detected.

The kit is 100% specific.

EurobioPlex Flu A/Flu B/SARS-CoV-2

Diagnostic sensibility and specificity

The validation of the performance was performed on:

538 negative samples	
505 SARS-CoV-2 negative samples	Characterised with EBX-041 kit, CE marked
33 positive respiratory samples	RSV-A/B, Parainfluenzae, Rhinovirus, Metapneumovirus, CoV NL63

305 positive samples	
140 SARS-CoV-2 positive samples	Characterized by EurobioPlex EBX-041 CE-IVD kit
59 Influenza B positive samples	Characterized with the Seegene Allplex SARS-Cov-2 Flu A / Flu B / RSV Assay
106 Influenza A positive samples	Characterized with the Seegene Allplex SARS-Cov-2 Flu A / Flu B / RSV Assay

For this study, sample extraction was performed on Maelstrom 9600 (TANBead) with the OptiPure Viral Auto Plate kit (ref: W665A46), and RT-PCR on CFX96 (Bio-Rad).

Overall performances are:

		EBX-043		
		Flu A		
		POSITIVE	NEGATIVE	TOTAL
Pretested	POSITIVE	106	0	106
	NEGATIVE	2	707	709
		Flu B		
		POSITIVE	NEGATIVE	TOTAL
Pretested	POSITIVE	59	0	59
	NEGATIVE	3	753	756
		SARS-CoV-2		
		POSITIVE	NEGATIVE	TOTAL
Pretested	POSITIVE	140	0	140
	NEGATIVE	0	505	505

EurobioPlex Flu A/Flu B/SARS-CoV-2

	Sensitivity %	Specificity %	Concordance %
FLU A	> 98%	> 99%	> 99%
FLU B	95.16%	> 99%	> 99%
SARS-CoV-2	> 99%	> 99%	> 99%

SPECIFICITIES OF THE REAL-TIME PCR T-COR 8®-IVD INSTRUMENT

T-COR 8®-IVD is a real-time PCR instrument with 8 independently programmable wells that can be used, patient by patient, independently of each other, in terms of the assays performed, the thermal programs, and the timing of the PCR run.

All 8 wells can also be used simultaneously with the same assay.

Note: Mix well after placing the sample in the tube by pipetting back and forth, avoiding the formation of bubbles and making sure that the volume of liquid is well located at the bottom of the tube. Close each tube tightly after each deposit of control or sample to avoid contamination.

Controls

On T-COR 8®-IVD, the positive and negative controls should be tested at least the first time a new kit is used.

Thereafter, the endogenous control allows to ensure that the extraction has been performed correctly, and that within the tested sample, there is no inhibition of the PCR, and that the PCR is working correctly.

Fluorophores

There are 4 fluorophore combinations available on the T-COR 8®-IVD.

For all EBX, including EBX-043, the FAM/HEX/Texas Red/Cy5 combination is the one chosen. If some channels are not used in some kits, do not take this into account when analyzing the results.

Use of the barcodes

1- Select Menu > New Analysis/New Run

EurobioPlex Flu A/Flu B/SARS-CoV-2

- 2- Select Barcode,
- 3- Scan the desired barcode on the right side of the instrument:
 - either for the positive control (Barcode EBX-043 Pos Ctrl),
 - or for the negative control (Barcode EBX-043 Neg Ctrl),
 - or for a sample (Barcode EBX-043)

The instrument automatically selects one of the 8 available wells. Follow the instructions given by the instrument.

- 4- Lift the flap of the well x selected by the instrument, and check that the blue LED light is on, then select "Yes".
- 5- Place the tube in the corresponding well and fold the valve back, then select " Next ".
- 6- If you wish to name it (optional), select " Sample x ", and name it. Select " Accept ".
- 7- Select " Next ".
- 8- To add/test another control or sample, select " Add well " and start again at point 3-.
- 9- Select " Start Run ".

Note: the nomination of a sample or run cannot be modified or added after the end of the run. It must be done before or during the run.

Presentation and Interpretation of Results

Ct results and amplification curves are available from the SmartCTTM Values / SmartCTTM Table (Ct values assigned on each channel), and Ct vs. run graphs on each channel, both of which can be viewed in real time on the machine.

An analysis for the negative and positive controls, as well as for the samples, is generated automatically. This is available in "Interpretations", at the end of the run, in the "View" window.

For the positive and negative control, the following results are possible:

"Neg Ctrl Fail": Not valid
"Neg Ctrl Valid": Valid
"Pos Ctrl Fail": Not valid
"Pos Ctrl Valid": Valid

Determination of status for samples:



"Détecté(e-s) / Detected" : Positive → green box
"Not detected / Non détecté" : Negative → red box
"Non Valide / Invalid" : Invalid result -> retest -> yellow box

Example of a presentation of automatic interpretation of results


EurobioPlex Flu A/Flu B/SARS-CoV-2

Summary									
	Well	Sample	Assay	FAM HEX TxR Cy5			Call	Note	
	1	Sample 1	EBX-043	12.1		17.3	Detected • FLUA	FLUB, SARS-CoV-2 Negative	
	2	Sample 2	EBX-043		18.8	19.0	Detected • FLUB	FLUA, SARS-CoV-2 Negative	
	3	Sample 3	EBX-043			21.2 18.4	Detected • SARS-CoV-2	FLUA, FLUB Negative	
	4	Ctlr POS	EBX-043 Pos Ctrl	24.1	22.2	26.7 24.4	Detected • Pos Ctrl	Pos Ctrl Valid	
	5	Ctrl POS	EBX-043 Pos Ctrl	20.1			Invalid	Pos Ctrl Fail	
	6	Ctrl NEG	EBX-043 Neg Ctrl				Detected • Neg Ctrl	Neg Ctrl Valid	
	7	Sample 4	EBX-043				Invalid	Retest	
	8	Sample 5	EBX-043			18.1	Not Detected	FLUA, FLUB, SARS-CoV-2 Negative	

Example of amplification curves


EurobioPlex Flu A/Flu B/SARS-CoV-2

EurobioPlex Flu A/Flu B/SARS-CoV-2

Barcodes on T-COR8®-IVD for EBX-043


EBX-043 Pos Ctrl
Positive Control

EBX-043 Neg Ctrl
Water = Negative Control (CN-H2O)

EBX-043

14. Bibliography

Xiaojing Wu, Ying Cai, Xu Huang, Xin Yu, Li Zhao, Fan Wang, Quanguo Li, Sichao Gu, Teng Xu, Yongjun Li, Binghuai Lu, Qingyuan Zhan, Co-infection with SARS-CoV-2 and Influenza A Virus in Patient with Pneumonia, China, Emerg Infect Dis. 2020 Jun;26(6):1324-1326.

Yan Li , Jiangshan Wang , Chunting Wang , Qiwen Yang , Yingchun Xu , Jun Xu , Yi Li Xuezhong Yu , Huadong Zhu , Jihai Liu Characteristics of respiratory virus infection during the outbreak of 2019 novel coronavirus in Beijing. Int J Infect Dis. 2020 Jul;96:266-269.

Chan JF, Kok KH, Zhu Z, Chu H *et al.* Genomic characterization of the 2019 novel human-pathogenic coronavirus isolated from a patient with atypical pneumonia after visiting Wuhan. Emerg microbes infect. 2020 Dec; 9(1):221-236.

Cohen J, Normile D. New SRAS-like virus in China triggers alarm, Science. 2020 Jan 17;367(6475):234-235.

Drosten C, Muth D, Corman VM, Hussain R *et al.* An observational, laboratory-based study of outbreaks of middle east respiratory syndrome coronavirus in Jeddah and Riyadh, kingdom of Saudi Arabia, 2014. Clin infect dis. 2015 Feb 1; 60(3):369-77.

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EurobioPlex Flu A/Flu B/SARS-CoV-2

15. Quality control

In accordance with Eurobio's ISO EN 13485 certified quality management system, each batch of EurobioPlex Flu A/Flu B/SARS-CoV-2 is tested according to predefined specifications to ensure consistent product quality.

16. Waste disposal

Eliminate all waste in accordance with the legislation on DASRI.

17. Incident report

Serious incidents related to the system must be reported to Eurobio Scientific.

18. Technical assistance

For assistance with our products, please contact our technical support.

Eurobio Scientific's customer service can be contacted by e-mail at adv@eurobio-scientific.com or by phone at +33 (0)1.69.79.64.80.

Change tracking

Date	Version	Revision
2023/02/17	v2.00	Additional clinical study on influenza A and B samples. Validation of EBX-043 on a new instrument: MIC®. Study on the new FLU A H3N2 Danish variant.
2022/05/03	v1.00	Creation

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





EBX-043

14. Literatur

Xiaojing Wu, Ying Cai, Xu Huang, Xin Yu, Li Zhao, Fan Wang, Quanguo Li, Sichao Gu, Teng Xu, Yongjun Li, Binghuai Lu, Qingyuan Zhan, Co-infection with SARS-CoV-2 and Influenza A Virus in Patient with Pneumonia, China, *Emerg Infect Dis.* 2020 Jun;26(6):1324-1326.

Yan Li, Jiangshan Wang, Chunting Wang, Qiwen Yang, Yingchun Xu, Jun Xu, Yi Li, Xuezhong Yu, Huadong Zhu, Jihai Liu Characteristics of respiratory virus infection during the outbreak of 2019 novel coronavirus in Beijing. *Int J Infect Dis.* 2020 Jul;96:266-269.

Chan JF, Kok KH, Zhu Z, Chu H *et al.* Genomic characterization of the 2019 novel human-pathogenic coronavirus isolated from a patient with atypical pneumonia after visiting Wuhan. *Emerg microbes infect.* 2020 Dec; 9(1):221-236.

Cohen J, Normile D. New SRAS-like virus in China triggers alarm, *Science.* 2020 Jan 17;367(6475):234-235.

Drosten C, Muth D, Corman VM, Hussain R *et al.* An observational, laboratory-based study of outbreaks of middle east respiratory syndrome coronavirus in Jeddah and Riyadh, kingdom of Saudi Arabia, 2014. *Clin infect dis.* 2015 Feb 1; 60(3):369-77.

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