EurobioPlex SARS-CoV-2 Fast-SVT Screening and VOC typing (Omicron BA.1, Omicron BA.2, Omicron BA.4/BA.5, Delta)

RT-PCR EN TEMPS REEL

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14. Bibliographie

Cherian, Sarah, Varsha Potdar, Santosh Jadhav, Pragya Yadav, Nivedita Gupta, Mousmi Das, Partha Rakshit, et al. 2021. « Convergent Evolution of SARS-CoV-2 Spike Mutations, L452R, E484Q and P681R, in the Second Wave of COVID-19 in Maharashtra, India ». Preprint. *Molecular Biology*. <https://doi.org/10.1101/2021.04.22.440932>.

Hu, Ben, Hua Guo, Peng Zhou, et Zheng-Li Shi. 2021. « Characteristics of SARS-CoV-2 and COVID-19 ». *Nature Reviews Microbiology* 19 (3): 141-54. <https://doi.org/10.1038/s41579-020-00459-7>.

Lauring, Adam S., et Emma B. Hodcroft. 2021. « Genetic Variants of SARS-CoV-2—What Do They Mean? » *JAMA* 325 (6): 529. <https://doi.org/10.1001/jama.2020.27124>.

Plante, Jessica A., Brooke M. Mitchell, Kenneth S. Plante, Kari Debbink, Scott C. Weaver, et Vineet D. Menachery. 2021. « The Variant Gambit: COVID-19's next Move ». *Cell Host & Microbe* 29 (4): 508-15. <https://doi.org/10.1016/j.chom.2021.02.020>.

Seyed Hosseini, Elahe, Narjes Riahi Kashani, Hossein Nikzad, Javid Azadbakht, Hassan Hassani Bafrani, et Hamed Haddad Kashani. 2020. « The Novel Coronavirus Disease-2019 (COVID-19): Mechanism of Action, Detection and Recent Therapeutic Strategies ». *Virology* 551 (décembre): 1-9. <https://doi.org/10.1016/j.virol.2020.08.011>.

Srivastava, Surabhi, Sofia Banu, Priya Singh, Divya Tej Sowpati, et Rakesh K. Mishra. 2021. « SARS-CoV-2 Genomics: An Indian Perspective on Sequencing Viral Variants ». *Journal of Biosciences* 46 (1): 22. <https://doi.org/10.1007/s12038-021-00145-7>.

European Centre for Disease Prevention and Control : <https://www.ecdc.europa.eu/en/covid->

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
Reference EBX-048 Version 3.00 – 13/06/2022

Validated on:

- CFX96™ Real Time PCR detection system (Bio-Rad) with analysis on CFX Manager version 3.1 (Bio-Rad)

Instructions for use

Available on www.eurobio-scientific.com

Fiche technique EurobioPlex SARS-CoV-2 Fast-SVT

Table of contents

1. Introduction.....	29
2. Purpose of the system.....	31
3. Symbols	31
4. Principle of detection	32
5. Content of the kit	33
6. Storage.....	34
7. Materials required not provided.....	34
8. Real-time PCR instrument	34
9. Cautions and note	34
10. Procedure	36
11. Validation of the experiment	40
12. Data analysis and interpretation.....	43
13. Performance analysis	45
14. Bibliography.....	53
15. Quality control.....	53
16. Waste disposal	54
17. Incident report	52
18. Technical assistance	52

Fiche technique EurobioPlex SARS-CoV-2 Fast-SVT

1. Introduction

The SARS-CoV-2 virus emerged at the end of 2019 in the city of Wuhan in China. It is classified as a member of the Coronaviridae family and the Sarbecovirus subgroup. Its genome consists of 29903 ribonucleic acid (RNA) bases. SARS-CoV-2 is the 7th identified coronavirus infecting humans, after human coronaviruses (HCoV), SARS-CoV (coronavirus causing severe acute respiratory syndrome) and MERS-CoV (coronavirus inducing Middle East respiratory syndrome).

The sequences of SARS-CoV-2 show similarities to those of the betacoronaviruses found in bats. SARS-CoV-2 is genetically distinct from other human coronaviruses.

The clinical profile is variable, ranging from flu-like symptoms, fever, cough, difficulty breathing, to pneumonia, to a severe respiratory syndrome that can be fatal. The lethality rate of the virus is estimated to be between 2-3%. SARS-CoV-2 is highly contagious with an estimated 173 million confirmed cases worldwide since early June 2021.

Several variants and mutations of SARS-CoV-2 have demonstrated an impact on public health (increased transmissibility, severity of infection or immune escape). They emerged at the end of 2020 in distinct geographical areas and are designated as variants of concern (VOC):

- Alpha variant from the United Kingdom (UK): lineage B.1.1.7 with mutation N501Y.
- Beta variant from South Africa (SA): lineage B.1.351 with mutation E484K.
- Gamma variant from Brazil (BR): lineage P.1 with mutation E484K.
- Delta variant from India: lineage B.1.617.2 with mutation L452R.
- Kappa variant from India: lineage B.1.617.1 with the mutations L452R and E484Q.
- The L452R mutation is also found in other lineages, including A.27, the origin of which is unclear, and C.36, detected in Egypt. An L452Q variation of lineage C.37 has also been detected in Peru.
- A variant called Omicron appeared in South Africa in November 2021 carrying up to 61 mutations in several regions of the viral genome including the K417N mutation.
 - o Omicron now features multiple BA.1, BA.1.1, BA.2, BA.3 sublineages.
 - o BA.2 appeared in Israel in January 2022. We are interested in BA.2 because it supplanted BA.1 in Denmark, it is progressing in the United Kingdom and now in France. It presents in addition to the K417N mutation, the V213G mutation which is absent in the BA.1 sub-lineage.
 - o Sublineages BA.4 and BA.5 began to be identified in January 2022 in South Africa, before spreading to several countries around the world (Portugal, France, ...)

The interest of the differentiation of the variants of Covid-19 is therefore to focus on the Delta variant, as well as the Omicron variant and its BA.x and BA.2 sub-lines on which the currently available vaccines could not have little or no effect, since these are mutations responsible for immune escape. The new variant identification strategy aims to identify immune escape

Fiche technique EurobioPlex SARS-CoV-2 Fast-SVT

mutations as a priority, rather than distinguishing variants. This strategy will allow faster management of patients with variants with this type of mutation.

2. Purpose of the system

The EurobioPlex SARS-CoV-2 Fast-SVT test is a real-time polymerase chain reaction (RT-PCR) assay which, in a single test, identifies the prevalence of SARS-CoV-2 and distinguishes between the K417N (Omicron variant BA.x and BA.2 screening), V213G (BA.2 sublineage of Omicron variant screening) and L452R (Delta variant screening) mutations in a qualitative manner. The presence of the V213G makes it possible to differentiate the BA.2 sub-lineage of the Omicron variant. The presence of this 3 mutations indicate the BA.4 or BA.5 sub-lineages.

The RNA extract is the starting material for the EurobioPlex SARS-CoV-2 Fast-SVT kit. It is up to the user to use extraction methods adapted to the tested samples. The kit has been validated on nasopharyngeal swabs.

The EurobioPlex SARS-CoV-2 Fast-SVT Kit is exclusively intended for use for the detection of SARS-CoV-2 and mutations of interest to differentiate the Delta variant and the BA.x , BA.2, BA.4/BA.5 sublineages of the Omicron variant.

The EurobioPlex SARS-CoV-2 Fast-SVT test must be used by qualified medical laboratory personnel. It is for single use only and should not be recycled after use.

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Manufacturer

CE marked product

In Vitro Diagnostic Medical Device

Instructions for use

Store away from light

Do not use if packaging is damaged

Caution

Fiche technique EurobioPlex SARS-CoV-2 Fast-SVT

4. Principle of detection

The EurobioPlex SARS-CoV-2 Fast-SVT kit is a qualitative ribonucleic acid (RNA) amplification test that allows the analysis of:

- The presence of two targets (RdRp1/RdRp2) of SARS-CoV-2 on the RdRp gene, in the same reading channel.
- The presence of the L452R (Delta variant screening), K417N (Omicron variant screening BA.1 and BA.2) and V213G (BA.2 sublineage of Omicron variant screening) mutations and distinguish them in two different detection channels. The presence of the 3 mutations indicate a BA.4/BA.5 sub-lineage.
- The presence of a human housekeeping gene used as a control for the quality of the sample, for RNA extraction and RT-PCR inhibition. The latter allows to determine possible variations that may occur during the RNA extraction and real-time RT-PCR amplification steps. It ensures that a negative result cannot be due to poor RNA extraction and/or the presence of too many RT-PCR inhibitors.

The RNA extract is the starting material for the EurobioPlex SARS-CoV-2 Fast-SVT kit. The RNA is extracted from a patient's nasopharyngeal swab.

SARS-CoV-2 RNA is detected using specific probes for each mutation or gene labeled in FAM, HEX, Texas Red and Quasar 705. The endogenous internal control is detected using a Cy5-labeled probe (Table 1). During elongation of the amplification product, the probes emit specific fluorescence following their hydrolysis. Measurement of the fluorescence intensity in real time, emitted and measured individually by an optical system during PCR, is relative to the accumulation of amplification products.

Table 1: Detection of targets by fluorophores

Targets	Fluorophore	Excitation	Emission
Target 1 / V213G mutation	FAM	495 nm	515 nm
Target 2 / K417N mutation	HEX	535 nm	555 nm
Target 3 / RdRP1/RdRP2 gene	Texas Red	585 nm	605 nm
Target 4 / L452R mutation	Quasar 705	681 nm	698 nm
Endogenous internal control	Cy5	650 nm	670 nm

Fiche technique EurobioPlex SARS-CoV-2 Fast-SVT

Recommended equivalent channels on different PCR instruments:

- **FAM channel** (CFX96, QuantStudio 5, QuantStudio 6), Green channel (Rotor-Gene Q),
- **HEX channel** (CFX96), VIC channel (QuantStudio 5, QuantStudio 6), Yellow channel (Rotor-Gene Q),
- **Texas Red channel** (CFX96, QuantStudio 6), Orange channel (Rotor-Gene Q), Rox channel (QuantStudio 5),
- **Cy5 channel** (CFX96, QuantStudio 5, QuantStudio 6), Red channel (Rotor-Gene Q),
- **Quasar 705 channel** (CFX96, QuantStudio 5, QuantStudio 6), Purple channel (Rotor-Gene Q).

5. Content of the kit

The EurobioPlex SARS-CoV-2 Fast-SVT real-time RT-PCR kit is ready to use and contains the reagents and enzymes necessary for the detection of the virus (Table 2).

Fluorescence is emitted and measured individually by an optical system during the PCR process. Detection of the amplified fragments is performed by a fluorometer using the channels listed in Table 1.

Table 2: Content of the kit

Cap color	Content of the kit	50 reactions	100 reactions	200 reactions	600 reactions	Reconstitution
Red	Enzymes	195 µL	450 µL	825 µL	3 x 825 µL	Ready to use
Transparent	Oligomix	156 µL	360 µL	660 µL	2 x 990 µL	Ready to use
Yellow	Positive Control (CP)	80 µL	160 µL	320 µL	320 µL	Ready to use
Blue	Water = Negative Control (CN-H2O)	1000 µL	1000 µL	2 x 1000 µL	6 x 1000 µL	Ready to use

Oligomix: contains primers and probes for the 4 targets and for the endogenous control

Fiche technique EurobioPlex SARS-CoV-2 Fast-SVT

6. Storage

All reagents should be stored between -15°C and -22°C.

All reagents can be used until the expiration date indicated on the kit label.

The sensitivity of the assay may decrease if the kit components undergo multiple freeze/thaw cycles. The kit can be used after the initial opening for up to 5 freeze/thaw cycles.

7. Materials required not provided

- ◇ Biological cabinet
- ◇ Real time PCR instrument
- ◇ Centrifuge for microtubes
- ◇ Vortex
- ◇ Plates / tubes for real-time PCR reaction
- ◇ Micropipettes
- ◇ DNase-free and RNase-free filter tips for micropipettes
- ◇ Sterile microtubes
- ◇ Gloves (talc-free)

8. Real-time PCR instrument

The EurobioPlex SARS-CoV-2 Fast-SVT kit has been developed and validated for use with the following real-time PCR instruments:

- CFX96™ Real Time PCR detection system (Bio-Rad) with analysis on CFX Manager version 3.1 (Bio-Rad)

9. Cautions and note

Read these instructions carefully before starting the procedure.

- ◇ This experimentation should be performed by medical laboratory technicians.
- ◇ The local and national biosafety regulations in place for the detection of SARS-CoV-2 must be followed strictly at all times, especially in laboratories and with laboratory equipment in agreement.

Fiche technique EurobioPlex SARS-CoV-2 Fast-SVT

- ◇ Ensure that instruments have been installed, calibrated and maintained in accordance with the manufacturer's recommendations.
- ◇ It is the responsibility of the user to use equipment other than that validated, and in such cases performance is not guaranteed.
- ◇ Clinical samples should be considered as potentially infectious material and should be prepared in a laminar flow cabinet.
- ◇ This experiment should be performed according to good laboratory practice.
- ◇ Do not use this kit after the stated expiration date.
- ◇ The kit is shipped under dry ice, and the kit components should arrive frozen. If one or more of the components arrive thawed, or if the tubes have been damaged in transit, contact Eurobio Scientific.
- ◇ Avoid freeze/thaw cycles of reagents as this may lead to a decrease in assay sensitivity.
- ◇ Once the reagents have been thawed, centrifuge the tubes briefly before use.
- ◇ The use of ice or an ice pack is recommended in case of long delays due to, for example, a large number of samples to be processed or high temperatures.
- ◇ It is recommended to define three separate work areas: 1) RNA isolation, 2) Preparation of the reaction mixture, and 3) Amplification/Detection of the amplified products.
- ◇ The fluorescent probes in the oligomix are light sensitive. Any prolonged exposure of the oligomix to light should be limited to the technical time required to prepare the PCR plate.
- ◇ It is recommended that the positive control be opened and handled separately from the biological samples to be tested and the other kit components to avoid cross-contamination.
- ◇ Distinct lab coat and gloves (talc-free) should be worn in each work area.
- ◇ Pipettes, reagents and other working materials should not be moved between these areas.
- ◇ Special care should be taken to maintain the purity of reagents and reaction mixtures.
- ◇ The internal endogenous control detects a detectable cell target in all samples of human origin but not in the CN-H2O negative control provided in the kit. The absence of a signal in the negative control prevents cross-contamination from occurring.
- ◇ Appropriate RNA preparation/extraction methods for quality RNA production and RT-PCR application should be used, especially to avoid contamination with RNAs.
- ◇ Use filter tips for micropipettes, RNase-free and DNase-free.
- ◇ Do not pipette with the mouth. Do not eat, drink or smoke in the laboratory.
- ◇ Avoid aérosols.

Fiche technique EurobioPlex SARS-CoV-2 Fast-SVT

10.Procedure

a. Samples collection

- ◇ Collect samples in sterile tubes.
- ◇ It is the responsibility of the user to control their own sample collection, transport, storage and extraction conditions to ensure that RNA extraction using appropriate systems produces RNA of high quality.
- ! It is recommended that samples are extracted immediately, or stored according to the sample storage recommendations before extraction (Table 3).

Table 3: Storage recommendations before use

Recommendations for maximum storage of samples before extraction	
Ambient temperature	2 h
4°C	72 h
-20°C (preferably -80°C)	Long-term storage

Caution	
	◇ The user may refer to the recommendations issued by the World Health Organization or the French National Authority for Health for the proper storage of samples. ◇ Extracted RNA should be stored at -80°C. ◇ The transport of clinical samples is subject to local regulations for the transport of infectious agents.

b. RNA extraction

It is the responsibility of the user to ensure that the nucleic acid extraction system used is compatible with real-time RT-PCR technology. We recommend using virus RNA extraction methods suitable for respiratory specimens, and to refer to the instructions of the supplier of the extraction kit used.

In the EBX-048 kit, the endogenous internal control on the Cy5 channel is already present in the clinical sample. Properly detected after amplification, the endogenous internal control is used to ensure the quality of the sample and the extraction. The performance described in this document was obtained with the following extraction technique: Maelstrom 9600 (TanBEAD) and the OptiPure Viral Auto Plate kit (ref: W665A46).

c. Real-time RT-PCR

General note:

- The positive control contains high concentrations of matrix. Handling must be done carefully to avoid contamination.
- To confirm that the RT-PCR is working correctly, it is necessary to test the positive control, as well as the negative control (water supplied = CN-H₂O) (see II-2/6 of the real-time RT-PCR protocol).

Fiche technique EurobioPlex SARS-CoV-2 Fast-SVT

Diagram of the procedure:

1 - PREPARATION OF MASTERMIX

Number of reactions	N+3
Enzymes	(N+3) x 3.75 µl
Water	(N+3) x 3.25 µl
Oligomix	(N+3) x 3 µl
Total Volume Mastermix	(N+3) x 10 µl

2 - PREPARATION OF REACTIONS

Sample

10 µl Mastermix
+
5µl RNA sample

Positive Control

10 µl Mastermix
+
5µl CP

Negative Control

10 µl Mastermix
+
5µl Molecular biology water (CN-H₂O)

3 – REAL TIME PCR INSTRUMENT

Program	Tempera-	Dura-	Cy-	
Reverse Tran-	45°C	5 min	1	-
Denaturation	98°C	20 sec	1	-
Amplification	98°C	3 sec	40	-
	58°C	10 sec		Acquisition of fluo-

Fiche technique EurobioPlex SARS-CoV-2 Fast-SVT

d. Detailed procedure

- 1) Ensure that the reagents are completely thawed. Homogenise the enzyme tubes, and vortex the Oligomix and CP for about 15 seconds, then centrifuge briefly.
- 2) Prepare the Mastermix as below. N is the number of reactions (including positive and negative controls. Plan to prepare enough reagents for at least N+3 reactions.

Number of reactions	N+3
Enzymes	(N+3) x 3.75 µl
Water	(N+3) x 3.25 µl
Oligomix	(N+3) x 3 µl
Total Volume Mastermix	(N+3) x 10 µl

**For smaller series (≤10): preparing for N+2 is sufficient.*

- 3) Homogenise the Mastermix prepared in 2) and centrifuge briefly
- 4) Dispense 10 µL of Mastermix into separate microplate tubes/wells using a micropipette and filter tips.
- 5) Add 5 µL of extracted RNA sample.
- 6) In parallel perform the following controls:
 - Positive Control:
 - 10 µL of Mastermix + 5 µL of CP.
 - Negative Control:
 - 10 µL of Mastermix + 5 µL supplied water (CN-H2O).
- 7) Immediately close the tubes, or plate with an adhesive film to avoid contamination.
- 8) Centrifuge briefly to collect the reaction mixture at the bottom of the tubes or microplate wells.
- 9) Run the following program on the real-time PCR instrument.

Program	Temperat	Durati	Cycle(s)	
Reverse	45°C	5 min	1	-
Denaturation	98°C	20 sec	1	-
Amplification	98°C	3 sec	40	-
	58°C	10 sec		Acquisition of fluo-


Note 1: On the CFX96™ (Bio-Rad), start the run from version 1.6 or later, then analyse with version 3.1 (see § Validation of the experiment). **Please ensure that white PCR plates are used for appropriate reading in all channels.**

Fiche technique EurobioPlex SARS-CoV-2 Fast-SVT


11. Validation of the experiment

Post-acquisition data analysis on a CFX96 PCR instrument (Bio-Rad) must be performed using version 3.1 of the CFX Manager software (Bio-Rad). In order to change to this version from a run started on an earlier version, please follow the procedure below: At the end of the run, the data file with the extension .pcrd must be opened and processed with CFX Manager (Bio-Rad) version 3.1.

If the run was started with CFX Manager v1.6 for example, to open a data file with CFX Manager v3.1, click on the CFX Manager v3.1 icon. The home screen appears.


Click on "File" and select "Open" then "Data File".


Select the file you need to analyse and click on "Open".

Fiche technique EurobioPlex SARS-CoV-2 Fast-SVT

The "Drift correction" option must be applied in the "Settings" tab as shown in the image below: click on the "Settings" tab, then on "Baseline Setting" and on "Apply Fluorescence Drift Correction".


To optimise the run analysis, select the "Log Scale" checkbox for each channel analysed. Then place the threshold bar above the background noise corresponding to the middle of the exponential phase. The 5 channels of interest have high RFUs and differ according to the channel considered, so the "Log Scale" option allows for a better reading and analysis of the run.

Once this step has been completed, the analysis can proceed for the 5 channels considered **FAM, HEX, Texas Red, CY5 and Quasar 705**.

Fiche technique EurobioPlex SARS-CoV-2 Fast-SVT

To validate the assay, the Ct values for the controls must be the following (Table 4). Outside of these values, the experiment cannot be validated.

Table 4: Run validation

Positive Control	
FAM	Ct ≤ 30
HEX	Ct ≤ 30
Texas Red	Ct ≤ 30
Cy5	Ct ≤ 30
Quasar 705	Ct ≤ 30
Negative Control	
FAM	Ct not determined
HEX	Ct not determined
Texas Red	Ct not determined
Cy5	Ct not determined
Quasar 705	Ct not determined

Fiche technique EurobioPlex SARS-CoV-2 Fast-SVT

12.Data analysis and interpretation

For clinical samples, the following results are possible:

* Ct threshold for sample: + Positive → positive Ct < 40

For SARS-CoV-2 targets and mutations in samples:

Analyse the results as described below.

SARS-CoV-2 (Texas Red), V213G mutation (FAM), K417N mutation (HEX), L452R mutation (Quasar 705), Endogenous internal control (Cy5) : Ct < 40

Results are analyzed in the **FAM**, **HEX**, **Texas Red**, **Quasar 705** and **Cy5** channels (Table 5).

Table 5: Detection of SARS-CoV-2 and mutations of interest

Test validity					
FAM	HEX	Texas Red	Quasar 705	Cy5	
-	-	+/-	+	+/-	SARS-CoV-2 positive and presence of L452R mutation (Delta suspected)
-	+	+/-	-	+/-	SARS-CoV-2 positive and presence of K417N mutation (Omicron suspected BA.x or BA.2)
-	-	+	-	+/-	SARS-CoV-2 positive without K417N, V213G and L452R mutations (Omicron, Beta et Delta excluded)
+	+	+/-	-	+/-	SARS-CoV-2 positive and presence of K417N and V213G mutations (BA.2 sublineage of Omicron suspected)
+	+	+/-	+		SARS-CoV-2 positive and presence of L452R, K417N and V213G mutations (BA.4/BA.5 sublineage of Omicron suspected)
-	-	-	-	+	SARS-CoV-2 negative

Fiche technique EurobioPlex SARS-CoV-2 Fast-SVT

-	-	-	-	-	NI ²
---	---	---	---	---	-----------------

¹+/- : if a signal is detected in one of the FAM/HEX/Quasar 705 channels but not in the Texas Red channel, the sample is considered to be SARS-CoV-2 positive because one or more mutations of interest are revealed.

²NI : not interpretable because of RT-PCR inhibition or extraction problems: no conclusion can be given. It is then recommended to collect a new sample and/or repeat the extraction and/or dilute the sample 5 times.

Typing for V213G, K417N and L452R mutations is presumptive. Only genomic sequencing provides a definitive result.

Limitations on use and interpretation:

- ❖ The EurobioPlex SARS-CoV-2 Fast-SVT kit is used for first-line diagnostic purposes
- ❖ All samples should be treated as potentially infected with SARS-CoV-2, and local biosafety regulations should be carefully followed.
- ❖ Interpretation of results should consider the possibility of false negatives and false positives.

False negatives may be due to:

- Inadequate collection of samples, or incorrect storage,
- Samples outside the viremia phase,
- Inappropriate extraction conditions or use of non-validated PCR instruments,
- Experimentation not respecting all elements of these instructions for use.

False positives may be due to:

- Contamination due to mishandling of high positive samples, the positive control, or the PCR amplification products,
- Failure to follow the procedure described in these instructions for use, especially to avoid sources of contamination.
- ❖ All results should be interpreted by medical professionals in the context of the patient's clinical situation, history and symptoms.
- ❖ This assay does not exclude the presence of pathogens other than SARS-CoV-2.
- ❖ A negative result of this assay does not absolutely exclude a possible SARS-CoV-2 infection.

Fiche technique EurobioPlex SARS-CoV-2 Fast-SVT

13. Performance analysis

A dilution range was performed using a plasmid mixture ranging from 10^5 to 1 copy/ μ L. This dilution range was used to determine the detection limit (at 100% detection) of the EBX-048 kit.

In addition, a detection limit was determined on the WHO 1st International Standard for SARS-CoV-2 RNA (NIBSC code 21/146).

Detection limitation/analytical sensitivity on CP:

SARS-CoV-2 : 5 copies/ μ L

K417N mutation : 5 copies/ μ L

L452R mutation : 5 copies/ μ L

V213G mutation : 10 copies/ μ L

IC Cy5 : 50 copies/ μ L

Detection limitation/analytical sensitivity on WHO 1st International Standard for SARS-CoV-2 RNA (cut-off 95%):

SARS-CoV-2: 3.197 copies/ μ L

Signal variability on FAM, HEX, Texas Red, Cy5 and Quasar 705 channels

- Intra-experimental variability:

Mean CV %				
V213G	K417N	L452R	SARS-CoV-2	Endogenous internal control
0,97%	0,51%	1.05%	0.54%	0,62%

CV: coefficient of variation

- Inter-batch variability:

Mean CV %

Fiche technique EurobioPlex SARS-CoV-2 Fast-SVT

V213G	K417N	L452R	SARS-CoV-2	Endogenous internal control
1.13%	1.35%	0.74%	2.92%	2.09%

CV: coefficient of variation

Diagnostic specificity

This validation included:

- 23 samples negative for all known coronaviruses, including respiratory panel samples characterised for other viruses to test for cross-reactions.
- 1 positive samples for CoV NL63.

Number of sample	Positive for	Coronavirus	Status SARS-CoV-2 EBX-048
2	Adenovirus	Negative	Negative
1	Parainfluenzae type 4	Negative	Negative
1	Rhinovirus	Negative	Negative
1	Legionella pneumoniae	Negative	Negative
2	Mycoplasma tuberculosis	Negative	Negative
6	Enterovirus	Negative	Negative
3	Streptococcus pneumoniae	Negative	Negative
2	Staph.epiderm/oralis/mitis/sanguis	Negative	Negative

Fiche technique EurobioPlex SARS-CoV-2 Fast-SVT

1	Candida albicans	Negative	Negative
2	Pseudomonas aeruginosa	Negative	Negative
1	CoV NL63	Positive	Negative
2	VRS	Negative	Negative

In silico analysis was used and demonstrated the specificity of the primers and probes for SARS-CoV-2. For example, on the sequences of MERS and SARS viruses, no significant homology of SARS-CoV-2 primers and probes that could amplify these viruses was found.

No aspecificity or cross reaction was detected.

The kit is 100% specific.

Diagnostic sensibility

The validation of the performance was performed on:

543 negative samples	
519 SARS-CoV-2 negative samples	Characterised with EurobioPlex EBX-047 CE-IVD kit
24 positive respiratory samples	RSV-A/B, Parainfluenzae, Rhinovirus, Metapneumovirus, CoV NL63, CoV OC43, CoV 229E, CoV KU1

391 positive samples	
122 omicron BA.2 sublineage variant samples with V213F mutations	Characterized by EurobioPlex EBX-047 CE-IVD kit and sequencing
61 Omicron variant BA.1	Characterized by EurobioPlex EBX-047 CE-IVD kit
4 Alpha variant or WT *	Characterized by EurobioPlex EBX-046 CE-IVD kit
160 Delta variant samples with the L452R mutation	Characterized by EurobioPlex EBX-047 CE-IVD kit

Fiche technique EurobioPlex SARS-CoV-2 Fast-SVT

44 omicron BA.4/BA.5 sublineage variant samples with L452R, K417N and V213F mutations

Characterized by sequencing

*None of EBX-047 or EBX-048 cannot differentiate a Alpha variant or a Wild-Type

For this study, sample extraction was performed on Maelstrom 9600 (TANBead) with the OptiPure Viral Auto Plate kit (ref: W665A46), and RT-PCR on CFX96 (Bio-Rad).

		EBX-048	
		SARS-CoV-2 Positive	SARS-CoV-2 Negative
Characterized by EurobioPlex EBX-046 CE-IVD kit	SARS-CoV-2 Positive	391	0
	SARS-CoV-2 Negative	0	543

Overall sensitivity: > 99% (391/391)

Overall specificity : > 99% (543/543)

Overall concordance: > 99% (934/934)

SARS-CoV-2 negative samples pre-tested with the EBX-047 CE-IVD kit: Specificity for SARS-CoV-2 for mutations of interest			
		SARS-CoV-2 (RdRp1/RdRp2 targets)	
SARS-CoV-2 negative	N=543	(543/543)	
Specificity to mutations		100%	

Fiche technique EurobioPlex SARS-CoV-2 Fast-SVT

SARS-CoV-2 positive samples pre-tested with the EBX-047 CE-IVD kit:

Sensitivity for SARS-CoV-2 and specificity for mutations of interest

N=391 (+) positif / (-) négatif		SARS-CoV-2 (RdRp1/RdRp2 targets)	Mutations		
			L452R	K417N	V213G
SARS-CoV-2 samples without L452R, K417N and V213G mutations	N=4	+ (4/4)	-	-	-
Sensitivity to SARS-CoV-2		100%	N/A		
Specificity to mutations		N/A	100%		
Omicron BA.2 variants samples with the K417N and the V213G mutation	N=122	+ (122/122)	-	+ (122/122)	+ (122/122)
Sensitivity to SARS-CoV-2		100%	N/A		
Specificity to mutations		N/A	N/A	100%	100%
Delta variants samples with the L452R mutation	N=160	+ (160/160)	+ (160/160)	-	-
Sensitivity to SARS-CoV-2		100%	N/A		
Specificity to mutations		N/A	100%	N/A	N/A

Fiche technique EurobioPlex SARS-CoV-2 Fast-SVT

Omicron BA.x variants samples with the K417N mutation	N=61	N=61	-	+ (61/61)	-
Sensitivity to SARS-CoV-2		100%	N/A		
Specificity to mutations		N/A	N/A	100%	N/A
Omicron BA.4/BA.5 variants samples with V213G, L452R and K417N mutations	N=44	N=44	+ (44/44)	+ (44/44)	+ (44/44)
Sensitivity to SARS-CoV-2		100%	N/A		
Specificity to mutations		N/A	100%	100%	100%

Note: The EBX-048 test positively detected a previously sequenced BA.2 virus strain provided by the RSS for Respiratory Infection Viruses (HCL – Lyon).

K417N sensitivity: > 99% (227/227)
K417N specificity: > 99% (164/164)
K417N concordance: > 99% (391/391)

L452R sensitivity: > 99% (204/204)
L452R specificity: > 99% (187/187)
L452R concordance: > 99% (391/391)

V213G sensitivity: > 99% (166/166)
V213G specificity: > 99% (225/225)
V213G concordance: > 99% (391/391)

Fiche technique EurobioPlex SARS-CoV-2 Fast-SVT

14. Bibliography

- Cherian, Sarah, Varsha Potdar, Santosh Jadhav, Pragya Yadav, Nivedita Gupta, Mousmi Das, Partha Rakshit, et al. 2021. « Convergent Evolution of SARS-CoV-2 Spike Mutations, L452R, E484Q and P681R, in the Second Wave of COVID-19 in Maharashtra, India ». Preprint. Molecular Biology. <https://doi.org/10.1101/2021.04.22.440932>.
- Hu, Ben, Hua Guo, Peng Zhou, et Zheng-Li Shi. 2021. « Characteristics of SARS-CoV-2 and COVID-19 ». *Nature Reviews Microbiology* 19 (3): 141-54. <https://doi.org/10.1038/s41579-020-00459-7>.
- Lauring, Adam S., et Emma B. Hodcroft. 2021. « Genetic Variants of SARS-CoV-2—What Do They Mean? » *JAMA* 325 (6): 529. <https://doi.org/10.1001/jama.2020.27124>.
- Plante, Jessica A., Brooke M. Mitchell, Kenneth S. Plante, Kari Debbink, Scott C. Weaver, et Vineet D. Menachery. 2021. « The Variant Gambit: COVID-19's next Move ». *Cell Host & Microbe* 29 (4): 508-15. <https://doi.org/10.1016/j.chom.2021.02.020>.
- Seyed Hosseini, Elahe, Narjes Riahi Kashani, Hossein Nikzad, Javid Azadbakht, Hassan Hassani Bafrani, et Hamed Haddad Kashani. 2020. « The Novel Coronavirus Disease-2019 (COVID-19): Mechanism of Action, Detection and Recent Therapeutic Strategies ». *Virology* 551 (décembre): 1-9. <https://doi.org/10.1016/j.virol.2020.08.011>.
- Srivastava, Surabhi, Sofia Banu, Priya Singh, Divya Tej Sowpati, et Rakesh K. Mishra. 2021. « SARS-CoV-2 Genomics: An Indian Perspective on Sequencing Viral Variants ». *Journal of Biosciences* 46 (1): 22. <https://doi.org/10.1007/s12038-021-00145-7>.
- European Centre for Disease Prevention and Control :
<https://www.ecdc.europa.eu/en/covid-19/variants-concern>
- Haut Conseil de la Santé Publique :
<https://www.hcsp.fr/Explore.cgi/AvisRapportsDomaine?clefr=963>
<https://doi.org/10.1101/2021.08.04.455140>

15. Quality control

In accordance with Eurobio's ISO EN 13485 certified quality management system, each batch of EurobioPlex SARS-CoV-2 Fast SVT is tested according to predefined specifications to ensure consistent product quality.

16. Waste disposal

Eliminate all waste in accordance with the legislation on DASRI.

Fiche technique EurobioPlex SARS-CoV-2 Fast-SVT

17. Incident report

Serious incidents related to the system must be reported to Eurobio Scientific.

18. Technical assistance

For assistance with our products, please contact our technical support.

Eurobio Scientific's customer service can be contacted by e-mail at adv@eurobio-scientific.com or by phone at +33 (0)1.69.79.64.80

7, avenue de Scandinavie
91940 Les Ulis Cedex
FRANCE