Instructions For Use EurobioPlex Chikungunya/Dengue/Leptospire

eurobio
SCIENTIFIC

EurobioPlex *Chikungunya/Dengue/Leptospire* **REAL TIME RT-PCR**

For real time **qualitative** RT-PCR

EBX-050-25 (for 25 reactions)

EBX-050-50 (for 50 reactions)

EBX-050-100 (for 100 reactions)

EBX-050-200 (for 200 reactions)

EBX-050-600 (for 600 reactions)

25/50/100/200/600 reactions

EBX-050 V 1.02 – 16/12/2024

Validated by:

- CFX96™ Real Time PCR detection system (Bio-Rad) based on analysis by CFX Manager version 3.1 (Bio-Rad)

Instructions For Use

Available on request at info@eurobio-scientific.com

As soon as it is operational, this *Résumé des Caractéristiques de Sécurité and Performance* (RCSP – Summary of Safety and Performance Characteristics) will be made available by the notified certifying organisation on EUDAMED. This Summary is also available on request.

Instructions For Use EurobioPlex Chikungunya/Dengue/Leptospire

Table of contents

ENGLISH.....	1
FRANCAIS.....	24
ESPANOL.....	48
ITALIANO	68
1. General Information.....	3
2. Test procedure	5
3. Symbols	6
4. Principles of Operation.....	7
5. Kit Components	8
6. Conservation and storage	8
7. Required materials not provided	8
8. Real time PCR instrument	9
9. Important preliminary instructions and precautions	9
10. Protocol	10
10.1 Sample collection	10
10.2 Nucleic acid extraction	11
10.3 Performance of real time RT-PCR procedures	11
10.4 Detailed protocol.....	13
11. Experimental validation.....	14
12. Data analysis and interpretation.....	16
13. Analysis of performance.....	17
14. Quantity control	22
15. Bibliography.....	22
16. Elimination of waste	22
17. Declaration of incident.....	22
18. Technical assistance	23

Instructions For Use EurobioPlex Chikungunya/Dengue/Leptospire

1. General Information

Chikungunya and Dengue Fever viruses are Arthropod-Borne viruses (aka Arbovirus), that is, they are transmitted by arthropods such as mosquitoes, ticks, sandflies, etc. The Dengue virus belongs to the *Flaviviridae* family and the Chikungunya virus to the *Togaviridae* family. These viruses are transmitted by various species of mosquito of the *Aedes* genus, among others *aegypti* (principally in tropical and subtropical America and in the south of the United States) and *albopictus* (principally found in Asia), *furcifer*, *africanus*, *taylori* and *luteophaeus* which are globally prevalent, but especially so in the tropics. To date, the two first-mentioned species have been identified as the major vectors of epidemics, this being due to their adaptability to zones of human settlement. France's southern regions are more particularly affected.

There are 4 serotypes for Dengue Fever (DENV 1-4), whose genome is a single strand RNA. Virus replication and transcription take place in the cytoplasm of the host cell. Initial clinical signs are fever associated with muscular and joint pain. In the majority of cases, there are no complications, but in 20% of cases Dengue symptoms may be severe to the point of life-threatening, particularly when secondary infections occur involving haemorrhagic shock, since infection by one strain does not protect the host from infection by another strain.

The Chikungunya (CHIKV) virus is a small, spherical, enveloped virus whose genome is a single strand RNA of approximate size 12 Kbp. The viraemic phase of Chikungunya infection lasts from 3 to 10 days. Major clinical signs are the rapid onset of high fever (> 40 °C) lasting from 24 to 48 hours, associated with severe often multiple joint pain and myalgia, particularly affecting the extremities, such as wrists, ankles or fingers, where joint swelling may be pronounced, associated with maculopapular or erythematous skin rash and subacute oedema. CHIKV infections may evolve into a chronic condition. The symptoms are similar to those of malaria or dengue, and not infrequently two infections may be present simultaneously. In a proportion varying from 10% to 70 % of cases, clinical relapse may be observable several months after the acute phase. There is no specific treatment for chikungunya.

Leptospire bacteria are often categorised as belonging to the *Leptospira interrogans* species and are the cause of the condition known as leptospirosis, which is considered to be an ubiquitous zoonosis. The number of reported cases has significantly risen in recent times. As a zoonosis, Leptospirosis affects animals such as rat or pig living in proximity to and infecting humans. Leptospire is a pathogen for humans.

Leptospirosis is prevalent in a number of parts of the world, such as the islands of La Réunion and Wallis and Futuna, as well as in the Caribbean, in French Guyana, in the tropics, and in South-East Asia, Cambodia, Vietnam and Laos. The leptospire has even been detected in the Poitou Charente region of France, and elsewhere in the country's south and south-west, in humid environments close to pig raising farms.

Leptospire is a spiral bacteria whose motility depends on their endoflagellae. The extremely wide range of non-specific symptoms of leptospirosis and the diversity of the organs affected make clinical diagnosis difficult, since the condition can be mistaken for a severe bout of influenza. Symptoms may

Instructions For Use EurobioPlex Chikungunya/Dengue/Leptospire

begin with diffuse or localised aches and pains, which, if not diagnosed in good time, may lead to impaired speech and mental faculties. In some cases, symptoms include jaundice and/or nephritis. Recovery is possible thanks to kidney dialysis and treatment by antibiotics, provided the condition remains moderately or less severe, in 5 or 6 weeks (although bacteria may still be found in urine several weeks after patients have become asymptomatic). Leptospirosis may develop into a severe form causing violent and even life-threatening haemorrhaging or severe renal dysfunction. There are numerous serovars (for example *Leptospira icterohaemorrhagiae*, *Leptospira canicola*, *Leptospira pomona*) which present no homogenous antigenic signatures. This makes it difficult to design effective vaccines.

The World Health Organisation recommends that the biological confirmation of leptospirosis takes place by means, either of isolation of the bacteria involved, or of the identification of nucleic acids in biological samples, or again, in the clinical context, by reliance on positive serological indicators taking into consideration epidemiological factors. The usefulness of bacteriological culture is extremely limited, given the long bacterial growth time. Early leptospirosis detection is possible in real time by PCR, generally using patient blood plasma or serum. Optimum time for taking blood samples is after 8 days of illness. There is no interest in real time PCR blood testing after 10 days from the onset of symptoms (end of bacterial period). The Haute Autorité de Santé (HAS – High Authority for Healthcare in France) has indicated that simple (i.e., NON-real time) PCR testing is of no interest for the specific biological diagnosis of leptospirosis. The HAS places the greatest value on real time PCR testing for the early specific biological diagnosis of leptospirosis. Testing requires blood samples taken within 8 days of onset of fever, preferably before the administration of antibiotics. But because of the normally low levels of leptospires in the blood, negative real time PCR results need to be confirmed by further serological investigation.

Instructions For Use EurobioPlex Chikungunya/Dengue/Leptospire

2. Test procedure

EurobioPlex Chikungunya/Dengue/Leptospire testing relies on inverse transcription amplification followed by real time polymerase chain reaction (RT-PCR), designed for the qualitative detection of the presence or absence of the Chikungunya and Dengue viruses and of the Leptospire in an extract of nucleic acids.

The EurobioPlex Chikungunya/Dengue/Leptospire test procedure is indicated for the diagnosis of suspected infection in humans, or alternatively, as an adjunct to serological diagnosis both determined and indeterminate. Extract of nucleic acids is the basic material required for the EurobioPlex Chikungunya/Dengue/Leptospire procedure. Nucleic acids are extracted from patient blood plasma, serum or urine.

This amplification system has been validated on samples of plasma, serum and urine, and by means of endogenous internal controls.

The EurobioPlex Chikungunya/Dengue/Leptospire diagnostic procedure relies on *in vitro* analysis and must be undertaken by qualified analytical medical laboratory personnel. The test package must not be recycled for further use. It is not intended for chikungunya / dengue/ leptospire blood or organ bank screening.

Instructions For Use EurobioPlex Chikungunya/Dengue/Leptospire

3. Symbols

Reference catalogue

Batch code

Maximum Temperature

Use-by date

Content for n number of tests

Manufacturer

Date of Manufacture

CE-marked product

Dispositif Médical de Diagnostic *in vitro*

Consult User Instructions

Beware!

Do not use if packaging is damaged

Store away from sunlight

Instructions For Use EurobioPlex Chikungunya/Dengue/Leptospire

4. Principles of Operation

The EurobioPlex Chikungunya/Dengue/Leptospire kit provides a test procedure based on the amplification of the ribonucleic acid (RNA) of the Chikungunya and Dengue viruses and on the amplification of the deoxyribonucleic acid (DNA) coding for ribosomal RNA and ribosomal RNA of the Leptospire, including the endogenous internal control of RT-PCR extraction and inhibition, using RT-PCR enabled real time amplification. Testing takes place of the nucleic acids extracted from a given sample by means of a single reaction in a single well. The kit's primers and probes have the capability of detecting (screening) the RNA of the 4 serotypes of Dengue without typing.

The RT-PCR endogenous internal control process ensures that a negative result cannot be caused by poor extraction and/or by excess of PCR inhibitors.

Chikungunya virus RNA is detected by an FAM marked probe. Dengue virus is detected by a HEX marked probe. The ribosomal RNA and DNA coding for the ribosomal RNA of the leptospire are detected by a Texas Red marked probe. The endogenous internal control of RT-PCR extraction and inhibition is detected by a CY5 marked probe. All probes emit specific fluorescent spectra.

The apparatus is not automated. Users may under their responsibility use other devices than those validated, in which case performance is not guaranteed.

Table 1: Target detection by fluorophore

Target	Fluorophore	Excitation	Emission
CHIKUNGUNYA	FAM	495 nm	515 nm
DENGUE	HEX	535 nm	555 nm
LEPTOSPIRE	TEXAS RED	585 nm	605 nm
Endogenous internal control	Cy5	650 nm	670 nm

Equivalent channels on different PCR instruments:

- FAM channel (ABI Systems, SmartCycler II, Mx Systems, Chromo4/CFX96, T-COR 8®-IVD), Channel 510 (LC 480), Green Channel (RotorGene)
- HEX channel (Chromo 4/CFX96, Mx Systems, T-COR 8®-IVD), Canal VIC (ABI Systems), Alexa532 Channel (SmartCycler II), 580 Channel (LC 480), Yellow Channel (RotorGene)
- Texas Red Channel (ABI Systems, SmartCycler II, Chromo 4/CFX96™, Mx Systems, T-COR 8™-IVD), LC Red 610 (LC480®), Orange Channel (RotorGene)
- Canal Cy5 (ABI Systems, Mx Systems, Chromo4/CFX96, T-COR 8®-IVD), Alexa647 Channel (SmartCycler II), 660 Channel (LC 480), Red Channel (RotorGene)

Note : On LC480 instrument II : Apply colour compensation for wavelengths FAM-HEX/VIC-Cy5 (465-510,533-580,618-660).

Instructions For Use EurobioPlex Chikungunya/Dengue/Leptospire

5. Kit Components

The Chikungunya/Dengue/Leptospire real time RT-PCR test package is delivered ready for use as a specific detector of the afore-named pathogens.

Fluorescence is emitted in the course of the PCR process and individually measured by an optical system. Amplified fragment detection is by a fluorimeter using the channels shown in Table 1.

The kit contains the reagents and enzymes required for the amplification of the viral RNA of Chikungunya and Dengue and for the amplification of the DNA coding for ribosomal RNA and ribosomal RNA of the Leptospire; as well as the reagents and enzymes required for the endogenous control of RT-PCR extraction and inhibition (Table 2).

Table 2 : Kit Components

Colour of Cap	Kit Content	EBX-050 (25 tests)	EBX-050 (50 tests)	EBX-050 (100 tests)	EBX-050 (200 tests)	EBX-050 (600 tests)	Reconstitution
Red	Enzymes	225 µl	450 µl	900 µl	2x900 µl	4x1350 µl	Ready for use
Transparent	Oligomix	80 µL	160 µL	320 µL	2X320 µL	3x640 µL	Ready for use
Blue	Water = Negative Control (NC-H ₂ O)	1,5 ml	1,5 ml	1,5 ml	1,5 ml	1,5 ml	Ready for use
Yellow	Positive Control (PC) Chikungunya/ Dengue/ Leptospire/ endogenous CI	75 µl	150 µl	300 µl	2 x 300 µl	3x300 µL	Ready for use

6. Conservation and storage

All reagents must be stored at temperatures in the range of -15°C to -22°C.

All reagents may be used prior to the use-by date shown on label.

Take care to avoid more than five repeated cycles of freezing/thawing (>5x) as this may impair analytical sensitivity.

7. Required materials not provided

- ◇ Biological extractor hood
- ◇ Real time PCR machine
- ◇ Centrifuge for microtubes

Instructions For Use EurobioPlex Chikungunya/Dengue/Leptospire

- ◇ Vortex
- ◇ Plates / tubes for PCR reactions in real time
- ◇ Micropipettes
- ◇ DNase-free tips and RNase-free filter tips for micropipettes
- ◇ Sterile microtubes
- ◇ Talc-free gloves

8. Real time PCR instrument

The EurobioPlex kit has been developed and validated for use with the Bio-Rad CFX96™ Real Time PCR detection system with analysis by Bio-Rad CFX Manager version 3.1 applying drift correction, between 15 and 45 cycles (cf. §11.).

9. Important preliminary instructions and precautions

Carefully read instructions before starting protocol

- ◇ Kit is delivered in a dry ice package, and its constituents in frozen condition. If any part arrives unfrozen, or if tubes have been damaged in transport, the technician must not use the kit.
- ◇ This experimental procedure must be carried out by medical biology analytical laboratory technicians.
- ◇ Ensure that instruments have been installed, calibrated and properly maintained in accordance with manufacturer's recommendations.
- ◇ Clinical samples must be considered as potentially infectious materials and must be prepared under a laminar flow extractor hood.
- ◇ Experimentation must be performed in accordance with good laboratory practice.
- ◇ Do not use this kit after its use-by date.
- ◇ Avoid cycles of reagent thawing and refreezing, as this may impair test sensitivity.
- ◇ Avoid prolonged exposure of reagents to light, and limit exposure to the time required for preparation of the PCR plate, failing which performance cannot be guaranteed.
- ◇ Once the reagents are unfrozen, briefly centrifuge tubes before use.
- ◇ For very heavy viral loads ($Ct < 15$), bell curve shaped amplification may be observed. This is a positive signal. From 15 cycles onward, a classic exponential amplification curve is presented.
- ◇ It is recommended to separate out three distinct zones of work: 1) Isolation of nucleic acids, 2) Preparation of reaction mixture, and 3) Amplification / Detection of amplified products.
- ◇ In each zone of work, wear separate sets of lab coats and talc-free gloves.
- ◇ Pipettes, reagents and other work materials must not circulate from one zone to another.
- ◇ Special attention is required to protect the purity of reagents and reaction mixtures.
- ◇ Appropriate methods of nucleic acid preparation / extraction must be adopted to ensure high quality RNA/DNA production and RT-PCR operation, and in particular to avoid any risk of contamination of RNases/DNases.

Instructions For Use EurobioPlex Chikungunya/Dengue/Leptospire

- ◇ Use filter tips for micropipettes, RNase-free and DNase-free.
- ◇ No use of mouth for pipetting. No food, drink or smoking to be allowed in the laboratory.
- ◇ Avoid aerosols.
- ◇ The kit is not designed for once only use. When reused, recommendations must be followed, notably in regard to authorized thawing/refreezing cycles and precautions to prevent contamination of reagents.
- ◇ Kit constituents must not be used separately, i.e. with other reagents, or with reagents from other batches.
- ◇ The package is not automated. Users make use of other procedures than those validated under their sole responsibility, in which case performance is not guaranteed.
- ◇ Reagent temperatures after thawing must be carefully managed (+COC/+COC or room temperature) in order not to adversely affect performance.

10. Protocol

10.1 Sample collection

- ◇ Use sterile tubes for sample collection. When nucleic acid is extracted from serum, collect samples in dry tubes. When extracting from plasma, collect samples in EDTA tubes.
- ◇ Lipemic plasma may induce PCR inhibition and prevent pathogen detection

Use of heparin as anti-coagulant is prohibited.

Users must properly manage and control their own conditions of sample collection, transport, storage and extraction, and must use suitable means to ensure high quality nucleic acid extraction.

- ◇ It is recommended that samples be extracted immediately, or if kept, then in accordance with the recommendations for the storage of samples prior to extraction (Table 3).
- ◇ The bibliographical recommendations of section “15. Bibliography” provide indicative data on the stability of samples and RNAs.

Table 3: Recommendations for the storage of blood prior to the preparation of plasma or serum, and for the storage of plasma, serum and urine before extraction

Recommendations for whole blood storage (before plasma and serum preparation)	
<24h	
Recommendations for maximum sample storage times of urine, plasma or serum before extraction	
Room temperature	2h
+2°C/+8°C	5 days
<-70°C (preferable to -20°C)	Long term storage (>5 days and no more than 2 months at -COC)

Instructions For Use EurobioPlex Chikungunya/Dengue/Leptospire

Attention	
	◇ Users may take into consideration the sample storage recommendations of the World Health Organisation or of the Haute Autorité de Santé (High Authority for Healthcare in France) ◇ Nucleic acids after extraction must be stored at <-70°C to ensure their stability. If RNAs are stored at <-70°C for more than a year, or if they are stored at -20°C, the Ct readings obtained may increase. Re-extraction is advised if a biological sample has been stored for more than one year. It is advisable to limit to 3 the number of times RNAs are unfrozen. ◇ The transportation of clinical samples is governed by local regulations on the transportation of infectious materials.

10.2 Nucleic acid extraction

Users of the real time RT-PCR technology are responsible for ensuring its compatibility with the nucleic acid extraction system used. In regard to this package, we recommend methods of nucleic acid extraction from samples of plasma, serum or urine, and adherence to the instructions of the supplier of the extraction kit used.

Performance levels, including the interfering pathogens described in §13 of this document, were obtained using the following extraction techniques: Qiagen columns (QIAamp Ultrasens Virus kit Ref : 53704), or EZ1 extractor (Qiagen) or Maelstrom TANbeads (Taiwan Advanced Nanotech).

10.3 Performance of real time RT-PCR procedures

General comment:

Positive controls and the endogenous internal controls of RT-PCR extraction and inhibition high levels of matrix concentration. Great care must be taken when handling to avoid contamination. PCR operational controls involve testing Positive Control (PC) and Negative Control (water supplied = NC-H₂O) (see II-2/6) of RT-PCR protocol in real time.

Instructions For Use EurobioPlex Chikungunya/Dengue/Leptospire

Outline of procedure:

1 - REACTION MIX / MASTERMIX PREPARATION

Number of reactions	N+3
Enzymes	(N+3) x 7.5 µl
Oligomix	(N+3) x 2.5 µl
Total volume Mastermix	(N+3) x 10µl

2 - PREPARATION OF REACTIONS

Sample

10 µl Mastermix
+
5µl extracted sample

Positive control (PC)

10 µl de Mastermix
+
5µl CP

Negative control (NC-H2O)

10 µl de Mastermix
+
5µl Biological molecular water (NC-H2O)

3 – REAL TIME PCR INSTRUMENT

Programme	Temperature	Duration	Cycle(s)	
Reverse Transcription	45 °C	5 min	1	-
Denaturation	98°C	20 sec	1	-
Amplification	98°C	3 sec	45	-
	58°C	45 sec		Acquisition of fluorescence

See 11. Experimental validation for detailed Settings with drift correction and analysis between 15 and 45 cycles.

Instructions For Use EurobioPlex Chikungunya/Dengue/Leptospire

10.4 Detailed protocol

- 1) Homogenise the Enzyme tube and vortex the Oligomix before centrifuging
- 2) Prepare Mastermix as shown below, N being the number of reaction(s), prepare Mastermix quantity for minimum of N+3 reactions.

Number of reactions	N+3
Enzymes	(N+3) x 7.5 µl
Oligomix	(N+3) x 2.5 µl
Total volume Mastermix	(N+3) x 10µl

- 3) Homogenise Mastermix as prepared in 2) above and briefly centrifuge.
 - 4) Using micropipette and filter tips distribute 10 µL Mastermix into each microplate to set up real time PCR.
 - 5) Add 5 µL extracted sample.
 - 6) Perform following controls in parallel:
 - Positive Control: 10 µl Mastermix + 5 µl PC positive control
 - Negative Control: 10µl Mastermix + 5µl water supplied (NC-H2O)
 - 7) Prevent contamination by immediately sealing off under adhesive film or isolate under transparent stoppers.
 - 8) Briefly collect reaction mixture at bottom of tubes or microplate wells.
 - 9) Initiate and follow programme shown below on the real time PCR instrument.
- Approximate duration: 1h10

Programme	Temperature	Duration	Cycle(s)	
Reverse Transcription	45 °C	5 min	1	-
Denaturation	98°C	20 sec	1	-
Amplification	98°C	3 sec	45	-
	58°C	45 sec		Acquisition of fluorescence


Note: On Bio-Rad CFX96, initiate run using release 1.6 or later, followed by analysis on release 3.1 (see **11 Experimental validation for detailed Settings with drift correction and analysis between 15 and 45 cycles**).

Instructions For Use EurobioPlex Chikungunya/Dengue/Leptospire


11. Experimental validation

Analysis of post-acquisition data on Bio-Rad PCR CFX96 requires running release 3.1 of Bio-Rad CFX Manager software. To switch to this release from an earlier-versioned run, following procedure should be followed: at end of run, first open data file with suffix < .pcrd > to process by means of Bio-Rad CFX Manager release 3.1.

If running CFX Manager v1.6 software, for example, click on CFX Manager v3.1 icon to open data file running CFX Manager v3.1 software. Welcome screen appears.

- Click on File, select Open then Data File.


Instructions For Use EurobioPlex Chikungunya/Dengue/Leptospire

- Select file it is wished to analyse and click on Open.

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

Cycles to Analyze

In Settings tab -> Drift Correction option must be applied as shown below. Click on Settings tab, then on Baseline Setting and on Apply Fluorescence Drift Correction. -> Under Cycles to Analyze, analysis takes place between 15 and 45 cycles, as shown below.

- **Analysis:** using Log Scale, position thresholds between low and middle of PC positive control exponential curves, as shown below. On completion of this stage, analysis may begin on the FAM, HEX, Texas Red and Cy5 channels.

Instructions For Use EurobioPlex Chikungunya/Dengue/Leptospire

For dosage validation, control Ct values must be as shown in Table 4 below. Outlier values, if apparent, invalidate experimentation.

Table 4: Positive and Negative Control specifications to be validated

Positive Control	
Canal FAM	Ct ≤ 28
Canal HEX	Ct ≤ 28
Canal Texas Red	Ct ≤ 28
Canal CY5	Ct ≤ 30
Negative Control	
Canal FAM	Indeterminate Ct
Canal HEX	Indeterminate Ct
Canal Texas Red	Indeterminate Ct
Canal CY5	Indeterminate Ct or > 35

12. Data analysis and interpretation

As regards endogenous extraction and inhibition control of RT-PCR samples, evaluation of the correct operation of the RT-PCR reaction may take place on the Cy5 channel measuring endogenous control. A number of different situations may arise:

1/ Positive outcome from test of RT-PCR internal endogenous control of extraction and inhibition, when $Ct \leq 35$. In this case, the nucleic acids have been correctly extracted, there are no RT-PCR inhibitors and the result can be validated.

2/ Negative/indeterminate outcome from test of RT-PCR internal endogenous control of extraction and inhibition, when $Ct > 35$. In this case, the nucleic acids have not been correctly extracted, or reverse transcription is impaired, or the presence of PCR inhibitors has adversely affected the PCR reaction. This situation may arise with lipemic plasma. Repeat of the extraction process is recommended, or if not, dilution of the sample, except if a specific signal appears in the 3 FAM, HEX or Texas Red channels.

In regard to clinical samples, the outcomes shown in Table 5 are possible.

Instructions For Use EurobioPlex Chikungunya/Dengue/Leptospire

Table 5 : Analysis and interpretation of clinical sample results

Positivity threshold for FAM/HEX/Texas Red: Ct ≤ 45

Chikungunya	Dengue	Leptospire	Endogenous internal control	Test/Interpretation validity
FAM	HEX	Texas Red	CY5	
+	-	-	Ct ≤ 35	YES : Presence of Chikungunya nucleic acids, Negative for Dengue and Leptospire
-	+	-	Ct ≤ 35	YES : Presence of Dengue nucleic acids, Negative for Chikungunya and Leptospire
-	-	+	Ct ≤ 35	YES : Presence of Leptospire nucleic acids; Negative for Chikungunya and Dengue
+	+	-	Ct ≤ 35	YES : Presence of Chikungunya and Dengue nucleic acids, Negative for Leptospire
-	+	+	Ct ≤ 35	YES : Presence Dengue and Leptospire nucleic acids, Negative for Chikungunya
+	-	+	Ct ≤ 35	YES : Presence of Chikungunya and Leptospire nucleic acids, Negative for Dengue
+	+	+	Ct ≤ 35	YES : Presence of 3 targets
-	-	-	Ct ≤ 35	YES : Absence of 3 targets
-/+	-/+	- /+	Ct ND or- > 35	Yes for positive targets. NI for all negative samples

NI: Not Interpretable either because of RT-PCR inhibition or because of extraction problem: no definitive conclusion is possible for the 3 targets. **The recommendation is then to take a new sample and/or repeat extraction and/or dilute sample by factor of 5.**

ND: Non Determinate

13. Analysis of performance

Limits of detection/analytical sensitivity

EBX-050 sensitivity for Chikungunya was determined by the International Standard 1st World Health Organization International Standard for Chikungunya virus RNA for Nucleic acid amplification techniques (NAT)-based assays PEI code : 11785/16, extracted on Maelstrom TANbeads8, and by Amplirun Chikungunya virus RNA control (Vircell) on 3 batches of EurobioPlex kits.

Instructions For Use EurobioPlex Chikungunya/Dengue/Leptospire

EBX-050 sensitivity for the 4 serotypes of Dengue was determined by quantified Amplirun Dengue 1, Dengue 2, Dengue 3 and Dengue 4 virus RNA control (Vircell) on 3 batches of EurobioPlex kits.

EBX-050 for leptospire was determined by titrated bacterial culture *Leptospira interrogans* serovar *Canicola* extracted by Maelstrom TANbeads8 on 3 batches of EurobioPlex kits.

To determine LOD_{95%}, Probit analysis was performed: dilution within given ranges performed by a plasmidic mixture of the 3 targets, tested on 24 replicata for each concentration.

Table 6: Analytical Sensitivity

	LOD _{100%} for PC	LOD _{95%} for PC		LOD _{100%} for RNA
Chikungunya	50 copies/μl	8.8 copies/μl	Chikungunya by WHO Standard	7.5 IU/μl eluate
			Chikungunya by RNAcontrol	55 copies/μl eluate
Dengue	5 copies/μl	3 copies/μl	Dengue 1 by RNAcontrol	3 copies/μl eluate
			Dengue 2 by RNAcontrol	10 copies/μl eluate
			Dengue 3 by RNAcontrol	30 copies/μl eluate
			Dengue 4 by RNAcontrol	30 copies/μl eluate
Leptospire	50 copies/μl	7.8 copies/μl	Leptospire by titrated bacterial culture	500 bacteria/μl eluate (95% for 250 bacteria/μl eluate)


Tests were performed on Bio-Rad CFX96.

Linearity

From 5 to 10⁵ copies/μl of PC for the 3 targets.

	PCR % efficacy	Slope	R ²
Chikungunya/FAM	97.29	-3.39	0.998
Dengue/HEX	102.92	-3.25	0.999
Leptospire/Texas Red	96.98	-3.40	0.999

Instructions For Use EurobioPlex Chikungunya/Dengue/Leptospire

Detection of different Leptospire strains

The following strains were tested and detected with EBX-050.

- *Leptospira inadai* serovar *Lyme* strain 10
- *Leptospira borgpetersenii* serovar *Javanica veldrat*
- *Leptospira weilli* serovar *Celladoni*
- *Leptospira wolbachii* serovar *Codice*
- *Leptospira genomospecies5/yanagawae* serovar *Saopaulo*
- *Leptospira interrogans* serovar *Icterohaemorrhagiae*
- *Leptospira alexanderi* serovar *Ranarum*
- *Leptospira kirschneri* serovar *Cynopteri* strain 3522C
- *Leptospira fainei* serovar *Hurstbridge* strain BUT6
- *Leptospira santarosai* serovar *Shermani* strain 1342K
- *Leptospira noguchii*

Variability of RT-PCR signal for Chikungunya (FAM), Dengue (HEX) and Leptospire (Texas Red) channels

EBX-050 reproducibility was determined on quantified RNA Amplirun for Chikungunya, Dengue 2, Vircell virus RNA control, and for nucleic acids extracted from *Leptospira interrogans* bacteria serovar *Canicola*, on the 3 batches of EurobioPlex kits.

Instructions For Use EurobioPlex Chikungunya/Dengue/Leptospire

Table 7 : Repeatability and reproducibility

	Coefficient of variation (CV %) intra-run	Coefficient of variation (CV %) inter-run	Coefficient of variation (CV %) inter-lots
Chikungunya	0.89	0,45	1.14
Dengue serotype 1	0.99	3.00	1.79
Dengue serotype 2	1.93	1.70	1.12
Dengue serotype 3	2.64	2.20	1.17
Dengue serotype 4	3.14	1.67	1.21
Leptospire	0.06	1.44	1.09

Tests were performed on Bio-Rad CFX96.

Diagnostic performance

Tests were performed on positive or negative pre-characterised samples: 265 serum, 21 plasma-EDTA, 28 urine samples, 6 EDTA anti-coagulated whole blood. Numbers of positives and negatives are detailed in Table 8.

For this study, sample extraction was performed on Maelstrom TANbeads (Taiwan Advanced Nanotech). Tests were performed on Bio-Rad CFX96.

Table 8 : Table of diagnostic performance contingency

EBX-050 Chikungunya				
		POS	NEG	TOTAL
Characterised pre Chikungunya test	POS	8	0	8
	NEG	0	141	141
	TOTAL	8	141	149
EBX-050 Dengue				
		POS	NEG	TOTAL
Characterised pre Dengue test	POS	144	3	147
	NEG	0	83	83
	TOTAL	144	86	230
EBX-050 Leptospire				
		POS	NEG	TOTAL
Characterised pre Leptospire test	POS	49	3	52
	NEG	0	60	60
	TOTAL	49	63	112

Instructions For Use EurobioPlex Chikungunya/Dengue/Leptospire

Diagnostic sensitivity and specificity performance was therefore:

Table 9 : Diagnostic Sensitivity and Diagnostic Performance

	Chikungunya	Dengue	Leptospire
Diagnostic Sensitivity (%)	>99%	98% [confidence interval (CI) at 95%: 94.6%-99.9%]	94 % [confidence interval (CI) at 95%: 87.9%-100%]
Diagnostic Specificity (%)	>99%	>99%	>99%

Absence of cross reaction with other pathogens

The chikungunya and leptospire system was experimentally tested and found non-cross-reactive with the following pathogens:

- HEV
- Parechovirus type 1
- Norovirus GI
- Norovirus GII
- HHV6
- HSV1
- HSV2
- VZV
- EBV
- Zikavirus
- Mayaro virus
- *Legionella*
- *Coxiella*
- *Bartonella henselae*
- *Bartonella quintana*
- *Clostridium difficile*

Cibles de multimarqueur RP 1	Cibles de multimarqueur RP 2
Influenza A H3N2 (Brisbane/10/07)	Influenza A H1 (Nouvelle Calédonie/20/99)
Influenza A H1N1 (NY/02/2009)	Influenza B (Florida/02/06)
Rhinovirus (Type 1A)	VRS (Type A)
Adénovirus (Type 3)	Parainfluenza (Type 2)
Parainfluenza (Type 1)	Parainfluenza (Type 3)
Parainfluenza (Type 4)	Coronavirus (HKU-1 recombinant)
Métapneumovirus (Pérou 6-2003)**	Coronavirus (OC43)
C. pneumoniae (CWL-029)	Coronavirus (NL63)
M. pneumoniae (M129)	Coronavirus (229E)
Virus Cocksackie (Type A1)	Bordetella pertussis (A639)

The Dengue system was experimentally tested and found to be non-cross-reactive with chikungunya, leptospire, HEV and zikavirus.

An overall *in silico* analysis was performed at time of design and indicated absence of potential cross reactions.

Specific *in silico* analysis demonstrated an absence of potential cross reaction for the following viruses in particular Usutu, West Nile, Japanese encephalitis, Tick encephalitis, Yellow fever, Saint Louis encephalitis, Hepatitis A, Hepatitis C.

Instructions For Use EurobioPlex Chikungunya/Dengue/Leptospire

14. Quantity control

In accordance with the Eurobio system of Quality Management, certified ISO EN 13485, each EurobioPlex Chikungunya /Dengue/Leptospire batch is tested compliant to predefined specifications in order to assure constant product quality.

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16. Elimination of waste

Waste management and elimination compliant to DASRI legislation in France

17. Declaration of incident

All serious incidents (if any) arising in connection with this EUROBIO procedure is notified to EUROBIO SCIENTIFIC and to the competent authorities of the European Union Member State in which user and/or patient are established.

Instructions For Use EurobioPlex Chikungunya/Dengue/Leptospire

18. Technical assistance

For technical assistance regarding EUROBIO products, please contact our technical support.

EUROBIO SCIENTIFIC customer service can be contacted by e-mail sent to adv@eurobio-scientific.com or by telephone on +33 (0)1.69.79.64.80.

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15. Bibliografia

- Haute Autorité de Santé: Cryopréservation de tissus, cellules et liquides biologiques issus du soin, Recommandations de bonnes pratiques, HAS / Service des bonnes pratiques professionnelles / Septembre 2009; www.has-santé.fr.
- Organizzazione Mondiale della Sanità. <https://iris.who.int/handle/10665/334254> Test diagnostici per la SARS-CoV-2: guida provvisoria/interim guidance, 11 settembre 2020. Licenza: CC BY-NC-SA 3.0 IGO
- Fleige Simone, Pfaffl Michael W. Molecular Aspects of Medicine 27 (2006) 126–139. RNA integrity and the effect on the real-time qRT-PCR performance (review).
- Alfaro-Núñez Alonzo, Stephanie Crone, Shila Mortensen, Maiken Worsøe Rosenstjerne, Anders Fomsgaard, Ellinor Marving, Sofie Holdflod Nielsen, Michelle Grace Pinto Jørgensen, Charlotta Polacek, Arieh S. Cohen, Claus Nielsen, Transbound Emerg Dis. 2022; 69:189–194 SARS-CoV-2 RNA stability in dry swabs for longer storage and transport at different temperatures
- Yu Keke, Jie Xing, Jie Zhang, Ruiying Zhao, Ye Zhang, Lanxiang Zhao, Effect of multiple cycles of freeze-thawing on the RNA quality of lung cancer tissues Cell Tissue Bank 2017 Sep;18(3):433-44

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