Instructions for Use EurobioPlex DermaScreen

EurobioPlex

DermaScreen

Real-Time PCR

For **qualitative** real-time PCR

REF

EBX-073-25 /EBX-073-50/EBX-073-100/EBX-073-200/EBX-073-400/EBX-073-600

25/50/100/200/400/600 reactions

Reference EBX-073 Version 1.06 from 04/03/2025

For use and validated with :

- CFX96™ Real Time PCR detection system (Bio-Rad) with analysis on the CFX Maestro versions 1.1 and 2.3 (Bio-Rad) software or with the CFX Manager version 3.1 (Bio-Rad) software
- CFX Opus 96 Real Time PCR detection system (Bio-Rad) with analysis on the CFX Maestro versions 1.1 and 2.3 (Bio-Rad) software or with the CFX Manager version 3.1 (Bio-Rad) software
- QuantStudio 5/7 Pro (Thermofisher) with analysis on the Design and Analysis version 2.7.0 (Thermofisher) software
- LC480 II (Roche) with analysis on the LightCycler® 480 SW version 1.5.1 (Roche) software
- LC PRO (Roche) with analysis on the LightCycler PRO version 1.0 (Roche) software
- MIC (Bio Molecular Systems) with analysis on the mic PCR version 2.12.6 (Bio Molecular Systems) software

Instructions for Use

Available on request from info@eurobio-scientific.com

The Summary of Safety and Performance Characteristics (SPPC) will be made available by the Notified Body on EUDAMED once it is operational. It can also be obtained on request.

Instructions for Use EurobioPlex DermaScreen

ENGLISH.....	1
FRANÇAIS.....	22

Table of contents

1. Introduction.....	3
2. Intended use.....	3
3. Symbols.....	4
4. Principle.....	5
5. Content of the kit.....	6
6. Storage.....	6
7. Materials required but not provided.....	6
8. Real-time PCR instrument.....	7
9. Cautions and notes.....	7
10. Procedure.....	8
11. Validation of the PCR run.....	11
12. Data analysis and interpretation.....	15
13. Performances analysis.....	17
14. Bibliography.....	20
15. Quality Control.....	20
16. Waste disposal.....	20
17. Incident report.....	21
18. Technical assistance.....	21

Instructions for Use EurobioPlex DermaScreen

1. Introduction

Dermatophytes are common filamentous fungi that cause superficial mycoses by invading the keratin of the nails, hair and skin of humans and animals. The result is a variety of clinical symptoms, including onychomycosis, a very common nail condition. Onychomycosis is mainly caused by *Trichophyton rubrum* and *Trichophyton interdigitale* and rarely by other dermatophytes and fungi.

The identification of dermatophytes in daily practice is important, as an accurate diagnosis enables the selection of the correct (oral) antifungal treatment and avoids (long-term) treatment side effects. A cautious approach is particularly necessary in the elderly, children and pregnant women. In cases where systemic treatment is required, terbinafine is preferable, except during pregnancy and breast-feeding.

Diagnosis of onychomycosis often relies on culture of these dermatophytes and direct microscopic examination. However, identification of specific dermatophyte species is highly complex and dermatophytes are slow-growing and often contaminated by other fungi. The use of real-time polymerase chain reaction (PCR) can overcome these problems by enabling rapid identification of the most common dermatophyte infections in nails, enabling accurate and rapid treatment of onychomycosis.

2. Intended use

The EurobioPlex DermaScreen (EBX-073) is a qualitative, real-time multiplex PCR for the detection of clinically relevant dermatophytes and the detection of *Trichophyton rubrum* and *Trichophyton interdigitale* from extracted DNA. Extracted DNA is the starting material used for detection.

EBX-073 has been validated on the following type of sample:

- Nails

The EurobioPlex DermaScreen test aids in the diagnosis of dermatophyte infections of the nail. The diagnosis must always be combined with other diagnostic and clinical data.

Identification of dermatophytes in daily practice is important as an accurate diagnosis is essential for selecting the correct antifungal treatment and avoid side-effects of long-term treatment. A cautious approach is especially needed in pregnant women for whom several antifungals – topical (miconazole) and oral (terbinafine) – should be avoided during pregnancy and lactation, due to the potential toxicity.

The EurobioPlex DermaScreen test is an *in vitro* medical diagnostic device, it should be used by qualified medical biology analysis laboratory personnel. It should not be recycled after use.

Instructions for Use EurobioPlex DermaScreen

3. Symbols

	Reference
	Lot number
	Limit of storage temperature
	Expiry date
	Content sufficient for « N » reactions
	Manufacturer
	CE marked product
	<i>In Vitro</i> Diagnosis Medical Device
	Instructions for use
	Do not use if packaging is damaged
	Keep away from light
	Attention

Instructions for Use EurobioPlex DermaScreen

4. Principle

The targets detected by the EurobioPlex DermaScreen kit are described in Table 1.

The kit contains an oligomix containing primers and probes specific for the detection of dermatophyte species, which must be used in combination with the polymerase enzyme to enable amplification of the specific DNA. Amplified DNA is visualized by hydrolysis of the fluorescently-labelled probes in the oligomix and provides a specific semi-quantitative signal on a real-time PCR instrument. Detection of the amplified fragments is performed by fluorescence detection in the PCR instrument using the channels indicated in Table 2.

An inhibition control (IC) is provided to monitor variations that may occur during DNA extraction from nail samples using the DermaExtract kit (EBE-077). This allows confirmation that a negative sample is a true negative or the result of ineffective DNA extraction or PCR failure.

A positive control (PC) is also included in the kit and can be used in every PCR run. Signals should be observed in all four detection channels: FAM, HEX, Texas Red and Cy5.

The system is not automated. It is the user's responsibility to use equipment other than those validated (front page of this IFU), and in this case, performance is not guaranteed.

Table 1 : Target detection by fluorophores

Target	Fluorophore	Excitation	Emission
<i>Pan-dermatophytes (ITS)</i>	FAM	495 nm	515 nm
<i>Trichophyton rubrum/ Trichophyton soudanense (ITS)</i>	HEX	535 nm	555 nm
<i>Trichophyton interdigitale/Trichophyton mentagrophytes (ITS)</i>	Texas Red	585 nm	605 nm
Inhibition control	Cy5	650 nm	670 nm

ITS : Internal Transcribed Spacer DNA

Table 2 : Equivalent channels of other thermocyclers

		Channels			
CFX 96/ CFX Opus/LC Pro/ QuantStudio 6*		FAM	HEX	Texas Red	Cy5
QuantStudio 5/7		FAM	VIC	ROX	Cy5
LC480 II	Excitation (nm)	465	533	533	618
	Emission (nm)	510	580	610	660
MIC/Rotor-Gene Q*		Green	Yellow	Orange	Red

*These PCR instruments were not initially validated and require customers validation.

Instructions for Use EurobioPlex DermaScreen

5. Content of the kit

The EurobioPlex DermaScreen real-time PCR kit is ready to use and contains the necessary reagents and enzymes for the detection of these fungi (Table 3).

Table 3 : Content of the kit

Cap colour	Kit contents	25 reactions	50 reactions	100 reactions	200 reactions	400 reactions	600 reactions	Reconstitution
Red	Enzymes	275 µL	550 µL	1100 µL	2x 1100 µL	4x 1100 µL	6x 1100 µL	Ready to use
Transparent	Oligomix*	150 µL	275 µL	550 µL	1100 µL	2x 1100 µL	3x 1100 µL	Ready to use
Yellow	PC (Positive Control)	150 µL	275 µL	550 µL	1100 µL	2x 1100 µL	3x 1100 µL	Ready to use
Blue	NC-H2O (molecular biology grade water)	1000 µL	1000 µL	1000 µL	1000 µL	2x 1000 µL	2x 1000 µL	Ready to use
White	IC-PCR (Armored DNA) (inhibition control)	290 µL	575 µL	1150 µL	2 x 1150 µL	4 x 1150 µL	6x 1150 µL	Ready to use

*Oligomix: contains the primers and probes for the 3 targets and for the inhibition control

6. Storage

All the reagents must be stored between -15°C and -22°C.

All reagents can be used up to the expiration date indicated on the kit label.

Sensitivity may be reduced if kit components undergo multiple freeze/thaw cycles. The kit may be used for a maximum of 5 freeze/thaw cycles.

7. Materials required but not provided

- Biological cabinet/ PCR-/UV-cabinet
- Calibrated Real-time PCR instrument
- Centrifuge for microtubes
- Vortex
- 96-wells plates and appropriate seals / tubes for real-time PCR reaction (MIC PCR)
- Micropipettes
- DNase-free and RNase-free filter tips for micropipettes
- Sterile microtubes
- Gloves (talc-free)

Instructions for Use EurobioPlex DermaScreen

8. Real-time PCR instrument

The EurobioPlex DermaScreen has been developed and validated for use on the following real-time PCR instruments:

- CFX96™ Real Time PCR detection system (Bio-Rad) with analysis on the CFX Maestro versions 1.1 and 2.3 (Bio-Rad) software and with the CFX Manager version 3.1 (Bio-Rad)
- CFX Opus 96 Real Time PCR detection system (Bio-Rad) with analysis on the CFX Maestro versions 1.1 and 2.3 (Bio-Rad) software and with the CFX Manager version 3.1 (Bio-Rad)
- QuantStudio 5/7 Pro (Thermofisher) with analysis on the Design and Analysis version 2.7.0 (Thermofisher) software
- LC480 II (Roche) with analysis on the LightCycler 480 II version 1.5.1 (Roche) software
- LC PRO (Roche) with analysis on the LightCycler Pro version 1.0 (Roche) software
- MIC (Bio Molecular Systems) with analysis on the MIC PCR version 2.12.6 (Bio Molecular Systems) software

9. Cautions and notes

Read these instructions carefully before starting the procedure

- This experiment must be carried out by qualified laboratory employees.
- When conducting dermatophyte screening in laboratories, adherence to both local and national biosafety regulations is critical to ensure safety and compliance.
- Ensure that instruments have been installed, calibrated and maintained in accordance with the manufacturer's recommendations.
- It is the user's responsibility to use equipment other than that validated, in which case performance is not guaranteed.
- Clinical samples should be considered as potentially infectious material and should be prepared under a laminar flow hood.
- Experiments must be carried out in accordance with good laboratory practice.
- Do not use this kit after the expiration date, indicated on the kit label.
- The kit is shipped on dry ice, and kit components must arrive frozen. If one or more components arrive thawed, or if tubes have been damaged in transit, contact Eurobio Scientific.
- Avoid multiple freeze/thaw cycles of reagents, as this may lead to a decrease in test sensitivity.
- Once reagents have thawed, centrifuge tubes briefly before use.

Instructions for Use EurobioPlex DermaScreen

- The use of ice or an ice rack is recommended during preparation of the reactions in the event of long delays due, for example, to a large number of samples to be processed or high temperatures.
- Do not combine reagents of different DermaScreen lot numbers in an experiment.
- It is recommended to define three distinct work areas:
 1. DNA extraction or adding DNA targets to the reaction mixture
 2. Preparation of the reaction mixture
 3. Amplification/detection of the amplified products
- The fluorescent probes in the oligomix are sensitive to light. Prolonged exposure of the oligomix to light should be minimized and limited to the time necessary for preparation of the reagents.
- Pipettes, reagents and other working materials must not circulate between these areas.
- It is recommended to open and handle the positive control separately from the biological samples to be tested, and other kit components, to avoid cross-contamination.
- Wear separate labcoats and gloves (talc-free) in each work area.
- Particular care must be taken to maintain the purity of reagents and reaction mixtures.
- The EBE-077 extraction kit is recommended for DNA extraction.
- It is the user's responsibility to use an alternative extraction method, in which case performance is not guaranteed.
- Use DNase-free filter tips for micropipettes, RNase-free and DNase-free.
- Do not pipette by mouth. Do not eat, drink or smoke in the laboratory.
- Do not open PCR plates or tubes after PCR amplification to avoid aerosols.
- The system is not automated. It is the user's responsibility to use equipment other than those validated, and in this case performance is not guaranteed.
- All samples must be treated as potentially infected with the fungi of interest and local biosafety regulations must be followed.

10. Procedure

10.1 Sample collection

- Transfer small nail parts to sterile 1.5 mL tubes with safety lock before proceeding with DNA extraction. The internal control (IC) from the extraction kit must be used in the DNA-extraction procedure.
- It is the responsibility of the user to collect, transport and store samples before starting with the DNA-extraction of the nails
- The nail samples can be stored at room temperature.

Instructions for Use EurobioPlex DermaScreen

Caution	
	<ul style="list-style-type: none">◇ The user may refer to the recommendations issued by the World Health Organization or the Haute Autorité de Santé for the proper storage of samples.◇ Extracted DNA should be stored at -20°C up to 24 months.◇ The transport of clinical samples is subject to local regulations for the transport of infectious agents.

10.2 DNA Extraction

The DermaXtract kit (EBE-077) is recommended for extraction of nail samples. The exogenous inhibition control must be added at the time of extraction, according to the instructions in the EBE-077 manual. The inhibition control (armored DNA) is a DNA sequence in a protein coat and is used to differentiate true negative samples from negative samples resulting from inefficient DNA extraction or PCR failure. Add 10 µl per nail sample when using the EBE-077 extraction procedure.

It is the user's responsibility to use an alternative extraction method, in which case performance is not guaranteed.

For high-throughput extraction of samples, automation can be used. We recommend using the Qiacube (Qiagen) extraction system which offers identical performances as the DermaXtract kit (tested on 55 *T. rubrum* samples, 10 *T. interdigitale* samples and 14 samples negative for dermatophytes).

10.3 Real-time PCR

General remarks :

- Ensure reagents are completely thawed.
- Briefly vortex the tubes containing the reagents, then centrifuge briefly.
- To check that the qPCR is working correctly, it is necessary to include the positive control, as well as the negative control (water supplied = NC-H₂O) in each run.
- The positive control is a synthetic double-stranded DNA containing the regions targeted by the oligomix. No extraction is required. It should be handled with care to avoid contamination.
- Ensure that the correct (skirted/non-skirted) qPCR plates are used for each type of instrument. White plates must be used in the CFX instruments (CFX 96, CFX Opus).
- An overview of the PCR preparation is provided in the figure below

Instructions for Use EurobioPlex DermaScreen

Procedure overview :

1 – PREPARATION OF THE MASTERMIX (AREA 1)

Number of reactions	N+3
Enzyme mix	(N+3) x 10 µL
Oligomix	(N+3) x 5 µL
Total volume Mastermix	(N+3) x 15 µL

2 - PREPARATION OF THE REACTIONS (AREA 2)

Sample

15 µL Mastermix
+
5 µL DNA sample

Positive Control

15 µL Mastermix
+
5 µL PC

Negative Control

15 µL Mastermix
+
5 µL molecular biology graded water (NC-H₂O)

- Seal the PCR plate with a suitable adhesive film before transferring it to the real-time PCR instrument.
- Briefly centrifuge the plate before transferring it to the real-time PCR instrument.

3 – REAL TIME PCR INSTRUMENT (AREA 3)

Program	Temperature	Duration	Cycle(s)	
Denaturation	95°C	2 mins	1	-
Amplification	95°C	10 secs	44	-
	62°C	30 secs		Fluorescence acquisition*

*Detection in all channels must be selected

Instructions for Use EurobioPlex DermaScreen

11. Validation of the PCR run

11.1 Validation of the PCR run on CFX 96 and CFX Opus (Bio-Rad)

Post-acquisition data analysis on a CFX96 or CFX Opus (Bio-Rad) PCR instrument must be carried out using CFX Manager software version 3.1 (Bio-Rad) or CFX Maestro software version 1.1 or version 2.3 (Bio-Rad).

Run analysis must be performed on the data file with the .pcrd suffix. On the CFX 96, this is the format automatically saved at the end of the run. On the CFX Opus, the .zpcr file automatically generated at the end of the run must first be opened using the CFX Manager or CFX Maestro software, then saved in .pcrd format for analysis.

The “**Baseline Subtracted Curve Fit**” option must be applied in the “Settings” tab: click on the “Settings” tab, then on “Baseline Setting” and on “Baseline Subtracted Curve Fit”.

It is recommended to perform the analysis on a linear scale.

Ct values are determined using the regression Ct determination mode: click on the “Settings” tab, then on “Cq determination mode” and on “**Regression**”. This enables automatic determination of Ct values.

For the run to be valid, the Ct values obtained for the positive control must correspond to the values shown in Table 4. **The negative control must be negative in all detection channels.**

Table 4: Validation of the run on CFX96 and CFX Opus

Positive Control	
FAM	Ct ≤ 32
HEX	Ct ≤ 34
Texas Red	Ct ≤ 33
Cy5	Ct ≤ 31


11.2 Validation of the PCR run on QuantStudio 5/QuantStudio 7 (ThermoFisher)

Post-acquisition data analysis should be performed on the data file with the .eds suffix, which is automatically generated at the end of the run, using the Design and Analysis software version 2.7.0.

The analysis must be performed on each channel independently by linear scaling on the amplification curves: in the “Quality Check” tab, select the “Amplification plot” option from the drop-down menu. Select the settings icon at the top right of the amplification curves to select “Linear” for the **Y scale** parameter.

Instructions for Use EurobioPlex DermaScreen

The threshold for each fluorophore should then be set at the inflection point of the curve. The inflection point is the point at which a change in the direction of curvature occurs. It can be defined by tracing a straight line on the linear phase: the inflection point is the point where this straight line crosses the curve. An example is shown in the following figure.

For the run to be valid, the Ct values obtained for the positive control must correspond to the values shown in the following table (Table 5). **The negative control must be negative in all detection channels.**

Table 5: Validation of the run on QuantStudio5/QuantStudio7

Positive control	
FAM	Ct ≤ 33
HEX	Ct ≤ 36
Texas Red	Ct ≤ 34
Cy5	Ct ≤ 33

11.3 Validation of the PCR run on LC480 II (Roche)

On the LC480II, a cross-talk can be observed between the 533-580/533-610 channels and between the 533-610/618-660 channels. To prevent this from happening, a 4-channel (FAM, HEX, Texas Red and Cy5 - or equivalent) colour compensation kit can be used.

Post-acquisition data analysis should be performed on the data file with the .ixo suffix, which is automatically generated at the end of the run, using LightCycler 480 software version 1.5.1.

To generate the analysis, select the “Analysis” tab. Depending on whether a colour compensation kit has been used or not, either select the “Abs/2nd derivative max” option (without colour compensation) or the “Color Comp” option (with colour compensation). Each channel must then be analysed independently and Ct values are calculated automatically by clicking on the “Calculate” button.

Instructions for Use EurobioPlex DermaScreen

For the run to be valid, the Ct values obtained for the positive control must correspond to the values shown in the following table (Table 6). **The negative control must be negative in all detection channels.**

Table 6: Validation of the run on LC480 II

Positive control	
FAM	Ct ≤ 33
HEX	Ct ≤ 34
Texas Red	Ct ≤ 34
Cy5	Ct ≤ 33

Whether a colour compensation kit has been used or not does not change the analytical sensitivity in each of the specific channels and the Ct values remain comparable. However, data interpretation depends on whether a colour compensation kit has been used or not and details regarding data interpretation can be found in section 12. Data analysis and interpretation.

11.4 Validation of the PCR run on LC Pro (Roche)

Data analysis on the LC PRO (Roche) is carried out using the LightCycler® PRO (v1.0) software. Post-acquisition data can be imported into the software by creating a new project and then selecting the runtime file stored on USB or on the server. In the plate configuration section, sample IDs and sample roles can be indicated, including unknowns, positive matrix control (PTC) and no matrix control (NTC). Note that only one PTC and one NTC can be included per analysis.

Perform an analysis on a single plate, select qualitative analysis and define the channels (FAM, HEX, Texas Red and Cy5). In the "Define parameters" section, target information can be provided and color compensation can be selected for all detection channels. Inhibition control can be selected for the Cy5 detection channel. The 'Positive target confirmation without internal control' option can be selected for the FAM, HEX and Texas Red detection channels. In the 'Edit sample information' section, you can delete unused wells.

For the run to be valid, the Ct values obtained for the positive control must correspond to the values shown in the following table (Table 7). **The negative control must be negative in all detection channels.**

Table 7: Validation of the run on LC Pro

Positive Control	
FAM	Ct ≤ 34
HEX	Ct ≤ 36
Texas Red	Ct ≤ 34


Instructions for Use EurobioPlex DermaScreen

Cy5	Ct ≤ 34
-----	---------

11.5 Validation of the PCR run on the MIC (Biomolecular systems)

Post-acquisition data analysis should be performed on the data file with the .micrun suffix, which is automatically generated at the end of the run, using micpcr software version 2.12.6.

The analysis must be performed independently on each channel. In the “Analysis” section, select the following parameters:

For the run to be valid, the Ct values obtained for the positive control must correspond to the values shown in the following table (Table 8). **The negative control must be negative in all detection channels.**

Table 8: Validation of the run on the MIC

Positive Control	
FAM	Ct ≤ 32
HEX	Ct ≤ 34
Texas Red	Ct ≤ 33
Cy5	Ct ≤ 32

Instructions for Use EurobioPlex DermaScreen

12. Data analysis and interpretation

Data analysis and interpretation is based on the signal measurements recorded for each channel of the DermaScreen kit, as detailed in Table 9, which lists the corresponding dermatophyte species and controls.

Table 9: Data interpretation


Data interpretation EBX-073 DermaScreen (CFX96, CFX Opus, QuantStudio5/7, LC480 with colour compensation, LC PRO, MIC)				
FAM	HEX	TexasRed	Cy5	Result
Pos	Neg	Neg	Pos/Neg	Positive for dermatophytes species (1)
Pos*	Neg	Pos	Pos/Neg	<i>Trichophyton interdigitale</i> / <i>Trichophyton mentagrophytes</i> (2)
Pos*	Pos	Neg	Pos/Neg	<i>Trichophyton rubrum</i> / <i>Trichophyton soudanense</i> (3)
Neg	Neg	Neg	Pos	Negative for dermatophytes species (5)
Neg	Neg	Neg	Neg	Invalid
Data interpretation EBX-073 DermaScreen (LC480 II without colour compensation)				
FAM	HEX	TexasRed	Cy5	Result
Pos	Neg	Neg	Pos/Neg	Positive for dermatophytes species (1)
Pos*	Neg	Pos	Pos/Neg	<i>Trichophyton interdigitale</i> / <i>Trichophyton mentagrophytes</i> (2)
Pos*	Pos	Pos	Pos/Neg	<i>Trichophyton rubrum</i> / <i>Trichophyton soudanense</i> (3) (4)
Neg	Neg	Neg	Pos	Negative for dermatophytes species (5)
Neg	Neg	Neg	Neg	Invalid

*A negative FAM/panderma signal with a HEX/*T. rubrum* or Texas Red/*T.interdigitale* signal can occur but these are very low dermatophyte concentrations.

- (1) The FAM probe detects all dermatophyte species, including all clinically relevant *Trichophyton* spp, *Microsporum* spp, *Epidermophyton* spp, *Nannizzia* spp and *Paraphyton* spp. If only the FAM probe is positive, it will be reported as “dermatophyte positive” but not as *T. interdigitale* or *T. rubrum*.
- (2) The Texas Red probe detects also *T. mentagrophytes* (all subtypes), *T. indotineae*, *T. schoenleinii* and *T. quinckeanum*. Last two species result in amplification with limited fluorescence. The clinical relevance in nails of other species compared to *T. interdigitale* species has not been established, particularly in Europe.
- (3) The HEX probe detects also *T. soudanense*, this dermatophyte is however observed predominantly in hair samples.
- (4) Spectral cross-talk is observed between the detection channels on the LC480 II PCR instrument. As a result, if no colour compensation kit is used, signals from the 533-580 detection channel (*T. rubrum*) will also be observed in the 533-610 detection channel (*T. interdigitale*). Also signals in 533-610 (*T. interdigitale*) can be observed in the 618-660 detection channel (Inhibition control). When signals are present in the 533-580 as well as the 533-610 detection channel, this must be interpreted as *T. rubrum* only (considering a double infection as very rare). Important note: always compare amplification curves of samples to the PC signal for correct data analysis.
- (5) Inhibition control signals are only visible when the inhibition control is used during DNA extraction. The inhibition control is outperformed by high concentrations of dermatophyte DNA, resulting in a limited or absent Cy5 signal. The result is still valid.

Instructions for Use EurobioPlex DermaScreen

An example of a CFX96 result is shown in the figure below. The sample is positive for *T. rubrum*. The green curve indicates a dermatophyte signal (FAM), the blue curve is a *T. rubrum* signal (HEX), the purple curve is a suppressed inhibition control signal (Cy5) and the *T. interdigitale* signal (TexasRed) is negative.

Limits of use and interpretation:

- Interpretation of results must take into account the possibility of false negatives and false positives.

False negatives may be due to :

- Inadequate sample collection or storage,
- Unsuitable extraction conditions or use of non-validated PCR instruments,
- Performance that does not comply with all the instructions in this manual.

False positives may be due to :

- Contamination due to mishandling of high-positive samples, positive control or PCR amplification products,
- Failure to follow the procedures described in this manual, in particular to avoid sources of contamination.

- All results must be interpreted by medical personnel in the context of the patient's clinical situation, history and symptoms.
- This test does not exclude the presence of pathogens other than those detected by this kit.
- Cross-reaction may occur (see Table 14), with other rare dermatophytes species, or almost exclusively present in other types of biological samples like hairs or skin.

Instructions for Use EurobioPlex DermaScreen

- A negative result of this test does not absolutely exclude a possible infection with the fungi detected in this kit.

13. Performances analysis

Analytical sensitivity

A Probit analysis to determine the 95% cut-off per target was carried out using a dilution series of synthetic target DNA. The results for each of the kit's specific targets (*T. rubrum* and *T. interdigitale*) and per instrument are presented in the tables below:

Table 10: Analytical sensitivity (95% cut-off) in copies/μL

Analytical sensitivity (95% cut-off) in copies/μL						
Specific targets	CFX 96 (Bio-Rad)	CFX Opus (Bio-Rad)	QuantStudio 7 (ThermoFisher)	LC 480 II (Roche)	LC Pro (Roche)	MIC (Biomolecular systems)
<i>T. rubrum</i> (HEX)	0,97	1,19	0,84	0.51	0,88	1.53
<i>T. interdigitale</i> (Texas Red)	1	8,83	1,32	13	1,03	12.5

Signal variability on specific HEX and Texas Red channels

- Intra-experimental variability : CV (Coefficient of Variation) intra-lot

Table 11: CV intra-lot

CV intra-lot (%)						
Specific targets	CFX 96 (Bio-Rad)	CFX Opus (Bio-Rad)	QuantStudio 7 (ThermoFisher)	LC 480 II (Roche)	LC Pro (Roche)	MIC (Biomolecular systems)
<i>T. rubrum</i> (HEX)	3,46	3,23	1,65	1,54	0,91	1,05
<i>T. interdigitale</i> (Texas Red)	5,69	2,11	0,76	3,01	1,29	0,93

- Inter-lots variability : CV inter-lot

Table 12: CV inter-lot

CV inter-lot (%)						
Specific targets	CFX 96 (Bio-Rad)	CFX Opus (Bio-Rad)	QuantStudio 7 (ThermoFisher)	LC 480 II (Roche)	LC Pro (Roche)	MIC (Biomolecular systems)
<i>T. rubrum</i> (HEX)	3,04	2,29	1,76	1,95	1,08	2,2
<i>T. interdigitale</i> (Texas Red)	5,75	1,44	1,42	2,14	1,16	2,74

Instructions for Use EurobioPlex DermaScreen

Clinical sensitivity and specificity

The clinical specificity and sensitivity was determined based on the testing of clinical nail samples in a routine setting. All samples were first extracted with the EurobioPlex DermaXtract kit (EBE-077) before analysed with the EurobioPlex DermaScreen real-time PCR using the CFX96 and the LC480.

Table 13: Clinical sensitivity and specificity

	Molecular reference methods (n)	DermaScreen (n)	Sensitivity (%) (95% CI*)	Specificity (%) (95% CI*)
<i>T. rubrum</i>	121	120	99,2 (98,224-100,208)	
<i>T. interdigitale</i>	26	25	96,2 (95,312-97,235)	
Negative	68	62		91,2 (90,356-92,179)

*CI: confidence interval

Analytical specificity

The analytical specificity was determined by testing the confirmed genomic DNA of 17 dermatophyte strains but also fungal and bacterial strains which can be present on nails or in the environment. All strains were tested on the CFX96 and listed in the table below. Additionally, all primer and probe sequences were compared to sequences published in Genbank to check for possible homologies which might result in unwanted signals.

Table 14: Analytical specificity

Species	Result DermaScreen		
	Panderna (FAM)	<i>T. rubrum</i> (HEX)	<i>T. interdigitale</i> (TexasRed)
<i>T. mentagrophytes</i> Thailand*	Positive	ND	Positive
<i>T. indotineae</i> *	Positive	ND	Positive
<i>T. quinckeanum</i> **	Positive	ND	Positive
<i>T. schoenleinii</i> **	Positive	ND	Positive
<i>T. interdigitale</i>	Positive	ND	Positive
<i>T. mentagrophytes</i> *	Positive	ND	Positive
<i>T. rubrum</i>	Positive	Positive	ND
<i>T. tonsurans</i>	Positive	ND	ND
<i>T. violaceum</i>	Positive	ND	ND
<i>T. soudanense</i> ***	Positive	Positive	ND
<i>T. verrucosum</i>	Positive	ND	ND
<i>P. cookei</i>	Positive	ND	ND
<i>E. floccosum</i>	Positive	ND	ND

Instructions for Use EurobioPlex DermaScreen

<i>M. canis</i>	Positive	ND	ND
<i>M. audouinii</i>	Positive	ND	ND
<i>T. benhamiae</i>	Positive	ND	ND
<i>N. gypsea</i>	Positive	ND	ND
<i>C. albicans</i>	ND	ND	ND
<i>C. parapsilosis</i>	ND	ND	ND
<i>M. guilliermondii</i>	ND	ND	ND
<i>F. solani</i>	ND	ND	ND
<i>A. fumigatus</i>	ND	ND	ND
<i>R. oryzae</i>	ND	ND	ND
<i>S. brevicaulis</i>	ND	ND	ND
<i>T. marneffeii</i>	ND	ND	ND
<i>S. aureus</i>	ND	ND	ND
<i>S. epidermidis</i>	ND	ND	ND
Humane DNA	ND	ND	ND

**T. mentagrophytes* species including the subtype originating from Thailand but also the *T. indotinaea* result in a signal in the TexasRed detection channel. These closely related species are not considered typical nail pathogens

***T. quinckeanum* and *T. schoenleinii* are also detected in TexasRed with the *T. interdigitale* probe but the signals are much lower. Moreover, note that both species are almost never nail pathogens and thus not clinically relevant to EBX-073 which is validated only on nails.

****T. soudanense* is also detected in HEX with the *T. rubrum* probe but this is rarely a nail pathogen compared to *T. rubrum*, particularly in Europe.

Interfering substances

In total four interfering substances likely to be found in the samples were tested for their possible effect on kit efficacy: clear nail varnish, red nail varnish, miconazolnitate and terbinafine, both of which are commonly found in anti-fungal creams.

None of these 4 substances interfered with the extraction reaction or with the PCR.

Instructions for Use EurobioPlex DermaScreen

14. Bibliography

1. Irene Weitzman, Richard C. Summerbell, The Dermatophytes, *Clinical Microbiology Reviews*, Apr. 1995, p. 240-259 Vol. 8, No. 2
2. Dyanne P. Westerber, Michael J. Voyack, Onychomycosis: Current Trends in Diagnosis and Treatment, *Am Fam Physician*. 2013;88(11):762-770
3. Yvonne Gräser, Ditte M L Saunte, A Hundred Years of Diagnosing Superficial Fungal Infections: Where Do We Come From, Where Are We Now and Where Would We Like To Go?, *Acta Derm Venereol* 2020 Apr 20;100(9)
4. Pietro Nenoff, Dieter Reinel, Peter Mayser, Dietrich Abeck, Guntram Bezdold, Philipp P. Bosshard, Jochen Brasch, Georg Daeschlein, Isaak Effendy, Gabriele Ginter-Hanselmayer, Yvonne Gräser, Gudrun Hamm, Ulrich Hengge, Uta-Christina Hipler, Peter Höger, Alexandra Kargl, Annette Kolb-Mäurer, Constanze Krüger, Bartosz Malisiewicz, Johannes Mayer, Hagen Ott, Uwe Paasch, Martin Schaller, Silke Uhrlaß, Miriam Zidane *Guideline onychomycosis, JDDG: Journal der Deutschen Dermatologischen Gesellschaft: May, 2023, Volume 21, Issue 6, Pages: 571-710*
5. A.K. Gupta, M. Paquet, F. Simpson, A. Tavakkol, Terbinafine in the treatment of dermatophyte toenail onychomycosis: a meta-analysis of efficacy for continuous and intermittent regimens, *Volume 27, Issue 3, March 2013, Pages 267-272*
6. Subuhi Kaul, Savita Yadav, Sunil Dogra, Treatment of Dermatophytosis in Elderly, Children, and pregnant women, *Indian Dermatol Online J* 2017 ;8 :310-8
7. Pilmis B, Julien V, Sobel J, Lecuit M, Lotholary O, Charlier C, Antifungal drugs during pregnancy: An updated review, *J Antimicrob Chemother* 2015;70:14-22
8. Murase JE, Heller MM, Butler DC, Safety of dermatologic medications in pregnancy and lactation: Part I. Pregnancy, *J Am Acad Dermatol* 2014;70:401-159.
9. Patel VM, Schwartz RA, Lambert WC, Topical Antiviral and Antifungal Medications in Pregnancy: A Review of Safety Profiles, *J Eur Acad Dermatol Venereol* 2017;27
10. Christina Forsberg*, Lina Pettersson, Lina Boiso, The effect of freezing, thawing and long-term storage on forensic DNA extracts. *Forensic Science International: Genetics Supplement Series* 8(2022)77–78

15. Quality Control

In accordance with Eurobio's ISO EN 13485-2016:A21 certified quality management system, each batch of EurobioPlex DermaScreen is tested according to predefined specifications to guarantee consistent product quality.

16. Waste disposal

Dispose of all waste in accordance with DASRI legislation or to the local applicable regulations.

Instructions for Use EurobioPlex DermaScreen

17. Incident report

Any serious incident involving the reagents is reported to Eurobio Scientific and to the competent authority of the Member State in which the user and/or patient is established.

18. Technical assistance

For assistance with our products, please contact our technical support.

Eurobio Scientific's customer service can be contacted by e-mail at adv@eurobio-scientific.com or by phone at +33 (0)1.69.79.64.80

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14. Bibliographie

1. Irene Weitzman, Richard C. Summerbell, The Dermatophytes, Clinical Microbiology Reviews, Apr. 1995, p. 240-259 Vol. 8, No. 2
2. Dyanne P. Westerber, Michael J. Voyack, Onychomycosis: Current Trends in Diagnosis and Treatment, Am Fam Physician. 2013;88(11):762-770
3. Yvonne Gräser, Ditte M L Saunte, A Hundred Years of Diagnosing Superficial Fungal Infections: Where Do We Come From, Where Are We Now and Where Would We Like To Go?, Acta Derm Venereol 2020 Apr 20;100(9)
4. Pietro Nenoff, Dieter Reinel, Peter Mayser, Dietrich Abeck, Guntram Bezold, Philipp P. Bosshard, Jochen Brasch, Georg Daeschlein, Isaak Effendy, Gabriele Ginter-Hanselmayer, Yvonne Gräser, Gudrun Hamm, Ulrich Hengge, Uta-Christina Hipler, Peter Höger, Alexandra Kargl, Annette Kolb-Mäurer, Constanze Krüger, Bartosz Malisiewicz, Johannes Mayer, Hagen Ott, Uwe Paasch, Martin Schaller, Silke Uhrlaß, Miriam Zidane Guideline onychomycosis, JDDG: Journal der Deutschen Dermatologischen Gesellschaft: May, 2023, Volume 21, Issue 6, Pages: 571-710
5. A.K. Gupta, M. Paquet, F. Simpson, A. Tavakkol, Terbinafine in the treatment of dermatophyte toenail onychomycosis: a meta-analysis of efficacy for continuous and intermittent regimens, Volume 27, Issue 3, March 2013, Pages 267-272
6. Subuhi Kaul, Savita Yadav, Sunil Dogra, Treatment of Dermatophytosis in Elderly, Children, and pregnant women, Indian Dermatol Online J 2017 ;8 :310-8
7. Pilmis B, Julien V, Sobel J, Lecuit M, Lotholary O, Charlier C, Antifungal drugs during pregnancy: An updated review, J Antimicrob Chemother 2015;70:14-22
8. Murase JE, Heller MM, Butler DC, Safety of dermatologic medications in pregnancy and lactation: Part I. Pregnancy, J Am Acad Dermatol 2014;70:401-159.
9. Patel VM, Schwartz RA, Lambert WC, Topical Antiviral and Antifungal Medications in Pregnancy: A Review of Safety Profiles, J Eur Acad Dermatol Venereol 2017;27
10. Christina Forsberg*, Lina Pettersson, Lina Boiso, The effect of freezing, thawing and long-term storage on forensic DNA extracts. Forensic Science International: Genetics Supplement Series 8(2022)77-78

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