TQS
TETANOS QUICK STICK

References:
TE-15U (15 x TE-1U)
TE-50U (50 x TE-1U)

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11. BIBLIOGRAPHIE

1. Elkharrat D. et al, Impact of guidelines to alter antitetanus prophylaxis practices and reduce costs in the emergency department. Am J Ther. 1999; 6, (4): 203-9.
2. Thiebaut A. et al. Intérêt clinique et économique d'un test rapide de mi-ose en évidence de l'immunoprotection antitétanique. J. Pharm. Clin. 2003; 22, (1): 31-35.
3. Colombet I. et al. Diagnosis of Tetanus Immunization status: Multicenter Assessment of a Rapid Biological Test. Clin. & Diagn. Lab. Immunol. 2005; 12, (9): 1057-1062.
4. Stubbe M. et al. Improving tetanus prophylaxis in the emergency department: a prospective, double-blind cost-effectiveness study. Belgian Society of Emergency and Disaster Medicine. Emerg Med J. 2007; 24, (9): 648-53.
5. Galazka A.M. Les bases immunologiques de la vaccination: 3 Le tetanus. Organisation Mondiale de la Santé, Genève. 1993.
6. Tierney R. et al. The first international standard for anti-tetanus immunoglobulin, human: pharmaceutical evaluation and international collaborative study. Biologicals. Journal of the International Association of Biological Standardization, 1993; 21:67-75.
7. Borrow R. et al., WHO, The immunological basis for immunization series. Module 3 Tetanus updated. 2006.
8. Elkharrat D. et al. Evaluation of a bedside immunotest to predict individual anti-tetanus seroprotection: a prospective concordance study of 1018 adults in an emergency department. Emerg Med J. 2010; 27, (1): 36-42.

eurobio
SCIENTIFIC

7 Avenue de Scandinavie
ZA de Courtaboeuf
91 940 LES ULIS

Tél : + 33 (0)1 69 79 64 80
info@eurobio-scientific.com

ENGLISH

Immunochromatographic test for the detection of tetanus antitoxin antibodies in whole blood, plasma and serum

TE-15U (15 x TE-1U)
TE-50U (50 x TE-1U)

In vitro Diagnostic (Professional use only)

1. INTRODUCTION

1.1 Clinical context

Tetanus is an acute, serious and often lethal infection caused by a toxin produced by *Clostridium tetani*, the spores of which can often be found in dirt and animal faeces. Tetanus develops following contamination of a cut, and the germination of these spores produces a toxin that causes a neuromuscular disease. The most common clinical signs of tetanus are increased muscle contractions and generalised muscle spasms.

Tetanus is prevented by immunisation with tetanus toxoid, which gives significant and long-lasting immunity, as 91% of people are still immunised 10 years after being vaccinated. Tetanus has not been eradicated in industrialised countries, although it is rare, despite a significant vaccination programme.

For patients referred to hospital emergency departments for wounds, treatment is based on cleaning and disinfecting cuts and immunisation with a tetanus toxoid injection and/or tetanus antitoxin immunoglobulins. During the consultation, the immunisation status of the patient with regards to tetanus is best determined by the patient's vaccine record or, alternatively, and often in practice, by interviewing the patient. However, data regarding patients in France show that the vaccination record is not available in 90-95% of cases and various studies have shown that patients misunderstand or are unaware of their actual vaccination status¹. The difficulty in accessing this data leads a majority of emergency departments to resort to rapid immunoassays in order to determine the patient's tetanus vaccine status and choose the correct prophylaxis^{2,3,4}.

1.2 Intended use

Various studies have shown that a level of circulating antitoxin greater than 0.01 IU/mL of serum is protective. The WHO, in its recommendations for treatment of tetanus, specifies a factor of 10 for defining the positivity threshold. Patients with a level of antibodies greater than 0.1 IU/mL of anti-tetanus antibody serum are considered to be immunoprotected^{5,6,7}.

TQS TETANOS QUICK STICK detects antibody levels greater than 0.1 IU/mL in whole blood, plasma and serum.

A positive result on serum or plasma implies that the sample contains more than 0.1 IU/mL of anti-tetanus antibodies in the serum. A positive result on whole blood implies that the sample contains more than 0.1 IU/mL of anti-tetanus antibodies in whole blood, or 0.2 IU/mL in serum from the same donor.

The test should be used only by healthcare professionals in medical biological laboratories or in other types of structures, notably in emergency departments.

2. PRINCIPLE OF THE TEST

The rapid test TQS TETANOS QUICK STICK is a rapid immunochromatographic assay for the qualitative detection of tetanus antitoxin antibodies in whole blood, plasma and serum. TQS TETANOS QUICK STICK detects antibody levels greater than 0.1 IU/mL in whole blood, plasma and serum. The test result is obtained in 10 minutes.

The stick comprises:

- a sample and diluent dispensing well
- a reading window comprising 2 zones: a T zone where the test result appears and a C zone where the control line appears.

A sample of whole blood, plasma or serum is deposited into the stick dispensing well. The diluent is then added to the same well.

The diluent migrates, dragging along the gold conjugate complexed to the anti-tetanus immunoglobulins present in the sample.

These complexes react with the immobilised toxoid to form a coloured line in the T window.

If the amount of tetanus immunoglobulins present in the specimen is lower than the detection threshold of 0.1 IU/mL, no coloured line will appear in the T zone.

The excess gold conjugate binds to a control reagent immobilised in the C zone to form a pink line indicating that the test has been carried out correctly.

3. KIT CONTENTS

TE-15U (15 x TE-1U):

TE-1U:

The bag of the reference TE-1U contains everything needed to perform 1 test:

- 1 single-use stick and 1 pipette (20µL), individually packaged in an aluminium pouch with a desiccant bag
- 1 vial of diluent (0,5mL)
- 1 lancet
- 1 alcohol pad
- 1 summary of the test procedure
- 2 patient labels for result archiving and traceability

Each reference of TE-15U contains in addition:

- 1 package insert

TE-50U (50 x TE-1U):

TE-1U:

The bag of the reference TE-1U contains everything needed to perform 1 test:

- 1 single-use stick and 1 pipette (20µL), individually packaged in an aluminium pouch with a desiccant bag
- 1 vial of diluent (0,5mL)
- 1 lancet
- 1 alcohol pad
- 1 summary of the test procedure
- 2 patient labels for result archiving and traceability

Each reference of TE-50U contains in addition:

- 1 package insert

4. PRECAUTIONS FOR USE

1. For in vitro diagnostic use.
2. Healthcare professionals not working in a laboratory should be trained, especially in the use of the pipette, and supervised by a competent person in accordance with the institution's practices.
3. Do not use the test if the stick or aluminium pouch is damaged or if the pouch does not contain any desiccant bag.
4. Do not reuse a stick.
5. Handle all specimens as if they were potentially infectious. Dispose of all waste products in accordance with the legislation in force.
6. Wear a laboratory coat and disposable gloves when handling specimens.
7. Do not eat, drink or smoke in the room where specimens and kit reagents are handled.
8. The diluent contains 0.05% sodium azide as an antimicrobial agent. According to the REACH regulation, this solution is not considered dangerous at this concentration. However, it is recommended to avoid contact with the hands, eyes and nose when handling samples and to avoid its release into the environment.
9. Do not use the test if when opening the pouch the strip shows a pink color in the read window.

5. STORAGE AND STABILITY

The TQS TETANOS QUICK STICK kit should be stored between 4°C and 30°C. Do not freeze the kit. The kit can be used until the expiry date indicated on the packaging.

6. SPECIMEN COLLECTION

1. TQS TETANUS QUICK STICK can be performed on non-haemolysed blood, serum, plasma (citrate or EDTA) or re-calcified plasma.

2. Specimens should be tested immediately. If this is not possible, they can be refrigerated between 2°C and 8°C for a maximum of 3 days.

If testing within 3 days is not possible, serum and plasma specimens should be frozen at -20°C and brought to ambient temperature before the test is performed.

3. Specimens containing a precipitate may give inconsistent test results.

7. TEST PROCEDURE (see explanatory diagrams)

1. Open the stick's protective sachet by tearing along the split and place the stick on a level surface.
2. Label the test with the patient name and test number.
3. Collect the sample and deposit it into the dispensing well.

*** Whole blood**

Clean the patient's finger with soapy water or disinfectant. Use a single-use lancet to prick the patient's finger (see explanatory diagram), collect and dispense into the well a 20 µL sample of blood using the plastic pipette provided. .

*** Serum or plasma**

Collect and dispense 20 µL of the sample using the pipette supplied or a calibrated pipette into the well.

4. Add 3 drops of diluent immediately, holding the dropper vertically above the well and avoiding touching the bottom of the well with the tip of the dropper.

NB: if no liquid flux appears in the reading zone after 2 minutes, an additional drop of diluent can be added.

8. READING RESULTS

Read the results **no later than** 10 minutes after adding the 3 drops of diluent. A positive result can show up earlier than 10 minutes if the test is validated in the C zone. A delayed reading could give erroneous results.

1 Test validation

A pink line appears in the C zone, indicating that the test has been carried out correctly. The absence of a coloured line in the C zone means the test is invalid and it must be done again with a new stick.

2 Reading the test

- If no coloured line appears in the T zone: the test is negative.
- If a pink line is visible in the T zone, no matter how intense: the test is positive.
- Any positive reaction indicates that the level of antibodies in the specimen is greater than the test detection threshold of 0.1 IU/mL.

Note: In case of doubt regarding the presence or not of a coloured line in the T zone (e.g., for a patient whose antibodies level are close to the threshold), the test result is negative.

In no case should there be a comparison between the colour intensities of the lines T and C.

3 Automatic reading

The test can be read by an automatic reader validated to provide a positive or negative result depending on the appearance of the test line and the control line. The ZReader (ref: H-KP-001, Zentech) is validated for this use in incubation mode for 10 minutes or in kinetic mode (positive result can be displayed in fewer than 10 minutes if the test is validated in the C zone).

9. QUALITY CONTROL

The use of control reagents is recommended to ensure proper kit performance. The control kit (ref. TE2C / S-GZ-CZ), which includes positive and negative serum, may be used for this purpose.

10. TEST PERFORMANCE

A-Validation of the detection threshold

The detection threshold has been validated for whole blood, plasma and serum from 7 donors and on 3 different batches. Two-fold dilution steps of the WHO standard (NIBSC TE-3) in whole blood, plasma and serum matrices were carried out, encompassing the detection threshold of 0.1 IU/mL.

For each of the 7 donors, the IgG tetanus antitoxin titers from 3 supplemented matrices with the WHO standard were defined using a commercial ELISA test.

This study was able to confirm that the test detects antibody levels greater than 0.1 IU/mL in whole blood, plasma and serum.

B-Clinical validation for whole blood, plasma and serum

1. Internal studies

TQS TETANOS QUICK STICK was assessed using 865 samples from 289 donors (289 whole-blood samples, 287 plasma samples and 289 serum samples), and was compared with a commercial tetanus antitoxin ELISA kit. This study, carried out in 2015, enabled the following performance calculation by excluding discordant specimens around the threshold:

	Whole blood	Plasma	Serum
Sensitivity (PPV)	97.0% n=169 (100%)	97.2% n=216 (100%)	98.6% n=213 (100%)
Specificity (NPV)	100% n=75 (93.8%)	100% n=43 (87.8%)	100% n=44 (93.6%)

This performance was confirmed on a separate batch of sticks.

This study confirms previously obtained results:

	Whole blood (2008)	Serum (2005)
Sensitivity (PPV)	98.1% n=52 (100%)	96.4% n=160 (100%)
Specificity (NPV)	100% n=25 (96.2%)	100% n=34 (85%)

2. External studies

TQS TETANOS QUICK STICK has been the subject of several external studies in real-life conditions, which assessed its performance along with its clinical and medico-economic value in an emergency departments^{3,8}.

C-Interferences

The following blood components have not been shown to significantly interfere with the results of the test up to these concentrations: albumin (100 mg/mL), human haemoglobin (20 mg/mL), bilirubin (0.05 mg/mL) and cholesterol (2.75 mg/mL).

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Version tracking:

The areas highlighted enable to identify words, lines or paragraphs changed since the previous version

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(EN): Lancet manufactured by Sarstedt (Medical Device CE European Directive 93/42/EEC)

- Twist off the protective cap.
- Press the Safety-Lancet against the chosen and cleansed puncture site. Press the trigger button.
- Dispose of the Safety -Lancet immediately in a suitable sharps container. Gently massage the finger in the direction of the puncture site to obtain the required blood volume.

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
(EN) Do not use if cap has been previously removed

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




(en)

1. The finger should be cleaned properly. After the puncture of the finger with a single use safety lancet and press the finger in order to obtain a big drop of blood.
2. Hold the pipette bulb gently in a horizontal position and place the tip of the pipette into the blood. Do not squeeze the pipette bulb while sampling! Maintain this position until the flow of sample into the pipette has stopped ; capillary action will automatically draw the sample to the fill line and stop.
3. To expel the sample, hold the pipette vertically. When dispensing the blood sample, let the end of the pipette touch the sample target and squeeze the bulb to dispense the blood sample at the center of the sample well. If a sample won't expel, hold the pipette vertically and position a finger over the vent hole and squeeze the bulb.
4. Add 3 drops of diluent within 10 seconds, holding the flask vertically above the sample well and avoiding touching it with the dropper.
5. Read results at 10 minutes :
 - a) Negative result : one coloured band appears in control window (C)
 - b) Positive result : two coloured bands appear in test (T) and control windows (C)

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