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Medium for the fluidification of biological samples

REF

D0DIG00-AI

CE IVD

Version 6.01 from July 2022

Instructions for use

Available in www.eurobio-scientific.com

Instructions for use Digest-EUR ®

Table of contents

Table of contents.....	10
1. General information	11
2. Intended use.....	11
3. Symbols.....	12
4. Packaging.....	13
5. Principle of the test.....	13
6. Composition	13
7. Storage/stability	13
8. Material required but not provided	13
9. Cautions and recommendation for users.....	14
10. Sample collection	14
11. Techniques	14
12. Waste disposal	15
13. Incident Report	15
14. Technical assistance	15
15. Bibliography.....	16

Instructions for use Digest-EUR ®

1. General information

This reagent allows the fluidization of mucosal biological samples for culture for the detection and identification of pathogenic bacterial strains.

These products can also be used for research applications, and under appropriate circumstances, as primary components in subsequent production applications. It is the responsibility of the end user to qualify these products for their specific application. These products are not qualified or intended for any animal or human therapeutic use.

2. Intended use

Digest-EUR ® is a mucolytic agent, for the preparation of biological samples before their culture or diagnostic techniques (molecular biology).

Digest-EUR ® is a 10 X concentrated solution composed of dithiothreitol for rapid digestion and fluidification of biological samples. The very gentle action of Digest-EUR ® allows the liquefaction of samples without damage to the microbial flora or cells present.

This reagent has been designed for the recommended techniques in this technical sheet. Its use for other purposes must be validated by the user.

The Digest-EUR ® is a medical device for in vitro diagnosis use.

The Digest-EUR ® must be used by qualified medical laboratory personnel.

The Digest-EUR ® is suitable for use in In Vitro diagnostic procedures.

Instructions for use Digest-EUR ®

3. Symbols

	Reference
	Batch number
	Temperature limit
	Expiration date
	Keep out of direct sunlight
	Manufacturer
	Date of manufacture
	CE marked product
	In vitro diagnostic medical device
	Instructions for use
	Attention
	Do not use if packaging is damaged
	Pictogram of danger : corrosive

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Once ready, Digest-EUR ® shows a physiological solution which allows the sample liquefaction without damaging the microbial flora or present cells.

5. Principle of the test

- If used as described in this technical sheet, the reagent causes the oxidation (in disulfide) of thiol groups present in the sputum.
- The reagent is 10 times concentrated and needs to be diluted.:

6. Composition

The reagent contains dithiothreitol, which might be irritating for the skin, the eyes and the respiratory system .

7. Storage/stability

- The product must be stored at +15/+30°C before opening and is stable until the date printed on the vial's label.
- After opening, the reconstituted reagent, ready-to-use is stable for 48 hours at +2/+8°C.

8. Material required but not provided

- Pipette with sterile tip or disposable sterile pipette
- Sterile container for the dilution
- Centrifuge
- Conical tube for centrifuge
- Vortex mixer
- Sterile distilled water
- Sterile physiological water
- Benzalkonium 1 p. 3000 product n°CEZBZK00 if needed.

Instructions for use Digest-EUR ®

9. Cautions and recommendation for users

-  Do not use the product if the individual packaging is damaged.
- Upon receipt of the device, it is necessary to ensure the integrity of the packaging and containers that guarantee the quality of the reagents. Preliminary examination of the reagents must not reveal any turbidity, precipitate, particles or leaks.
- This reagent is intended for professional use in analytical laboratories by qualified personnel with expertise in microbiology.
- This reagent is exclusively designed for an in vitro diagnostic use.
- This product must be manipulated with all the precautions in compliance with the local regulation.
- The reagent contains dithiothreitol, which might be irritating for the skin, the eyes and the respiratory system .
- This reagent should be disposed of in accordance with local regulations.

10. Sample collection

- The sample (sputum, gastric fluid, sperma, joint aspiration, urine with mucus) is collected in a sterile container. It must be tested as soon as possible after collection.
- Products with animal origin must be considered as potentially infectious and must be manipulated with all the usual precautions.

11. Techniques

11.1 Solution reconstitution

Dilute the Digest-EUR ® 10 times concentrated with sterile distilled water, that is to say one vial of 10 ml for 90 ml of sterile distilled water. The obtained solution is ready-to-use. Prepare one bottle in once.

This ready-to-use solution is stable for 48 hours. It is possible to use a syringe to collect only fractions of 1 ml through the percutable cap.

11.2 Standard procedure

1. In a sterile conical tube for centrifuge, put 10 ml of sputum, for example, and the same quantity of diluted Digest-EUR®.
2. Homogenize with a vortex, incubate for 15 minutes at +15/+30°C.
3. Centrifuge 5 minutes at 2000 rpm and discard the supernatant

→ **For the flora and cells study:** dilute the centrifugation sediment in a sterile diluent (1 volume of centrifugation sediment for 10 volumes of diluent)

Instructions for use Digest-EUR ®

11.3 **Procedure for the culture of AFBS (Avid Fast Bacili)**

An additional decontamination with NaOH 1%, or even better with a sterile buffered solution of benzalkonium at 1 per 3000 is required.

Use the standard procedure and follow by decontamination:

1. Dilute the centrifugation sediment for 30 minutes at +15/+30°C with 20 ml of NaOH 1% or even better with 20 ml of sterile buffered solution of Benzalkonium, pH 7-7,2 1p.3000.
2. Centrifuge and discard the supernatant.
3. Wash the centrifugation sediment with distilled water. The sample is ready for culture.

Note: The use of Digest-EUR ® is particularly adapted for the membrane filtration of sputum, for example STERIFIL-MILLIPORE or SM SARTORIUS systems.

12. Waste disposal

Eliminate all waste in accordance with the local legislation.

13. Incident Report

Any serious incident occurring in connection with the device shall be notified to EUROBIO SCIENTIFIC and to the competent authority of the Member State in which the user and/or patient is established.

14. Technical assistance

For assistance with our products, please contact our technical support.

Eurobio Scientific customer service can be reached by e-mail at adv@eurobio-scientific.com or by phone at +33 (0)1.69.79.64.80.

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